

Catalog 2007/2008



Designed to work perfect

At-a-glance Table of contents

Mixing

Magnetic stirrers Overhead stirrers Shakers Laboratory kneader

10 – 29	
30 – 43	
44 – 58	
59	

Crushing

Dispersers Mills 62 – 85 86 – 91

Heating / Tempering

Heating baths Heating plates Termostats

Liquid Handling

Fixed volume pipettes 102 - 103

Distilling

Rotary evaporators

106 – 111

Mechanical accessories

Stands Fixing elements 114 – 115 116

Electronic

accessories

Temperature measuring	
instrument	117 – 118
Revolution counter	119
Vacuum controller	120
Vacuum pumps / valves	121

Laboratory reactors

Systems up to 2 liters 124 – 137

Rheology

instruments

138 – 139

Software

Laboratory software for control and data collection 142 – 151

Analytical line

Calorimeters	154 – 163
Decomposition system	164 – 165

General

IKA [®] world-wide	4 – 5
News	6 – 7
Questionaire	168
Information	169
Conditions of sale	
and delivery	170 – 171
HANDS for Children	172 – 173
Process Technology	174 – 175
Index	176 – 177
Notes	178

3

IKA[®] 4 Success through Creativity



IKA®-Werke Staufen, Germany

Success through Creativity

The IKA® group is a manufacturer of laboratory and analytical technology as well as mechanical engineering with world-wide locations. Over decades, we have built a close relationship to our international customers. This is based on our commitment to the highest quality and innovation of our products and our dedication to the development of solutions for complex applications. With an extraordinary research and development program, the company was able to achieve its leading world market position. IKA® has received many design and innovation awards for its innovative, leading edge products.

IKA® world-wide: www.ika.net

North America	Europe Middle East Africa
IKA® Works, Inc. 2635 North Chase Pkwy SE Wilmington NC 28405-7419 USA Phone: 800 733-3037 Phone: +1 910 452-7059 Fax: +1 910 452-7693 e-mail: usa@ika.net	IKA®-Werke GmbH & Co.KG Janke & Kunkel-Str. 10 D-79219 Staufen Phone: +49 7633 831-0 Fax: +49 7633 831-98 e-mail: sales@ika.de
Brasil	India
IKA® Works Inc. Av. das Américas, 15700, sala 235 Recreio dos Bandeirantes CEP 22790-701 Rio de Janeiro, RJ Brasil Phone: +55 21 2487-7743	IKA®-Werke GmbH & Co.KG Liaison Office India No. 31 (Old No. 264) 1st Floor, 10th Cross 1st "N" Block, Rajajinagar 560 010 Bangalore Phone: +91 80 41157736

Fax: +55 21 2487-7743 e-mail: fcabral@ika.net Fax: +91 80 41157735

e-mail: info@ikaindia.com

IKA[®] world-wide

The IKA® group



e-mail: sales@ika.com.my

IKA[®] News 2007

6 All news at-a-glance



ULTRA-TURRAX®

RCT basic safety control p. 10 – 12 Tube Drive

Magnetic stirrer with digital display

The new improved version of the laboratory bestseller with revolutionary technology and new functions:

new: Integrated temperature control new: Incl. PT 1000 temperature sensor (PT 1000.60)

p. 63

- World first: Universal disposable disperser system - Disperse, stir and grind using a single drive unit
- High level of user safety
- No possibility of cross-contamination - Hermetically sealable disposable sample tubes



MS 3 basic MS 3 digital p. 45

Small shakers

Compact, universal small shaker suitable for shaking tasks with all small vessels and microtiter plates.

- Attachment detection
- Continuous or touch operation

MS 3 digital additionally:

- Timer with countdown function



Electronic contact thermometers With optimised PID control and RESET function.

ETS-D6 additionally:

- With integrated pH measuring instrument (without pH electrode) and with multilingual user quide



HCT basic safety control

Heating plate with digital display With revolutionary technology and new functions: new: Integrated temperature control new: Incl. PT 1000 temperature sensor



(PT 1000.60)

p. 96

T 25 digital High-performance dispersing instrument Standard version with digital display

- For volumes from 1 - 2.000 ml (H₂O) - Three types of shaft bearings

IKA® RW 20

IKA[®] News 2007

All news at-a-glance



p. 66

Calorimeter system Recommended for learning and teaching purposes as well as for laboratories with low sampling volumes



p. 38

RW 20 digital



- COLOR SQUID
- BIG SQUID
- lab disc



lab dancer

Economic, compact test tube shaker with touch function. Its compact and clever design makes it an indispensable tool for every laboratory. Page 45







NEW!

RCT basic safety control IKAMAG® The laboratory bestseller: New with digital display for speed and temperature

ETS-D6

Electronic contact thermometer with PID control and RESET function, incl stainless steel temperature sensor H 62.51, page 24 Ident. No. 3378100

H 44 Boss head clamp, page 27 Ident. No. 2437700

H 38

Holding rod for fastening ETS-D5 or ETS-D6 with H 44 to the support rod H 16 V, page 27 Ident. No. 3547700

H 16 V

Support rod for all magnetic stirrers with M 10 threaded bushing, page 27 ldent. No. 1545100 _

RCT basic *safety control* IKAMAG[®] Magnetic stirrer with digital display, incl. protection cover H 100, page 12 Ident. No. 3380000 _



IKA[®] Mixing Magnetic stirrers with heating

12 Magnetic stirrers with heating



Ident. No. 3380000 230 V 50/60 Hz 3380001 115 V 50/60 Hz





included with unit Ident No. 3516800

RCT basic safety control IKAMAG®

The new improved version of the laboratory bestseller with revolutionary technology and new functions:

- new: Integrated temperature control new: Incl. PT 1000 temperature sensor (PT 1000.60)
- new: Exact temperature and speed setting via digital display, even when switched off
- new: Set safety temperature limit displayed digitally
- **new:** Hot Top indicator >> hot surface warning to prevent burns!
- new: Digital error code display
- With adjustable safety circuit of heating plate temperature (50 - 360 °C)
- Safety magnetic stirrer with heating, suitable for unsupervised operation
- Bushing according to DIN 12878 for connecting a contact thermometer, e.g. ETS-D5, enables precise temperature control
- High level of safety thanks to improved heat control technology
- Enclosed assembly (IP 42) guarantees long service life
- Highly polished aluminium heating plate for optimum heat transfer
- High magnetic adhesion
- Incl. protection cover H 100

Accessories (page):

Electronic contact thermometers (24): ETS-D5, ETS-D6, IKAFLON®- Stirring bars (29), TRIKA®-Stirring bars (29), RS 1 Set of stirring bars (29), RSE Stirring bar remover (29), Bath attachments (28): H 15, H 28, Oil bath attachments (28): H 29, H 30

Technical data	
Stirring quantity (H ₂ O)	20
Motor rating input	9,5 W
Motor rating output	3,5 W
Speed display	digital
Speed range	50 – 1.200 rpm
Max. magnetic bar (L x Ø)	80 x 10 mm
Heating function	
Heat output	600 W
Heating rate (1 H ₂ O)	6,5 K/min
Temperature range	RT – 310 °C
Setting accuracy	± 1 K
Temperature undulation without temper	rature sensor ± 2 K
Adjustable safety circuit	50 – 360 °C
Digital temperature	
limit display	50 – 360 °C
Control accuracy	PT 1000 / ± 1 K
with sensor	ETS-D5 / ± 0,5 K
	ETS-D6 / ± 0,2 K
Heating plate	
Material	aluminium alloy
Dimensions	Ø 135 mm
General data	
Dimensions (W x D x H)	160 x 270 x 85 mm
Weight	2,5 kg
Permissible ambient temperature	5 – 40 °C
Permissible relative humidity	80 %
Protection class acc. to DIN EN 60529	IP 42

Technical data	
Stirring quantity (H ₂ O)	20
Motor rating input	12 W
Motor rating output	5 W
Speed display	scale
Speed range	0 – 1.500 rpm
Max. magnetic bar (L x Ø)	80 x 10 mm
Heating function	
Heat output	600 W
Heating rate (1 H ₂ O)	7 K/mir
Temperature range	RT – 340 °C
Setting accuracy	± 10 ł
Adjustable safety circuit	50 – 380 °C
Sensor for temperature in medium	ETS-D5, ETS-D6
Control accuracy with sensor	ETS-D5 ± 0,5 ł
	ETS-D6 ± 0,2 H
Heating plate	
Material RET basic s. c.	stainless stee
RET basic C s. c.	stainless steel white coated
Dimensions	Ø 135 mm
General data	
Dimensions (W x D x H)	160 x 280 x 97 mn
Weight	2,4 k
Permissible ambient temperature	5 – 40 °C
Permissible relative humidity	80 %
Protection class acc. to DIN EN 6052	29 IP 42

RET basic safety control IKAMAG® RET basic C safety control IKAMAG®

Magnetic stirrer with new heating and stirring technology for fast, accurate temperature control. The RET basic (with stainless steel surface) and the RET basic C (with white coated surface, chemical resistant) are exemplified by: - Extremely fast heating times - Very broad temperature range (RT - 340 °C) - Adjustable safety circuit (50 - 380 °C) - Electronic speed control - Speed range from 0 - 1.500 rpm - Bushing according to DIN 12878 for the use of contact thermometers

- - (e.g. ETS-D5 or ETS-D6)
- Incl. protection cover H 99

```
Accessories (page):
Oil bath attachments (28): H 29, H 30
```

IKA[®] Mixing Magnetic stirrers with heating

Electronic contact thermometers (24): ETS-D5, ETS-D6, IKAFLON®- Stirring bars (29), TRIKA®-Stirring bars (29), RS 1 Set of stirring bars (29), RSE Stirring bar remover (29), R 380 Stand support (27), Bath attachments (28): H 15, H 28,



RET basic safe	ty control
ldent. No.	
3188800	230 V 50/60 Hz
3188801	115 V 50/60 Hz



RET basic C safety control Ident. No. 3197600 230 V 50/60 Hz 3197601 115 V 50/60 Hz





RET control-visc safety control		
ldent. No.		
3364000	230 V 50/60 Hz	
3364001	115 V 50/60 Hz	



RET control-visc C safety control		
ldent. No.		
3364100	230 V 50/60 Hz	
3364101	115 V 50/60 Hz	

RET control-visc safety control IKAMAG® RET control-visc C safety control IKAMAG®

Safety magnetic stirrer with heating, suitable for unsupervised operation.

- Option: 1 sensor for medium temperature (PT 100) or 2 separate temperature sensors for heat transfer fluid and medium (PT 1000) available (automatic identification)
- 2 adjustable safety circuits
- Stirring bar crack detection - With stainless steel surface or white coated surface (chemical resistant)
- Setting acc. medium temperature: 0,5 K - HOT warning display indicating presence of any
- residual heat when the unit is switched off - Easy-to-read backlit LCD display
- Actual medium temperature resolution displayed: 0,5 K (RT to 100 °C): 1 K (from 100 °C upwards)
- Fuzzy control and microprocessor technology guarantee maximum control accuracy
- PC control via RS 232 interface, with optional safety function
- Software labworldsoft® is available to control and document all measuring values via PC
- HOT visual warning for hot heating plate - 3 modes of operation, e.g. stirring and heating functions can be secured against inadvertent changes of set parameters
- Viscosity trend display
- Enclosed assembly (IP 42) guarantees long service life
- Incl. protection cover H 99

Accessories (page):

Temperature sensors (26): PT 100.50, PT 100.51, PT 100.52, PT 1000.50, PT 1000.51, IKAFLON®-Stirring bars (29), TRIKA®- Stirring bars (29), RS 1 Set of stirring bars (29), RSE Stirring bar remover (29), R 380 Stand support (27), Bath attachments (28): H 15, H 28, labworldsoft® (143), Oil bath attachments (28): H 29, H 30, AM 1 Analog module (120)

Stirring quantity (H ₂ O)	20
Motor rating input	12 W
Motor rating output	5 W
Speed display	digital
Speed range	0 – 1.500 rpm
Max. magnetic bar (L x Ø)	80 x 10 mm
Heating function	
Heat output	600 W
Heating rate (1 H ₂ O)	7 K/min
Temperature range	RT – 340 °C
Setting accuracy	0,5 K (bis 100 °C)
	1 K (ab 100 °C)
Adjustable safety circuit	50 – 350 °C
Digital temperature limit display	50 – 350 °C
Sensor for temperature	1 x PT 100
in medium	or 2 x PT 1000
Control accuracy with sensor	± 0,2 K
Heating plate	
Material RET control-visc s. c.	stainless steel
RET control-visc C s. c.	stainless steel
	white coated
Dimensions	Ø 135 mm
General data	
Dimensions (W \times D \times H)	160 x 280 x 97 mm
Weight	2,8 kg
Permissible ambient temperature	5 – 40 °C
Permissible relative humidity	80 %
Protection class acc. to DIN EN 60529	IP 42
Interface	RS 232 / analog

Technical data	
Stirring quantity (H ₂ O)	10
Motor rating input	15 V
Motor rating output	2 V
Speed display	scale (0 - 6
Speed range	100 – 2.000 rpn
Max. magnetic bar (L x Ø)	40 x 8 mn
Heating function	
Heat output	400 V
Heating rate (1 H ₂ O in H15)	3 K/mi
Temperature range	RT – 320 °C
Heating plate	
Material	stainless steel (AISI 304
Dimensions	Ø = 125 mn
General data	
Dimensions (W x D x H)	168 x 220 x 105 mn
Weight	2,4 k
Permissible ambient temperature	5 – 40 °(
Permissible relative humidity	80 %
Protection class acc. to DIN EN 60529	IP 2

Technical data		
Stirring quantity (H ₂ O)		15
Motor rating input		15 W
Motor rating output		2 W
Speed display		scale (0 – 6) / digital
Speed range		100 – 2000 rpm
Max. magnetic bar (L x Ø)		40 x 8 mm
Heating function		
Heat output		500 W
Heating rate (1 H ₂ O)		4,5 K/min
Temperature range		RT – 320 °C
Adjustable safety circuit		100 – 400 °C
Sensor for temperature in medium		ETS-D5, ETS-D6
Control accuracy with sensor	ETS-D5	± 0,5 K
	ETS-D6	± 0,2 K
Heating plate		
Heating plate, stainless steel		white coated
Dimensions		130 x 130 mm
General data		
Dimensions (W x D x H)	1	168 x 220 x 105 mm
Weight		2,4 kg
Permissible ambient temperature		5 – 40 °C
Permissible relative humidity		80 %

Protection class acc. to DIN EN 60529

RH basic 2 IKAMAG[®]

Economic magnetic stirrer with stainless steel heating plate.

- Fixed safety circuit 400 °C

- Soft-start stirring motor

Accessories (page): IKAFLON®- Stirring bars (29), TRIKA®- Stirring bars (29), RSE Stirring bar remover (29), Bath attachments (28): H 15, H 28

Universal magnetic stirrers with heating and bushing according to DIN 12878 for connecting an electronic temperature controller, e.g. ETS-D5, ETS-D6. RH digital KT/C safety control complete with digital display for set and actual temperature and actual speed.

- resistance
- Heat output 500 W
- Adjustable safety circuit for heating
- plate temperature
- Soft-start stirring motor

IP 21

Accessories (page):

Electronic contact thermometers: ETS-D5, ETS-D6 (24), IKAFLON®- Stirring bars (29), TRIKA®-Stirring bars (29), RSE Stirring bar remover (29), Bath attachments (28): H 15, H 28

IKA[®] Mixing Magnetic stirrers with heating



Ident. No.	
3339000	230 V 50/60 Hz
3339001	115 V 50/60 Hz

RH basic KT/C safety control IKAMAG® RH digital KT/C safety control IKAMAG®

- Heating plate with excellent chemical

- Long life cycle due to foil heating and solid-state switching for heat control

- Safety feature: in the event of motor failure the heating switches off automatically



RH basic KT/C safety control Ident. No. 230 V 50/60 Hz 3207100 115 V 50/60 Hz 3207101



RH digital KT/C safety control

Ident. No.	
3207000	230 V 50/60 Hz
3207001	115 V 50/60 Hz

16 Multi-position magnetic stirrers with heating



ldent. No.	
2930300	230 V 50/60 Hz
2930301	115 V 50/60 Hz

RT 5 power IKAMAG®

The RT 5 power is a high-performance multiposition magnetic stirrer with 5 stirring positions and integrated temperature control plate. Precise temperature distribution on the heating plate makes it possible to perform series experiments. - Simultaneously operating stirrers

- Absolute consistency over sample conditions for the individual samples

Accessories (page):

IKAFLON®- Stirring bars (29), TRIKA®- Stirring bars (29), RSE Stirring bar remover (29)

Technical data	
Stirring positions	5
Max. stirring quantity per stirrer (H ₂ O)	0,4
Distance between stirring places	90 mm
Motor rating input	7,2 W
Motor rating output	1,8 W
Speed display	scale (1 - 10)
Speed range	0 – 1.100 rpm
Deviation for individual stirring positions	5 %
Max. magnetic bar (L x Ø)	30 x 8 mm
Heating function	
Heat output	175 W
Temperature range (surface)	RT – 120 °C
Max. temperature medium (dep. on vessel)	70 °C
Heat control	scale (1 - 10)
Temperature consistancy in the medium	± 2 K
Heating plate	
Material	silicone
Dimensions	120 x 450 mm
General data	
Dimensions (W x D x H)	138 x 552 x 65 mm
Weight	3 kg
Permissible ambient temperature	5 – 40 °C
Permissible relative humidity	80 %
Protection class acc. to DIN EN 60529	IP 42

257	

Ident. No.	
2930500	230 V 50/60 Hz
2930501	115 V 50/60 Hz



Ident. No.	
2930700	230 V 50/60 Hz
2930701	115 V 50/60 Hz

RT 10 power	IKAMAG®
-------------	---------

Same features as RT 5 power, but with 10 stirring positions.

Accessories (page):

IKAFLON®- Stirring bars (29), TRIKA®- Stirring bars (29), RSE Stirring bar remover (29)

RT 15 power IKAMAG®

Same features as RT 5 power, but with 15 stirring positions.

Accessories (page):

IKAFLON®- Stirring bars (29), TRIKA®- Stirring bars (29), RSE Stirring bar remover (29)

echnical data	
Stirring positions	10
Notor rating input	14,4 W
Notor rating output	3,6 W
leating function	
leat output	375 W
leating plate	
Dimensions	180 x 450 mm
General data	
Dimensions (W \times D \times H)	198 x 552 x 65 mm
Veight	4,2 kg

Technical data	
Stirring positions	15
Motor rating input	21,6 W
Motor rating output	5,4 W
Heating function	
Heat output	580 W
Heating plate	
Dimensions	270 x 450 mm
General data	
Dimensions (W \times D \times H)	288 x 552 x 65 mm
Weight	6 kg

5 0,4 l 90 mm 7,2 W 1,8 W scale (1 – 10) 0 – 1.100 rpm 5 %
5 0,4 l 90 mm 7,2 W 1,8 W scale (1 – 10) 0 – 1.100 rpm 5 %
0,4 l 90 mm 7,2 W 1,8 W scale (1 – 10) 0 – 1.100 rpm 5 %
90 mm 7,2 W 1,8 W scale (1 – 10) 0 – 1.100 rpm 5 %
7,2 W 1,8 W scale (1 – 10) 0 – 1.100 rpm 5 %
1,8 W scale (1 – 10) 0 – 1.100 rpm 5 %
scale (1 – 10) 0 – 1.100 rpm 5 %
0 – 1.100 rpm 5 %
5 %
30 x 8 mm
ss steel (AISI 304)
120 x 450 mm
122 x 552 x 65 mm
2,3 kg
5 – 40 °C
80 %
00 /0

Technical data	
Stirring positions	10
Motor rating input	14,4 W
Motor rating output	3,6 W
Set-up plate	
Dimensions	180 x 450 mm
General data	
Dimensions (W x D x H)	182 x 552 x 65 mm
Weight	3,2 kg

of liquids.

- Optimum use of laboratory space - Including removable PUR cover

Accessories (page): IKAFLON®- Stirring bars (29), TRIKA®- Stirring bars (29), RSE Stirring bar remover (29)

R0 10 power IKAMAG®

Same features as RO 5 power, but with 10 stirring positions.

Accessories (page): IKAFLON®- Stirring bars (29), TRIKA®- Stirring bars (29), RSE Stirring bar remover (29)

Technical data	
Stirring positions	15
Motor rating input	21,6 W
Motor rating output	5,4 W
Set-up plate	
Dimensions	270 x 450 mm
General data	
Dimensions (W \times D \times H)	272 x 552 x 65 mm
Weight	4,7 kg

RO 15 power IKAMAG®

Same features as RO 5 power, but with 15 stirring positions.

Accessories (page): IKAFLON®- Stirring bars (29), TRIKA®- Stirring bars (29), RSE Stirring bar remover (29)

IKA[®] Mixing

Multi-position magnetic stirrers without heating

17

R0 5 power IKAMAG®

Multi-position magnetic stirrer with 5 stirring positions, without heating. The stainless steel surface covers the unit allowing easy cleaning and providing protection against the penetration



Ident. No.	
2930200	230 V 50/60 Hz
2930201	115 V 50/60 Hz



Ident. No.		
2930400	230 V	50/60 Hz
2930401	115 V	50/60 Hz



Ident. No. 2930600 230 V 50/60 Hz 2930601 115 V 50/60 Hz

IKA[®] Mixing Magnetic stirrers without heating



Mini MR standard IKAMAG®

Magnetic stirrer without heating, for stirring quantities up to 800 ml (H₂O). - Infinitely variable speed from 0 - 1.500 rpm

- White set-up plate suitable for observing color reactions

Accessories (page): IKAFLON®- Stirring bars (29), TRIKA®- Stirring bars (29), RSE Stirring bar remover (29)

ldent. No.	
0040000	

2812000 230 V 50/60 Hz 2812001 115 V 50/60 Hz



Ident. No. Motif 3277700 white 100 - 240 V 50/60 Hz

KMO 2 basic IKAMAG®

Small, powerful magnetic stirrer without heating. - Strong magnetic field

- Strong magnetic fieldMotor with optoelectronic speed control
- Infinitely variable speed from
- 0 1.100 rpm
- Stainless steel casing facilitates cleaning and sterilization
- Includes M 10 thread for H 16 V support rod

Accessories (page):

IKAFLON®- Stirring bars (29), TRIKA®- Stirring bars (29), RS 1 Set of stirring bars (29), RSE Stirring bar remover (29), H 16 V Support rod (27)

lab disc IKAMAG®

Ultra-flat compact magnetic stirrer, guaranteed with modern magnet coil technology. Wear-free drive with no moving parts. The lab disc can reverse direction of rotation automatically every 30 seconds to ensure better mixing.

- Explosion hazard zone 2 (see techn. data)
- Rotation direction is reversible
- High IP protection class (IP 65)
- Wear-free
- Set-up plate and casing made from chemically resistant materials
- Slip-proof, safe stand

Accessories (page):

IKAFLON®- Stirring bars (29), TRIKA®- Stirring bars (29), RSE Stirring bar remover (29)

Technical data	
Stirring quantity (H ₂ O)	0,8
Motor rating input	4 W
Motor rating output	0,5 W
Speed range	0 – 1.500 rpm
Max. magnetic bar (L x Ø)	30 x 8 mm
Set-up plate	
Material	polyester
Dimensions	115 x 115 mm
General data	
Dimensions (W x D x H)	115 x 130 x 40 mm
Weight	0,23 kg
Permissible ambient temperature	5 – 40 °C
Permissible relative humidity	80 %
Protection class acc. to DIN EN 60529	IP 42

Technical data	
Stirring quantity (H ₂ O)	5 I
Motor rating input	14 W
Motor rating output	4 W
Speed display	scale
Speed range	0 – 1.100 rpm
Max. magnetic bar (L x Ø)	50 x 8 mm
Set-up plate	
Material	stainless steel (AISI 304)
Dimensions	140 x 120 mm
General data	
Dimensions (W x D x H)	140 x 200 x 75 mm
Weight	1,4 kg
Permissible ambient temperature	5 – 40 °C
Permissible relative humidity	80 %
Protection class acc. to DIN EN 60529	IP 21

Technical data	
Stirring quantity (H ₂ O)	800 ml
Motor rating input	5 W
Motor rating output	3 W
Speed range	15 – 1.500 rpm
Reversion of rotation direction (switchab	every 30 s
Max. magnetic bar (L x Ø)	30 x 8 mm
Set-up plate	
Material	polyester
Dimensions	Ø 90 mm
General data	
Dimensions (W \times D \times H)	114 x 161 x 12 mm
Weight	0,3 kg
Permissible ambient temperature	5 – 40 °C
Permissible relative humidity	80 %
Protection class acc. to DIN EN 60529	IP 65
Ex-proof	EX II 3 G EE x nC IIB T6

lab disc IKAMAG[®] – the ultra-flat magnetic stirrer with new motifs

Height only

	ldent. No.	Motif		Ident. No.	Motif
1	3277700	white	4	3578600	brown
2	3362300	Sunny Side Up	5	3578200	green
3	3578000	orange	6	3578400	pink



3 orange

4 brown

IKA[®] Mixing Magnetic stirrers without heating





1 white





2 Sunny Side Up





COLOR SQUID IKAMAG[®] – the compact magnetic stirrer with attractive designs

Technical data on page 22.



	Ident. No.	Motif		Ident. No.	Motif
1	2669917	Number one	6	3576600	bighead
2	2669936	Mighty Mouse	7	3576900	balloon
3	2669922	Wildcat	8	3576000	curly
4	2669913	Zip	9	3576300	sunny
5	2669910	Harry's notes			



BIG SQUID IKAMAG® the magnetic stirrer with the extra large set-up plate

Technical data on page 22.

	Ident. No.	Motif		Ident. No.	Motif
1	3050009	White	6	3573000	pineapple
2	3050002	Froggy	7	3573200	strawberry
3	3050001	Star	8	3573600	orange
4	3050005	IKAmäleon	9	3573400	kiwi
5	3050004	Ocean			



2 Mighty Mouse



3 Wildcat



4 Zip



5 Harry's notes







8 curly



9 sunny



2 Froggy





6 pineapple

7 strawberry

Number one

IKA[®] Mixing Magnetic stirrers without heating 21







Ident. No. Motif 2669917 number one 100 – 240 V 50/60 Hz (29), RSE Stirring bar remover (29)



Ident. No. Motif 3050009 White 100 – 240 V 50/60 Hz

COLOR SQUID IKAMAG®

Small magnetic stirrers without heating, in a variety of motifs.

- Recyclable materials
- Very good chemical resistance due to glass top and synthetic bottom made of Hytrel®
- Electronically controlled motor with infinitely variable speed from 0 - 1.500 rpm

Accessories (page):

IKAFLON®- Stirring bars (29), TRIKA®- Stirring bars

Motifs on page 20.

Technical data	
Stirring quantity (H ₂ O)	0,8
Motor rating input	2 W
Motor rating output	1 W
Speed display	none
Speed range	0 – 1.500 rpm
Max. magnetic bar (L x Ø)	30 x 8 mm
Set-up plate	
Material	glass
Dimensions	Ø 110 mm
General data	
Dimensions (W \times D \times H)	130 x 135 x 50 mm
Weight	0,48 kg
Permissible ambient temperature	5 – 40 °C
Permissible relative humidity	80 %
Protection class acc. to DIN EN 60529	IP 54

Technical data	
Stirring quantity (H ₂ O)	1,5
Motor rating input	2 W
Motor rating output	1 W
Speed display	none
Speed range	0 – 1.500 rpm
Max. magnetic bar (L x Ø)	30 x 8 mm
Set-up plate	
Material	glass
Dimensions	Ø 160 mm
General data	
Dimensions (W \times D \times H)	180 x 180 x 48 mm
Weight	1 kg
Permissible ambient temperature	5 – 40 °C
Permissible relative humidity	80 %
Protection class acc. to DIN EN 60529	IP 54

Technical data	
Stirring quantity (H ₂ O)	20
Motor rating input / output	9,5 W / 3,5 W
Speed display	digital
Speed range	0 – 1.500 rpm
Max. magnetic bar (L x Ø)	80 x 10 mm
Set-up plate	
Material	stainless steel
Dimensions	Ø 135 mm
General data	
Dimensions (W x D x H)	160 x 270 x 85 mm
Weight	2,4 kg
Permissible ambient temperature	5 – 40 °C
Permissible relative humidity	80 %
Protection class acc. to DIN EN 60529	IP 42

Technical data

Stirring quantity (H₂O)

Motor rating input

Motor rating output

Max. magnetic bar (L x Ø)

Permissible ambient temperature

Permissible relative humidity Protection class acc. to DIN EN 60529

Speed display

Speed range

Set-up plate

Dimensions

General data Dimensions (W x D x H)

Technical data

Stirring quantity (H₂O)

Material

Weight

REO basic IKAMAG® REO basic C IKAMAG®

- Digital speed display

- Infinitely variable speed

- Incl. protection cover H 101

request.

Accessories (page): (29), RS 1 Set of stirring bars (29), RSE Stirring bar remover (29)

Midi MR 1 digital IKAMAG®

- Flat, sturdy stainless steel casing
- Non-locking motor digital 0 – 1.000 rpm

50 I

70 W

19 W

9 kg 5 – 40 °C

80 %

IP 21

150 I

80 x 10 mm

stainless steel (AISI 304)

360 x 430 x 110 mm

- Infinitely variable speed - Digital LED speed display
- Timer (0 56 min) or continuous operation
- For stirring quantities up to 50 liters (H₂O) 350 x 350 mm

Accessories (page):

Maxi MR 1 digital IKAMAG®

Same features as Midi MR 1 digital.

Accessories (page):

Motor rating input	80 W
Motor rating output	35 W
Speed display	digital
Speed range	0 – 600 rpm
Max. magnetic bar (L x Ø)	155 x 27 mm
Set-up plate	
Material	stainless steel (AISI 304)
Dimensions	500 x 500 mm
General data	
Dimensions (W \times D \times H)	505 x 585 x 110 mm
Weight	16 kg
Permissible ambient temperature	5 – 40 °C
Permissible relative humidity	80 %
Protection class acc. to DIN EN 60529	IP 21

Magnetic stirrer in five motifs. Same features as COLOR SQUID, but with a larger set-up plate:

- Diameter: 160 mm
- Stirring quantity (H2O) up to 1,5 l

Accessories (page):

IKAFLON®- Stirring bars (29), TRIKA®-Stirring bars (29), RSE Stirring bar remover (29)

Motifs on page 21.

IKA[®] Mixing Magnetic stirrers without heating 23

- Classic magnetic stirrer without heating.
- Non-locking, electronically controlled motor
- Constant speed even during changes in load
- REO basic C with white coated set-up plate on
- REO basic IKAMAG® available 3. quarter 2007.
- IKAFLON®- Stirring bars (29), TRIKA®- Stirring bars



REO basic I	KAMAG®	
ldent. No.		

3384200	230 V	50/60 Hz
3384201	115 V	50/60 Hz

- Powerful magnetic stirrer without heating.
- IKAFLON®- Stirring bars (29), TRIKA®- Stirring bars (29), RSE Stirring bar remover (29)



ldent. No.	
2621900	230 V 50/60 Hz
2621901	115 V 50/60 Hz

- For stirring quantities up to 150 I (H₂O)
- IKAFLON®- Stirring bars (29), TRIKA®- Stirring bars (29), RSE Stirring bar remover (29)



24 Magnetic stirrers accessories



Electronic Contact Thermometers ETS-D5 and ETS-D6

Ensures perfect temperature control without overshooting the set temperature, even in the case of quick heating. With optimised PID control and RESET function, incl. stainless steel sensor H 62.51. For all magnetic stirrers with contact thermometer bushing according to DIN 12878, class 2 (e.g. IKA®, Heidolph and Corning with adapter AD-C1, Ident. No. 3414000, please order separately).

ETS-D6 additionally:

- With integrated pH measuring instrument (without pH electrode)
- Large, graphic LCD display with multilingual user guide

Patented: 3 modes of operation guarantee optimum adjustment to your working method: Operating mode A

Suitable for work with varying parameters (from -50 °C to 450 °C). Safety temperature adjustable.

Operating mode B Suitable for series operation under uniform conditions.

Operating mode C Suitable for unsupervised operation.

All values are taken from the memory. This ensures perfect protection against inadvertent improper adjustment.

Accessories ETS-D5 and ETS-D6 (page): Sensor (25): H 62.51, H 66.51, H 70 Extension cable (25), H 52 Power pack set (25), H 16 V Support rod (27), H 44 Boss head clamp (27), H 38 Holding rod (27)

Temperature	
Temperature measuring rang	e -50 - 450 °C
Resolution	0,1 K
Measuring accuracy	± 0,2 K + Sensor tolerance PT 1000
	DIN IEC 751 class A
Setting accuracy	0,1 K
Control deviation	± 0,5 K
General data	
Supply voltage	8 – 16 VDC
Power consumption	10 mA (at 9 V)
Max. ON time	100 %
Plug	6 pin DIN 45322
Connection	DIN 12878 class 2
Dimensions (W x D x H)	82 x 22 x 83 mm
	(without sensor)
Weight	0,2 kg
Permissible ambient tempera	ature 0-60 °C
Permissible relative humidity	80 %
Protection class acc. to DIN	EN 60529 IP 54

Temperature	
Temperature measuring range	-50 – 450 °C
Resolution	0,01 K
Measuring accuracy ± 0,05 K +	Sensor tolerance PT 1000
	DIN IEC 751 class A
Setting accuracy	0,1 K
Control deviation	± 0,2 K
pH measurement	
Measuring range	0 – 14 pH
Accuracy	± 0,1 pH
Resolution	± 0,01 pH
pH connection	BNC-bushing
General data	
Supply voltage	8 – 16 VDC
Power consumption	15 mA (at 9 V)
Max. ON tim	100 %
Plug	6 pin DIN 45322
Connection	DIN 12878 class 2
Dimensions (W \times D \times H)	96 x 45 x 98 mm
	(without sensor)
Weight	0,2 kg
Permissible ambient temperature	0 - 60 °C
Permissible relative humidity	80 %
Protection class acc. to DIN EN 6052	29 IP 54

General data	
Depth of immersion	230 mm
Diameter	3 mm
Length	260 mm

		H 66.51 Staini
General data		glass-coated
Depth of immersion	230 mm	J
Diameter	6 mm	For use with ETS
Length	260 mm	with aggressive r
		solutions.

General data

General data

Analog output

General data

Diameter

Depth of immersion

Length

H 70 Extension cable

ETS-D6.

1 m

10 mV/K

To separate the casing from the sensor. The casing with the electronics may thus be kept away from dangerous vapor released by the medium (for use with ETS-D5 and ETS-D6).

H 52 Power pack set

The power pack set is required in order to operate older magnetic stirrer models (prior to 1990) with ETS-D5 and ETS-D6. If you have any questions, please contact our service department. In addition, the power pack set features an analog output to document signals on a recorder.

Accessories (page): labworldsoft® (143), DC 2 DATACONTROL (147), AK 2.1 Analog cable (148)

PC documentation is also possible in combination with DC 2 DATACONTROL and labworldsoft®.

PT 1000.60 Temperature sensor

Made of stainless steel, for use with 230 mm RCT basic safety control. 3 mm

ETS-D6 Ident. No. 3378100

NEW!

IKA[®] Mixing

Magnetic stirrers accessories

H 62.51 Stainless steel sensor

Spare sensor for use with ETS-D5 and

Ident. No. 2735451

H 66.51 Stainless steel sensor

For use with ETS-D5 and ETS-D6, for work with aggressive media such as acid and alkaline

> Ident. No. 2735551



2735600



Ident. No. 8010600 230 V 50/60 Hz 8010601 115 V 50/60 Hz



3516800

25

26 Magnetic stirrers accessories



Ident. No. 2601900



Ident, No. 2600300



Ident. No. 2847700



Ident. No. 3367600



Ident, No. 3377700

PT 100.50 Temperature sensor

For use with RET control-visc safety control and RET control-visc C safety control.

Accessories (page):

alkaline solutions.

H 16 V Support rod (27), H 44 Boss head clamp (27), H 38 Holding rod (27)

PT 100.51 Temperature sensor

General data	
Depth of immersion	230 mm
Diameter	8 mm

230 mm

230 mm

3 mm

3 mm

Accessories (page): H 16 V Support rod (27), H 44 Boss head clamp (27), H 38 Holding rod (27)

visc safety control and RET control-visc C safety

For use with RET control-visc safety control and

RET control-visc C safety control, glass-coated

for work with aggressive media such as acid and

General data	
Depth of immersion	60 mm
Diameter	3 mm



H 16 V Support rod

stirrers with M 10 threaded bushing.

Accessories (page):

H 16.1 Extension

R 380 Stand support

For fitting along the multifunction strips of the magnetic stirrers RET basic safety control, RET basic C safety control, RET control-visc safety control, RET control-visc C safety control and REO basic C. It allows the support rod H 16 V to be fixed at any given position. This makes it possible to use several support rods.

Accessories (page): H 16 V Support rod (27)

H 44 Boss head clamp

For fastening the holding rod H 38 (p. 27) to the support rod H 16 V (p. 27).

H 38 Holding rod

For fastening ETS-D5 or ETS-D6 with H 44 (p. 27) to the support rod H 16 V (p. 27).

General data	
Depth of immersion	230 mm
Diameter	3 mm

PT 100.52 Temperature sensor	
	General data
Made of stainless steel, for use with RET control-	Depth of immersion

General data

Diameter

Depth of immersion

General data

Diameter

Depth of immersion

control.			

Accessories (page): H 16 V Support rod (27), H 44 Boss head clamp (27), H 38 Holding rod (27)

PT 1000.50 Temperature sensor

2 separate steel sensors for heat transfer fluid and medium. Ideal for the magnetic stirrers RET control-visc safety control and RET controlvisc C safety control.

Accessories (page):

H 16 V Support rod (27), H 44 Boss head clamp (27), H 38 Holding rod (27)

PT 1000.51 Temperature sensor

As per PT 1000.50, but glass-coated for work with corrosive media such as acids and lyes.

Accessories (page): H 16 V Support rod (27), H 44 Boss head clamp (27), H 38 Holding rod (27)

IKA[®] Mixing Magnetic stirrers accessories

Stainless steel support rod for all magnetic

R 380 Stand support (27), H 16.1 Extension (27), H 44 Boss head clamp (27), H 38 Holding rod (27)

For work with bath attachment over 180 mm Ø.



Ident, No 5000500

Ident. No.

1545100



Ident. No. 2636700



Ident. No. 2437700



3547700

27

28 Magnetic stirrers accessories



H 99 Protective cover H 100 Protective cover H 101 Protective cover

Resistant to most acids, alkaline solutions and organic solvents. The protective cover is included with the magnetic stirrer.

Ident. No.
0551300

Ident. No. 2167400



H 28 Bath attachment

H 15 Bath attachment

tempering 0,5 and 1 l flasks.

Stainless steel bath attachment, suitable as a sand bath basin.

Stainless steel bath attachment, suitable for

H 29 Oil bath attachment H 30 Oil bath attachment

The oil bath attachments H 29 and H 30 can be used as oil baths together with an IKAMAG[®] magnetic stirrer with heating or with an IKATHERM[®] heating plate having a diameter of 135 mm.

- Positioning border prevents sliding on the heating plate
- Safety grips protect you from burns caused by hot oil
- The bath attachment is made of aluminum. This ensures good heat transfer and quick heating-up of the
- tempering medium
- Easy cleaning
- The bath attachments can only be used as an oil bath

General data	
Material	silicone
Max. temperature	135 °C
Protective cover	
H 99	for RET series and RCT basic s.c. (3378500)
H 100	for RCT basic s.c. (3380000)
H 101	for REO basic

General data	
Inner diameter	140 mm
Height	125 mm
Volume	1,5

General data	
Inner diameter	140 mm
Height	70 mm
Volume	11
Max. temperature	350 °C

General data		
Inner diameter	H 29	136 – 180 mm
	H 30	136 – 190 mm
Height	H 29	81 mm
	H 30	110 mm
Volume	H 29	11
	H 30	1,5

General data	
For heating plate diameter	135 mm
Material	aluminium
Diameter enlarged to	200 mm

H 12 / 135 Supporting plate

200 mm.

Ident. No.	Description	Length	Ø
1572000	IKAFLON® 10*	10 mm	6 mm
1572100	IKAFLON® 15*	15 mm	6 mm
1572200	IKAFLON® 20*	20 mm	8 mm
1572300	IKAFLON® 25*	25 mm	8 mm
1572400	IKAFLON® 30*	30 mm	8 mm
1572500	IKAFLON® 40*	40 mm	8 mm
1572600	IKAFLON® 50*	50 mm	8 mm
1572800	IKAFLON® 80*	80 mm	10 mm
0793300	IKAFLON [®] 110	110 mm	27 mm
1129000	IKAFLON® 155	155 mm	27 mm

Description	Length
TRIKA [®] 25*	25 mm
TRIKA [®] 40*	40 mm
	Description TRIKA® 25* TRIKA® 40*

Triangular, PTFE-coated, especially suited for stirring liquids which have a low solids content and where sedimentation is not desired.

RS 1 Set of magnetic stirring bars

Consisting of the IKAFLON® and TRIKA® Magnetic stirring bars marked with *, see above.

RSE Stirring bar remover

PTFE-coated.

		H 11
Ident. No.		
1091500	Euro plug	Spare
1091700	USA plug	
2410700	UK plug	
1091600	CH plug	

H 11 Mains cable

IKA[®] Mixing Magnetic stirrers accessories

29

For increasing the heating plate to a diameter of



Ident. No. 0771700



Round, PTFE-coated.

IKAFLON[®] Magnetic stirring bars

TRIKA[®] Magnetic stirring bars



Ident. No. 1358600

For all stirring bars up to 80 mm in length,

Ident. No. 1293100



30 IKA[®] Mixing Overhead stirrers



EUROSTAR power control-visc Powerful, digital laboratory stirrer

EUROSTAR power control-visc

Stirrer for quantities up to 40 l, with RS 232 interface, page 35 Ident. No. 2600000

R 271 Boss head clamp, page 116 Ident. No. 2664000

R 2723 Telescopic stand, page 115 Ident. No. 1412100 ____

R 1331 Anchor stirrer, page 40 Ident. No. 2022400

RH 5 Strap clamp, page 116 Ident. No. 3159000

With labworld*soft*® you can network up to 64 laboratory devices and control these from a PC, see page 143

IKA[®] Mixing Electronic overhead stirrers 31







RW 11 basic "Lab egg"

Small-sized stirrer available in four attractive colors.

- Glass housing resistant to chemicals
- Max. stirring quantity 2 I (H₂O)
- Incl. paddle stirrer R 1001 and extension arm

Accessories (page):

R 103 Stand (114), R 1001 Spare paddle stirrer (42), R 1002 Screw-type stirrer (42)

Technical data	
Stirring quantity (H ₂ O)	2
Max. viscosity	100 mPas
Motor rating input	8 W
Motor rating output	1 W
Output at stirring shaft	1 W
Max. ON-time	100 %
Max. torque (plug-in coupling)	0,8 Ncm
Speed range	0 – 2.000 rpm
Speed display	none
Plug-in coupling Ø	4 mm
Support holder Ø	integrated (10 mm)
General data	
Dimensions (W x D x H)	86 x 175 x 89 mm
Weight	0,39 kg
Permissible ambient temperature	5 – 40 °C
Permissible relative humidity	80 %
Protection class acc. to DIN EN 60529	IP 42

4 transparent



echnical data	
tirring quantity (H ₂ O)	81
lax. viscosity	10.000 mPas
lotor rating input	28,5 W
lotor rating output	17 W
utput at stirring shaft	17 W
1ax. ON-time	100 %
lax. torque at chuck	8 Ncm
peed range	100 – 2.000 rpm
peed display	scale
huck range	0,5 – 8 mm
iameter / length of extension arm	13 / 160 mm

Diameter / length of extension arm General data

Dimensions	
without extension arm (W \times D \times H)	70 x 176 x 197 mm
Weight	2 kg
Permissible ambient temperature	5 – 40 °C
Permissible relative humidity	80 %
Protection class acc. to DIN EN 60529	IP 20

RW 14 basic

viscosity changes.

- Constant speed thanks to electronic control - Infinitely adjustable speed - Very smooth thanks to direct drive

- Easy to operate

- Slimline

Accessories (page):

Stands (114): R 1825, R 1826, 1827, R 182 Boss head clamp (116), RH 3 Strap clamp (116), R 301 Stirring shaft protection (42), R 301.1 Support holder (42), Stirring elements (40 / 41): e.g. R 1342, FK 1 Flexible coupling (42)

EUROSTAR digital

Stirrer for quantities up to 20 l, page 34 Ident. No. 2482000 _

R 182 Boss head clamp, page 116 Ident. No. 2657700

RW 16 basic Stirrer for quantities up to 10 l, page 34 Ident. No. 2572100 _

R 1373 Paddle stirrer, page 40 Ident. No. 0757600

R 1330 Anchor stirrer, page 40 Ident. No. 2022300

RH 3 Strap clamp, page 116 Ident. No. 3008600

R 1825

Plate stand, page 114 Ident. No. 3160000

IKA[®] Mixing Electronic overhead stirrers 33

Quiet, economical laboratory stirrer with electronic infinitely adjustable speed. For stirring substances of low to medium viscosity. The laboratory stirrer is suitable for reproducibly setting the speed or processing media with substantial temporary

- Non-locking, overload capabilities





34 Electronic overhead stirrers



Ident, No.

2482000	230 V	50/60 Hz
2482001	115 V	50/60 Hz

RW 16 basic

Laboratory stirrer for simple stirring tasks of up to 10 liters (H₂O) with ideal speed range from 40 - 1.200 rpm. Especially suitable for schools, universities and inspection laboratories.

- Infinitely adjustable without gear shifting
- Slim casing
- Quiet operation - Safety circuit
- Non-locking, overload capabilities

Accessories (page):

Stands (114): R 1825, R 1826, R 1827, R 182 Boss head clamp (116), FK 1 Flexible coupling (42), RH 3 Strap clamp (116), DZM control. o Revolution counter (119), R 301 Stirring shaft protection (42), Stirring elements (40 / 41): e.g. R 1342, R 1330, R 1373

EUROSTAR digi	tal
---------------	-----

Laboratory stirrer that can be used up to the "medium viscosity" range.

- Constant speed by microprocessor control - Digital display presents set and actual speed
- Infinitely adjustable without gear shifting
- Slim casing
- Quiet operation
- Safety circuit
- Non-locking, overload capabilities
- Push-through agitator shafts
- Enhanced safety as a result of smooth start

Accessories (page):

Stands (114): R 1825, R 1826, R 1827, R 182 Boss head clamp (116), FK 1 Flexible coupling (42), RH 3 Strap clamp (116), R 301 Stirring shaft protection (42), Stirring elements (40 / 41): e.g. R 1342, R 1330, R 1373

Stirring quantity (H ₂ O)	10
Max. viscosity	10.000 mPas
Motor rating input	75 W
Motor rating output	55 W
Output at stirring shaft	53 W
Max. ON-time	100 %
Max. torque at chuck	40 Ncm
Speed range	40 – 1.200 rpm
Speed display	scale (1 - 10)
Chuck range	0,5 – 10 mm
Hollow shaft, inner diameter	11 mm
Diameter / length of extension arm	13 mm / 160 mm
General data	
Dimensions (W x D x H)	80 x 190 x 222 mm
Weight	2,8 kg
Permissible ambient temperature	5 – 40 °C
Permissible relative humidity	80 %
Protection class acc. to DIN EN 60529	IP 42

Protection class acc. to DIN EN 60529

Technical data

Technical data	
Stirring quantity (H ₂ O)	20
Max. viscosity	10.000 mPas
Motor rating input	75 W
Motor rating output	55 W
Output at stirring shaft	53 W
Max. ON-time	100 %
Max. torque at chuck	30 Ncm
Speed range	50 – 2.000 rpm
Speed display	digital
Chuck range	0,5 – 10 mm
Hollow shaft, inner diameter	11 mm
Diameter / length of extension arm	13 mm / 160 mm
General data	
Dimensions (W x D x H)	80 x 190 x 222 mm
Weight	2,8 kg
Permissible ambient temperature	5 – 40 °C
Permissible relative humidity	80 %
Protection class acc. to DIN EN 60529	IP 42

Technical data	
Stirring quantity (H ₂ O)	40 I
Max. viscosity	50.000 mPas
Motor rating input	130 W
Motor rating output	110 W
Output at stirring shaft	105 W
Max. ON-time	100 %
Max. torque at chuck	60 Ncm
Speed range	50 – 2.000 rpm
Speed display	scale
Chuck range	0,5 – 10 mm
Hollow shaft, inner diameter	11 mm
Diameter / length of extension arm	16 mm / 200 mm
General data	
Dimensions (W x D x H)	80 x 190 x 253 mm
Weight	3,8 kg
Permissible ambient temperature	5 – 40 °C
Permissible relative humidity	80 %
Protection class acc. to DIN EN 60529	IP 42

Interface

Technical data

Max. viscosity

Max. ON-time

Speed range

Speed display

Chuck range

General data

Weight

Interface

Motor rating input

Motor rating output

Output at stirring shaft

Max. torque at chuck

Hollow shaft, inner diameter

Dimensions (W x D x H)

Torque measurement

Diameter / length of extension arm

Permissible ambient temperature

Protection class acc. to DIN EN 60529

Permissible relative humidity

Stirring quantity (H₂O)

EUROSTAR power basic

"high viscosity" range.

- Slim casing - Quiet operation
- Safety circuit
- Non-locking, overload capabilities
- Push-through agitator shafts
- Analog recording of speed parameters
- is possible

Accessories (page): (40 / 41): e.g. R 1345, R 1375

EUROSTAR power control-visc

EUROSTAR power basic, additionally: and documented via a PC.

- Digital display presents rated / actual speed - Integrated torque trend display for viscosity
- Analog interface for recording speed and torque
- RS 232 interface

Accessories (page):

100 % 60 Ncm 50 – 2.000 rpm digital 0,5 – 10 mm control 11 mm 16 mm / 200 mm

analog

40 I

130 W

110 W

105 W

50.000 mPas

80 x 190 x 253 mm

3,8 kg

80 %

IP 42

trend

5-40 °C

RS 232 / analog

IKA[®] Mixing Electronic overhead stirrers

- Powerful laboratory stirrer for tasks up to the
- Constant speed by microprocessor control - Infinitely adjustable without gear shifting
- Enhanced safety as a result of smooth start

Stands (114 / 115): R 2722, R 2723, R 271 Boss head clamp (116), FK 1 Flexible coupling (42), RH 5 Strap clamp (116), VK 600 control (139), DZM control.o Revolution counter (119), R 301 Stirring shaft protection (42), Stirring elements



Powerful, digital laboratory stirrer for tasks up to the "high viscosity" range. Same features as labworldsoft[®] software is available to allow speed and torgue parameters to be controlled, regulated

Stands (114 / 115): R 2722, R 2723, R 271 Boss head clamp (116), FK 1 Flexible coupling (42), RH 5 Strap clamp (116), VK 600 control Torque measurement instrument (139), R 301 Stirring shaft protection (42), Stirring elements (40 / 41): e.g. R 1345, R 1375, labworldsoft® (143), RC 1 Remote control (43), AM 1 Analog module (120)



IKA[®] Mixing Electronic overhead stirrers







EUROSTAR power control-visc 6000

High-performance digital laboratory stirrer for tasks up to the "medium viscosity" range. Same features as EUROSTAR power control-visc (page 35), additionally:

- Speed range up to 6.000 rpm
- Agitator elements are not push-through
- Cone seat for precision shaft, incl. with delivery (stirring elements can be screw connected, please order separately, see page 42)
- Analog output of speed and torque

Accessories (page):

Stands (114 / 115): R 2722, R 2723, R 271 Boss head clamp (116), RH 5 Strap clamp (116), R 301 Stirring shaft protection (42), R 1402 Dissolver (42), R 1405 Propeller (42), R 1401 Propeller (42), labworldsoft[®] (143), RC 1 Remote control (43), AM 1 Analog module (120)

Technical data	
Stirring quantity (H ₂ O)	20
Max. viscosity	10.000 mPas
Motor rating input	130 W
Motor rating output	110 W
Output at stirring shaft	95 W
Max. ON-time	100 %
Max. torque at chuck	15 Ncm
Speed range	150 – 6.000 rpm
Speed display	digital
Diameter / length of extension arm	16 mm / 220 mm
General data	
Dimensions (W x D x H)	80 x 190 x 317 mm
Weight	4,8 kg
Permissible ambient temperature	5 – 40 °C
Permissible relative humidity	80 %
Protection class acc. to DIN EN 60529	IP 42
Interface	RS 232 / analog
Torque measurement	trend

EUROSTAR power control-visc P1

Powerful, digital laboratory stirrer for tasks up to the "high viscosity" range.

- Constant speed by microprocessor control
- Infinitely adjustable without gear shifting
- Slim casing
- Quiet operation
- Safety circuit
- Non-locking, overload capabilities
- Push-through agitator shafts
- Enhanced safety as a result of smooth start
- Digital display presents rated- and actual-speed
- Integrated torque trend display for viscosity control
- Analog interface for recording speed and torque
- RS 232 interface
- Software labworldsoft® is available to control and document all measuring values via PC

Accessories (page):

Stands (114 / 115): R 2722, R 2723, R 271 Boss head clamp (116), RH 5 Strap clamp (116), VK 600 control Torque measurement instrument (139), R 301 Stirring shaft protection (42), Stirring elements (40 / 41): e.g. R 1331, R 1312, labworldsoft® (143), RC 1 Remote control (43), AM 1 Analog module (120)

60 I
70.000 mPas
153 W
134 W
126 W
100 %
100 Ncm
50 – 1.200 rpm
digital
0,5 – 10 mm
11 mm
16 mm / 200 mm
80 x 190 x 253 mm
4 kg
5 – 40 °C
80 %
IP 42
RS 232 / analog
trend

EUROSTAR power control-visc P4

transmission reduction	4-1010
Stirring quantity (H ₂ O)	40
Max. viscosity	100.000 mPas
Motor rating input	153 W
Motor rating output	134 W
Output at stirring shaft	126 W
Max. ON-time	100 %
Max. torque at chuck	200 Ncm
Speed range	14 – 530 rpm
Speed display	digital
Chuck range	0,5 – 10 mm
Hollow shaft, inner diameter	11 mm
Diameter / length of extension arm	16 mm / 200 mm
General data	
Dimensions (W x D x H)	80 x 190 x 330 mm
Weight	4,9 kg
Permissible ambient temperature	5 – 40 °C
Permissible relative humidity	80 %
Protection class acc. to DIN EN 60529	IP 42
Interface	RS 232 / analog
Torque measurement	trend
EUROSTAR power control-visc P7	

Transmission reduction 7-fold Stirring quantity (H₂O) 40 I Max. viscosity 150.000 mPas Motor rating input 153 W Motor rating output 134 W Output at stirring shaft 126 W Max. ON-time 100 % Max. torgue at chuck 380 Ncm Speed range 8 – 290 rpm Speed display digital Chuck range 0,5 – 10 mm Hollow shaft, inner diameter 11 mm Diameter / length of extension arm 16 mm / 200 mm General data Dimensions (W x D x H) 80 x 190 x 330 mm Weight 4,9 kg Permissible ambient temperature 5 – 40 °C Permissible relative humidity 80 % Protection class acc. to DIN EN 60529 IP 42 Interface RS 232 / analog Torque measurement trend

EUROSTAR power control-visc P4 EUROSTAR power control-visc P7

Powerful, digital laboratory stirrer for tasks up to the "high viscosity" range.

- Slim casing
- Quiet operation
- Safety circuit
- Non-locking, overload capabilities
- Push-through agitator shafts
- Integrated torque trend display for viscosity
- control
- RS 232 interface

P4 with 4-fold transmission reduction and P7 with 7-fold transmission reduction; without pushthrough agitator shafts.

Accessories (page): Stands (114 / 115): R 2722, R 2723, R 271 Boss head clamp (116), RH 5 Strap clamp (116), VK 600 control Torque measurement instrument (139), R 301 Stirring shaft protection (42), Stirring elements (40 / 41): e.g. R 1331, R 1312, labworldsoft® (143), RC 1 Remote control (43), AM 1 Analog module (120)

IKA[®] Mixing Electronic overhead stirrers

- Constant speed by microprocessor control - Infinitely adjustable without gear shifting

- Enhanced safety as a result of smooth start

- Digital display presents rated- and actual-speed

- Analog interface for recording speed and torque

- Software labworldsoft® is available to control and document all measuring values via PC



2850700	230 V 50/60 Hz
2850701	115 V 50/60 Hz

Mechanical overhead stirrers 38





Ident, No. 2760000 230 V 50/60 Hz 2760001 115 V 50/60 Hz

RW 20 digital

Overhead stirrer with digital display. Technical improvements on the trusted RW 20 series designs. new: With digital display new: Robust, slimline, ergonomic design

With constant power-drive

- Two speed ranges for universal use from 60 - 2.000 rpm
- Push-through agitator shafts (only when stationary)

Accessories (page):

Stands (114): R 1825, R 1826, R 1827, FK 1 Flexible coupling (42), RH 3 Strap clamp (116), VK 600 control Torque measurement instrument (139), R 301 Stirring shaft protection (42), Stirring elements (40 / 41): e.g. R 1342, R 1381, VK 60/01 Adapter (139)

Fechnical data	
Stirring quantity (H ₂ O)	20
Max. viscosity	10.000 mPas
Notor rating input	70 W
Notor rating output	35 W
Output at stirring shaft	26 W
Max. ON-time	100 %
Max. torque at chuck	150 Ncm
Speed range I (per 50 Hz)	60 – 500 rpm
Speed range II (per 50 Hz)	240 – 2.000 rpm
Speed display	digital
Chuck range	0,5 – 10 mm
Diameter / length of extension arm	13 mm / 160 mm
General data	
Dimensions (W x D x H)	88 x 212 x 294 mm
Weight	3,1 kg
Permissible ambient temperature	5 – 40 °C
Permissible relative humidity	80 %
Protection class acc. to DIN EN 60529	IP 20

Technical data	
Stirring quantity (H ₂ O)	80
Max. viscosity	50.000 mPas
Motor rating input	270 W
Motor rating output	180 W
Output at stirring shaft	135 W
Max. ON-time	100 %
Max. torque at chuck	
per 60 rpm	1.515 Ncm
per 100 rpm	911 Ncm
per 1.000 rpm	91 Ncm
Speed range I (per 50 Hz)	40 – 370 rpm
Speed range II (per 50 Hz)	120 – 1.400 rpm
Speed range I (per 60 Hz)	48 – 444 rpm
Speed range II (per 60 Hz)	144 – 1.680 rpm
Speed display	scale
Chuck range	1 – 10 mm
Diameter / length of extension arm	16 mm / 160 mm
General data	
Dimensions (W x D x H)	140 x 279 x 468 mm
Weight	9,3 kg
Permissible ambient temperature	5 – 40 °C
Permissible relative humidity	80 %

RW 28 D

d pilot plant stations.

Push-through agitator shafts

cessories (page): pport holder (42)

RW 28 basic

Powerful, mechanically controlled stirrer. Suitable for quantities up to 80 I (H₂O) for use in laboratories and pilot plant stations.

- Two selectable speed ranges for high viscosity (range I) or intensive mixing (range II)
- Push-through agitator elements

Accessories (page):

Stands (114 / 115): R 2722, R 2723, R 271 Boss head clamp (116), FK 1 Flexible coupling (42), RH 5 Strap clamp (116), R 301 Stirring shaft protection (42), Stirring elements (40 / 41): e.g. R 1345, R 1300, R 301.1 Support holder (42)

echnical data	
tirring quantity (H ₂ O)	80
flax. viscosity	50.000 mPas
Notor rating input	220 W
Notor rating output	90 W
Output at stirring shaft	90 W
/lax. ON-time	100 %
/lax. torque at chuck	
er 60 rpm	1.144 Ncm
er 100 rpm	900 Ncm
er 1.000 rpm	86 Ncm
ipeed range I (per 50 Hz)	60 – 400 rpm
peed range II (per 50 Hz)	240 – 1.400 rpm
peed range I (per 60 Hz)	72 – 480 rpm
peed range II (per 60 Hz)	288 – 1.680 rpm
peed display	scale
Chuck range	1 – 10 mm
lollow shaft, inner diameter	10,5 mm
Diameter / length of extension arm	16 mm / 145 mm
General data	
)imensions (W x D x H)	123 x 252 x 364 mm
Veight	7,4 kg
ermissible ambient temperature	5 – 40 °C
Permissible relative humidity	80 %
rotection class acc. to DIN EN 60529	IP 42

Technical data	
Stirring quantity (H ₂ O)	200
Max. viscosity	100.000 mPas
Motor rating input	513 W
Motor rating output	370 W
Output at stirring shaft	300 W
Max. ON-time	100 %
Max. torque at chuck	
per 60 rpm	4.642 Ncm
per 100 rpm	3.000 Ncm
per 1.000 rpm	285 Ncm
Speed range I (per 50 Hz)	57 – 275 rpm
Speed range II (per 50 Hz)	275 – 1.300 rpm
Speed range I (per 60 Hz)	69 – 330 rpm
Speed range II (per 60 Hz)	330 – 1.560 rpm
Speed display	scale
Chuck range	3 – 16 mm
Hollow shaft, inner diameter	13 mm
Fixing	flange
General data	
Dimensions (W x D x H)	145 x 340 x 445 mm
Weight	15 kg
Permissible ambient temperature	5 – 40 °C
Permissible relative humidity	80 %
Protection class acc. to DIN EN 60529	IP 54

Protection class acc. to DIN EN 60529

RW 47 D

IP 54

The most powerful IKA® stirrer for laboratories, pilot plant stations and small-scale production. - For stirring tasks up to 200 I (H₂O) - Two speed ranges for highly viscous media and intensive mixing - Cables with plugs not included in delivery

Accessories (page):

Safety switch (43), Fixing devices (43): SI 472, SI 474

IKA[®] Mixing Mechanical overhead stirrers 39

- werful, mechanically controlled stirrer with AC otor and high IP protection class. Suitable for antities up to 80 I (H₂O) for use in laboratories
- wo selectable speed ranges for high viscosity range I) or intensive mixing (range II)
- Cables with plugs not included in delivery

ands (114 / 115): R 2722, R 2723, R 271 Boss ad clamp (116), RH 5 Strap clamp (116), 301 Stirring shaft protection (42), Stirring ements (40 / 41): e.g. R 1345, R 1301, R 301.1



3297000 3 x 400 V 50 Hz 3297006 3 x 230 V 60 Hz

R 472 Floor stand (115), R 474 Telescopic stand (115), R 302 Shaft protection (43), Stirring elements (40 / 41): e.g. R 2305, R 2311, SI 400



Ident. No.		
1602000	3 x 400 V	50 Hz
1602010	3 x 230 V	60 Hz

40 Stirring elements (stainless steel AISI 316L)



Propeller stirrer, 4-bladed

Standard stirring element. For drawing the material to be mixed from the top to the bottom. Local shearing forces. Generates axial flow in the vessel. Used at medium to high speeds

Propeller stirrer, 3-bladed

Flow-efficient design. For drawing the material to be mixed from the top and the bottom. Minimum shearing forces. Used at medium to high speeds.

Turbine stirrer

For drawing the material to be mixed from above. Generates axial flow in the vessel. Minimum danger of injury when contact is made with vessel. Minimum shearing forces. Used at medium to high speeds.

Dissolver stirrer

Radial flow, for drawing the material to be mixed from the top and the bottom. High turbulence, high shearing forces. Particle reduction. Used at medium to high speeds.

Centrifugal stirrer

Two-bladed, blades open with increasing speed. For stirring in round vessels with narrow necks. Effect is similar to that of a 4-bladed propeller stirrer. Medium to high speeds required.

Paddle stirrer

Tangential flow, minimum turbulence, good heat exchange, gentle treatment of product. Used at low to medium speeds.

Anchor stirrer

Tangential flow, high shearing rate at edges, minimum deposits on the vessel wall. Used at low speeds. Polymer reactions, even distribution of high mineral contents in liquids. The ideal stirrer for medium to highly viscous fluids.

Kneading stirrer

Tangential flow with oscillating compacting between the kneading surfaces. Minimum deposits on vessel. Used at low speeds.

	Ident. No.	Stirrer-Ø	Shaft Ø	Shaft length	Max.	RW 14 basic	EUROSTAR	EUROSTAR power
					speed	RW 16 basic	digital	power control-visc
Propeller stirrer, 4-bladed								
R 1342	0741000	50 mm	8 mm	350 mm	2.000 rpm	•	•	•
R 1345	0741300	100 mm	8 mm	540 mm	800 rpm			•
R 2305	0739300	150 mm	13 mm	550 mm	1.300 rpm			
R 2302	0739000	150 mm	13 mm	800 mm	600 rpm			
Propeller stirrer, 3-bladed								
R 1381	1296000	45 mm	8 mm	350 mm	2.000 rpm	•	•	•
R 1382	1295900	55 mm	8 mm	350 mm	2.000 rpm	•	•	•
R 1385	0477700	140 mm	10 mm	550 mm	800 rpm			
R 1388	0477800	140 mm	10 mm	800 mm	400 rpm			
R 1389 (PTFE-coated)	2343600	75 mm	8 mm	350 mm	800 rpm	•	•	•
Turbine stirrer								
R 1311	2332900	30 mm	8 mm	350 mm	2.000 rpm	•	•	•
R 1312	2333000	50 mm	8 mm	350 mm	2.000 rpm	•	•	•
R 1313	2333100	70 mm	10 mm	400 mm	800 rpm			•

Dissolver stirrer								
R 1300	0513500	80 mm	8 mm	350 mm	2.000 rpm		•	•
R 1302	2387900	100 mm	10 mm	350 mm	1.000 rpm			•
R 1303	2746700	42 mm	8 mm	350 mm	2.000 rpm	•	•	•

Centrifugal stirrer								
R 1352	0756900	60 / 15 mm	8 mm	350 mm	2.000 rpm	•	•	•
R 1355	1132700	100 / 24 mm	8 mm	550 mm	800 rpm			•

Paddle stirrer								
R 1373	0757600	70 mm	8 mm	350 mm	1.000 rpm		•	•
R 1375	0757700	70 mm	8 mm	550 mm	800 rpm			•
R 1376	0757800	150 mm	10 mm	550 mm	800 rpm			
R 2311	0739500	150 mm	13 mm	800 mm	600 rpm			
Anchor stirrer								
R 1330	2022300	45 mm	8 mm	350 mm	1.000 rpm	•	•	•
R 1331	2022400	90 mm	8 mm	350 mm	1.000 rpm			•
R 1332 (PTFE-coated)	2343700	60 mm	8 mm	350 mm	800 rpm	•	•	•
R 1333	2747400	150 mm	10 mm	550 mm	800 rpm			

Kneading stirrer								
R 1335	2022500	45 mm	8 mm	350 mm	2.000 rpm	•	•	•

IKA[®] Mixing

Stirring elements (stainless steel AISI 316L) 41

basic / / P1	EUROSTAR power control-visc P4 / P7	RW 20 digital	RW 28 basic / RW 28 D	RW 47 D
	•			
	•		•	•
				•
				•
	•	•		
	•	•		
	•		•	•
	•		•	•
	•	•		
	•	•		
	•	•		
	•		•	
		•	•	
			•	•
		•		
		•		
			•	•
	•	•	•	
	•		•	
	•		•	•
				•
	-	_		
	•	•		
	•		•	
	•	•	•	
	•		•	•

٠

٠

Overhead stirrers accessories 42



R 1001 Paddle stirrer	
	Spare for use with RW 11 basic.

R 1002 Screw-type stirrer

For use with RW 11 basic.

Ident. No.	
0527400	

ldent. No.	
0527500	



Ident No 1242900 R 1401 1243300 R 1402 1289800 R 1405







ł	1401	Propeller	

R 1402 Dissolver

R 1405 Propeller

structural variances.

EUROSTAR series.

For use with EUROS

General data	
Shaft length	160 mm
Shaft Ø	4 mm
Stirrer Ø	34 mm

General data	
Shaft length	140 mm
Shaft Ø	4 mm
Stirrer Ø	12 mm

R 1401 Propeller	
Working range	1 – 30 I
Rotor diameter	55 mm
R 1402 Dissolver	
Working range	1 – 30 I
Rotor diameter	42 mm
R 1405 Propeller	
Working range	0,25 - 30
Rotor diameter	45 mm

A second s	Statement of the local division in which the local division in the local division in the local division in the

FK 1 Flexible coupling	
	General data
Required for stirring tasks using glass stirring	Clamping range
rods. The flexible coupling compensates for any	Max. torque

6 – 10 mm
10 Ncm

R 301 Stirring shaft protection		
	General data	
Prevents potential injuries at the rotating shafts and	Length adjustment	190 – 310 mn
stirring elements. Can be directly attached to the	Material	plexiglas
stirring motors RW 16 basic, RW 20 digital and the		

General data	
Length	275 mm
Material	plexiglass

General data	
Dimensions (W \times D \times H)	139 x 99 x 250 mm
Material	macrolon

Prevents potential injuries due to the rotating shafts and stirring elements. Can be directly attached to the stirrer RW 47 D.

SI 400 Safety switch Dimensions end switch (B x T x H) 84 x 19 x 16 mm Dimensions switch contact (B x T x H) 73 x 10 x 19 mm 1 normally closed contact

Casing materia Protection class according to DIN EN 60529 Operating temperature max. 250 VAC / 2A Voltage / current

General data

Contact

plastic (ABS) IP 67 -10-65 °C

SI 474. The stirring unit RW 47 can only be switched on through the SI 400, when the agitator is adjusted in the mixing vessel to the user designated height. The power of the RW 47 automatically shuts off if the stirring unit is lifted off the designated height. Also suitable for dispersing instrument T 65 D ULTRA-TURRAX®.

Accessories (page): Fixing devices (43): SI 472, SI 474

SI 472 Fixing device

stand R 472.

SI 474 Fixing device

ttach the safety switch SI 400 to the telescostand R 474 and to the telescopic stand T 653 (for T 65 D ULTRA-TURRAX®).

RC 1 Remote control

Remote control to operate the Eurostar power control-visc (also P1, P4 and P7) over a 10 m cable. - Provides problem-free control of stirrers

- even under load
- overload status

R 301.1 Support holder

For fixing the stirring shaft protection R 301 to the stand.

R 301 Stirring shaft protection

Accessories (page): Boss head clamp (116): R 182, R 270

General data	
ength	275 mm
Naterial	plexiglass

General data	
Dimensions	80 x 80 mm

General data Dimensions (W x D x H)

95 x 83 x 20 mm	To at
	pic s
	/ 6

General data	
Power supply	Two 1.5 V batteries
	(included with unit)
Max. cable length	10 m
Power consumption remote	
Off-state	ca. 7 μA
On-state	ca. 7 mA
Dimensions (W x D x H)	65 x 140 x 30 mm
Weight (incl. battery)	0,3 kg

	VVOLKI
STAR power control-visc 6000.	Rotor
	R 140
	Worki
	Rotor

ropeller	
range	1 – 30 I
meter	55 mm
issolver	
range	1 – 30 I



	Working range
	Rotor diameter
	R 1402 Dissolver
	Working range
).	Rotor diameter
	R 1405 Propeller
	Working range

IKA[®] Mixing Overhead stirrers accessories

R 302 Stirring shaft protection



The SI 400 consists of an end switch (normally closed contact / switch) and a magnetic switch contact (actuator) which is mounted on the floor stand R 472 with the fixing device SI 472 and on the telescopic stand R 474 with the fixing device

To attach the safety switch SI 400 to the floor

Ident. No. 3294800



Ident. No. 3264400

- Displays actual speed, target speed and



IKA[®] Mixing Shakers

44



NEW!

MS 3 digital Orbital shaker: Timer with countdown function

Technical data			
Shaking movement orbita		orbital	
Orbital diamet	ter		4,5 mm
Max. permitte	d shaking weight (incl	. attachment)	0,5 kg
Motor rating i	nput		10 W
Motor rating o	output		8 W
Permissible O	N time		100 %
Infinitely adjust	stable speed range		0 – 3.000 rpm
Speed display	r		scale
Timer	MS 3 basic		no
	MS 3 digital		yes
Time setting	MS 3 basic		-
	MS 3 digital		1 s – 999 min
Operating	MS 3 basic		Continuous /
mode			touch operation
	MS 3 digital	Timer and c	ontinous mode,
			touch operation
General data			
Dimensions (W x D x H) 148 x 205 x 63		x 205 x 63 mm	
Weight 2,9 k		2,9 kg	
Protection class acc. to DIN EN 60529 IP 2		IP 21	

Technical data

Orbital diameter

Motor rating input

Speed (fixed)

General data

Attachment

Bottom

Weight

Material Casing

Dimensions (Ø x H)

Permissible ambient temperature Permissible relative humidity

Protection class acc. to DIN EN 60529

Motor rating output

Shaking movement

Shaken quantity (1 test tube)

MS 3 basic MS 3 digital

Compact, universal small shaker suitable microtiter plates

- Wide range of attachments - Attachment detection

- Continuous or touch operation

(with standard attachment) - Two operating modes: limited to 1.300 rpm. Mode B (without attachment detection) attachments.

- Sturdy zinc die cast casing

MS 3 digital additionally: Timer with countdown function

Accessories (page):

lab dancer

Economic, compact test tube shaker with touch function.

- 1,2 W
 - centrifuge tubes, Eppendorf vessels
 - Excellent mixing action
 - are made from inert plastic
- zinc, coated 100 x 70 mm
- 5 40 °C

orbital 4,5 mm max. 50 ml

0,8 W

PP

TPU

80 %

IP 40

2.800 rpm

- indispensable tool for every laboratory - The lab dancer can be used with all small vessels
- of up to 30 mm in diameter, e.g. test tubes,
- The upper casing and the test tube surface
- Secure stand thanks to coated zinc die cast base - Includes light 12 V power pack set 0,55 kg



- for shaking tasks with all small vessels and
- Mode A (safe mode with attachment detection) The maximum speed of 3.000 rpm is only reached with the standard attachment in touch mode. When using other attachments the speed is
- A speed of 3.000 rpm is possible with all
- Stable in all speed ranges
- Attachments (50): MS 1.31, MS 1.32, MS 1.33
- Included with delivery (page):
- MS 3.1 Standard attachment (50), MS 3.3 Universal attachment (50), MS 1.21 One-hand insert (50), MS 3 digital complete with MS 3.4 Microtiter attachment and MS 1.32 Test tube insert (50)



MS 3 basic

Ident. No.	
3617000	230 V 50/60 Hz
3617001	115 V 50/60 Hz





MS 3 digital Ident. No.

3319000	230 V	50/60 Hz
3319001	115 V	50/60 Hz

NEW!

- Its compact and clever design makes it an



Ident. No. 3365000 100 - 240 V 50/60 Hz

Orbital shakers 46



Ident. No.	
3340000	230 V 50/60 Hz
3340001	115 V 50/60 Hz

VORTEX Genius 3

New vortex shaker suitable for short-time operation (touch function), activated by pressing shaker attachment or continuous operation.

Technical data

Orbital diameter

Speed display

Speed setting

General data

Weight

Shaking movement

Motor rating input

Motor rating output

Permissible ON time

Dimensions (W x D x H)

Infinitely adjustable speed range*

Permissible ambient temperature

Protection class acc. to DIN EN 60529

Permissible relative humidity

- Wide speed range, infinitely adjustable - Different applications thanks to 3 inter-
- changeable attachments and 7 inserts (e.g. Eppendorf tubes, microtiter plates, Erlenmeyer
- flasks 250 ml etc.), please order separately Attachments securely click onto appliance in any position
- Special strap (VG 3.36, page 51) ensures
- easy handling of round/Erlenmeyer flasks - Sturdy metal zinc die cast casing
- Compact design
- Short-time operation activated by pressing
- attachment (touch function) - Stable at high speeds thanks to special feet (silicon base with ultra high vibration damping)
- Eccentric with ball bearings
- Suitable for continuous operation with low self heating thanks to self ventilation of motor

Accessories (page):

Attachments (51): VG 3.1, VG 3.2, VG 3.3 Inserts (51): VG 3.31, VG 3.32, VG 3.33, VG 3.34, VG 3.35, VG 3.36, VG 3.37



Ident. No. 2819000 230 V 50/60 Hz 2819001 115 V 50/60 Hz

VXR basic Vibrax®

Optoelectronically controlled small shaker with a very wide speed range.

- Suitable for continuous operation
- New design and improved drive system
- Circular shaking motions
- Slow speeds are well maintained
- Attachments are interchangeable

Accessories (page):

Attachments (52 / 53): VX 1, VX 2, VX 2E, VX 7, VX 8, VX 8.1, VX 11, VX 11.1, VX 11.2, VX 11.3, VX 11.4

Technical data	
Shaking movement	orbital
Orbital diameter	4 mm
Max. shaking weight (with attachment)	2 kg
Motor rating input	35 W
Motor rating output	13,2 W
Permissible ON time	100 %
Speed range	0 – 2.200 rpm
Speed display	scale
General data	
Dimensions (W x D x H)	157 x 247 x 130 mm
Weight (without attachment)	6,1 kg
Permissible ambient temperature	5 – 40 °C
Permissible relative humidity	80 %
Protection class acc. to DIN EN 60529	IP 21

Technical data Shaking movement orbital Orbital diameter 3 mm 2 or 4 microtiter plates Max. shaking weight Motor rating input 35 W Motor rating output 13,2 W 100 % l time 0 – 1.100 rpm scale ∞ / 1 – 99 min digital x D x H) 185 x 320 x 105 mm 2,7 kg

Special shaker for shaking two or four microtiter plates.

- Electronic speed control
- Digital timer

Technical data		
Shaking movem	ent	orbital
Orbital diameter		4 mm
Max. shaking w	eight (with attachment)	2 kg
Motor rating inp	ut	45 W
Motor rating out	put	10 W
Permissible ON	time	100 %
Speed range		80 – 800 rpm
Speed display	KS 130 basic	LED line
	KS 130 control	digital
Timer	KS 130 basic	∞ / 5 – 50 min
	KS 130 control	∞ / 0 – 9 h 59 min
General data		
Dimensions (W	x D x H)	270 x 316 x 98 mm
Weight	KS 130 basic	8,8 kg
	KS 130 control	9,8 kg
Permissible amb	pient temperature	5 – 50 °C

KS 130 basic KS 130 control

80 %

IP 21

RS 232 / analog

Small, guiet shaker ensures long life with ideal swivel motion, for a maximum shaking weight of 2 kg.

- Electronic adjustment of speed and timer
- sizes of vessels
 - separately

KS 130 control additionally with:

- modes
- 80 % - Electronic time switching clock: IP 21

 - documented with labworldsoft® - Special version with reverse rotating
 - direction on request

Accessories (page): Attachments (54): AS 130.1, AS 130.2, AS 130.3, AS 130.4 KS 130 control additionally: labworldsoft® (143), PC 1.5 Cable (148)

500 – 2.500 rpm	Permissible ON
scale 0 – 6	Speed range
knob, front	Speed display
	Timer
127 x 149 x 136 mm	Timer display
4,5 kg	General data
5 – 40 °C	Dimensions (W

orbital

4 mm

58 W

10 W

100 %

80 %

IP 21

*depending on attachment and loading

Weight Permissible ambient temperature 5 – 40 °C Permissible relative humidity

Protection class acc. to DIN EN 60529

Permissible relative humidity

Interface

Protection class acc. to DIN EN 60529

KS 130 control



MTS 2/4 digital microtiter shaker

- Alarm to indicate set time has expired - Incl. attachment (without microtiter plate)



ldent. No.		
3208000	230 V	50/60 Hz
3208001	115 V	50/60 Hz

- LED display for speed and time adjustment - Wide range of attachment combinations makes it possible to use almost all shapes and

- Attachments are not included, please order

- A digital display makes it possible to read the speed, timer function and operating

0 - 9 h 59 min or continuous operation (∞) - With integrated end point positioning (for automated robot-controlled sampling) - All functions can be controlled and



KS 130 basic Ident. No. 230 V 50/60 Hz 2980000 115 V 50/60 Hz 2980001



KS 130 control Ident. No. 230 V 50/60 Hz 2980100 115 V 50/60 Hz 2980101



Orbital shakers 48



KS 260 basic Ident. No. 2980200 230 V 50/60 Hz 2980201 115 V 50/60 Hz



KS 260 control Ident. No. 2980300 230 V 50/60 Hz 2980301 115 V 50/60 Hz



-	June -		Call	
-	Non In			
		N AS		
	L	-0-		
		**		

230 V 50/60 Hz
115 V 50/60 Hz

KS 260 basic KS 260 control

Compact, flat shaker with ideal swivel motion, for a maximum shaking weight of 7,5 kg.

- Electronic adjustment of speed and timer - LED display for speed and time adjustment
- Wide range of attachment combinations makes it possible to use almost all shapes and
- sizes of vessels - Attachments are not included, please order
- separately

KS 260 control additionally with:

- Digital display makes it possible to read the speed, timer function and operating modes
- Electronic time switching clock:
- 0 9 h 59 min or continuous operation (∞) With integrated end point positioning
- (for automated robot-controlled sampling) All functions can be controlled and
- documented with lab Special version with direction on request

Accessories (page):

Attachments (54 / 55): AS 260.1, AS 260.2, AS 260.3 KS 260 control additionally: labworldsoft® (143), PC 1.5 Cable (148)

controlled and	
oworldsoft®	
reverse rotating	

KS 501	digital
--------	---------

of up to 15 kg.

- Infinitely variable speed control of 0 300 rpm - Digital display
- Ideal for vessels with a volume of more than 250 ml, e.g. round flasks, Erlenmeyer flasks, culture flasks and culture bottles
- Guaranteed continuous operation (∞) even under extreme loads
- Includes timer
- Attachments are not included, please order separately

Acc Atta AS

Technical data		
Shaking movem	ent	orbital
Orbital diameter		10 mm
Max. shaking w	eight (with attachment)	7,5 kg
Motor rating inp	ut	45 W
Motor rating out	tput	10 W
Permissible ON	time	100 %
Infinitely adjusta	ble speed range	
	KS 260 basic	20 – 500 rpm
	KS 260 control	10 – 500 rpm
Speed display	KS 260 basic	LED line
	KS 260 control	digital
Timer	KS 260 basic	∞ / 5 – 50 min
	KS 260 control	∞ / 9 h 59 min
Timer display	KS 260 control	digital
General data		
Dimensions (W	x D x H)	360 x 420 x 98 mm
Weight	KS 260 basic	8,5 kg
	KS 260 control	8,8 kg
Permissible ambient temperature		5 – 50 °C
Permissible rela	tive humidity	80 %
Protection class acc. to DIN EN 60529		IP 21
Interface	KS 260 control	RS 232 / analog

Technical data		
Shaking movem	ent	reciprocating
Orbital diameter		20 mm
Max. shaking we	eight (with attachment)	7,5 kg
Motor rating inp	ut	45 W
Motor rating out	put	10 W
Permissible ON	time	100 %
Infinitely adjusta	ble speed range	
	HS 260 basic	20 – 300 rpm
	HS 260 control	10 – 300 rpm
Speed display	HS 260 basic	LED line
	HS 260 control	digital
Timer	HS 260 basic	∞ / 5 – 50 min
	HS 260 control	∞ / 9 h 59 min
Timer display	HS 260 control	digital
General data		
Dimensions (W	x D x H)	360 x 420 x 100 mm
Weight	HS 260 basic	8,5 kg
	HS 260 control	8,8 kg
Permissible amb	pient temperature	5 – 50 °C
Permissible rela	tive humidity	80 %
Protection class	acc. to DIN EN 60529	IP 21
Interface	HS 260 control	RS 232 / analog

HS 260 basic HS 260 control

- Wide range of attachment combinations makes it possible to use almost all shapes and sizes of vessels
- please order separately
 - HS 260 control additionally:
 - Digital display makes it possible to read the speed, timer function and operating mode - Electronic time switching clock:
 - 0 9 h 59 min or continuous operation (∞) - With integrated endpoint positioning

 - (for automated robot-controlled sampling) - All functions can be controlled and

Accessories (page): AS 260.3, AS 260.5 labworldsoft® (143)

Low profile laboratory shaker with a pleasant design, large mounting surface and load capacity

- Orbital diameter 30 mm

achments (55 / 56): AS 501.1, AS 501.4	ļ,
501.5	

echnical data	
Shaking movement	orbital
Drbital diameter	30 mm
Nax. shaking weight (with attachment)	15 kg
Notor rating input	70 W
Notor rating output	19 W
Permissible ON time	100 %
nfinitely adjustable speed range	0 – 300 rpm
speed display	digital
imer	∞ / 0 – 56 min
General data	
Dimensions (W x D x H)	505 x 585 x 120 mm
Veight	26 kg
Permissible ambient temperature	5 – 50 °C
Permissible relative humidity	80 %
Protection class acc. to DIN EN 60529	IP 21

Technical data	
Shaking movement	reciprocating
Orbital diameter	30 mm
Max. shaking weight (with attachment)	15 kg
Motor rating input	70 W
Motor rating output	19 W
Permissible ON time	100 %
Infinitely adjustable speed range	0 – 300 rpm
Speed display	digital
Timer	∞ / 0 – 56 min
General data	
Dimensions (W x D x H)	505 x 585 x 120 mm
Weight	26 kg
Permissible ambient temperature	5 – 50 °C
Permissible relative humidity	80 %
Protection class acc. to DIN EN 60529	IP 21

IS 501 digital

- f up to 15 kg.
- Digital display
- Orbital diameter 30 mm
- - under extreme loads Includes timer
- please order separately

Accessories (page): Attachments (55 / 56): AS 501.1, AS 501.2, AS 501.3, AS 501.4, AS 501.5, AS 501.6

IKA[®] Mixing Horizontal shakers

Compact, flat shaker with ideal swivel motion, for a maximum shaking weight of 7,5 kg. - Electronic adjustment of speed and timer - LED display for speed and time adjustment

- Attachments are not included in delivery,

- documented with labworldsoft® software
- Attachments (54 / 55): AS 260.1, AS 260.2,
- HS 260 control additionally: PC 1.5 Cable (148),



HS 260 basic

Ident. No.	
3066600	230 V 50/60 Hz
3066601	115 V 50/60 Hz



HS 260 control

Ident. No.		
3066700	230 V	50/60 Hz
3066701	115 V	50/60 Hz

ow profile laboratory shaker with a pleasant lesign, large mounting surface and load capacity

Infinitely variable speed control of 0 - 300 rpm

Ideal for all lying vessels, e.g. separating funnels Guaranteed continuous operation (∞) even

Attachments are not included in delivery,



2527001	115 V 50/60 Hz
2527000	230 V 50/60 Hz
Ident. No.	

IKA[®] Mixing 50 Shakers accessories (MS 3 basic / digital)



3426300





Ident. No. 3426400













For test tubes and small vessels up to Ø 50 mm, included with the minishakers MS 3 basic and MS 3 digital.

MS 3.3 Universal attachment

For various foam inserts, included with the minishakers MS 3 basic and MS 3 digital.

MS 3.4 Microtiter attachment

For use with a microtiter plate, included with the minishakers MS 3 basic and MS 3 digital.

MS 1.21 One-hand insert

For inserting into the universal attachment, included with the minishakers MS 3 basic and MS 3 digital.



For inserting into universal attachment, for 14 test tubes Ø 10 mm, material: ethylvinyl-acetate.

MS 1.32 Test tube insert

For inserting into the universal attachment, for 6 test tubes Ø 12 mm. Material: ethylvinyl-acetate.

MS 1.33 Test tube insert

For inserting into the universal attachment, for 4 test tubes Ø 16 mm. Material: ethylvinyl-acetate. Included with the minishaker MS 3 digital.

MS 1.34 Test tube insert

For inserting into the universal attachment, for any number of bore holes. Material: ethylvinyl-acetate.

VG 3.1 Standard a Standard attachment vessels (continuous / included with delivery

VG 3.2 One-hand One-hand attachment rubber insert (continu

VG 3.3 Universal Universal attachment insert (continuous / to

VG 3.31 Test tube For 54 Eppendorf tub

VG 3.32 Test tube For 18 reagent glasse operation).

VG 3.33 Test tube For 12 reagent glasse operation).

VG 3.34 Test tube For 8 reagent glasses operation).

VG 3.35 Test tube For 8 reagent glasses operation).

VG 3.36 Erlenmey For 1 Erlenmeyer / ro 250 ml (continuous op

VG 3.37 Microtite For 1 standard microt operation).

IKA[®] Mixing Shakers accessories (VORTEX Genius 3) 51

attachment		
for reagent glasses / small		1 mil
touch operation),	Ident. No.	
1.	3341200	and the second second
attachment		
t. 88 mm. round. with	Ident, No.	
ious / touch operation).	3342300	Starses .
attachment		
, 150 mm, with rubber	Ident. No.	
ouch operation).	3342400	
attachmont*		2 2 2 2 2 2 3 3 2 3
es (continuous operation)		1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1
	Ident, No.	
	3344300	- Anno
1		
attachment*		
es, 10 mm (continuous		and the second s
	Ident No	
	3343900	
attachment*		
es, 12 mm (continuous		The second s
	11 · N	SP
	Ident. No.	
	3344000	
attachment*		
s, 16 mm (continuous		and the second s
		- F
	Ident. No.	
	3344100	Charles and the second s
attachment*		
s, 20 mm (continuous		
	Ident. No.	
	3344200	HACCOM STATE
er flask attachment*		
und flask from 100 to		
peration).		
	Ident. No.	
	3342100	and the second second
r nlate attachment*		
titer plate (continuous		
		and the second s

Ident, No. 3344400

IKA[®] Mixing 52 Shakers accessories (VXR basic)



VX 1 One-hand	attachment

For shaking single, non-fixed vessels of 1 - 250 ml.

	Sec. of Sec.
ldent. No.	
0568900	

VX 2 Test tube attachment

For up to 36 test tubes or centrifugal tubes with a diameter of 16 mm.

General data		

130 x 135 x 40 mm

140 x 145 x 115 mm

macrolon

300 g

160 g

Gaparal data		VX 11.1
Bore holes (number)	70	Attachme
Hole Ø	10 mm	



VX 2E "Eppendorf" attachment

For intensive mixing of up to 64 "Eppendorf" tubes (1.5 ml).

VX 7 Dish attachment

Petri dishes, ect.

For careful mixing of culture bottles,

General data	
Dimensions (W \times D \times H)	210 x 210 x 65 mm
Weight	790 g

General data	
Bore holes (number)	41
Hole Ø	12 mm

VX 11.2 Test tube insert

Attachment for test tubes.



0953300

Ident. No.

3375400



VX 8 Universal attachment

For rapid and secure clamping, e.g. 2 Erlenmeyer flasks up to 500 ml.

General data	
Dimensions (W x D x H)	265 x 136 x 60 mm
Clamping range	25 – 135 mm
Min. height of vessel	80 mm
Weight	760 g

General data	
Bore holes (number)	32
Hole Ø	16 mm

VX 11.3 Test tube insert

Attachment for test tubes.

		VX 11.4
General data		
Bore holes (number)	18	Attachme
Hole Ø	20 mm	

VX 8.1 Clamping roll (without fig.)

Spare for use with VX 8 universal attachment.

General data

Weight

Dimensions (W x D x H)

General data

Weight

Material

Weight

Dimensions (W x D x H)

Dimensions (W x D x H)

Gene	ral d	ata							
D.:			-			005	100	~~	

(x D x H)	265 x 136 x 60 mm
9	25 – 135 mm
vessel	80 mm
	760 a

410 x 210 x 40 mm

740 g

IKA[®] Mixing Shakers accessories (VXR basic)

VX 11 Basic holder

Attachment for test tube inserts.



Test tube insert

ent for Eppendorf tubes or test tubes.



Ident. No. 3659000





Ident. No. 3659100



Ident. No. 3659200



ent for test tubes.



Ident. No. 3659300

53

IKA[®] Mixing 54 Shakers accessories (KS 130 und HS / KS 260)



Ident, No. 8017300



3115000



3120000





AS 130.1 Universal attachment

For use with various types of vessels by means of universal, infinitely variable clamping rolls.

Included with delivery (page): 1 x AS 1.30 Basic holder (57), 3 x AS 1.31 Clamping roll (57), 6 x AS 1.5 Fastening screw (57)

AS 130.2 Fixing clip attachment

For processing round flasks, measuring flasks and Erlenmeyer flasks. Please order fixing clips separately.

Accessories (page): Fixing clips (58): AS 2.1, AS 2.2, AS 2.3, AS 2.4, AS 2.5

AS 130.3 Dish attachment

AS 130.4 Test tube attachment

For intensive shaking, e.g. small tubes,

test tubes, cuvettes, centrifuge tubes.

AS 260.1 Universal attachment

Included with delivery (page):

1 x AS 1.60 Basic holder (57),

4 x AS 1.61 Clamping roll (57),

8 x AS 1.5 Fastening screw (57)

vessels.

For various types of vessels. Infinitely variable

clamping rolls allow universal adaptation to the

For smooth shaking operations in the low viscosity range, e.g. Petri dishes or culture bottles. With integrated slip-resistant foil (PP).

General data	
Dimensions (W x D x H)	420 x 270 x 32 mm
Set-up plate	220 x 340 mm
Veight	370 q

325 x 234 x 88 mm

230 x 230 x 24 mm

20 x AS 2.1 (25 ml)

12 x AS 2.2 (50 ml)

12 x AS 2.3 (100 ml)

4 x AS 2.4 (250 ml)

4 x AS 2.5 (500 ml)

650 g

220 x 220 mm

850 g

General data

Set-up plate

General data

Capacity:

Weight

Dimensions (W x D x H)

Number of fixing clips (volume)

Weight

Dimensions (W x D x H)

General data	
Dimensions (W x D x H)	220 x 230 x 95 mm
Capacity	64
Vessel Ø	10 – 16 mm
Min. height of vessel	80 mm
Weight	670 g

General data	
Dimensions (W \times D \times H)	425 x 335 x 135 mm
Set-up plate	320 x 320 mm
Weight	1.600 g

General data	
Dimensions (W x D x H)	330 x 330 x 24 mm
Capacity:	
Number of fixing clips (volume)	56 x AS 2.1 (25 ml)
	23 x AS 2.2 (50 ml)
	23 x AS 2.3 (100 ml)
	11 x AS 2.4 (250 ml)
	9 x AS 2.5 (500 ml)
Weight	1.290 g

General data	
Dimensions (W x D x H)	410 x 370 x 32 mm
Set-up plate	320 x 320 mm
Weight	460 g

S 260.2 Fixing clip attachment

or shaking flasks, Erlenmeyer flasks and bottles vith a round crosssection (without fixing clips).

ccessories (page): ixing clips (58): AS 2.1, AS 2.2, AS 2.3, AS 2.4, AS 2.5

For smooth movement for cell cultures, nutrient media in Petri dishes, culture bottles and vessels with a low center of gravity. With integrated slipresistant foil (PP).

General data	
Dimensions (W x D x H)	334 x 425 x 145 mm
Capacity: (number of separating	6 x 50 ml
funnels per volume, pear-shaped)	5 x 100 ml
	3 x 250 ml
	3 x 500 ml
Weight	1.550 g

480 x 500 x 120 mm

420 x 420 mm

4.000 g

2.640 g

General data

Set-up plate Weight

General data

Weight

Dimensions (W x D x H)

Capacity: Number of fixing clips (volume)

Dimensions (W x D x H)

For shaking out, salting out, extracting, eluting, enriching. The 3 clamping rolls (included in delivery) are height-adjustable for adaption to different separating funnel sizes. The separating funnels are secured with O-rings (5 O-rings included).

AS 501.1 Universal attachment

For various types of vessels with a minimum volume of 50 ml. Ideally more than 250 ml. The clamping rolls may be adjusted to two levels.

Included with delivery (page): 1 x AS 1.10 Basic holder (57), 6 x AS 1.11 Clamping roll (57), 12 x AS 1.6 Fastening screw (57)

AS 501.4 Fixing clip attachment

475 x 460 x 95 mm 50 x AS 2.1 (25 ml) 48 x AS 2.2 (50 ml) 25 x AS 2.3 (100 ml) 16 x AS 2.4 (250 ml) 12 x AS 2.5 (500 ml)

Accessories (page): Fixing clips (58): AS 2.1, AS 2.2, AS 2.3, AS 2.4, AS 2.5

Ident No. 8017400

IKA[®] Mixing

Shakers accessories (HS / KS 260 und HS / KS 501) 55

Ident. No. 3115500

AS 260.3 Dish attachment

AS 260.5 Separating funnel attachment



Ident. No. 3120900

For shaking flasks, Erlenmeyer flasks and pear-shaped flasks (without fixing clips).





IKA[®] Mixing Shakers accessories (HS / KS 501) 56





For smoothly shaking dishes, but also for smooth mixing in vessels with a large, flat bottom (widenecked Erlenmeyer flasks and beakers). A plastic foil with mild adhesive prevents the vessel from slipping.

General data

Set-up plate

Weight

Weight

Dimensions (W x D x H)

	E	6
Ł	 	A

Ident. No 8000300



Ident. No. 8000400



Ident. No. 8000500

AS 501.2 Separating funr	nel attachment
--------------------------	----------------

For shaking out, eluting, extracting, gassing out, dissolving, enriching, etc. Adjustment for the clamping rolls is infinitely variable, the set-up height can be changed by means of clamping devices.

Included with delivery (page): 1 x AS 1.10, 6 x AS 1.11, 6 x AS 1.6 (57), 6 x AS 1.7 (58)

S 501.3 Separating funnel attachment	
--------------------------------------	--

Same features as AS 501.2.

Included with delivery (page): 1 x AS 1.10, 4 x AS 1.11, 4 x AS 1.6 (57), 4 x AS 1.7 (58)

AS 501.6 Separating funnel attachment

Same features as AS 501.2. This attachment will

hold 4 x 1.000 ml separating funnels.

Included with delivery (page): 1 x AS 1.10, 4 x AS 1.6 (57), 4 x AS 1.12,

8 x AS 1.13 (58)

General data	
Dimensions (W x D x H)	480 x 505 x 225 mm
Capacity: (number of separating funnels	
per volume, pear-shaped)	4 x 1.000 ml
Weight	5.500 q

Concernel data	
General data	
Dimensions (W x D x H)	480 x 505 x 190 mm
Capacity: (number of separating funnels	12 x 50 ml
per volume, pear-shaped)	10 x 100 ml
	6 x 250 ml

450 x 450 x 45 mm

420 x 420 mm

1.120 g

4.180 g

General data	
Dimensions (W x D x H)	252 x 234 x 88 mm

AS 1.30 Basic holder

Accessories (page): AS 1.31, AS 1.5 (57)

AS 1.60 Basic holder

For use with universal attachment AS 260.1. 348 x 335 x 135 mm

> Accessories (page): AS 1.61, AS 1.5 (57)

General data 480 x 480 x 120 mm Dimensions (W x D x H)

General data

Dimensions (W x D x H)

For use with universal attachment AS 501.1 and separating funnel attachments AS 501.2, AS 501.3 and AS 501.6.

Accessories (page): AS 1.11, AS 1.6 (57), AS 1.7, AS 1.8, AS 1.12,

		Clamping roll
AS 1.31		
For basic holder	AS 1.30	AS 1.31
Length	228 mm	
AS 1.61		AS 1.61
For basic holder	AS 1.60	
Length	335 mm	A C 1 11
AS 1.11		AS 1.11
For basic holder	AS 1.10	
Length AS 1.11	410 mm	

AS 1.5 Fastening screw

corresponding basic holder.

AS 1.6 Fastening screw

Two AS 1.6 clamping devices are required for fastening a clamping roll to the corresponding basic holder (for basic holder AS 1.10 only).

General data	
Dimensions (W x D x H)	480 x 505 x 190 mm
Capacity: (number of separating funnels	4 x 500 ml
per volume, pear-shaped)	3 x 1.000 ml
	2 x 2.000 ml
Weight	3.720 g

Α

IKA[®] Mixing Shakers accessories

For use with universal attachment AS 130.1.

AS 1.10 Basic holder





Ident. No.	
3030500	AS 1.31
3030501	AS 1.61
2339800	AS 1.11

Fastening screw for the universal attachments AS 130.1, AS 260.1 and the separating funnel attachment AS 260.5. Two AS 1.5 fastening screws are required for fastening a clamping roll onto the



Ident. No. 1268400

Ident. No.

2979400



Shakers accessories 58



Ident. No. 1269200



Ident. No. 1268900





Ident. No 2597000

3 AS 2.3 Fixing clip



	Ident. No.	
1	1234300	AS 2.1
2	1234400	AS 2.2
3	1234500	AS 2.3
4	1234600	AS 2.4
5	1234700	AS 2.5

AS 1.7 Clamping device

Two AS 1.6 and two AS 1.7 clamping devices are required for fastening two clamping rolls one above the other (for clamping separating funnels). For basic holder AS 1.10 only.

AS 1.8 Supporting clamping device

Two AS 1.6 clamping devices and two AS 1.8 supporting clamping devices are required if a clamping roll is to be attached at a higher position (e.g. for fixing a vessel which has a higher point of gravity). For basic holder AS 1.10 only.

AS 1.12 Supporting bar

For attaching two AS 1.13 ground section holders for fixing 1.000 ml separating funnels. For basic holder AS 1.10 only.

Accessories (page): AS 1.13 (58)

AS 2.1 Fixing clip

2 AS 2.2 Fixing clip

4 AS 2.4 Fixing clip

5 AS 2.5 Fixing clip

AS 1.13 Ground section holder

For attaching separating funnels with ground opening NS 29 (2x AS 1.13 necessary per separating funnel). For basic holder AS 1.10 only.

General data		
For flask volume	AS 2.1	25 ml
	AS 2.2	50 ml
	AS 2.3	100 ml
	AS 2.4	200 ml / 250 ml
	AS 2.5	500 ml

437 mm

General data

Length

General data	
Shape of kneading blades	duplex
Trough	
Useful volume min. / max.	100 / 300 ml
Total volume	600 ml
Attainable vacuum	50 mbar
Trough base for heating up to	210 °C
Bore hole for accomodating	yes
temperature measuring sensor PT 100.27	
Materials in contact with medium	stainless steel
	(AISI 316 Cb)
Drive	
Motor rating input	320 W
Motor rating output	180 W
Motor principle	asynchron
Motor protection	thermo contact
Nominal torque	48 Nm
Speed of front kneading blade	35 rpm
Speed of back kneading blade	18 rpm
Safety device	cover contact
Conoral data	

General da Dimensions (W x D x H) 660 x 250 x 380 mm Weight 27 kg 5 – 40 °C Permissible ambient temperature Permissible relative humidity 80 % Protection class acc. to DIN EN 60529 IP 54

D-T 06 D \®-High-performance oratory kneader

processing non-flowable, highly viscous media. orm mixing is based on intensive processing neans of wide-bladed kneading elements. The ading medium is moved within the trough both zontally and vertically. Additional media quanmay be added during the kneading operation.

e double-walled kneading chamber ows cooling or heating of the product e product temperature may be measured

- ough can easily be removed
- leading blades can easily be removed ort kneading time
- - nd trough wall ensures efficient wipe-off

Accessories (page):

HKD 06.2 Plunger

Double-walled plunger seat, for heating and cooling, which presses down the kneading material by means of spring-action. This considerably improves heat conduction in the measuring kneader. The central supply opening can be closed with a plunger. The HKD 06.2 plunger is not suitable for vacuum operation.

HKD 06.10 Kneading blade

Special lined front kneading blade for viscous, elastic products to prevent dead zones. Alternative to the standard kneading blade.

IKA[®] Mixing Laboratory kneader

- ectly behind the kneading blades
- e narrrow gap between the kneading blades
- andard version equipped for vacuum operation
- Trough cover with inspection glass and safety screen

HKD 06.2 Plunger (59), HKD 06.10 Kneading blade (59), DTM 12 Digital temperature measuring device (118), CC3-308B vpc Circulation thermostat (97), VC 2 IKAVAC[®] Vacuum controller (120)

Ident. No. 2936000





Ident. No. 3134800



1911800	3 x 400 V	50 Hz
1911803	3 x 230 V	60 Hz



M 23 Star-shaped cutter

Used to crush fibrous substances such as paper and vegetation, but also for plastics and material with a low specific weight. Page 89

Dispersers Mills 62 - 85 86 - 91

Servi

Alm Bie

Crushing



62 Dispersers (batch operation)



NEW!

ULTRA-TURRAX[®] Tube Drive World first! The revolutionary disposable disperser

Technical data Rating input 20 W 17 W Rating output Speed range, infinitely adjustable 300 – 6.000 rpm Timer 1 – 59 s (300 – 6.000 rpm) 1 – 29 min (300 – 4.000 rpm) Display time digital Dimensions (W x D x H) 100 x 160 x 40 mm Weight 0,75 kg Protection class acc. to DIN EN 60529 IP 20

Application areas: Human medicine, pathology, veterinary medicine, animal hygiene institutes, clinical diagnosis research, foodstuffs testing laboratories, diagnostic laboratories, toxicology, medical research, pharmaceutical research, biological research, tumour biology, immunology, chemistry, cosmetics

ULTRA-TURRAX[®] Tube Drive

World first: Universal disposable disperser system with hermetically sealable disposable sample tubes. Protection and security for: Infectious sample materials, toxic substances, high-odour substances.

- Disperse, stir and grind using a single drive unit - No possibility of cross-contamination - Hermetically sealable disposable sample tubes - γ-sterilized tubes on request

- High level of user safety

- Quick and simple

- Hygienic and clean - Covers with pierceable membranes available

- on request - Volumes from 2–15 ml
- - Anti-locking function

 - Chemical-resistant plastic
 - Patent pending

ULTRA-TURRAX[®] Workstation

Included with delivery (page): 1 x ULTRA-TURRAX[®] Tube Drive (63), 2 x ST-20 Tube with stirring device (64), 2 x DT-20 Tube with rotor-stator element (64), 2 x BMT-20 G / S Tube for grinding with glass balls (G) or with stainless steel balls (S) (64), 1 x removal hook for removal the rotor-stator unit, power supply



IKA[®] Crushing Dispersers (batch operation)

- Suitable for individual use and use in series - Increases safety due to low voltage (24 V)





Ident. No. 3646000 100 - 240 V 50/60 Hz





Ident. No. 3645000 100 - 240 V 50/60 Hz

64 Dispersers (batch operation)



ldent. No.		
3599500	ST-20	
3681500	ST-20 y-sterilized	



Ident. No.	
3599400	DT-20
3681600	DT-20 γ-sterilized



Ident. No. 3599600 glass balls 3553700 stainless steel balls 3681700 stain. steel, y-sterilized 3681800 glass, y-sterilized

NEW!

ST-20 Tube with stirring device

Suitable	for:
- Mixing	

- Stirring

- Extractions
- Preparation of soil sample suspensions
- Solubility testing for medications
- Volumes from 2 15 ml
- 25 items per pack

DT-20 Tube with rotor-stator element

Suitable for: - Dispersion

- Homogenisation
- Suspensions
- Pharmocokinetics
- Metabolism studies
- Diagnosis
- Volumes from 5 15 ml
- 25 items per pack

BMT-20 G / S Tube for grinding with glass balls (G) or with stainless steel balls (S)

Suitable for:

- Dry milling of dry and brittle samples (e.g. kaolin, gypsum, coloured pigments, tablets)
- Cell maceration
- Processing of materials mixed with fluids
- Volumes from 2 15 ml
- 25 items per pack

y-sterilized tubes

All sample tubes available in y-sterilized versions.



Technical data	
Motor rating input	125 W
Motor rating output	75 W
Volume range (H ₂ O)	0,5 – 100 ml
Max. viscosity	5.000 mPas
Speed adjustment	stepless
Speed range	8.000 – 30.000 rpm

Speed range	8.000 – 30.000 rpm
Speed stability	< 6 %
Speed display	scale
Noise without	
dispersing element	65 dB (A)
Overload protection	yes
Permitted ON-time (ON / OFF)	max. 10 min /
	min. 5 min

General data

Dimensions (W x D x H) 45 x 60 x 180 mm Weight 5 – 40 °C Permissible ambient temperature Permissible relative humidity Protection class acc. to DIN EN 60529

T 10 basic

Competitively priced dispersing instrument for mic finish.

- the dispersing elements easy
- Immense speed stability with various media
- light weight and ergonomic form
- -
- 0,4 kg be dismounted without tools 80 % -IP 30
 - analysis
 - spare seals and clamp R 200

Accessories (page):

R 200 Clamp (116), R 104 Stand (114), H 44 Boss head clamp (116), Dispersing elements (72): S 10 N – 5 G, S 10 N – 8 G, S 10 N - 10 G, Plastic dispersing elements (74): S 10 D - 7 G - KS - 65, S 10 D - 7 G - KS - 110

IKA[®] Crushing Dispersers (batch operation)

volumes of 0,5 to 100 ml. A wide speed range allows you to work at high circumferential speeds even with small rotor diameters. Perfect ergono-

Quick-release coupling makes changing

thanks to high performance 125 Watt drive - Ideal for manual operation thanks to its - Extremely mobile thanks to direct mains operation (no transformer required) Stainless steel dispersing elements (5 mm, 8 mm and 10 mm diameter) can be cleaned quickly and easily as they can Plastic disposable dispersing elements in two sizes, particularly suitable for PCR

- Included with delivery: empty storage case (for drive, clamp, dispersing elements) and

Ident, No. 3420000 230 V 50/60 Hz 3420001 115 V 50/60 Hz



66 Dispersers (batch operation)



Ident. No. 3565000 230 V 50/60 Hz 3565001 115 V 50/60 Hz

T 18 basic ULTRA-TURRAX®

Competitively priced dispersing instrument for volumes of 1 to 1.500 ml (H_2O). A wide speed range allows you to work at high circumferential speeds.

- Electronic speed control
- Electronic overload protection

- Quick release button for dispersing element - As standard, the T 18 is equipped with a

connection for a revolution counter Dispersing elements not included with delivery.

Accessories (page):

Dispersing instruments (70), Stands (114): R 1825, R 1826, R 1827, R 182 Boss head clamp (116), DZM control.o Revolution counter (119), RH 3 Strap clamp (116)

echnical data	
Notor rating input	500 W
Notor rating output	300 W
olume range (H ₂ O)	1 – 1.500 ml
flax. viscosity	5.000 mPas
peed adjustment	stepless
peed range (under load)	3.500 – 24.000 rpm
peed display	scale
loise without	
ispersing element	73 dB (A)
Overload protection	yes
Diameter / length of extension arm	13 mm / 175 mm
General data	
Dimensions (W x D x H)	65 x 80 x 240 mm
Veight	1,6 kg
ermissible ambient temperature	5 – 40 °C
ermissible relative humidity	80 %
rotection class acc. to DIN EN 60529	IP 20

T 25 digital ULTRA-TURRAX®

Dispersing instrument for quantities up to approx. 2.000 ml, page 66 Ident. No. 3565000

T 18 basic ULTRA-TURRAX[®]

Dispersing instrument for quantities up to approx. 1.500 ml, page 66 Ident. No. 3561000

R 182

Boss head clamp, page 116 Ident. No. 2657700

S 18 N – 19 G

Dispersing element for quantities between 10 – 1.500 ml, page 70 Ident. No. L004640

High-performance dispersing instrument for volumes from 1 - 2.000 ml (H_2O). The spectrum of applications ranges from homogenizing waste water samples to the use in laboratory reactors, to dispersion tasks under vacuum / pressure and sample preparation in medical diagnostics.

- Three types of shaft bearings
- Standard version with digital display

T 25 digital ULTRA-TURRAX®

- Rotor-Stator configurations have thirty years of proven, guaranteed comparability of test results - Wide range of dispersing elements (not included
- with delivery, page 70 / 71

Accessories (page):

Dispersing instruments (70 / 71), Stands (114): R 1825, R 1826, R 1827, R 182 Boss head clamp (116), RH 3 Strap clamp (116)

Technical data	
Notor rating input	500 W
Notor rating output	300 W
/olume range (H ₂ O)	1 – 2.000 ml
Max. viscosity	5.000 mPas
Speed adjustment	stepless
Speed range (under load)	3.400 – 24.000 rpm
Speed display	digital
loise without	
lispersing element	73 dB (A)
Overload protection	yes
Diameter / length of extension arm	13 mm / 175 mm
General data	
Dimensions (W \times D \times H)	65 x 80 x 240 mm
Veight	1,6 kg
Permissible ambient temperature	5 – 40 °C
Permissible relative humidity	80 %
Protection class acc. to DIN EN 60529	IP 20

S 25 N – 18 G Dispersing element for quantities between 10 – 1.500 ml, page 71 Ident. No. 0593400

RH 3

Strap clamp, page 116 Ident. No. 3008600



R 1827 Plate stand, page 114 Ident. No. 3160200



IKA[®] Crushing Dispersers (batch operation) 67

68 Dispersers (batch operation)



Ident. No.	
2953100	230 V 50/60 Hz
2953101	115 V 50/60 Hz

T 50 basic ULTRA-TURRAX®

- High-performance dispersing instrument for volumes from $0,25 - 30 \mid (H_2O)$
- Three types of shaft bearings
- Several rotor-stator configurations - Agitator shaft R 50 allows the use of the T 50 basic as a "high-speed stirrer" (not included in delivery, page 76)
- Infinitely variable speed control, for continuous operation
- · Reproducible operations due to constant speed even with viscosity changes
- Large selection of dispersing elements - Plug-in connectors facilitate exchange
- of dispersing elements
- Electronic safety circuit and smooth start - As standard, the T 50 basic is equipped with a connection for the revolution counter
- Wide range of dispersing elements (not included in delivery, page 72 / 73)

Accessories (page):

Dispersing elements (72 / 73), Special dispersing elements (76), Stands (114 / 115): R 2722, R 2723, R 271 Boss head clamp (116), DZM control.o Revolution counter (119), RH 5 Strap clamp (116)

T 50 basic ULTRA-TURRAX®

Dispersing instrument for quantities up to approx. 30 l, page 68 Ident. No. 2953100

R 271

Boss head clamp, page 116 Ident. No. 2664000

S 50 N – G 45 G

Dispersing element for coarse crushing, page 72 Ident. No. 8003000

RH 5

Strap clamp, page 116 Ident. No. 3159000

R 2723

Telescopic stand, page 115 ldent. No. 1412100

S 50 N – G 45 F

Dispersing element for subsequent fine crushing, page 73 Ident. No. 8003900

Technical data	
Motor rating input	1.100 W
Motor rating output	700 W
Volume range (H ₂ O)	0,25 – 30 l
Max. viscosity	5.000 mPas
Speed adjustment	stepless
Speed range	4.000 – 10.000 rpm
Speed stability	1 %
Speed display	scale
Noise without	
dispersing element	72 dB (A)
Diameter / length of extension arm	16 mm / 220 mm
Overload protection	yes
General data	
Dimensions (W x D x H)	125 x 120 x 367 mm
Weight	6 kg
Permissible ambient temperature	5 – 40 °C
Permissible relative humidity	80 %
Protection class acc. to DIN EN 60529	IP 20



Technical data

Motor rating input	1.800 W
Motor rating output	1.500 W
Volume range (H ₂ O)	2 – 50 l
Max. viscosity	5.000 mPas
Speed, fixed	7.200 rpm
Speed stability	5 %
Noise without	
dispersing element	75 dB (A)
Overload protection	yes
General data	
Dimensions (W x D x H)	190 x 580 x 380 mm
Weight	28 kg
Permissible ambient temperature	5 – 40 °C
Permissible relative humidity	80 %
Protection class acc. to DIN EN 60529	IP 54

T 65 D ULTRA-TURRAX®

guantities from 2 - 50 I (H₂O).

- mPas delivery)
 - of dispersing elements
 - Speed controller on request - Dispersing instruments for the
 - technology catalogs - Cables and plugs not included with delivery

Accessories (page):

Nomenclature dispersing elements

The variety of media to be processed also requires a variety of rotor-stator configurations and seals. In many cases it is neccessary to use subsequently two dispersing elements, for pre-crushing and fine crushing. The plug-in connectors facilitate the exchange of the dispersing elements.



For dispersing	Dispersing element	With seal or	Generator	With outer diameter (mm)	Degree of fineness
instrument	Shaft / Agitator shaft	bearing type*	or element**		achieved***
T 10	S 10	Ν	-	5/8/10	G
T 18	S 18	Ν	-	10/19	G
T 25	S 25	N / KR / KV / NK	-	8 / 10 / 18 / 19 / 25	G / F
T 50	S / R 50	N / KV / KR / KG – HH	G/W	45 / 65 / 80	G / M / F
T 65	S 65	KG – HH	G	65	G / M / F

* N = PTFE bearing, KR = Ball bearing with FKM- seal, KV = Ball bearing with vacuum-tight sliding-ring seal with silicon carbide seal rings, NK = PTFE bearing with additional ball bearing without seal, KG - HH = Ball bearing with sliding-ring seals of hard metal allow with FFPM seal rings ** G = proved configuration, W = special element

***G = coarse, M = medium, F = fine

IKA[®] Crushing Dispersers (batch operation)





IKA[®] Crushing Dispersing elements













Dispersing elements T 18 basic, T 25 digital

Dispersing element		S 18 N – 10 G	S 18 N – 19 G	S 25 N – 8 G	S 25 N – 10 G	S 25 N – 10 G – VS	S 25 N – 18 G	S 25 KR – 18 G	S 25 KV – 18 G
Ident. No.		L004639	L004640	1024200	0594000	1899000	0593400	0560300	2348000
Fig.		without fig.	without fig.	1	2	without fig.	3	without fig.	without fig.
Suitable for dispersing instrument		T 18 basic	T 18 basic	T 25 digital	T 25 digital	T 25 digital	T 25 digital	T 25 digital	T 25 digital
Working range		1 – 100 ml	10 – 1.500 ml	1 – 50 ml	1 – 100 ml	1 – 100 ml	10 – 1.500 ml	10 – 1.500 ml	10 – 1.500 ml
Stator diameter		10 mm	19 mm	8 mm	10 mm	10 mm	18 mm	18 mm	18 mm
Rotor diameter		7,5 mm	12,7 mm	6,1 mm	7,5 mm	7,5 mm	12,7 mm	12,7 mm	12,7 mm
Gap between rotor and stator		0,35 mm	0,4 mm	0,25 mm	0,35 mm	0,35 mm	0,3 mm	0,3 mm	0,3 mm
Circumferential speed		9,4 m/s	15,9 m/s	7,7 m/s	9,4 m/s	9,4 m/s	15,9 m/s	15,9 m/s	15,9 m/s
Min. / max. immersion depth		25 / 70 mm	35 / 170 mm	27 / 85 mm	22 / 85 mm	22 / 85 mm	40 / 165 mm	40 / 185 mm	40 / 225 mm
Shaft length		108 mm	204 mm	108 mm	105 mm	105 mm	194 mm	194 mm	270 mm
Materials in contact with medium		PTFE, AISI 316L	PTFE, AISI 316L	PTFE, AISI 316L	PTFE, AISI 316L	PTFE, AISI 316L	PTFE, AISI 316L	FKM, AISI 316L	FFPM / SIC, AISI 316L
pH range		2 – 13	2 – 13	2 – 13	2 – 13	2 – 13	2 – 13	2 – 13	2 – 13
Suitable for solvents		yes	yes	yes	yes	yes	yes	no	yes
Suitable for abrasive substances		yes	yes	yes	yes	yes	yes	no	no
Max. temperature		180 °C	180 °C	180 °C	180 °C	180 °C	180 °C	80 °C	220 °C
Sterilization methods		all methods	all methods	all methods	all methods	all methods	all methods	wet chemical	wet chemical
Min. vacuum		-	-	-	_	_	_	50 mbar	1 mbar
Max. pressure		-	-	-	-	_	_	-	6 bar
Ultimate fineness, suspensions		10 – 50 μm	10 – 50 μm	10 – 50 μm	10 – 50 μm	10 – 50 μm	10 – 50 μm	10 – 50 μm	10 – 50 μm
Ultimate fineness, emulsions		1 – 10 µm	1 – 10 µm	1 – 10 µm	1 – 10 µm	1 – 10 µm	1 – 10 µm	1 – 10 µm	1 – 10 µm
Dispersing element	S 25 NK – 19 G	S 25 N – 25 G	S 25 KR – 25 G	S 25 KV – 25 G	S 25 N – 25 F	S 25 KR – 25 F	S 25 KV – 25 F	S 25 KV – 25 G – IL	S 25 KV – 25 F – IL
Ident. No.	2494700	1713300	1713400	2466900	1713800	1713900	2404000	2563000	2830200
Fig.	without fig.	4	without fig.	without fig.	5	without fig.	without fig.	without fig.	without fig.
Suitable for dispersing instrument	T 25 digital	T 25 digital	T 25 digital	T 25 digital	T 25 digital	T 25 digital	T 25 digital	T 25 digital	T 25 digital
Working range	25 – 1.500 ml	50 – 2.000 ml	50 – 2.000 ml	50 – 2.000 ml	100 – 2.000 ml	100 – 2.000 ml	100 – 2.000 ml	Inline	Inline
Stator diameter	19 mm	25 mm	25 mm	25 mm	25 mm	25 mm	25 mm	25 mm	25 mm
Rotor diameter	12,7 mm	17 mm	17 mm	17 mm	18 mm	18 mm	18 mm	17 mm	18 mm
Gap between rotor and stator	0,3 mm	0,5 mm	0,5 mm	0,5 mm	0,5 mm	0,5 mm	0,5 mm	0,5 mm	0,5 mm
Circumferential speed	15,9 m/s	21,4 m/s	21,4 m/s	21,4 m/s	22,6 m/s	22,6 m/s	22,6 m/s	21,4 m/s	22,6 m/s
Min. / max. immersion depth	40 / 165 mm	40 / 165 mm	40 / 185 mm	40 / 225 mm	40 / 165 mm	40 / 185 mm	40 / 225 mm	40 / 85 mm	40 / 85 mm
Shaft length	194 mm	194 mm	194 mm	270 mm	194 mm	194 mm	270 mm	110 mm	110 mm
Materials in contact with medium	PTFE, AISI 316L	PTFE, AISI 316L	FKM, AISI 316L	FFPM / SIC, AISI 316L	PTFE, AISI 316L	FKM, AISI 316L	FFPM / SIC, AISI 316L	FFPM / SIC, AISI 316L	FFPM / SIC, AISI 316L
pH range	2 – 13	2 – 13	2 – 13	2 – 13	2 – 13	2 – 13	2 – 13	2 – 13	2 – 13
Suitable for solvents	yes	yes	no	yes	yes	no	yes	yes	yes
Suitable for abrasive substances	yes	yes	no	no	yes	no	no	no	no
Max. temperature	120 °C	180 °C	80 °C	220 °C	180 °C	80 °C	220 °C	220 °C	220 °C
Sterilization methods	wet chemical	all methods	wet chemical	wet chemical	all methods	wet chemical	wet chemical	wet chemical	wet chemical
Min. vacuum	-	-	50 mbar	1 mbar	-	50 mbar	1 mbar	1 mbar	1 mbar
Max. pressure	-	-	-	6 bar	-	-	6 bar	6 bar	6 bar
Ultimate fineness, suspensions	10 – 50 μm	15 – 50 μm	15 – 50 μm	15 – 50 μm	5 – 25 µm	5 – 25 µm	5 – 25 µm	15 – 50 μm	5 – 25 µm
Ultimate fineness, emulsions	1 – 10 um	1 – 5 um	1 – 5 um	1 – 5 um	1 – 10 µm	1 – 5 µm			

Stator diameter	19 mm	25 mm	25 mm	25 mm	25 mm	25 mm	
Rotor diameter	12,7 mm	17 mm	17 mm	17 mm	18 mm	18 mm	
Gap between rotor and stator	0,3 mm	0,5 mm	0,5 mm	0,5 mm	0,5 mm	0,5 mm	
Circumferential speed	15,9 m/s	21,4 m/s	21,4 m/s	21,4 m/s	22,6 m/s	22,6 m/s	
Min. / max. immersion depth	40 / 165 mm	40 / 165 mm	40 / 185 mm	40 / 225 mm	40 / 165 mm	40 / 185 mm	
Shaft length	194 mm	194 mm	194 mm	270 mm	194 mm	194 mm	
Materials in contact with medium	PTFE, AISI 316L	PTFE, AISI 316L	FKM, AISI 316L	FFPM / SIC, AISI 316L	PTFE, AISI 316L	FKM, AISI 316L	
pH range	2 – 13	2 – 13	2 – 13	2 – 13	2 – 13	2 – 13	
Suitable for solvents	yes	yes	no	yes	yes	no	
Suitable for abrasive substances	yes	yes	no	no	yes	no	
Max. temperature	120 °C	180 °C	80 °C	220 °C	180 °C	80 °C	
Sterilization methods	wet chemical	all methods	wet chemical	wet chemical	all methods	wet chemical	
Min. vacuum	_	-	50 mbar	1 mbar	-	50 mbar	
Max. pressure	_	-	_	6 bar	-	-	
Ultimate fineness, suspensions	10 – 50 µm	15 – 50 μm	15 – 50 μm	15 – 50 μm	5 – 25 µm	5 – 25 µm	
Illtimate finances, emulsions	1 10 μm	1 10.000	1 10 um	1 10 um	1 5	1 5 4 50	_

SW 18 Slab rotor

Additional rotor for dispersing elements: S 25 N – 18 G S 25 KR – 18 G S 25 KV – 18 G

12,8 mm
0,35 mm
16,1 m/s
stainl. steel AISI 316L
viscous, fibrous tissue

Ident. No.	
8011900	

IKA[®] Crushing Dispersing elements

For nomenclature see page 69
72 Dispersing elements









Dispersing elements T 50 basic

For nomenclature see page 69

Dispersing element	S 10 N – 5 G	S 10 N – 8 G	S 10 N – 10 G
Ident. No.	3304000	3305500	3370100
Fig.	1	2	3
Suitable for dispersing instrument	T 10 basic	T 10 basic	T 10 basic
Working range	0,5 – 10 ml	1 – 50 ml	1 – 100 ml
Stator diameter	5 mm	8 mm	10 mm
Rotor diameter	3,8 mm	6,1 mm	7,6 mm
Gap between rotor and stator	0,1 mm	0,25 mm	0,2 mm
Min. / max. immersion depth	20 / 75 mm	20 / 95 mm	20 / 100 mm
Shaft length	92 mm	115 mm	115 mm
Materials in contact with medium	PTFE, AISI 316L	PTFE, AISI 316L	PTFE, AISI 316L
pH range	2 – 13	2 – 13	2 – 13
Suitable for solvents	yes	yes	yes
Suitable for abrasive substances	yes	yes	yes
Max. temperature	180 °C	180 °C	180 °C
Sterilization methods	all methods	all methods	all methods
Min. vacuum	_	-	-
Max. pressure	-	-	-
Ultimate fineness, suspensions	5 – 25 µm	5 – 25 µm	5 – 25 µm
Ultimate fineness, emulsions	1 – 10 um	1 – 10 um	1 – 10 µm

Dispersing elements T 65 D

Dispersing element	S 65 KG – HH – G 65 G	S 65 KG – HH – G 65 M	S 65 KG – HH – G 65 F
Ident. No.	8005500	8005700	8005900
Fig.	1	2	3
Suitable for dispersing instrument	T 65 D	T 65 D	T 65 D
Working range	2 – 50 l	2 – 40 l	2 – 30 l
Stator diameter	65 mm	65 mm	65 mm
Rotor diameter	58 mm	58 mm	58 mm
Circumferential speed	21,9 m/s	21,9 m/s	21,9 m/s
Min. / max. immersion depth	90 / 450 mm	80 / 450 mm	80 / 450 mm
Shaft length	520 mm	510 mm	500 mm
Materials in contact with medium	FFPM / SIC, AISI 316L	FFPM / SIC, AISI 316L	FFPM / SIC, AISI 316L
pH range	2 – 13	2 – 13	2 – 13
Suitable for solvents	yes	yes	yes
Suitable for abrasive substances	no	no	no
Max. temperature	180 °C	180 °C	180 °C
Sterilization methods	wet chemical	wet chemical	wet chemical
Min. vacuum	1 mbar	1 mbar	1 mbar
Max. pressure	6 bar	6 bar	6 bar
Ultimate fineness, suspensions	25 – 75 µm	20 – 50 μm	5 – 20 μm
Ultimate fineness, emulsions	5 – 25 µm	5 – 15 µm	1 – 10 µm









Dispersing element	S 50 N – G 45 G	S 50 KR – G 45 G	S 50 N – G 45 M	S 50 KR – G 45 M	S 50 N – G 45 F	S 50 KR – G 45 F	S 50 KV – G 45 G – IL
Ident. No.	8003000	8003100	8003300	8003400	8003900	8004000	8015800
Fig.	1	without fig.	2	without fig.	3	without fig.	without fig.
Suitable for dispersing instrument	T 50 basic	T 50 basic	T 50 basic	T 50 basic	T 50 basic	T 50 basic	T 50 basic
Working range	0,5 – 20 l	0,5 – 20	0,5 – 15	0,5 – 15 l	0,25 – 10	0,25 – 10	Inline
Stator diameter	45 mm	45 mm	45 mm	45 mm	45 mm	45 mm	45 mm
Rotor diameter	36 mm	36 mm	40,5 mm	40,5 mm	40 mm	40 mm	36 mm
Circumferential speed	18,8 m/s	18,8 m/s	21,2 m/s	21,2 m/s	20,9 m/s	20,9 m/s	18,8 m/s
Min. / max. immersion depth	70 / 250 mm	70 / 260 mm	70 / 250 mm	70 / 260 mm	70 / 250 mm	70 / 260 mm	70 mm
Shaft length	300 mm	300 mm	290 mm	290 mm	290 mm	290 mm	105 mm
Materials in contact with medium	PTFE, AISI 316L	FKM, AISI 316L	PTFE, AISI 316L	FKM, AISI 316L	PTFE, AISI 316L	FKM, AISI 316L	FFPM / SIC, AISI 316L
pH range	2 – 13	2 – 13	2 – 13	2 – 13	2 – 13	2 – 13	2 – 13
Suitable for solvents	yes	no	yes	no	yes	no	yes
Suitable for abrasive substances	yes	no	yes	no	yes	no	no
Max. temperature	180 °C	80 °C	180 °C	80 °C	180 °C	80 °C	220 °C
Sterilization methods	all methods	wet chemical	all methods	wet chemical	all methods	wet chemical	wet chemical
Min. vacuum	-	100 mbar	-	100 mbar	-	100 mbar	1 mbar
Max. pressure	-	-	-	-	-	-	6 bar
Ultimate fineness, suspensions	40 – 100 μm	40 – 100 µm	25 – 50 μm	25 – 50 μm	10 – 30 µm	10 – 30 µm	40 – 100 μm
Ultimate fineness, emulsions	10 – 30 µm	10 – 30 µm	5 – 20 µm	5 – 20 µm	1 – 10 µm	1 – 10 µm	10 – 30 µm

S 50 N - Special length shafts also available in 430 mm (order label S 50 N 1)

IKA[®] Crushing

Dispersing elements



For nomenclature see page 69

For nomenclature see page 69

74 Dispersing elements

Plastic dispersing elements

Plastic dispersing elements are ideal for those applications where absolutely no cross-contamination is permitted. They are disposable and can be thrown away after a single use. The element is disposable and designed for one-way use. However, it can be re-used several times in applications where this is permitted. If you decide to re-use the element, make sure that you follow the cleaning instructions carefully. Example use: homogenizing tissue samples.

For disperser	Dispersing element	Seals	Diameter stator	Degree of fineness	Material
	shaft		(mm)	achieved	
T 10	S 10	D = without seal	7	G = coarse	KS = plastic
T 18	S 18	D = without seal	10 / 14	G = coarse	KS = plastic
T 25	S 25	D = without seal	10 / 14	G = coarse	KS = plastic



S 10 D – 7G – KS – 65	
Ident. No.	
3433212	12 pcs.
3433225	25 pcs.



S 10 D – 7G – KS – 110 ldent. No. 3433312 12 pcs. 3433325 25 pcs.

Plastic dispersing elements for T 10 basic

Dispersing element	S 10 D – 7 G – KS – 65	S 10 D – 7 G – KS – 110
Ident. No. [Packing unit]	3433212 [12 pcs.]	3433312 [12 pcs.]
	3433225 [25 pcs.]	3433325 [25 pcs.]
Suitable for dispersing instrument	T 10 basic	T 10 basic
Working range	1 – 20 ml	1 – 40 ml
Stator diameter	7 mm	7 mm
Rotor diameter	4,8 mm	4,8 mm
Min. / max. immersion depth	20 / 50 mm	20 / 90 mm
Shaft length	65 mm	110 mm
Materials in contact with medium	Polycarbonate (PC)	Polycarbonate (PC)
	Polysulfon (PSU)	Polysulfon (PSU)
Max. temperature	100 °C	100 °C
Sterilization methods	yes, autoclavable	yes, autoclavable

Plastic materials used approved by FDA.

Plastic dispersing elements for T 18 basic

Dispersing element	S 18 D – 10 G – KS
Ident. No. [Packing unit]	3452000 [5 pcs.*]
	3452400 [10 pcs.*]
Suitable for dispersing instrument	T 18 basic
Working range	10 – 100 ml
Stator diameter	10 mm
Rotor diameter	6,75 mm
Min. / max. immersion depth	15 / 85 mm
Shaft length	150 mm
Materials in contact with medium	Polycarbonate (PC)
	Polyetheretherketon (PEEK)
Max. temperature	100 °C
Sterilization methods	yes, autoclavable

Plastic materials used approved by FDA.

* incl. 1 Disposable tube

Plastic dispersing elements for T 25 digital

Dispersing element	S 25 D – 10 G – KS
Ident. No. [Packing unit]	3451800 [5 pcs.*]
	3452200 [10 pcs.*]
Suitable for dispersing instrument	T 25 digital
Working range	10 – 100 ml
Stator diameter	10 mm
Rotor diameter	6,75 mm
Min. / max. immersion depth	15 / 85 mm
Shaft length	150 mm
Materials in contact with medium	Polycarbonate (PC)
	Polyetheretherketon (PEEK)
Max. temperature	100 °C
Sterilization methods	yes, autoclavable

PP

Plastic materials used approved by FDA. * incl. 1 Disposable tube

General data Material

Disposable tube S 18 / 25-ET50

50 ml for attaching onto plastic tools from S 18 D and S 25 D series. Allows dispersing in a closed system (splash guard).

IKA[®] Crushing

Dispersing elements

S 18 D – 10 G – KS		
Ident. No.		
3452000	5 pcs.*	
3452400	10 pcs.*	

S 18 D – 14 G – KS
3451900 [5 pcs.*]
3452300 [10 pcs.*]
T 18 basic
10 – 500 ml
14 mm
9,5 mm
15 / 85 mm
150 mm
Polycarbonate (PC)
Polyetheretherketon (PEEK)
100 °C
yes, autoclavable

S 18 D – 14 G – KS	
Ident. No.	
3451900	5 pcs.*
3452300	10 pcs.*

100 - 100 MTT

				-
			-	
			1000	
		1000	21	
1	No. of Concession, Name	-		
1.2.21				
100				

S 25 D – 10 G – KS	
Ident. No.	
3451800	5 pcs.
3452200	10 pcs.

S 25 D – 14 G – KS
3451700 [5 pcs.*]
3452100 [10 pcs.*]
T 25 digital
10 – 500 ml
14 mm
9,5 mm
15 / 85 mm
150 mm
Polycarbonate (PC)
Polyetheretherketon (PEEK)
100 °C
yes, autoclavable



S 25 D – 14 G – KS	
Ident. No.	
3451700	5 pcs.
3452100	10 pcs.



Ident. No. 3452500

IKA[®] Crushing 76 Special dispersing elements (T 50 basic)



ldent. No 1689300

Ident No 1289800





R 1402 Dissolver



S 50 ... - W 80 SMK

R 50 "high speed" stirring shaft

With the stirring shaft R 50, the T 50 basic is

guickly converted into a high speed stirrer. 700 W

and 10.000 rpm are provided for rapid mixing, dis-

solving, and disagglomerating pigment agglome-

rates. The conical shaft is supported by means of

ball bearings, the mixing elements have a screw

connection. For operational safety a protective

cage is fitted around the mixing element.

Dispersing elements (76): R 1405, R 1402

Included with delivery (page):

R 1402 Dissolver (76)

Accessories (page):

R 1405 Propeller

For shortening mixing a vertical flow and the high circumferential speed up to 10.000 rpm ensure intensive mixing. The head is used for adding gases or liquids, for lumpfree suspension of difficult to dissolve powders or for dissolving sedimented, already hardened substances.

ldent. No.	
8006300	S 50 N – W 80 SMK
8006400	S 50 KR – W 80 SMK



S 50 N – W 65 SK Cutting head

To crush large pieces (up to 50 mm) of fibrous materials, such as vegetation, vegetables and fruit.

General data	
Min. / max. immersion depth	80 / 350 mm
Working range	1 – 10
Generator diameter	65 mm
Available seals	S 50 N

Technical data	
Flow rate (H ₂ O)	11,6 l/min
Speed range	6.500 – 24.000 rpm
Materials in contact with medium	stainl. steel (AISI 316L)
	FFPM
Max. operating temperature	180 °C
Dimensions (W x D x H)	450 x 100 x 120 mm
Weight	3.8 ka

Chamber volume

Min. vacuum

Technical data

Flow rate (H₂O)

Materials in contact with medium

Permissible ambient temperature

Protection class acc. to DIN EN 60529

Permissible relative humidity

Max. operating temperature

Dimensions (W x D x H)

Chamber volume

Min. vacuum

Max. pressure

Speed range

Weight

00 – 24.000 rpm	laboratory.
teel (AISI 316L)	 Simple, compact ar
FFPM	- Sterilizable, autocla
180 °C	 Table-top or stand-stand-stand-stand-stand-standard
(100 x 120 mm	space requirement

- Easy disassembly 26 ml 1 mbar 6 bar rmissible ambient temperature
 - 5 40 °C 80 °C IP 20
 - emulsyfying and desagglomeration - For vacuum or pressurized operation (up to 6 bar) - If the DK 25.11 is used, air induction is also prevented in batch operation
 - Not self-priming - A pump can be integrated between intake nozzle and vessel. As a result, viscous fluids
 - can be processed
 - cyclical continuous operation

Included with delivery (page): T 25 digital (66), AD 25 Mounting (78), DK 25.11 Flow chamber (78), S 25 KV - 25 G - IL Dispersing element (71)

Accessories (page): dispersing element S 25 KV – 25 F – IL (71)

UTL 50 basic Inline ULTRA-TURRAX®

For circulation or flow-through processings in the laboratory or pilot plant stations. - Stand-supported device, low space requirement

24 l/min

FFPM

180 °C

1 mbar

6 bar

80 °C

IP 21

5 – 40 °C

4.000 – 10.000 rpm

stainl. steel (AISI 316L)

- 130 x 150 x 500 mm 6,1 kg
 - prevented in batch operation
 - 94 ml

Additional features as UTL 25 digital inline.

Included with delivery (page): T 50 basic (68), DK 50.11 Flow chamber (78), S 50 KV – G 45 G – IL Dispersing element (73)

Accessories (page): R 2723 Telescopic stand (115), R 271 Boss head clamp (116)

	General data
	Working range
	Rotor diameter
C Jet mixer head	
	General data
and dissolving times. The	Min. / max. imm

General data	
Min. / max. immersion depth	140 / 350 mm
Working range	1 – 50 I
Generator diameter	80 mm
Available seals	S 50 N
	S 50 KR

General data	
Working range	0,25 - 10
Rotor diameter	45 mm

ieneral data	
Vorking range	1 – 30
otor diameter	42 mr

General data	
Working range	
Rotor diameter	

General data

Working range

Materia

Immersion depth

Max. circumferential speed

Max. permissible rotor diameter

General data	
Working range	1 – 3
Potor diamotor	42 n

General data Working range Rotor diameter		
General data Working range Rotor diameter		
Working range Rotor diameter	General data	
Rotor diameter	Working range	
	Rotor diameter	

-		

General data	
Working range	1 – 30
Rotor diameter	42 m

Max. pressure
Permissible ambient temperature
Permissible relative humidity
Protection class acc. to DIN EN 60529

180 mm

50 mm

0,25 - 30 |

15,7 – 23 m/s

stainl. steel (AISI 316L)

IKA[®] Crushing Dispersers (inline operation)

UTL 25 digital Inline ULTRA-TURRAX®

For circulation or flow-through processings in the

- act and sturdy modular design utoclave-compatible tand-supported device, low
- Large delivery capacity of 4,4 to 11,6 l/min with open outlet (the mounting of a valve can reduce the flow rate)
- For air-free, sterile, and inline suspension,

- Not suitable for continuous operation or

- Large flow rate of 24 l/min with open outlet (the mounting of a valve reduces the delivery capacity) - For vacuum or pressurized operation to 6 bar - If the DK 50.11 is used, air induction is also

- Not suitable for continuous operation or cyclical continuous operation



Ident. No. 8014400 230 V 50/60 Hz 8014401 115 V 50/60 Hz



Example application



230 V 50/60 Hz 115 V 50/60 Hz

8015900

8015901

78 Accessories Dispersers



Ident. No. 2562500



DK 25.11 Flow chamber

For S 25 KV - 25 ... - IL dispersing elements. Allows inline operation mode, see UTL 25 digital, page 77.

Batch operation (see fig.):

DK 25.11 is mounted around the dispersing element. The DK 25.11 must be at a lower elevation than the surface of the liquid during operation. With this operating mode, no air is drawn in as a result of turbulence in the vessel.

AD 25

Mounting sup

General data	
Chamber volume	94 ml
Vacuum	1 mbar
Pressure	6 bar

Technical data	
Power	1,5 kW
Speed range (up to 40 m/s)	3.160 – 13.750 rpm
Flow rate	approx. 300 – 700 l/h
Circumferential speed	23,5 (9,4 - 41) m/s
Materials in contact with medium	stainl. steel AISI 316L/316T
Seal material	
Shaft sealing ring	PTFE-compound
Elastomer chamber	standard FPM
Dimensions (W x D x H)	450 x 250 x 350 mm
Generator	4M

LABOR-PILOT 2000/4

Small, multi-functional table-top dispersing system with the possibility of upscaling to production size.

- Belt drive

- Generator 4M (medium) - PTFE-shaft sealing ring Can be extended to a high-speed dispersing machine by means of the LABOR-PILOT-Controller.

Special features:

- capable

PROCESS-PILOT 2000/4

The PROCESS-PILOT is equipped with a double-acting mechanical seal complete with the necessary pressure locking system. This allows, in addition to the other modules, the use of the CMS module for easy and dust-free incorporation of powders into liquids in batch operation.

Q1Q1VMG-BQ1VMG standard FPM

- Operates under pressure/vacuum
- Works at elevated temperatures
- 450 x 250 x 930 mm - Low maintenance mechanical seal 4M

- Suitable for dry-running DN 25
- (ATEX 95) DN 15

- Three-phase motor
- Belt drive
- Single stage dispersing chamber UTL
- combination Q1Q1VMG-BQ1VMG
- Generator 4M (medium)

- O-rings made of FPM

neral data	
amber volume	94 ml
cuum	1 mbar
	0.1

For S 50 KV - G 45 ... - IL dispersing elements.

Allows operation in inline mode, see UTL 50

If used in batch mode: DK 50.11 is mounted

around the dispersing element. Additional fea-

DK 50.11 Flow chamber

basic, page 77.

tures as DK 25.11.

General data	
Chamber volume	94 ml
/acuum	1 mbar
Pressure	6 bar

port for DK 25.11		

General data

Vacuum

Pressure

Chamber volume

	Technical data	
	Power	2,2 kW
	Speed range (up to	40 m/s) 3.160 – 13.750 rpm
	Standard speed	7.900 rpm
	Flow rate	approx. 300 – 700 l/h
94 ml	Circumferential spe	ed 23,5 (9,4 - 41) m/s
1 mbar	Materials in contact	with medium stainl. steel AISI 316L/316Ti

Seal material

Generator

Inlet connector

Outlet connector

Elastomer chamber

Dimensions (W x D x H)

(including pressure locking)

Double acting mechanical seal

26 ml

1 mbar

6 bar

IKA[®] Crushing Dispersers (LABOR-PILOT / PROCESS-PILOT)

Basic unit with module ULTRA-TURRAX® UTL: - Three-phase motor with on/off-switch

- Single stage dispersing chamber UTL

- Multiple use for mixing, dispersing or wet milling - Modular design - the basic unit can be converted easily by using various mixing heads to reach circumferential speeds up to 40 m/s - CIP/SIP (clean-in-place / sterilization-in-place)

- All parts in contact with the medium are made of stainless steel AISI 316 L/316 Ti - The control system allows a product-specific adaptation to rheological characteristics



Ident, No T056762 3 x 220 - 240 V 50/60 Hz T055396 3 x 380 - 420 V 50/60 Hz S097950 for FC-operation

The advantages of the PROCESS-PILOT:

- Available in Ex-protected design acc. to 94/9/EG

Basic unit with module ULTRA-TURRAX® UTL:

- Double acting mechanical seal, material - All parts in contact with the medium are made of stainless steel AISI 316 L/316 Ti



Ident, No. T058102

for FC-operation

Other generators available.

IKA[®] Crushing Modules (LABOR-PILOT / PROCESS-PILOT) 80











Module colloid mill MK	
Ident. No.	
T054917	LABOR-PILOT
T058583	PROCESS-PILOT
Module cone mill MKO	
Ident. No.	
T061069	LABOR-PILOT
T061674	PROCESS-PILOT

LABOR-PILOT-Controller PROCESS-PILOT-Controller

The LABOR-PILOT and the PROCESS-PILOT are controlled via the Controller, complete with:

- Frequency converter
- Speed indicator
- Temperature indicator - Timer
- Adapter DN 15 with temperature sensor PT 100 at the product outlet

The controller enables flexible setting and measurement of processing parameters.

General data	
Power	
LABOR-PILOT-Controller	1,5 kW
PROCESS-PILOT-Controller	4 kW
Frequency range	20 – 87 Hz
Dimensions (W \times D \times H)	
LABOR-PILOT-Controller	220 x 340 x 380 mm
PROCESS-PILOT-Controller	220 x 340 x 450 mm

LABOR-PILOT-Controller

Ident. No. T054916 3 x 220 - 240 V 50/60 Hz T055171 3 x 380 - 420 V 50/60 Hz PROCESS-PIL OT-Controller Ident, No. T058761 3 x 380 - 420 V 50/60 Hz

Module DISPAX-REACTOR® DR

Extension for a three-stage dispersing machine with two additional generators, 2G (coarse) and 6F (fine).

Can be extended to the high-speed dispersing machine DRS by means of the LABOR-PILOT Controller.

Material module: stainl. steel (AISI 316L/316Ti) / FPM. Other generators available.

Module colloid mill MK / Module cone mill MKO

Extension for getting a colloid mill/cone mill:

Tools:

- Milling tool for MK: spiral geared, stainless steel (AISI 316L/316Ti)
- Milling tool for MKO: WCCO-coated (tungsten carbide-cobalt)

Material module stainl. steel (AISI 316L/316Ti) / FPM

General data	
unction mode	3-step inline
Generator	Standard 2G – 4M – 6F
c	other configurations available as optional
low rate 2G – 4M – 6F	up to 350 l/h
peed range (20 – 87 Hz) 3.160 – 13.750 rpm
itandard speed (50 Hz)	7.900 rpm
nlet connector	DN 25
Jutlet connector	DN 15

General data	
Function mode	two counteropposing cones
	for setting milling grade
Flow rate	approx. 30 – 300 l/h
Speed range	3.160 – 13.750 rpm
Standard speed	7.900 rpm
Inlet connector	DN 25
Outlet connector	DN 15

General data Function mode continuous "solid/liquid" incorporation Flow rate approx. 50 – 200 l/h Speed range 3.160 – 13.750 rpm Standard speed 7.900 rpm Inlet connector (solid) Inlet connector (liquid) Outlet connector DN 15

Generator

General data

Funnel connector

3-Way stop cock

Ball valve connector

Product temperature

Process temperature

Circulating operation

Funnel

Two-way ball valve

RW 28 basic, page 38

PTFE-wipers for anchor stirrer

Sealing set for circulation bend

R 271 Boss head clamp, page 116

Circulation bend with intermediate piece

Telescopic stand

Anchor stirrer

3-Way stop cock

Blind nut DN 15

Hose nozzle DN 15 K

Hose nozzle K / M DN 25

Hose nozzle DN 15 K

Ident. No.

T056458

T056468

T055330

2760000

T058564

T058714

2664000

T056482

T057796

T056730

T056484

T056486

T056486

Continuous operation Ident. No. T057251

Capacity

Module MHD

Extension for continuous incorporation and dispersion of powdery materials in liquids.

Tools:

- Injection unit DN 50

DN 15

DN 15

up to 160 °C

up to 3 bar

2P

- Feeding screw - Mixing blades

Material module: stainl. steel (AISI 316L/316Ti) / FPM

General data	
Function mode	"solid/liquid" incorporation
	in batch mode
Flow rate	up to 6.000 l/h
Standard speed	9.000 rpm
Inlet connector (solid)	DN 25
Inlet connector (liquid)	DN 25
Outlet connector	DN 25

Module CMS (for PROCESS-PILOT only!)

n	Extension for incor
le	dery materials into
′h	- Housing parts for

- Rotors - Adapter for housing and shaft

Material module: Stainless steel (AISI 316L 316Ti) / FPM

Funnel 8 I, double-jacketed, stainl. steel. 81 DN 25 - Two-way ball valve

- DN 25
 - Telescopic stand
 - RW 28 basic (page 38)
 - Anchor stirrer
 - PTFE-wipers for anchor stirrer

 - 3-Way stop cock

 - Blind nut DN 15

 - Hose nozzle DN 15

Accessories for continuous operation

Inlet side: Hose nozzle K / M DN 25 Outlet side: Hose nozzle DN 15 K

For more inform
www.ikaprocess

IKA[®] Crushing Modules (LABOR-PILOT / PROCESS-PILOT)





- Generator 2P (coarse with pump blades)

Ident, No T055142 LABOR-PILOT PROCESS-PILOT T058148

rporation and dispersion of powliquids in batch operation with: powder and liquid inlet





- R 271 Boss head clamp (page 116)

- Circulation bend with intermediate piece - Sealing set for circulation bend



nation go to: .com



One basic module – Seven configurations



Technical data Motor rating 900 W 3.000 rpm – 26.000 rpm Speed range (40 m/s) Power supply 220 – 240 V 50 – 60 Hz Frequency Process pressure up to 2,5 bar Product temperature in continuous operation up to 80 °C in short applications (up to 18/min) up to 120 °C Materials in contact with stainl. steel (AISI 316L and AISI 316Ti)

medium

Seal material

Shaft sealing ring

Elastomer chamber

PTFE-compound Standard FPM optional EPDM (FDA), FFPM Dimensions (W x D x H) 170 x 270 x 215 mm Weight, basic module 7 kg

magicLAB[®]

equipment.

- The magicLAB®

- processing
- series production parameters
- Modular construction
- Easy to upgrade, flexible upgrading

 - rential speeds up to 40 m/s
- panel (included as standard)

The magicLAB® basic module is fitted with the ULTRA-TURRAX[®] UTL (4M generator) and is designed for one-step, inline dispersion.

Included with delivery (basic module):

and carrying case.

Optional extension modules and related applications

The following extension modules are available as optionals for the basic module; for a description see pages 84 and 85.

Optional extension modules	Applications
DISPAX-REACTOR® DR	3-step dispersion
MK colloid mill	Wet and fine milling
MKO cone mill	
MHD	Continous powder incorporation
CMS	Intermittent powder incorporation
ULTRA-TURRAX® UTC	Single-step dispersion in batch mode
MICRO-PLANT	Dispersion in circulating mode – accessory for UTL

Module CMS

magicLAB[®] with module UTC

IKA[®] Crushing magicLAB®

Seamless transfer from laboratory to production

- Suitable for a wide range of applications - Designed for incorporation, dispersion, wet milling and incorporating powders into fluids in laboratory quantities - Ideal for continuous, circulating and batch

- Enables same working methods from prescription development right through to

- Makes scale-up easy thanks to constant work

- Simple heating and cooling of all modules

- Rotor-stator construction operating at circumfe-

- Flexible speed setting and temperature control from the magicLAB® controller and the operator

Basic module with inline dispersion element ULTRA-TURRAX[®] UTL with 4M generator, operator panel, controller, temperature probe, storage



magicLAB [®] with module UTL	
Ident. No.	
U078310	magicLAB®
U077727	230 V power supply cable (EU)
U077729	115 V power supply cable (USA)



Multifunction storage and carrying case



, DR, MK and MKO modules

magicLAB® 84





- Single-step, inline dispersion
- In cases where standard stirring and mixing are not enough (e.g. with instable material) or to save time
- Dispersing materials of different viscosity
- To speed up dissolving processes



Applications: Extension to modules with 2G and 6F generators

Module DISPAX-REACTOR[®] DR

- Three-step, inline dispersion
- When dispersion with ULTRA-TURRAX® is not sufficient
- For narrow dispersion spectra after just one pass
- For optimu of emulsio

Module N

Colloid mill with spiral-geared milling tool / Cone mill module with WCCO-coated milling tool. With two counteropposing cones with adjustable milling grade and friction. Applications:

- For wet and fine milling of hard and grainy material
- Making fine emulsions and pastes

0	0
	9

Module MHD

Applications:

- For continuous incorporation of powders in liquids, homogenising and dispersion in a single pass
- For continuous preparation of suspensions
- For rapid wetting of solids without lumps
- For solid fractions up to 80 % and for viscosities up to 50 Pas (depending on process)

General data	
Function mode	continuous "solid/liquid"
	incorporation
Generator	2P
Flow rate	approx. 60 l/h
Speed range	3.000 – 26.000 rpm
Standard speed	16.000 rpm

Generator	28
Flow rate	approx. 60 l/h
Speed range	3.000 – 26.000 rpm
Standard speed	16.000 rpm
Inlet connector (solid)	inner diam. 25 mm / outer diam. 35 mm
Inlet connector (liquid)	outer calibrated
	6 mm hose
Outlet connector	1/2" clamp

General data	
Function mode	"solid/liquid" incorpora
	in batch m
Flow rate	approx. 1.00
Speed range	3.000 - 16.000
Standard speed	11.000
Inlet connector (solid)	3⁄4" cl
Inlet connector (liquid)	3⁄4" cl
Outlet connector	3⁄4" cl

General data

Hopper connector

Product temperature

General data Function mode

Generator

Volume range

Speed range

Duty cycle

Three-way ball valve connector

Process pressure (optional)

Capacity

Module CMS

Applications:

- mode
- 00 l/h liquids) rpm) rpm clamp hopper

¾" clamp

1/2" clamp

- clamp in liquids clamp
 - processes

MICRO-PLANT

- With 1 litre, single-wall funnel 1litre (standard)
- With three-way ball valve 2 litres (optional)
 - With circulation bend

Optionals:

- up tp 120 °C up to 2,5 bar

Applications:

- modules

Module ULTRA-TURRAX[®] UTC

Applications:

same dispersing element as T25 up to 2.000 ml 10.000 – 24.000 rpm

Dispersing in batch processes

- 2 litres in a glass flask 25 % (= 15 min/h)
 - - in structure, viscosity, etc
 - Requires a stand.

Accessories for module UTC (page):

- Stand R 2722, R 2723 (114 / 115)
- Boss head clamp R 271 (116)

Further accessories on request.

Standard speed Inlet connector Outlet connector General data

General data

Function mode

Generator

Flow rate

Speed range

Function mode	inline, three-step
Generator	standard: 2G, 4M, 6F
other	configurations available as optional
Flow rate 2G – 4M – 6F	approx. 80 l/h
Speed range	3.000 – 26.000 rpm
Standard speed	16.000 rpm
Inlet connector	¾" clamp
Outlet connector	1/2" clamp

inline, single-step

optional: 2G, 2P, 6F

approx. 100 – 200 l/h

3.000 – 26.000 rpm

16.000 rpm

¾" clamp

1/2" clamp

standard: 4M (medium)

General data	
Function mode	two counteropposing cones
	for setting milling grade
Flow rate	approx. 200 l/h
Speed range	3.000 – 26.000 rpm
Standard speed	16.000 rpm
Inlet connector	34" clamp
Outlet connector	½" clamp

1K / Module MKO
um homogeneity and stability ns and suspensions
/ dispersion spectra after just one



magicLAB®



- The unit creates a vacuum at the solid inlet so that the powder is directly sucked in from the

- For dustless, lump-free incorporation of powders

- For enriching liquids with powders in circulating

- For solid fractions up to approx. 65%

- Heated, double-wall, 2 litre funnel - Cover and special extensions for working with overpressure or vacuum

- For dispersion in circulating processes. - As extension kit for UTL, DR, MK and MKO

- For prescription development and defining optimum incorporation and mixing module - For calculating process parameters such as speed, shear frequency, temperature and time



- For dispersion in batch processing - For formulating small product quantities up to - For rapid identification of rheological changes

- For estimating the shear stress required





A 11 basic Batchmill for two milling processes

Technical data

Motor rating input	300 M
Motor rating output	160 W
Speed	28.000 rpm (fixed
Useful volume	80 m
Duty cycle ON / OFF	5 min / 10 mir
Overload protection	yes
Circumferential speed	53 m/s
Max. granularity of task	10 mm
Dimensions (W x D x H)	85 x 85 x 240 mn
Weight	1,5 kg
Permissible ambient temperature	5 – 40 °C
Permissible relative humidity	80 %
Protection class acc. to DIN EN 60529	IP 43

A 11 basic Analytical mill

Batch mill for 2 different grinding procedures: Impact grinding of hard, brittle or non-elastic grinding materials with high-grade stainless steel beater. This beater can be used for a Mohs hardness up to 6 (incl. with delivery). Cutting grinding for pulverizing soft, fibrous materials with a cutting blade (not incl. with delivery).

- Moist and gluey materials can be pulverized by adding water - Grinding chamber made of Tefcel (ETFE, glass fiber-reinforced) with stainless steel inlet (AISI 316L), useful volume 80 ml (incl. with delivery). For embrittle ment of grinding materials with liquid nitrogen in the grinding chamber - Optionally, a 250 ml grinding chamber is

available (page 88)

Accessories (page): A 11.1 Spare beater (87), A 11.2 Cutting blade (87), A 11.3 Beater (87), A 11.4 Grinding chamber (88), A 11.5 Spare grinding chamber (88), A 11.6 Double beater (88), A 11.7 Funnel (88)

A 11.1 Spare beater

stainl. steel (AISI 420)

stainl. steel (AISI 440B)

stainl. steel (AISI 440B)

For pulverizing substances with a Mohs hardness up to 6. Included with the analytical mill A 11 basic.

A 11.2 Cutting blade

For pulverizing soft, fibrous grinding materials. Not included with the analytical mill A 11 basic.

General data Material

General data Material

General data

Material

A 11.3 Beater

For pulverizing substances with a Mohs hardness up to 9, coated with chromium carbide. Not included with the analytical mill A 11 basic.

IKA[®] Crushing Analytical mill accessories

Ident. No. 2905200

2900000

2900001

Ident. No.

2904600

Ident. No. 2983000





230 V 50/60 Hz

115 V 50/60 Hz



87

Accessories Analytical mill 88



A 11.4 Grinding chamber

Made of polycarbonate with stainless steel inlet. Not suitable for cooling with N₂, only applicable with double beater A 11.6. Not included with the analytical mill A 11 basic.

General data	
Useful volume	250 ml
Material	stainl. steel (AISI 316L)

General data

Useful volume

Material

Ident. No. 2904100

Ident, No.

2983100



A 11.5 Spare grinding chamber

Made of Tefcel (ETFE, glass fibre-reinforced) with stainless steel inlet. Excellent resistance to chemicals and low temperatures (- 200 °C). Included with the analytical mill A 11 basic.

ldent. No. 3302900	

A 11.6 Double beater

For use up to Mohs hardness 3. Only applicable with grinding chamber A 11.4. Not included with the analytical mill A 11 basic.

<	7

Ident. No. 3048700

General data	
Material	titanium, surface-hardened

80 ml

stainl. steel (AISI 316L)

A 11.7 Funnel

Prevents splashing by pouring in liquid nitrogen in the grinding chamber A 11.5. Not included with the analytical mill A 11 basic.

General data	
Material jacket	PTFE
Material sieve	stainl. steel (AISI 316L)

Technical data Motor rating input Motor rating output 20.000 rpm (fixed) Speed Circumferential speed Overload protection current limitation Useful volume Material grinding chamber stainl. steel (AISI 304) Material cover stainl. steel (AISI 304) Max. granularity of task max. 6 – 7 mm Duty cycle ON / OFF

(with cooling)	7 min / 10 min
Weight	6,6 kg
Dimensions (W x D x H)	170 x 170 x 350 mm
Permissible ambient temperature	5 – 40 °C
Permissible relative humidity	80 %
Protection class acc. to DIN EN 60529	IP 21

M 20 Universal mill

440 W brittle substances. 225 W - Double-walled grinding chamber 72 m/s hose adapters

250 ml - Two grinding chambers can be

Accessories (page):

stainl. steel (1.4122)

ungsten carbide (86,5 WC 13,5 Col)

Suitable for crushing materials up to Mohs hardness 5. Included with M 20.

M 22 Hard metal cutter

Made of tungsten carbide for hard materials up to Mohs hardness 9. Not included with M 20.

M 23 Star-shaped cutter

Used to crush fibrous substances such as paper and vegetation, but also for plastics and material with a low specific weight. Not included with M 20.

M 20.1 Grinding chamber

A second grinding chamber ensures effective processing. The grinding chambers can be placed on the drive alternately. One chamber is cleaned and filled while the other is being processed. Cutters are not included with M 20.

Accessories (page): M 21 Spare cutter (89), M 22 Hard metal cutter (89), M 23 Star-shaped cutter (89)

General data

General data

General data

Material

Material

stainl. steel (AISI 304) Material

IKA[®] Crushing Universal mill and accessories

Batch mill suitable for dry grinding of hard and

- can be cooled with water through two
- Removable grinding chamber, easy to clean alternately operated using one drive - M 21 blade incl. with delivery
- M 21 Spare cutter (89), M 22 Hard metal cutter (89), M 23 Star-shaped cutter (89), M 20.1 Grinding chamber (89)



Ident. No. 1603600 1603603



Ident. No. 0521800

Ident No 0328200



Ident. No. 1443400





Ident, No. 8006200

90 Microfine grinder and accessories



Ident. No. 2836000 230 V 50/60 Hz 2836001 115 V 50/60 Hz



ldent. No. 2870900



ldent. No. 2871000



MF 10 basic Microfine grinder drive

Continuously operating universal grinder. - Powerful drive

- Easy to clean working surface made of stainless steel
- Two different grinding heads can be
- attached to the drive
- Heads are easily changeable
- Grinding heads not incl. with delivery

Accessories (page):

MF 10.1 Cutting-grinding head (90), MF 10.2 Impact grinding head (90)

MF 10.1 Cutting-grinding head

For crushing fibrous substances such as paper and vegetation, but also for plastics and material with a low volume weight. Before being discharged, the ground material passes through a sieve. This sieve is interchangeable and available in different hole sizes (not incl. with delivery). The ground material can then be collected using an NS 29 standard ground vessel.

Accessories (page): MF Sieve (90)

MF 10.2 Impact grinding head

For crushing brittle, hard materials such as minerals, building materials up to Mohs hardness 6. Before being discharged, the ground material passes through a sieve. This sieve is interchangeable and available in different hole sizes (not incl. with delivery). The ground material can then be collected using an NS 29 standard ground vessel.

Accessories (page): MF Sieve (90)

MF Sieve

Interchangeable sieves for insertion into the grinding heads ensure maximum particle size filtering.

Fechnical data	
Notor rating input	1.000 W
Notor rating output	500 W
Speed range	3.000 – 6.500 rpm
Circumferential speed	
Cutting-grinding head	22,5 m/s
mpact grinding head	31,4 m/s
Vaterials in contact with medium	stainl. steel (AISI 316L)
Duty cycle ON / OFF	120 / 30 min
Overload protection	yes
Weight	9,7 kg
Dimensions (W x D x H)	320 x 300 x 380 mm
Permissible ambient temperature	5 – 40 °C
Permissible relative humidity	80 %
Protection class acc. to DIN EN 60529	IP 22

Technical data	
Circumferential speed	22,5 m/s
Max. granularity of task	max. 15 mm
Dimensions including MF 10 basic	320 x 300 x 560 mm
Weight incl. MF 10 basic	10,5 kg
Materials in contact with medium	stainl. steel
Grinding channel and cover	(AISI 304)
Blades	(AISI 440B)
Shaft, rotor, screws	(AISI 316L)

31,4 m/s

11 kg

max. 10 mm

stainl. steel

(AISI 304)

(AISI 304)

(AISI 316L)

320 x 300 x 450 mm

Technical data

Hammer beater

Shaft, rotor, screws

Circumferential speed

Max. granularity of task

Weight incl. MF 10 basic

Grinding channel and cover

Dimensions including MF 10 basic

Materials in contact with medium

MF 10 basic

Drive for inline microfine grinder. Grinding head and sieves not incl. with delivery, page 90 Ident. No. 2836000

MF 10.1

Cutting-grinding head, interchangeable with impact grinding head MF 10.2, page 90 Ident. No. 2870900

MF 10.2

Impact grinding head, interchangeable with cutting-grinding head MF 10.1, page 90 Ident. No. 2871000



MF 0.5

Sieve for insertion into cutting-grinding head MF 10.1 or impact grinding head MF 10.2, with hole size 0.5 mm, page 90 Ident. No. 2939000

MF 2.0

Sieve for insertion into cutting-grinding head MF 10.1 or impact grinding head MF 10.2, with hole size 2,0 mm, page 90 Ident. No. 2939400

General da	ata		
Material		stainl. stee	I (AISI 304)
Hole size (diameter)		
MF 0.25	0,25 mm	MF 2.0	2,0 mm
MF 0.5	0,5 mm	MF 3.0	3,0 mm
MF 1.0	1,0 mm	Wider holes	on request





HCT basic

Universal laboratory heating plate with revolutionary technology and new functions: new: Integrated temperature control new: Incl. PT 1000 temperature sensor (PT 1000.60) Page 96



Heating / Tempering



93

IKA[®] Heating / Tempering



Heating baths 94



HBR 4 digital Heating bath with digital display

Heating function Heat output Temperature range Setting tolerance Deviation Temperature display Safety class acc. to DIN 12877 Stirring function Stirring function General data Useful volume Material stainl. steel (AISI 304) Outer diameter Inner diameter

Outer height

Inner height

Heating function

Temperature range

Temperature display

Deviation (3 | H₂O, 90 °C)

Safety class acc. to DIN 12877

Setting tolerance

Stirring function

Stirring function

Speed range General data

Heat output

Permissible ambient temperature

Protection class acc. to DIN EN 60529

Permissible relative humidity

Weight

250 mm 160 mm 3,9 kg

HBR 4 digital Heating bath

Useful volume	4
Material	stainl. steel (AISI 304)
Outer diameter	340 mm
Inner diameter	200 mm
Outer height	250 mm
Inner height	160 mm
Weight	4,4 kg
Permissible ambient temperature	5 – 40 °C
Permissible relative humidity	80 %
Protection class acc. to DIN EN 60529	IP 21

The heating bath is characterized by the following features:

- Cylindrical bath shape
- ±5 K ±5 K - High-grade recyclable materials
 - The heating elements are situated
 - underneath the bath vessel
 - Either low viscosity oil (50 mPas) or water
 - can be used as the heat transfer fluid
 - Useful volume of 4 liters
 - Heat output 1.000 W - Infinitely adjustable safety temperature
- limiter acc. to DIN 12877 340 mm - Double jacket provides protection against burns 200 mm

- Accessories (page): H 240 Ring set (98)

2

yes 100 – 800 rpm

5-40 °C

80 %

IP 21

1.000 W RT – 225 °C

scale

2

no

4 |

- HBR 4 digital additionally: 1.000 W - Digital display presents rated, actual and
- RT 200 °C ±1K
 - Fuzzy logic control ±1K
 - digital
 - Integrated magnetic stirring drive to circulate the tempering fluid, which
 - contributes to improved heat distribution
 - The safety elements are checked
 - when the unit is switched on

Accessories (page):

(98), IKAFLON®-Stirring bars (29)

IKA[®] Heating / Tempering

Heating baths 95

HB 4 basic Heating bath



Ident. No. 2520000 230 V 50/60 Hz 2520001 115 V 50/60 Hz

- safety temperature as well as speed
- H 240 Ring set (98), H 159 Intermediate bottom



Ident. No.		
2602300	230 V	50/60 Hz
2602301	115 V	50/60 Hz

IKA[®] Heating / Tempering

96 Heating plates



Ident, No 2609900 230 V 50/60 Hz



Ident. No.	
3384100	230 V 50/60 Hz
3384101	115 V 50/60 Hz

NEW!

HP 30 digital IKATHERM®

Heating plate with a large set-up surface (300 x 300 mm).

- Adjustable heating capacities of 25 %, 50 % and 100 %
- Homogenous temperature distribution across entire plate
- Digital temperature display and electronic control
- Safety features to prevent skin contact - Surface is made of corrosive-resistant materials

Accessories (page):

HP 30.1 Temperature sensor (98)

HCT basic safety control IKATHERM®

Universal laboratory heating plate with revolutionary technology and new functions: new: Integrated temperature control new: Incl. PT 1000 temperature sensor

- (PT 1000.60) new: Exact temperature setting via digital display, even when switched off
- new: Set safety temperature limit displayed digitally
- new: Hot Top indicator >> hot surface warning to prevent burns!
- new: Digital error code display
- With adjustable safety circuit of heating plate temperature (50 - 360 °C)
- Heating plate, suitable for unsupervised operation - Bushing according to DIN 12878 for connecting
- a contact thermometer, e.g. ETS-D5, enables
- precise temperature control - High level of safety thanks to improved heat control technology
- Enclosed assembly (IP 42) guarantees long service life
- Highly polished aluminium heating plate for optimum heat transfer

Accessories (page):

Electronic contact thermometers (117): ETS-D5, ETS-D6, Bath attachments (28): H 15, H 28, Oil bath attachments (28): H 29, H 30

leating function	
Heat output	2.500 W
lemperature range	RT – 300 °C
Setting tolerance	±1K
Sensor for temperature in medium	PT 100
Control accuracy with temperature sensor	± 1 K
Heating plate	
Heating plate material	cast iron
Heating plate dimensions	300 x 300 mm
General data	
Dimensions (W x D x H)	320 x 400 x 140 mm
Weight	15 kg
Permissible ambient temperature	5 – 40 °C
Permissible relative humidity	80 %
Protection class acc. to DIN EN 60529	IP 21

Heating function	
Temperature display	digital
Heat output	600 W
Heating rate (1 H ₂ O in H 15)	6,5 K/min
Temperature range	RT – 310 °C
Setting tolerance	± 1 K
Temperature undulation	
without temperature sensor	± 2 K
Adjustable safety circuit	50 – 360 °C
Digital temperature	
limit display	50 – 360 °C
Control accuracy	PT 1000 / ± 1 K
with sensor	ETS-D5 / ± 0,5 K
	ETS-D6 / ± 0,2 K
Heating plate	
Material	alloy
Dimensions	Ø 135 mm
General data	
Dimensions (W x D x H)	160 x 270 x 85 mm
Weight	2,5 kg
Permissible ambient temperature	5 – 40 °C
Permissible relative humidity	80 %
Protection class acc. to DIN EN 60529	IP 42

Operating temperature range Operating temperature range 28 – 300 °C Min. temperature with refrigerator -20 °C Temperature stability at 70 °C 0,02 K Temperature adjustment / digital temperature indication digital 0,1 K Resolution of display Absolute accuracy setup for calibration PT 100 Temperature control internal External sensor PT 100 Analog Interface In/Out 4 – 20 mA alternative 1 – 5 V Safety classification FL Heating power 3 kW Pressure pump With 12 mm connection 33 l/min Delivery pressure (head) 0,7 bar With 12 mm connection 22 l/min Delivery suction pressure (head) 0,4 bar Pump connection M 16 x 1 General data 50 mm²/s Max. permissible kin. viscosity Bath volume 8.51 5,2 I Bath capacity with displacement rack 130 x 110 / 155 mm Width bath opening WxD/ bath depth Dimensions (W x D x H) 240 x 405 x 390 mm Height of bath opening 190 mm 18 kg Weight 230 V 1~ 50/60 Hz Power supply requirement Max. current 14,1 A

Fuse Permissible ambient temperature

Heating function 1.500 W Heat output 25 – 100 °C Temperature range Temperature display Temperature stability (70 °C) ± 0,12 K Adjustable temperature limitation 25 – 200 °C Max. pump pressure 0,08 bar Max. delivery rate 5 l/min General data Dimensions (W x D x H) 105 x 139 x 319 mm Weight Per

2,3 kg
5 – 40 °C
80 %
IP 31
yes

Heating circulator bath with housing, bath and all moistened parts are made of stainless steel. With cooling coil for water-cooling, pressure- and suction pump. With adjustable overtemperature protection according to DIN 12876.

Complete functions: With level protection and maximum and minimum setpoint for additional safety, external temperature sensor connection, external temperature control and temperature programmer (50 segments, may be split into 10 programs), interactive, contains a digital RS 232/ RS 485 interface as well as a (4...20mA) analogue interface for bidirectional communication. Plug & Play Technology - new generation of microprocessor controlled compatible control. Simple operation with a rotary knob and digital display, easy control, clear text, menu-driven, set point limiting, visually and acoustically alarm, mains failure automatic, programmable.

Accessories (page): LT 5.20 Hose (99), Hose adapters (99): LT 5.22, LT 5.23, LT 5.24, PC 2.3 Cable (99), labworldsoft® (143), Temperature sensors (99): PT 100.5, PT 100.7

16 A

scale

5 – 32 °C

For temperature control of liquids (NFL/I) up to 100 °C in open baths (min. bath depth 160 mm,

- min. usable depth 75 mm). - Complies with all safety requirements for
- electrically operated devices - Intended for supervised use
- only
- bath vessels

Accessories (page): Bath vessels (98)

IKA[®] Heating / Tempering Thermostats

CC3-308B vpc Circulation thermostat



Ident. No.	
3658800	230 V 50/60 Hz
3658801	115 V 50/60 Hz

NEW!

EH 4 basic Immersion thermostat

- For operation with non-flammable liquids

- With universal clamp, suitable for all standard



Ident. No.	
3164000	230 V 50/60 Hz
3164001	115 V 50/60 Hz

IKA[®] Heating / Tempering

Accessories heating baths / heating plates and thermostats 98



H 240 Ring set

To cover the heating baths HB 4 basic and HBR 4 digital. Prevents dust penetration, uncontrolled heat dissipation as well as water absorption and the formation of oil mist when working with oil.



H 159 Intermediate bottom

Allows vessels to be inserted in the heating bath HBR 4 digital without obstructing movement of the rotating magnetic bars.



HP 30.1 Temperature sensor

For heating plate HP 30 digital. For controlling the temperature in the medium.

General data	
Immersion depth	290 mm
Diameter	6 mm
Material	stainl. steel (AISI 304)

10

25 – 185 mm

General data

Number of rings

Diameter of opening (variable)

General data	
Material	metal
Length	1 m
Max. temperature	300 °C

LT 5.20 Hose

Coated metal hoses for circulation thermostat CC3-308B vpc. Package contains 2 hoses.

Accessories (page): Hose adapters (99): LT 5.22, LT 5.23, LT 5.24

General data R 1/8" x M 16 x 1 Dimensions adapter LT 5.22 Dimensions adapter LT 5.23 R 1/4" x M 16 x 1

For connection to the kneader HKD-T 06 D. LR 2000.3 and LR 2000.4

LT 5.24 Hose adapter R 1/8" x M 16 x 1

For connection to the reactor vessels LR 2000.1 and LR 2000.2

General data Length

General data

Dimensions adapter

PC 2.3 Cable

For connecting the circulation thermostat CC3-308B vpc control to a PC (9 pin interface)

3 m

Ident. No 1752900

R 350 Universal clamp

For clamping flask necks, condensers, etc. up to 11 cm diameter.





Ident. No. 3335000 EH 4.1 (5 l) 3335100 EH 4.2 (11 l) 3335200 EH 4.3 (18 I)

Polycarbonate bath vessels, suitable for use with the immersion thermostat EH 4 basic, up to 100 °C.

General data		
Material		polycarbonate
Volume without vessels		5, 11, 18
Outer dimensions (W x D x H)	EH 4.1	132 x 280 x 160 mm
	EH 4.2	350 x 313 x 168 mm
	EH 4.3	350 x 473 x 168 mm
Inner dimensions (W x D x H	EH 4.1	120 x 262 x 150 mm
	EH 4.2	302 x 295 x 150 mm
	EH 4.3	302 x 455 x 150 mm

General data	
Length	255 mm
Diameter	6 mm
Material	stainl. steel (AISI 316L)

PT 100.5

Temperature sensor for use with laboratory reactor systems LR 2000.

		PT 100.7
General data		
_ength	135 mm	Temperatu
Diameter	3 mm	laboratory
Vlaterial	stainl. steel (AISI 316L)	

IKA[®] Heating / Tempering

Accessories Thermostats 99

LT 5.22 Hose adapter LT 5.23 Hose adapter

LT 5.23: For connection to the reactor vessels



Ident. No. 2807000 LT 5.22 2235000



LT 5.23











Ident. No. 2611500



IKA®-PET

High-precision air cushion piston stroke pipette for 25 μ l, 50 μ l and 100 μ l for science, research and routine work in the field of liquid handling. **Page 102**

Fixed volume pipettes 102 - 103

Liquid Handling

IKA[®] Liquid Handling

102 Fixed volume pipettes and accessories



- pe

IKA®-PET for 25 µl, 50 µl and 100 µl

High-precision air cushion piston stroke pipette for science, research and routine work in the field of liquid handling.

- Single channel air cushion piston stroke pipette (fixed volume)
- · Pipette and tips with certificate of conformity acc. to DIN EN ISO 8655
- ika A - Maximum accuracy and precision of pipette and tip system
- High quality tip with ULR (Ultra low retention) surface
- High quality casing and flask material
- Minimum manpower required even for blow out and tip ejection
- Simple operation due to ergonomic design
- Easy to maintain thanks to the assembly tool provided
- Adjustable even for liquids with different density to water
- Includes software for measurement calibration and analysis

Model 25 µl	
Nominal volume	25 µl
Accuracy systematic	
error [e _s]	± 1,0 %
Precision random	
error [CV]	<u>≤</u> 0,3 %
Model 50 µl	
Nominal volume	50 µl
Accuracy systematic	
error [e _s]	± 0,7 %
Precision random	
error [CV]	<u>≤</u> 0,3 %
Model 100 µl	
Nominal volume	100 µl
Accuracy systematic	
error [e _s]	± 0,6 %
Precision random	
error [CV]	<u>≤</u> 0,2 %
General data	
Dimensions (W \times D \times H)	28 x 52 x 206 mm
Weight	0,075 kg

IKA®-PET Rack 1 and Rack 3
Pipette stands fo IKA®-pipettes.

IKA®-PET soft

- with fixed and variable volumes
- Testing in accordance with
- DIN EN ISO 8655 - Automatic transfer of measurements from scales, air pressure, humidity and temperature
- by RS 232 interface
- Automatic correction of environmental air
- pipette is possible
- manufacturers already stored
- Inventory of individual pipettes
- Date-controlled monitoring of calibration cycles - Measuring results saved in a database
- Log according to GLP guidelines with details of all individual measurements, average, accuracy and precision

ULR-pipette tips

Considerably improved inner surface, similar to Lotus Effect.

Description	Pack (pcs.)	sterile	Colour	Ident. No.
IKA®-TIP 200 µl – bag	1.000	no	natural	3330700
IKA®-TIP 200 μl – box	2 x 1.000	no	natural	3330701
IKA®-TIP 200 µl – 96-rack	10 × 96	no	natural	3330702
IKA®-TIP 200 µl – 96-rack	10 × 96	yes	natural	3330703
IKA®-TIP Filter 1 – 100 µl – bag	1.000	no	natural	3330800
IKA®-TIP Filter 1 – 100 µl – bag	1.000	yes	natural	3330801
IKA®-TIP Filter 1 – 200 µl – bag	1.000	no	natural	3330900
IKA®-TIP Filter 1 – 200 µl – bag	1.000	yes	natural	3330901
IKA®-TIP Filter 100 µl – 96-rack	10 × 96	yes	natural	3331000
IKA®-TIP Filter 200 µl – 96-rack	10 × 96	yes	natural	3331100

IKA[®] Liquid Handling

Accessories 103



Ident. No.	
3224000	Rack 1
3224001	Rack 3

Is a user-friendly program for the automatic recording, analysis and management of gravimetric measurements for pipette calibration. - Suitable for single and multi-channel pipettes

- pressure and temperature conditions
- Programm controlled recalibrating of the
- Master data for over 580 pipettes by other

- Compatible with all Microsoft® Win32 operating systems (Win®95/98/ME/NT/2000/XP/Vista etc.)



Ident. No. 3220300



Glass condenser

With space-saving vertical condenser. Included with delivery of set of glassware RV 06.2. Page 110





IKA[®] Distilling 106 Rotary evaporators



et of glassware	RV 06.1
ype of condenser	diagonal
ooling surface	1.200 cm ²
rive	RV 06-ML
lotor type	DC motor
lotor rating input / output	45 / 36 W
peed range	10 – 240 rpm
ead inclination, adjustable	± 10 °
ift	
troke min. / max.	74 / 150 mm
rive	motor
1ax. load	10 kg
lignment end stop on top	76 mm
ower limit switch	fixed
eating bath	HB 4 basic
emperature range	RT – 225 °C
eating output	1.000 W
etting accuracy	± 5 K
ontrol deviation	± 5 K
eneral data	
imensions (W x D x H)	840 x 390 x 880 mm
leight	18 kg
ermissible ambient temperature	5 – 40 °C
ermissible relative humidity	80 %
rotection class acc. to DIN EN 60529	IP 21

Set of glassware	RV 06.2
Type of condenser	vertical
Cooling surface	1.200 cm ²
Drive	RV 06-ML
Motor type	DC motor
Motor rating input / output	45 / 36 W
Speed range	10 – 240 rpm
Head inclination, adjustable	± 10 °
Lift	
Stroke min. / max.	74 / 150 mm
Drive	motor
Max. load	10 kg
Alignment end stop on top	76 mm
Lower limit switch	fixed
Heating bath	HB 4 basic
Temperature range	RT – 225 °C
Heating output	1.000 W
Setting accuracy	± 5 K
Control deviation	± 5 K
General data	
Dimensions (W x D x H)	640 x 390 x 1.130 mm
Weight	18 kg
Permissible ambient temperature	5 – 40 °C
Permissible relative humidity	80 %
Protection class acc. to DIN EN 60529	IP 21

RV 06-ML 1-B

adjustment.

the glassware is not moved

- Condenser geometry with 1.200 cm² cooling surface, higher yield and increased condensation power

- No chimney effect - Rodaviss screw joint allows bonded ground joints to be released easily,

```
Accessories (page):
controller (120)
```

RV 06-ML 2-B

Consisting of heating bath HB 4 basic, set of glassware RV 06.2 with vertical condenser (spacesaving), drive RV 06-ML with lift with electrical height adjustment. - Wearless drive with brushless DC motor

- the glassware is not moved - Condenser geometry with 1.200 cm² cooling surface, higher yield and increased
- condensation power - No chimney effect
- components
- Suitable for DIN EN 12697-3 (Asphalt test for hot asphalt)

Accessories (page): RV 06.1 Set of glassware (109), VC 2 Vacuum controller (120)

IKA[®] Distilling

Rotary evaporators 107



- Wearless drive with brushless DC motor
- Convenience and safety with motorized lift;
- removable screw connections facilitate
- cleaning of glass components

RV 06.2 Set of glassware (110), VC 2 Vacuum

ldent. No.		
8010000	230 V	50/60 Hz
8010001	115 V	50/60 Hz

- Convenience and safety with motorized lift;
- Rodaviss screw joint allows bonded ground joints to be released easily, removable screw connections facilitate cleaning of glass



Ident. No.	
8010100	230 V 50/60 Hz
8010101	115 V 50/60 Hz

IKA[®] Distilling

108 Rotary evaporators



Ident, No. 230 V 50/60 Hz 8017900 8017901 115 V 50/60 Hz



Consisting of heating bath HB 4 basic, set of glassware RV 06.1 with diagonal condensor, drive RV 05 basic, telescopic stand RV 05.3 and boss head clamp R 271.

- Easy and jolt-free raising and lowering of the rotary evaporator
- Telescopic stand tiltable to the side
- Rodaviss screw joint allows bonded ground joints to be released easily, removable screw connections facilitate cleaning of glass components
- Condenser geometry with 1.200 cm² cooling surface, higher yield and increased condensation power
- No chimney effect

Accessories (page):

RV 06.2 Set of glassware (110), VC 2 Vacuum controller (120)

Set of glassware	RV 06.1
Type of condenser	diagonal
Cooling surface	1.200 cm ²
Drive	RV 05 basic
Motor type	asynchronous
Motor rating input / output	133 / 65 W
Speed range	46 – 260 rpm
Head inclination, adjustable	any
Lift	RV 05.3
Stroke	190 mm
Max. load	10 kg
Swivel feature	90 °
Heating bath	HB 4 basic
Temperature range	RT – 225 °C
Heating output	1.000 W
Setting accuracy	± 5 K
Control deviation	± 5 K
General data	
Dimensions (W \times D \times H)	520 x 450 x 900 mm
Weight	12 kg
Permissible ambient temperature	5 – 40 °C
Permissible relative humidity	80 %
Protection class acc. to DIN EN 60529	IP 21

Drive	
Motor type	asynchronous
Motor rating input	133 W
Motor rating output	65 W
Speed range	46 – 260 rpm
Speed display	scale
General data	
Dimensions (W x D x H)	130 x 200 x 260 mm
Weight	4,5 kg
Permissible ambient temperature	5 – 40 °C
Permissible relative humidity	80 %
Protection class acc. to DIN EN 60529	IP 21

General data

Max. load

Stroke

Height

Diameter of support rod

Dimensions (W x D x H)

RV 05 basic Drive

drive output is transmitted directly to the vapor tube via a control gear with secondary torque coupling. bearing

heavy loads

Accessories (page): R 271 Boss head clamp (116)

RV 05.3 Telescopic stand

Raising is made easier by a jolt-free pneumatic spring.

34 mm

190 mm 710 – 900 mm

diagonal

1.200 cm²

580 x 450 x 900 mm

10 kg



ldent. No.	
8018000	230 V 50/60 Hz
8018001	115 V 50/60 Hz

RV 05 basic 2-B

Consisting of heating bath HB 4 basic, set of glassware RV 06.2 with vertical condenser (space-saving), drive RV 05 basic, telescopic stand RV 05.3 and boss head clamp R 271.

- Easy and jolt-free raising and lowering of the rotary evaporator
- Telescopic stand tiltable to the side - Rodaviss screw joint allows bonded ground
- joints to be released easily, removable screw connections facilitate cleaning of glass components
- Condenser geometry with 1.200 cm² cooling surface, higher yield and increased condensation power
- No chimney effect
- Suitable for DIN EN 12697-3 (Asphalt test for hot asphalt)

Accessories (page):

RV 06.1 Set of glassware (109), VC 2 Vacuum controller (120)

Set of glassware	RV 06.2
Type of condenser	vertical
Cooling surface	1.200 cm ²
Drive	RV 05 basic
Motor type	asynchronous
Motor rating input / output	133 / 65 W
Speed range	46 – 260 rpm
Head inclination, adjustable	any
Lift	RV 05.3
Stroke	190 mm
Max. load	10 kg
Swivel feature	90 °
Heating bath	HB 4 basic
Temperature range	RT – 225 °C
Heating output	1.000 W
Setting accuracy	± 5 K
Control deviation	± 5 K
General data	
Dimensions (W x D x H)	580 x 480 x 1.000 mm
Weight	12 kg
Permissible ambient temperature	5 – 40 °C
Permissible relative humidity	80 %
Protection class acc. to DIN EN 60529	IP 21

General data	
Type of condenser	
Cooling surface	

RV 06.1 Set of glassware

distillation tasks.

and without difficulty

Accessories (page):

IKA[®] Distilling Rotary evaporators accessories 109

- Rotary evaporator drive, newly designed. The
- Condenser motor supported by means of ball
- Constant operation is guaranteed even with

RV 05.3 Lift (109), Set of glassware (109 / 110): RV 06.1, RV 06.2, HB 4 basic Heating bath (95),



Ident. No. 230 V 50/60 Hz 3075000 115 V 50/60 Hz 3075001



Ident. No. 3154100



- Diagonally mounted condenser for all standard
- Can be assembled and disassembled quickly
- Included with delivery: one 1.000 ml evaporating flask and one 1.000 ml receiving flask

Evaporating flasks (111): RV 06.4, RV 06.5, RV 06.6, RV 06.7 Receiving flask (111), RV 06.11 Vapor tube (110), Seals (110): RV 06.13, RV 06.15



IKA[®] Distilling

110 Rotary evaporators accessories



RV 06.2 Set of glassware

With space-saving vertical condenser. The distributor part is fitted with a condensate blocker as well as a discharge channel which prevents the condensate from getting into contact with the seal. A PTFE inlet pipe facilitates constant infeed of the distillate. Included with delivery: one 1.000 ml evaporating flask and one 1.000 ml receiving flask.

Accessories (page):

Evaporating flasks (111): RV 06.4, RV 06.5, RV 06.6, RV 06.7 Receiving flasks (111), RV 06.11 Vapor tube (110), Seals (110): RV 06.13, RV 06.15

ldent. No.	
1958000	

RV 06.11 Vapor tube

For set of glassware RV 06.1 and RV 06.2.

General data	
Diameter	21,6 mm

vertical

1.200 cm²

General data

Cooling surface

Type of condenser

General data	
Volume	11
Material	borosilicate glass

Material



RV 06.13 Seal

For RV 06.11.





RV 06.15 Seal

For RV 06.11, resistant to solvents.

General data	
Material	PTFE

Evaporating flasks, NS 29

		RV 06.4
General data		
Volume	11	
Material	borosilicate glass	
		RV 06.5
General data		
Volume	2	
Material	borosilicate glass	
		RV 06.6
General data		
Volume	0,1	

RV 0	
1	
oato alaco	h

borosilicate glass

IKA[®] Distilling Rotary evaporators accessories



06.7 Receiving flask, KS 35



ETS-D5

Electronic contract termometer ensures perfect temperature control without overshooting the set temperature, even in the case of quick heating. With optimised PID control and RESET function, incl. stainless steel sensor H 62.51. Page 117

IKA' ETS-D5

Stands	114 – 115
Fixing elements	116
Temperature	
measuring instrument	117 – 118
Revolution counter	119
Vacuum controller	120
Vacuum pump / valve	121

Accessories



IKA[®] Mechanical accessories

114 Stands



R 103 Plate stand

R 104 Stand

Suitable for small instruments such as the overhead stirrer RW 11 basic.

Accessories (page): Boss head clamp H 44 (116) H 44 Boss head clamp (116)

Accessories (page): R 200 Clamp (116),

Small stand for T 10 basic.

With slip resistant foil.

Plate stands

R 1825

R 1826

R 1827

Accessories (page): Boss head clamp R 182 (116), RH 3 Strap clamp (116)

R 2722 H-Stand

Particularly stable stand with H-shape base which prevents the stand from tipping backwards. Provides optimum stability required for larger, heavier instruments and attachments, for example with rheological measurements using overhead stirrers. The stand has an adjustment screw which can be used to compensate for an uneven laboratory table surface.

Accessories (page): Boss head clamps (116): R 270, R 271, RH 5 Strap clamp (116)



R 2723 Telescopic stand RV 05.3 Telescopic stand T 653 Telescopic stand R 474 Telescopic stand R 472 Floor stand

Similar to R 2722, additional- Specially designed for the ly equipped with a pneumatic spring stand rod, which enables heavy instruments / attachments to be raised and spring. lowered smoothly without difficulty, e.g. with rheological measurements using overhead stirrers. The stand has an adjustment screw which can be used to compensate for an uneven laboratory table surface.

Accessories (page): Boss head clamps (116): R 270, R 271, RH 5 Strap clamp (116)

rotary evaporator drive RV 05. Raising is made easier by a jolt-free pneumatic

Accessories (page): Boss head clamp R 271 (116)

Specially designed for the dispersing instrument T 65 D. The stand is equipped with a pneumatic spring which enables effortless raising and lowering of the dispersion unit.

Descripton	R 103 Plate stand	R 104 Stand	R 1825	R 1826	R1827	R 2722 H-Stand	R 2723 Telescopic stand	RV 05.3 Telescopic stand	T 653 Telescopic stand	R 474 Telescopic stand	R 472 Floor stand
dent. No.	2972500	3386000	3160000	3160100	3160200	1412000	1412100	3154100	1608000	1643000	0738700
Diameter of support rod	10 mm	10 mm			16 mm	34 mm	34 mm	34 mm	48 mm	48 mm	-
Plate diameter	160 mm	-			-	-	_	_	_	_	-
Dimensions (W x D)	-	242 x 355 mm		20	0 x 316 mm	460 x 420 mm	460 x 420 mm	580 x 450 mm	460 x 530 mm	460 x 530 mm	80 x 80 mm
Heigt	360 mm	370 mm	560 mm	800 mm	1.000 mm	1.010 mm	620 – 1.010 mm	710 – 900 mm	1.200 mm	1.200 mm	2.020 mm
Vlax. load	1 kg	0,7 kg			5 kg	10 kg	10 kg	10 kg	-	-	-
Stroke	_	-			-	_	390 mm	190 mm	500 – 1.000 mm	500 – 1.000 mm	980 – 1.860 mm

IKA[®] Mechanical accessories Stands



Specially designed for the overhead stirrer RW 47 D; can be adapted for use with other instruments. The stand adapted for use with other is equipped with a pneumatic spring which enables effortless raising and lowering of the dispersion unit.

Accessories (page): SI 400 Safety switch (43), SI 474 Fixing device (43)

Mobile floor stand, specially designed for the overhead stirrer RW 47 D; can be instruments.

Accessories (page): SI 400 Safety switch (43), SI 472 Fixing device (43)

115

IKA[®] Mechanical accessories

116 Fixing elements



Ident. No 3372000



For fastening the T 10 basic to the stand R 104 (page 114) (included with delivery of T 10 basic).

RH 3 Strap clamp

For securing vessels against walls or for synchronized rotation during stirring or dispersing.

General data	
For stand diameter	8 – 16 mm
For vessel diameter	40 – 300 mm

General data

General data

General data Clamping range - stand

General data

General data

Diameter of extension arm

Length of extension arm

Material

Clamping range - stand

Clamping range - extension arm

Material

Material

Clamping range - stand

Clamping range - extension arm

Clamping range - extension arm

Material

Clamping range - stand

Clamping range - extension arm

10 – 11 mm

6 – 16 mm

6 – 16 mm

cast aluminium

25 – 36 mm

5 – 21 mm

34 mm

16 mm

8 mm

130 mm

cast aluminium

cast aluminium

cast aluminium

11 mm

Ident No 3159000

RH 5 Strap clamp

For securing vessels against walls or for synchronized rotation during stirring or dispersing, incl. boss head clamp R 270 (page 116).

General data	
For stand diameter	25 – 36 mm
For vessel diameter	40 – 300 mm

Temperature	
Temperature measuring range	-50 – 450 °C
Resolution	0,1 K
Measuring accuracy ±	0,2 K + Sensor tolerance PT 1000
	DIN IEC 751 class A
Setting accuracy	0,1 K
Control deviation	± 0,5 K
General data	
Supply voltage	8 – 16 VDC
Power consumption	10 mA (at 9 V)
Max. ON time	100 %
Plug	6 pin DIN 45322
Connection	DIN 12878 class 2
Dimensions (W x D x H)	82 x 22 x 83 mm
	(without sensor)
Weight	0,2 kg
Permissible ambient temperat	ure 0 - 60 °C
Permissible relative humidity	80 %
Protection class acc. to DIN EI	N 60529 IP 54

Electronic Contact Thermometers ETS-D5 and ETS-D6

Ensures perfect temperature control without overshooting the set temperature, even in the case of quick heating. With optimised PID control and RESET function, incl. stainless steel sensor H 62.51. For all magnetic stirrers with contact thermometer bushing according to DIN 12878, class 2 (e.g. IKA®, Heidolph and Corning with adapter AD-C1, Ident. No. 3414000, please order separately).

ETS-D6 additionally:

(without pH electrode) user guide

Patented: 3 modes of operation guarantee optimum adjustment to your working method. Operating mode A Suitable for work with varying parameters (from -50 °C to 450 °C). Safety temperature adjustable. Operating mode B Suitable for series operation under uniform conditions.

Operating mode C Suitable for unsupervised operation.

All values are taken from the memory. This ensures perfect protection against inadvertent improper adjustment.

H 38 Holding rod (27)

-50 – 450 °C

DIN IEC 751 class A

± 0,05 K + Sensor tolerance PT 1000

0,01 K

0,1 K

± 0,2 K

0 – 14 pH

Accuracy	± 0,1 pH
Resolution	± 0,01 pH
pH connection	BNC bushing
General data	
Supply voltage	8 – 16 VDC
Power consumption	15 mA (at 9 V)
Max. ON time	100 %
Plug	6 pin DIN 45322
Connection	DIN 12878 class 2
Dimensions (W x D x H)	96 x 45 x 98 mm
	(without sensor)
Weight	0,2 kg
Permissible ambient temperature	0 – 60 °C
Permissible relative humidity	80 %
Protection class acc. to DIN EN 60529	IP 54

Temperature

Measuring accuracy

Setting accuracy

Control deviation

pH measurement

Measuring range

Resolution

Temperature measuring range

IKA[®] Electronic accessories

ETS-D5

Ident. No.

3378000

NEW!

Temperature measuring instrument

IKA" ETS-D5

- With integrated pH measuring instrument - Large, graphic LCD display with multilingual

Accessories ETS-D5 and ETS-D6 (page): Sensors (25): H 62.51, H 66.51, H 70 Extension cable (25), H 52 Power pack set (25), H 16 V Support rod (27), H 44 Boss head clamp (27),

ETS-D6

Ident. No. 3378100



IKA[®] Electronic accessories

118 Temperature measuring instrument



Ident. No. 90 – 240 V 50/60 Hz 3113200







	Ident. No.	
1	3122100	PT 100.23
2	3122200	PT 100.24
3	3122300	PT 100.25
4	3122500	PT 100.27

Ident, No.

3127800

2616800



DTM 12 IKATRON[®] Digital temperature measuring instrument

For measuring temperatures between -200 °C up to +400 °C.

- LED display - Analog output (1 $^{\circ}C = 1mV$)
- Almemo interface for PC connection
- Sensor connection: Almemo

Accessories (page): Temperature sensors (118): PT 100.23, PT 100.24, PT 100.25, PT 100.27, DTM 12.10 Data cable (118), labworldsoft® (143)

Protective pipe, glass-coated. For use in acid and

E.g. for use with IKA® laboratory reactors in

combination with sensor receptacle

LR 2000.60 (page 134).

Temperature sensors

Standard sensor for a wide range of

PT 100.23

PT 100.24

PT 100.25

laboratory tasks.

alkaline solutions.







пт	100	77
	100	.27

With screw joint. Specially designed for IKA[®] laboratory kneader HKD-T 0.6 D.

DTM 12.10 Data cable, 9 pins (F)

Data cable with RS 232 interface to connect the DTM 12 with a PC.

PC 1.2 Adapter, 25 pins

9 pins (F) to 25 pins (F).

Measuring device	
Sensor	PT 100
Measuring range	-200 - 400 °C
Temperature display	digital
Resolution	0,01 K
General data	
Interface	Almemo, analog
Dimensions (W x D x H)	125 x 150 x 70 mm
Weight	1,1 kg
Permissible ambient temperature	0 – 50 °C
Permissible relative humidity	80 %
Protection class acc. to DIN EN 60529	IP 50

T 100.23	
laterial of protective pipe	stainl. steel (AISI 316L)
liameter	3 mm
ength	250 mm
leasuring range	-50 – 200 °C
esolution	0,01 K

PT 100.24	
Material of protective pipe	borosilicate glass
Diameter	8 mm
Length	250 mm
Measuring range	-50 – 200 °C
Resolution	0,01 K
PT 100.25	

Resolution

Material of protective pipe	stainl. steel (AISI 316L)
Diameter	6 mm
Length	255 mm
Measuring range	-50 - 400 °C
Resolution	0,1 K
PT 100.27	
Material of protective pipe	stainl. steel (AISI 316L)
Diameter	3 mm
Length	135 mm
Measuring range	-50 – 200 °C

0,01 K

Ge	eneral data	
Se	nsor	DZM-S.o
Sp	eed range	0 – 50.000 rpm
Ter	mperature co-efficient	0,005 %/°C
Me	easurement error of measured value	0,4 % = 1 Digit
An	alog output (0 – 4.000 rpm)	1 mV
	(> 4.000 rpm)	0,1 mV
Int	(> 4.000 rpm) rerface	0,1 mV RS 232
Int Dir	(> 4.000 rpm) terface mensions (W x D x H)	0,1 mV RS 232 70 x 180 x 75 mm
Int Dir We	(> 4.000 rpm) ereface mensions (W x D x H) eight	0,1 mV RS 232 70 x 180 x 75 mm 0,2 kg
Int Dir We Pe	(> 4.000 rpm) ereface mensions (W x D x H) eight rrmissible ambient temperature	0,1 mV RS 232 70 x 180 x 75 mm 0,2 kg 5 - 40 °C
Int Dir We Pe	(> 4.000 rpm) ereface mensions (W x D x H) eight rmissible ambient temperature rmissible relative humidity	0,1 mV RS 232 70 x 180 x 75 mm 0,2 kg 5 - 40 °C 80 %
Int Dir We Pe Pe	(> 4.000 rpm) erface mensions (W x D x H) eight rmissible ambient temperature rmissible relative humidity otection class acc. to DIN EN 60529	0,1 mV RS 232 70 x 180 x 75 mm 0,2 kg 5 - 40 °C 80 % IP 42

DZM control.o Revolution counter

Connection of an opto-electronic sensor enables measurement of the speed of rotating shafts from 0 - 50.000 rpm. The signals received are displayed on the monitor. This enables all IKA® overhead stirrers and dispersing instruments to be retrofitted with a speed display, thereby allowing reproducible work. A RS 232 interface allows the values to be recorded digital on a PC. An analog output signal for a recorder is also available. The monitor can be used as a table-top device or mounted on a stand rod.

Included with delivery (page): DZM-M Monitor (119), DZM-S.o Sensor optical (119), Power pack

Accessories (page): DZM-K Extension cable (119), labworldsoft® (143)

General data		
Length		

DZM-S.o Sensor optical

Spare sensor for DZM control.o.

1 m

DZM-M Monitor

Spare monitor for the revolution counter DZM control.o.

		DZM-K Extens
General data		
Length	1 m	Enables the sepa
Max. distance between monitor / sensor	2 m	sensor DZM-S.o.

IKA[®] Electronic accessories

Revolution counter (optical) 119



ldent. No.	
8014200	230 V 50/60 Hz
8014201	115 V 50/60 Hz





tension cable

eparation of monitor DZM.M and



Ident, No. 2808900

IKA[®] Electronic accessories

120 Vacuum controller and accessories magnetic / overhead stirrers



ldent. No.	
2300000	230 V 50/60 Hz
2300001	115 V 50/60 Hz

VC 2 IKAVAC® Vacuum controller

Used to create a controlled partial vacuum in laboratory applications. Typical tasks are the evacuation of desiccators, vacuum apparatus, etc. Solvent recovery rates of up to 99% are possible if rotary evaporators are used.

- Microprocessor-controlled
- Minimum solvent loss
- Considerable reductions in water costs
- Integrated air release valve
- Easy operation
- Space-saving stand-supported instrument

- The analog input values for speed and the analog

input for temperature and torque are converted

- Automatic setpoint correction - Clearly organized membrane keyboard
- Accessories (page):

to current

to normal signals

Accessories (page):

Analog cable (148): AK 2.3, AK 2.8

VC 1.1 Water jet pump (121)

-	

Ident. No. 2829300 230 V 50/60 Hz

Technical data Power input 14 W 1 – 1.200 mbar Control range Setting accuracy 1 mbar digital (LED) Display Dimensions (W x D x H) 150 x 57 x 85 mm Weight 1,0 kg Permissible ambient temperature 5 – 40 °C Permissible relative humidity 80 % Protection class acc. to DIN EN 60529 IP 50

Low water consumption.

VC 1.3 Magnetic solenoid valve

solenoid valve.

VC 2.4 Pump control

The pump control is required when using the chemical diaphragm pump MZ-2C or other electrical vacuum pumps, in conjuction with the vacuum controller VC 2. The pump is disconnected from the mains and then reconnected.

Advantage over VC 1.3: Due to the interruption of the pumps current lead, noise levels and energy costs are reduced.

Included with delivery: Magnetic solenoid valve, power pack





IKA[®] Electronic accessories

Vacuum pump / valve 121

VC 1.1 Water jet pump

With valves for water jet and cooling water. Automatic cooling water cut-off at end of distillation. Suitable for rotary evaporators.

> Ident. No. 1980700

In conjunction with the vacuum controller VC 2, the solenoid valve can be used to regulate an inhouse vacuum, the vacuum of uncontrolled water jet pumps or electrical vacuum pumps. The pump works constantly, the pipe is disconnected by the

Ident, No. 2163500











instruments

138 - 139



Anchor stirrer

With PTFE scarper or flow with bornigs, for all laboratory reactors.

Page 132

Flow breaker

Page 132

Laboratory reactors / Rheology

124 LR-2.ST Version 1

LR-2.ST Version 1 (without reactor vessel)



EUROSTAR power control-visc P7 Overhead stirrer, page 37 Ident. No. 2850700

LR 2000.80 Reactor cover, page 132 Ident. No. 2508200

LR 2000.11 Anchor stirrer with flow borings, page 132 Ident. No. 2509500

LR 2000.1 Double-walled reactor vessel, page 133

Laboratory reactor system LR-2.ST

The systems LR-2.ST and LR 2000 are modularly expandable laboratory reactors, designed and planned for reproducing and optimizing chemical reaction processes as well as mixing, dispersing and homogenization processes at laboratory scales.

- Some examples for these processes are: - Manufacturing of cremes, lotions, emulsions, and liposome preparations in the pharmaceutical and cosmetic sector
- Mixing of solids such as calcium carbonate, talc, titanium oxide, etc. into liquid polymers - Mixing of additives and solid polymer
- compounds into mineral oils - Grinding and disintegrating of solids and
- fibers in liquids and polymers

The cost efficient LR-2.ST laboratory reactors are available for vacuum applications.

The laboratory reactors of the series LR 2000 P (pressure) and LR 2000 V (vacuum) are especially designed for the use in the pharmaceutical and cosmetic sector.

The systems can be adapted individually to a wide range of different applications and specific requirements. IKA® laboratory devices, e.g. temperature measuring instruments, laboratory stirrers and dispersing instruments, pumps and thermostats can be combined and controlled via PC using labworldsoft[®]. The torque measuring instruments VK 600 control VISCOKLICK® or VM 600 basic allow for evaluation of rheological properties.

The IKA® laboratory reactors features among others are:

- Modularly expandable to accommodate groand joints)
- discharge valve

Technical data	
Min. volume (anchor stirrer)	500 ml
Min. volume (T 25 digital)	800 ml
Max. volume	2.000 ml
Max. temperature Kalrez	230 °C
Attainable vacuum	25 mbar
Max. viscosity	
Visco module VM 600	150.000 mPas
Speed range	
(EUROSTAR power control-visc P7)	8 – 290 rpm
Height of telescopic stand	620 – 1.010 mm
Dimensions (B x T x H)	460 x 420 x 1.240 mm
Materials in contact with medium	stainl. steel (AISI 316L)
	Kalrez (FFPM)
	borosilicate glass 3.3



IKA[®] Laboratory reactors

Modular and expandable

interchangeable instruments for various applications (3 x NS 29 and 2 x NS 14

- Single- and double-walled jacketed 2 liter vessels available made of borosilicate glass or stainless steel, with or without bottom

- Sealing materials (FFPM) resist solvents and temperatures for applications up to 230 °C

LR-2.ST System variants 126



LR-2.ST Version 1

[1] LR-2.ST

Basic package with reactor cover (sealing material: FFPM) consisting of:

- Stand system LR-2.ST
- LR-2.SI Safety disconnection
- EUROSTAR power control-visc P7
- LR 2000.11 Anchor stirrer
- with flow borings

[2] LR 2000.1

Double-walled reactor vessel, page 133 Ident. No. 2508300

Safety accessory for Version 1 and 2 (page): LR-2.SP Splinter protection (135)

LR-2.ST Version 2

[1] LR-2.ST

Basic package with reactor cover (sealing material: FFPM) consisting of:

- Stand system LR-2.ST
- LR-2.SI Safety disconnection
- EUROSTAR power control-visc P7
- LR 2000.11 Anchor stirrer with flow borings

[2] LR 2000.1

Double-walled reactor vessel, page 133 Ident. No. 2508300

[3] VM 600 basic Visco module, page 135

Ident. No. 8016600

LR-2.ST Version 3

[1] LR-2.ST

Basic package with reactor cover (sealing material: FFPM) consisting of:

- Stand system LR-2.ST
- LR-2.SI Safety disconnection
- EUROSTAR power control-visc P7
- LR 2000.11 Anchor stirrer with flow borings

[4] HBR 4 digital

Heating bath, page 95 Ident. No. 2602300

[5] LR 2.1

Single walled reactor vessel, page 133 Ident. No. 3070000

Configuration possibilities:

Basic package (page 124 – 126)

LR-2.ST Laboratory reactor system

consisting of: - LR-2.ST Stand system - LR-2.SI Safety disconnection - EUROSTAR power control-visc P7 - LR 2000.11 Anchor stirrer

ldent. No. 8016500

LR 2000.2

Hoses

Ident. No. 2509600

Ident. No. 2606700

CC3-308B vpc

Ident. No. 3658800

Circulation thermostat

Reactor vessels (page 133) and accessories (chapter Heating / Tempering)

LR 2000.1 Double-walled reactor vessel borosilicate glass Ident. No. 2508300

LT 5.24 Hose adapter (2 pieces required) Ident. No. 2578100

LT 5.20 Hoses Ident. No. 2606700 LT 5.24 Hose adapter (2 pieces required) Ident. No. 2578100 LT 5.20

CC3-308B vpc Circulation thermostat Ident. No. 3658800

Add-on units

VK 600 control Torque measurement instr., p. 139, Ident. No. 8015700

DTM 12 IKATRON® Digital temperature measuring instr., p. 118, Ident. No. 3113200

Software (page 142 – 146)

labworldsoft[®] PC software Ident. No. 2970000

Accessories (page 132)

LR 2000.10 Anchor stirrer with PTFE scraper Ident. No. 2508400

LR 2000.11 Anchor stirrer with flow borings Ident. No. 2509500

IKA[®] Laboratory reactors

LR-2.ST System variants



128 LR 2000 P System variants (pressure)

LR 2000 P System variant pressure



Please contact IKA® or your local dealer for a detailed quotation.

EUROSTAR power control-visc P7 Overhead stirrer, page 37, Ident. No. 2850700

T 25 digital Disperser can also be attached, page 66 Ident. No. 3565000

S 25 KV – 18 G Appropriate dispersing element, page 71 Ident. No. 2348000

LR 2000.40 Shaft receptacle, page 134, Ident. No. 2509200

LR 2000.85 Reactor cover, page 132, Ident. No. 2598100

LR 2000.11 Anchor stirrer with flow borings, page 132 Ident. No. 2509500

LR 2000.3 Reactor vessel, stainl. steel, page 133, Ident. No. 2509700

LR 2000.75 Stand for pressure variant, Ident. No. 2598000

Technical data	
Min. volume (anchor stirrer)	500 ml
Min. volume (T 25 digital)	800 ml
Max. volume	2.000 ml
Max. temperature FFPM	230 °C
Attainable pressure	6 bar
Max. viscosity	150.000 mPas
Speed range	
(EUROSTAR power control-visc P7)	8 – 290 rpm
Lift of telescopic stand	260 mm
Dimensions (W x D x H)	500 x 500 x 1.350 mm
Weight of basic device	30 kg
Materials in contact with medium	stainl. steel (AISI 316L)
	Kalrez (FFPM)

Configuration possibilities:

Basic components

EUROSTAR power control-visc P7 Overhead stirrer, p. 37 Ident. No. 2850700

LR 2000.75 Stand for pressure variants, p. 128 Ident. No. 2598000

Accessories (page 132)

LR 2000.10 Anchor stirrer with PTFE scraper Ident. No. 2508400

LR 2000.11 Anchor stirrer with flow borings Ident. No. 2509500

Reactor vessels (page 133) and accessories (chapter Heating / Tempering)

LR 2000.3 Double-walled reactor vessel stainl, steel Ident. No. 2509700

LT 5.23 Hose adapter (2 pieces required) Ident. No. 2235000

Double-walled reactor vessel with bottom outlet valve, stainless steel Ident. No. 3064900

LR 2000.4

LT 5.23 Hose adapter (2 pieces required) Ident. No. 2235000

CC3-308B vpc

Ident. No. 3658800

Circulation thermostat

LT 5.20 Hoses Ident. No. 2606700

CC3-308B vpc Circulation thermostat Ident. No. 3658800

LT 5.20

Hoses

Ident. No. 2606700

Add-on units

VK 600 control Torque measurement instrument, p. 139 Ident. No. 8015700

LR 2000 VK Attachment kit for LR 2000 P, p. 135 Ident. No. 2984600

Software (page 142 - 146)

DTM 12 IKATRON® Digital temperature measuring instrument, p. 118, Ident. No. 3113200

labworldsoft[®] PC software Ident. No. 2970000

Necessary components

IKA[®] Laboratory reactors

LR 2000 P System variants (pressure)



LR 2000 V System variants (vacuum) 130

LR 2000 V System variant vacuum



Please contact IKA® or your local dealer for a detailed quotation.

EUROSTAR power control-visc P7 Overhead stirrer, page 37, Ident. No. 2850700

T 25 digital Disperser, can also be attached, page 66, Ident. No. 3565000

S 25 KV – 18 G Appropriate dispersing element, page 71 Ident. No. 2348000

LR 2000.40 Shaft receptacle, page 134, Ident. No. 2509200

LR 2000.80 Reactor cover, page 132, Ident. No. 2508200

LR 2000.11 Anchor stirrer with flow borings, page 132 Ident. No. 2509500

LR 2000.1 Reactor vessel, page 133, Ident. No. 2508300

LR 2000.70 Stand for vacuum variant, Ident. No. 2509000

Technical data	
Min. volume (anchor stirrer)	500 ml
Min. volume (T 25 digital)	800 ml
Max. volume	2.000 ml
Max. temperature FFPM	230 °C
Attainable vacuum	25 mbar
Max. viscosity	150.000 mPas
Speed range	
(EUROSTAR power control-visc P7)	8 – 290 rpm
Lift of telescopic stand	260 mm
Dimensions (W x D x H)	500 x 500 x 1.350 mm
Weight of basic device	30 kg
Materials in contact with medium	stainl. steel (AISI 316L)
	Kalrez (FFPM)
	borosilicate glass 3.3

Configuration possibilities:

Basic components

EUROSTAR power control-visc P7 Overhead stirrer, p. 37 ldent. No. 2850700

LR 2000.70 Stand for vacuum variant, p. 130 Ident. No. 2509000

Accessories (page 132)

LR 2000.10 Anchor stirrer with PTFE scrapers Ident. No. 2508400

LR 2000.11 Anchor stirrer with flow borings Ident. No. 2509500

Reactor vessels (page 133) and accessories (chapter Heating / Tempering)

LR 2000.1 Double-walled reactor vessel borosilicate glass Ident. No. 2508300

Double-walled reactor vessel and bottom discharge valve, borosilicate glass Ident. No. 2509600

Ident. No. 2578100

Ident. No. 2606700

CC3-308B vpc

Ident. No. 3658800

Circulation thermostat

LR 2000.2

LT 5.24

LT 5.20

Hoses

LT 5.24 Hose adapter (2 pieces required) Ident. No. 2578100

LT 5.20 Hoses Ident. No. 2606700

CC3-308B vpc Circulation thermostat Ident. No. 3658800

Add-on units

VK 600 control Torque measurement instr., p. 139, Ident. No. 8015700

Digital temperature measuring instr., p. 118, Ident. No. 3113200

DTM 12 IKATRON®

LR 2000 VK Attachment kit for LR 2000 V, p. 135, Ident. No. 2984600

Software (page 142 – 146)

labworldsoft[®] PC software Ident. No. 2970000

IKA[®] Laboratory reactors

LR 2000 V System variants (vacuum)



T 25 digital ULTRA-TURRAX® Disperser, p. 66 Ident. No. 3565000

VC 2 IKAVAC® Vacuum controller, p. 120 Ident. No. 2300000

IKA® Laboratory reactors

132 Laboratory reactors accessories



Ident. No.	
2508200	LR 2000.80
2598100	LR 2000.85
2498900	LR 2000.54
2661200	LR 2000.57







LR 2000.80 Reactor cover For LR 2000 V (stand LR 2000.70).

For LR 2000 V (stand LR 2000.70). Incl. 3 x NS 29 and 2 x NS 14 / 23 groand joints.

Accessories (page): LR 2000.54 Sealing set (132)

LR 2000.85 Reactor cover (without fig.)

Accessories (page): LR 2000.57 Sealing set (132)

For LR 2000 P (stand LR 2000.75).

LR 2000.54 Sealing set (without fig.)

Spare, for LR 2000 V.

LR 2000.57 Sealing set (without fig.)

Spare, for LR 2000 P.

LR 2000.10 Anchor stirrer

With PTFE scraper, for all laboratory reactors.

LR 2000.11	Anchor stirrer

With flow borings, for all laboratory reactors.

LR 2000.20 Flow breaker

Only for LR 2000 V and LR-2.ST.

LR 2000.21 Flow breaker (without fig.)

Only for LR 2000 P in connection with LR 2000.40 (page 134).

General data	
Material of threaded seal	FFPM

General data	
Material of threaded seal	FFPM

stainl. steel (AISI 316L), PTFE

stainl. steel (AISI 316L)

General data Useful volume 2.000 ml Material stainl. steel (AISI 316L) Max. temperature 230 °C

General data

Useful volume

General data

Height

Max. temperature

Material

D 2000 1 D - - - + -

2.000 ml

230 °C

22 cm

borosilicate glass 3.3

Double-walled, with quick-action connectors, for LR-2.ST and LR 2000 V.

LR 2000.2 Reactor vessel (without fig.)

Double-walled, with quick-action connectors and bottom discharge valve, for LR-2.ST and LR 2000 V.

Accessories (page): LR 2000.53 Stand lower set (133), LT 5.24 Hose adapter (2 pieces required) (99), LT 5.20 Hose (99)

LR 2000.3 Reactor vessel

Double-walled for LR 2000 P (stand LR 2000.75).

LR 2000.4 Reactor vessel (without fig.)

Double-walled with bottom outlet valve, for LR 2000 P (stand LR 2000.75).

Accessories (page): LR 2000.53 Stand lower set (133), LT 5.23 Hose adapter (2 pieces required) (99), LT 5.20 Hose (99)

LR 2000.53 Stand lower set

To raise the laboratory reactor vessels LR 2000.2 and LR 2000.4. Only in connection with LR 2000.70 and LR 2000.75.

General data	
Material	stainl. steel (AISI 316L)
Installation length	180 mm

General data

General data

Material

Material

General data	
Material	stainl. steel (AISI 316L)
Installation length	180 mm

Ident. No.	
2571200	

IKA® Laboratory reactors

Laboratory reactors accessories

LR 2000.1

LR 2000.2

LR 2.1

Ident. No. 2508300

3070000

2509600

LR 2.1 Reactor vessel (without fig.)

Single-walled, for LR-2.ST.

LR 2000.1 Reactor vessel



Ident. No.	
2509700	LR 2000.3
3064900	LR 2000.4



ldent. No. 2509800 133

134 Laboratory reactors accessories



LR 2000.40 Shaft receptacle

To install the dispersing elements S 25 KV (page 71) and the flow breaker LR 2000.21 (page 132).

General data	
Material of seal	FFPM

Ident. No. 2509200



LR 2000.60 Sensor receptacle

To install the temperature sensors PT 100.25 (page 118) and PT 100.5 (page 99).

Ident. No.	
2509300	

LR 2000.30 Vacuum gauge

Only for LR 2000 V. Alternative to the vacuum controller VC 2 IKAVAC® (page 120).

General data	
Vaterial of seal	FFPM
Veasuring range	1 – 1.020 mbar
Measuring accuracy acc. to DIN 16005	class 1

FFPM

60 °C

General data

Material of seal

Max. temperature

Ident, No. 2509400



LR 2000.90 Drip funnel

For dosing, with ground joint NS 29. Only for LR-2.ST and LR 2000 V.

General data	
Volume	250 ml

Ident. No. 2277000

LR 2000.52 Tool set (without fig.)

Ident. No. 2508800

Spare. Included in the packages of the laboratory reactors.

LR 2000.VK Attachment kit (without fig.)

For LR 2000 V and LR 2000 P.

Accessories (page): Torque measurement instrument VK 600 control VISCOKLICK® (139)

VM 600 basic Visco module

LR-2.ST, consisting of adapter kit and

LR-2.SP Splinter protection (without fig.)

hot reactor vessel.

IKA[®] Laboratory reactors

Ident. No.

2984600

Laboratory reactors accessories

Torque measurement instrument for VK 600 control VISCOKLICK® (page 139).



Ident. No. 8016600

Prevents potential injuries caused by broken glass and burns as a result of accidentally touching the

Ident. No. 3326400

Optional components 136

Data processing: software, cable and adapters



IKA[®] Laboratory reactors

Optional components

Dispersing / Homogenizing

IKA[®] Rheology

138 Torque measurement instrument



EUROSTAR power control-visc Overhead stirrer, page 35 Ident. No. 2600000

Boss head clamp, page 116 Ident. No. 2664000

VK 600 control VISCOKLICK[®] Torque measurement instrument, page 139 Ident. No. 8015700

Telescopic stand, page 115 Ident. No. 1412100

Paddle stirrer, page 40 / 41 Ident. No. 0757800

Technical data	
Measuring range	0 – 600 Ncm
Display	digital
Flange-Ø	60 / 62 mm
Flange Height	≥ 10 mm
Linearity of Display:	
0 – 60 Ncm	± 0,5
60 – 600 Ncm	± 1,0
Reproducibility:	
Static	± 0,1 Ncm
Dynamic	± 0,5 Ncm

VK 600 control VISCOKLICK®

Rheological material properties such as viscosity, flow and deformation behavior are among the most important characteristics of any material: - They determine the application-technical manufacturing process of a product The structural composition of a material can be established from its viscosity behavior The sequence of chemical reactions can be

documented

The VK 600 control can be combined with all IKA® EUROSTAR overhead stirrers. The appropriate stirrer is "clicked" into the VK 600 control. During stirring, a force transducer determines a reaction force at the stirring shaft proportional to the torque.

- Simple assembly

- RS 232 interface and analog output

Accessories (page): VK 60 / 01 Adapter (139), labworldsoft® (143)

VK 60 / 01 Adapter

For adaption of IKA® overhead stirrer RW 20 digital.

IKA[®] Rheology Torque measurement instrument

Torque measurement instrument

- PC-controllable with labworldsoft® - Measuring system is overload-proof - Offset correction to eliminate errors

> Ident. No. 8015700 230 V 50/60 Hz 115 V 50/60 Hz 8015701



Ident. No. 2854100





$labworldsoft^{\scriptscriptstyle{(\!\!R)}}$

Eases life in the laboratory. With this laboratory software, you can network up to 64 laboratory devices simultaneously via one PC. That makes the automation of your laboratory experiments and processes possible. Page 142 / 143



Software

141

IKA[®] Software labworldsoft[®]

142



labworldsoft[®]

With this laboratory software, you can network up to 64 laboratory devices simultaneously via one PC.

labworldsoft[®]

With this laboratory software, you can network up to 64 laboratory devices simultaneously via one PC. That makes the automation of your laboratory experiments and processes possible.

Measurements and processes may be run independently from one another. This helps to avoid long waits and you increase your productivity. The communication between PC and laboratory device is performed via the serial interface RS 232 (COM1 or COM2).

With the help of plug-in cards and Ethernet RS 232 servers, up to 64 laboratory devices can be used simultaneously via one PC. All laboratory instruments can be controlled independently from each other and the measured values (speed, temperature, torque, pH, etc.) can be documented separately.

Hard- and software requirements:

Pentium 90 with at least 16 MB RAM, and a mouse. VGA display: monochrome with at least 16 levels of grey or color. Windows 95/98/2000/NT/ME/XP/Vista.

Accessories (page):

PCI 8.2 Plug-in card (147), PC 4.1 RS 232 Server (147), DC 2 DATACONTROL (147), DA 2 DATACONT-ROL (147), IO 2 DATACONTROL (148)

Networking, monitoring automated using labworldsoft®.

Controlling

Desired temperature and speed sequences can be precisely controlled by means of freely selectable ramp functions. The ramp functions can be graphically generated, stored, and then loaded again at any time.

Recording, evaluating

labworldsoft® enables a fast and easy recording of many physical parameters which are required in the laboratory, such as pH, conductivity, temperature, torque, weight, pump rates etc.

Exporting

Data recorded using labworldsoft® can be directly written to an Excel sheet or exported to any standard application at a later stage.

Storing / reproducing measured data

Do your test arrangements repeat themselves? With labworldsoft® all test arrangements can be stored. The stored data is available to reproduce the test, at a mouse click. The reproductability of tests is warranted within the scope of ISO 9000 and within GLP.

Documentation

For documentation purposes, all measuring results as well as the measurement flowcharts can be printed or plotted according to GLP, ISO and QA.

For more information and a download of your free trial version please visit: www.labworldsoft.com

IKA[®] Software

labworldsoft[®]

143

With labworldsoft[®] you can network up to 64 laboratory instruments simultaneously via one PC. From sample preparation to synthesis, all steps of research and development in the lab can be



Ident. No. 2970000

(**0**) labworld*soft* °

Manufactures with interface devices compatible to labworldsoft®:

- Ahlborn
- B. Braun Biotech
- Martin Christ
- Corning Inc.
- Ehret
- Eyela
- Fluid
- Fritsch
- Gerhardt
- GEL
- Harvard
- Heidolph
- Hermle
- Huber
- IKA®
- Ilmvac
- Infors
- Ismatec

- Kern
- KNE
- Knick
- Labovisco
- Lauda
- Metrohm
- Mettler-Toledo
- MLT
- PolyScience
- Sartorius
- Scaltec
- Sigma
- Telab
- Thermo Haake
- Thermo Neslab
- Troemner
- Vaccubrand
- yellowline

Interfaces to additional devices from other manufacturers will soon be available. Please ask for a current reference list.

- Julabo
IKA[®] Software

144 labworldsoft®



Presentation of results

The measuring results are directly displayed online or offline graphically with a selectable coordination system or numerically. Several numerical displays as well as four-channel displays are possible.

Storing a measuring configuration

The complete measuring configuration with all current parameters and the position of all opened windows can be stored. As a result, preconfigured flowcharts which are immediately ready for operation can be provided for the widest variety of tasks.

Fig. 1: Configuration example of a laboratory reactor with peripherals in operation. The speed of an overhead stirrer, the target temperature of a thermostat and a pump are controlled. Torque and temperature of the medium are recorded and are represented in a y/ t-graphic (fig. 2). By means of a IO 2 DATACONTROL, additional external sensors or valves are possible. Configuration example - Recording rheological data during the stirring process

labworldsoft®

Laboratory software for control and data collection, page 143 Ident. No. 2970000

EUROSTAR power control-visc

Stirrer, page 35 Ident. No. 2600000

R 270 Boss head clamp, page 116 Ident. No. 2657800

VK 600 control VISCOKLICK®

Torque measurement instrument, page 139 Ident. No. 8015700

PC 1.5 Cable, page 148 Ident. No. 2756000 _

R 1373 Paddle stirrer, page 40 Ident. No. 0757600

RH 5

Strap clamp for securing the vessel, incl. boss head clamp R 270, page 116 Ident. No. 3159000

R 2723 Telescopic stand, page 115 Ident. No. 1412100

PCI 8.2 Plug-in card

for mounting in the PC to control up to 8 instruments, page 147 Ident. No. 8017500





Figure 2: y/t-graphic: Shows torque and temperature changes in medium.



145

labworldsoft[®] 146

Configuration example - Controlling and recording temperature data during magnetic stirring with heating



Laboratory software for control and data collection,

Holding rod for casing of the PT 100.50 sensor, page 27

Boss head clamp, page 116

Temperature sensor for RET control-visc, page 26

Support rod for attachment to

Plug-in card for mounting in the PC to control up to 8 instruments, page 147

Protection cover included with the RET control-visc,

RET control-visc safety control Safety magnetic stirrer with RS 232 interface, page 14, incl. protection cover H 99, page 28

PCI 8.2 Plug-in card

For mounting in the PC to connect up to 8 instruments simultaneously. Plug-in cards for up to 64 instruments available on request.

PC 4.1 RS 232 Server

Up to 4 lab units can be controlled through the ethernet with the PC 4.1 RS 232 server. The server supports 4 RS 232 ports with a 10/100 mbps ethernet interface by TCP/IP. The server can be set-up through the ethernet and works as a transparent serial COM-Port without restrictions of platform and distance. Server for connection of up to 64 instruments available on request.

DC 2 DATACONTROL

For PC documentation of analog signals from up to 4 instruments.

Accessories (page): PC 1.5 Cable (148), PC 2.2 Adapter (148), AK 2.4 Analog cable (148)

DA 2 DATACONTROL

To convert digital signals into analog signals. In this manner, devices with analog control inputs (industrial controllers, temperature controllers) can be controlled using labworldsoft®. Connection box included.

Accessories (page): PC 1.5 Cable (148), PC 2.2 Adapter (148), Analog cable (148): AK 2.6, AK 2.7

General data 0 - 1 / 0 - 5 / 0 - 10 V Voltage input Current input 0 – 20 / 4 – 20 mA

0 – 20 / 4 – 20 mA

General data Voltage output 0 - 1 / 0 - 5 / 0 - 10 V Current output

IKA[®] Software labworldsoft[®] accessories



Ident. No. 3192000



Ident. No. 8015600 230 V 50/60 Hz 8015601 115 V 50/60 Hz



ident. No.	
8017200	230 V 50/60 Hz
8017201	115 V 50/60 Hz

IKA[®] Software

labworldsoft[®] accessories 148



ldent. No.	
3006000	230 V 50/60 Hz
3006001	115 V 50/60 Hz

IO 2 DATACONTROL

With the IO 2 DATACONTROL, the power switch IO 2.1 DATACONTROL and labworldsoft® any device without any interface (heaters, solenoid valves, etc.) can be turned on and off based on an event (a threshold value being exceeded, controller output, etc).

This opens up numerous control possibilities in connection with the PID, relay and trigger modules of labworldsoft[®]. In addition, using the 8 inputs on the IO 2 DATACONTROL, signals from switches etc. can be recorded by labworldsoft[®].

Accessories (page):

IO 2.1 Power switch (148), PC 1.5 Cable (148), PC 2.2 Adapter (148)

IO 2.1 DATACONTROL Power switch

Technical data	
Max. power of the	
connected devices	1,2 kW
Cable length	0,6 m

30 V / 1 A

0 – 24 V

Purpose	Digital control of up to 4 instruments and digital recording of measurement data via PC over ethernet	Digital control of up to instruments and digita recording of measurer data via PC with interfa plug-in card
IKA® instruments with interface	Magnetic stirrers RET control-visc <i>safety control</i> (p. 14 Overhead stirrers EUROSTAR power control-visc (p. 38 Shakers KS / HS control (p. 47 – 49) Revolution counter DZM control (p. 119) Torque measuring instrument VK 600 VISCOKLICK (S.	
Device interface		1!
Required components		[
	15 pin HD Sub-D (M)	
	Adapter PC 1.4	
	9 pin Sub-D (F)	
	9 pin Sub-D (M)	
	Cable PC 2.1	
	9 pin Sub-D (F)	
		25 pin Sub-D (M), 8 x
		Plug-in card PCI 8
		Installation in PC
Interface PC / recorder	9 pin Sub-D (M)	Installation in PC

PC 4.1 RS 232 Server with RJ 45 network connection

PC with unused slot

Ident, No. 3062000 230 V 50/60 Hz 3062001 115 V 50/60 Hz

Cable and Adapter (without fig.)

Cable	Length	Ident. No.
PC 1.1	3 m	2616700
PC 1.5	2,5 m	2756000
PC 2.1	5 m	2700700
PC 2.3	3 m	3036200
DTM 12.10	2,5 m	3127800
Adapter		
PC 1.2		2616800
PC 1.4		2755900
PC 2.2		2753200
PC 5.1		2621500
Analog cable		
AK 2.1	2,5 m	2734300
AK 2.2	2 m	2756100
AK 2.3	2 m	2801200
AK 2.4	2 m	2801300
AK 2.5	2 m	2845800
AK 2.6 (blue)	1,5 m	1719400
AK 2.7 (red)	1,5 m	1719300
AK 2.8	1,8 m	2907800

,	
	Technical data
	Technical data
	Max. power of the
	connected devices
	Cable length
	EURO connector (other connectors available on request)

Technical data

8 digital outputs

(relay contact)

8 digital inputs

(Voltage)

IKA[®] Software

Overview connection possibilities



IKA[®] Software

150 Overview connection possibilities



IKA[®] Software

Overview connection possibilities



C 14

The disposable crucible makes handling much easier because there is no longer any need for a quartz or stainless steel crucible. Optimises sample combustion. No

crucible to clean. Direct contact with ignition wire. No ignition thread required. Page 163



Calorimeters 154 – 163 Decomposition system 164 – 165

IKA® Analytical line

154 Calorimeters



C 5000 The calorimeter offers three user-selected operating modes.

Technical data		
Input power max.		120 W
Rated voltage		24 V DC, 5 A
Fuse		1 x 2.5 AT
Max. On-time	con	tinuous operation
Range of measurement		40.000 J
Measuring mode /	isoperibol	up to 17 min
Measuring time	dynamic	up to 8 min
	manuel (isoperibol)	up to 17 min
	time-controlled	up to 14 min
Reproducibility		
based on analysis of 1 g		
benzoic acid NBS 39i		0,1 % RSD
Operating oxygen pressure		30 bar
General data		
Dimensions (W x D x H)	400	0 x 400 x 400 mm
Weight		30 kg
Protection class		111
Interfaces	1	l x serial (RS 232)
	1 x pa	rallel (Centronics)
Ambient temperature	20 -	- 25 °C (constant)
Ambient humidity		80 %
Protection class according t	to DIN EN 60529	IP 21

C 200

New, compact low cost combustion calorimeter to determine calorific values of liquid and solid samples. Designed for teaching and training applications for use in schools, technical colleges, universities and practical training courses. Also suitable for laboratories with small numbers of analyses.

- In the manual mode (learning mode) the user triggers ignition and the end of measurement. The temperature changes are recorded at minute intervals. All calculations are manual.
- In the other operating modes ignition and calculation of calorific values are automatic. The calorific value is shown on the display. Acid correction of the calorific value and calculation of the heat values are performed manually.
- The C 5010 decomposition vessel can be used with a C 14 disposable crucible.
- The C 200 can also be operated with the "CalWin C 5040" calorimeter software. This enables control of up to eight C 200 measurement cells from a PC.

Functions:

- time-controlled

- GOST-certified
- Automatic sample ignition
- Compact modular design

Consisting of:

C 248 Oxygen station

IKA[®] Analytical line

Calorimeter C 200



Ident, No. 8802500 100 - 240 V 50/60 Hz

- Working methods: isoperibol, manual, dynamic,

- Validation according to DIN 51900, ISO 1928, ASTM D240, ASTM D4809, ASTM D5865, ASTM D1989, ASTM D5468, ASTM E711

- Highly operator maintenance friendly - Complies with all global voltages, from 100 - 240 V - Powered with a low operating voltage 24 V DC

Basic device C 200 incl. power pack and ignition adaptor, C 5010 Decompostion vessel standard,

NEW!

IKA[®] Analytical line

156 Calorimeters C 2000



C 2000 basic, C 2000 control, C 2000 basic high pressure und C 2000 control high pressure

The C 2000 basic and C 2000 control calorimeters are the tried-and-tested systems from IKA® for determining gross calorific values of liquid and solid samples.

A high level of automation with extremely simple handling characterizes these instruments. In addition to the isoperibolic measurement procedure (static jacket), a dynamic (reduced-time) working method is also available. Halogen resistant decomposition vessels of the C 5012 series for quantitative decomposition of sulfur and halogens in parallel to determining gross calorific values are available.

To provide the calorimeters with cooling water, they need to be connected to a thermostat like the KV 600 (page 131) or a firmly installed water connection.

The C 2000 basic is equipped with a very convenient console to operate the unit. The C 2000 control is delivered with the proven C 5040 Cal-Win calorimeter software in order to control the system via PC. Network connection and special configuration for data exchange with LIMS can be implemented at any time.

The C 2000 high pressure is a combination of the C 2000 basic / C 2000 control and the C 62 digestion container (up to 1200 bar operating pressure), see page 162.

C 2000 basic Version 1

Ident. No. Version 1 8801800 230 V 50/60 Hz 115 V 50/60 Hz 8801801 230 V 50/60 Hz Version 2 8801900 8801901 115 V 50/60 Hz high pressure 8802300 230 V 50/60 Hz 115 V 50/60 Hz 8802301

Consisting of:

C 2000 basic C 5010 Decomposition vessel, standard

C 2000 basic Version 2

Consisting of: C 2000 basic C 5012 Decomposition vessel, halogen resistant

C 2000 basic high pressure

Consisting of: C 2000 basic C 62 Decomposition vessel, high pressure C 60 Conversion set

Functions:

- Automatic water handling system includes tempe-

- ring, filling and emptying of calorimeter inner vessel
- Automatic oxygen filling of decomposition vessel - Automatic decomposition vessel identification
- Automatic sample ignition
- Validation according to DIN 51900,
- ASTM 240 D, ISO 1928, BSI etc.
- GOST-certified
- Working methods:

isoperibol, measurement time: approx. 22 min dynamic, measurement time: approx. 7 min

- Compact, integrated modular design for convenient operation
- Cooling water supply via thermostat (KV 600) or firmly installed water connection (C 25 pressure regulating valve recommended, page 162)
- Interface connections for each of the following: scale, printer, monitor and sample rack C 5020
- User-friendly software C 5040 CalWin for controlling the calorimeter and administration of measuring data (page 161)
- LIMS integration is possible
- Special halogen resistent vessel for quantitative decomposition of halogens and sulfur
- The decomposition vessel can be changed over to use disposable crucible C 14 burns during measuring (page 162/163)
- Up to 8 calorimeters can be controlled by a single PC, using a multi-serial plug-in card

Technical data		
Input power max.		1,8 kW
Power ON-time	continu	ous operation
Range of measurement		40.000 J
Reproducibility		
based on analysis of 1 g	isoperibol	0,05 % RSD
benzoic acid NBS 39i	dynamic	0,1 % RSD
Working modes / Start temperature	isoperibol	25 °C
	isoperibol	30 °C
	dynamic	25 °C
	dynamic	30 °C
Measurement time	isoperibol	up to 22 min
	dynamic	up to 7 min
Operating oxygen pressure		30 bar
Cooling medium		tap water
Min. flow rate		60 l/h
Operated with KV 600		
Pressure		0,3 bar
Temperature		
(depending on working mode)		18/25 °C
Operated at firmly installed water con	inection	
Pressure after C 25 pressure regulating va	alve	1 – 1,5 bar
Temperature		
(depending on working mode)		12 – 28 °C
Max. pressure at the tap		6 bar
General Data		
Dimensions (WxDxH)	440 x 4	50 x 500 mm
Weight		46 kg
Ambient temperature	20 – 25	°C (constant)
Ambient humidity		80 %
Protection class according to DIN EN 60	529	IP 21

C 2000 control Version 1

Consisting of: C 2000 control

C 2000 control Version 2

Consisting of: C 2000 control C 5012 Decomposition vessel, halogen resistant C 5040 CalWin, calorimeter software

C 2000 control high pressure

Consisting of: C 2000 control C 60 Conversion set

°C

C 2000 Ausbaugerät

Consisting of: without decomposition vessel), C 5041.10 Connection cable (for 8 x interface box)





IKA[®] Analytical line

Calorimeters C 2000

C 5010 Decomposition vessel, standard C 5040 CalWin, calorimeter software

C 62 Decomposition vessel, high pressure

A PC is required to operate the C 2000 control.

C 2000 control (without calorimeter software,

Ident. No. 8802200 230 V 50/60 Hz 8802201 115 V 50/60 Hz







Ident, No

8802000

8802001

8802100

8802101

8802401

high pressure 8802400

Version 1

Version 2

230 V 50/60 Hz

115 V 50/60 Hz

230 V 50/60 Hz

115 V 50/60 Hz

230 V 50/60 Hz

115 V 50/60 Hz

IKA[®] Analytical line

Calorimeters C 5000 158



C 5000 control und C 5000 duocontrol

The IKA® calorimeter C 5000 is the only calorimeter in the world that offers a free selection of 3 working methods: Thus it is possible to perform determinations of gross calorific values of liquid and solid samples in adiabatic (approx. 14 - 18 min), isoperibolic (approx. 22 min) and dynamic (reduced time: approx. 10 min) mode.

A high level of automation in addition to an extensive range of accessories leave nothing more to wish for.

Functions:

- Automatic water handling system includes tempering, filling and emptying of calorimeter inner vessel
- Automatic oxygen filling and degassing of the decomposition vessel
- Automatic decomposition vessel identification
- Validation according to DIN 51900, ASTM 240 D,
- ISO 1928, BSI etc.
- GOST-certified
- Interface connections for each of the following:
- scale, printer, monitor and sample rack C 5020 - User-friendly software C 5040 CalWin for controlling
- the calorimeter and administration of measuring data (page 161)
- LIMS integration is possible
- Special halogen resistent vessel for quantitative decomposition of halogens and sulfur
- The decomposition vessel can be changed over to use disposable crucible C 14 burns during measuring (page 162 / 163)
- The C 5000 control can be expanded to the duocontrol system with two measurement cells at any time

Technical data		
Input power max.		
(with one measuring cell)		1,3 kW
Power ON-time	continu	ous operation
Range of measurement		40.000 J
Reproducibility	adiaba	tic / isoperibol
based on analysis of 1 g		0,05 % RSD
benzoic acid NBS 39i	dynamic	0,1 % RSD
Working modes		adiabatic
		isoperibol
		dynamic
Measurement time	adiabatic	up to 15 min
	isoperibol	up to 22 min
	dynamic	up to 10 min
Operating oxygen pressure		30 bar
Cooling medium (C 5004)		tap water
Flow rate		18–42 l/h
Operated (C 5004) with KV 600		
Temperature		15 – 20 °C
Operated at firmly installed water co	onnection	
Min. / max. temperature		10 / 19 °C
Max. pressure at the tap		9 bar
General Data		
Dimensions (W x D x H)		
C 5000 control Package 1	740 x 3	380 x 400 mm
C 5000 control Package 2	560 x 3	380 x 400 mm
C 5000 duocontrol Package 3	1.440 x 3	380 x 400 mm
Weight Package 1		61 kg
Ambient temperature	20 – 25	°C (constant)
Ambient humidity		80 %
Protection class according to DIN EN 6	60529	IP 21

C 5000 control Package 2/10

Cooling water supply via thermostat KV 600 (page 147) or firmly installed water connection.

Consisting of: C 5000 Controller C 5003 Measurement cell C 5010 Decomposition vessel, standard C 5004 Heat exchanger

C 5000 control Package 2/12

Cooling water supply via thermostat KV 600 (page 147) or firmly installed water connection.

Consisting of:

C 5000 Controller C 5003 Measurement cell C 5004 Heat exchanger

C 5000 duocontrol Package 3/10

System with two measurement cells. nm nm

Consisting of:

C 5000 Controller nt) C 5003 Measurement cell (2 pieces) C 5002 Cooling system

C 5000 duocontrol Package 3/12

System with two measurement cells.

Consisting of:

C 5000 Controller C 5003 Measurement cell (2 pieces) (2 pieces) C 5002 Cooling system

C 5000 control Package 1/10

Package 1/10 8801000 230 V 50/60 Hz 8801001 115 V 50/60 Hz Package 1/12 8801500 230 V 50/60 Hz 115 V 50/60 Hz 8801501

Ident No

Consisting of: C 5000 Controller C 5003 Measurement cell C 5010 Decomposition vessel, standard C 5001 Cooling system

C 5000 control Package 1/12

Consisting of:

C 5000 Controller C 5003 Measurement cell C 5012 Decomposition vessel, halogen resistant C 5001 Cooling system

IKA[®] Analytical line

Calorimeters C 5000

C 5012 Decomposition vessel, halogen resistant

C 5010 Decomposition vessel, standard (2 pieces)

C 5012 Decomposition vessel, halogen resistant



	Ident. No.	
Package 2/10	8801200	230 V 50/60 Hz
	8801201	115 V 50/60 Hz
Package 2/12	8801600	230 V 50/60 Hz
	8801601	115 V 50/60 Hz



	Ident. No.	
Package 3/10	8801100	230 V 50/60 Hz
	8801101	115 V 50/60 Hz
Package 3/12	8801700	230 V 50/60 Hz
	8801701	115 V 50/60 Hz

IKA[®] Analytical line

Calorimeters C 7000 160

C 7000

The C 7000 is the first IKA® calorimeter with a completely dry system for measuring the gross calorific value of solid and liquid samples. The temperature is measured directly in the decomposition system. This results in measurement times in the range of 3 to 7 minutes (depending on the sample). The system can manage up to 8 different decomposition vessels using a code ring scheme.



Functions:

- High sample frequency
- Precise and reproducible determination of gross calorific values according to ISO 1928
- Reduction of routine work through automatic application flow
- Automatic decomposition vessel identification
- Interface connections for scale, printer and PC
- User-friendly software C 5040 CalWin for controlling the calorimeter and administration of measuring data (page 161)
- Special halogen resistent vessel for quantitative decomposition of halogens and sulfur
- The decomposition vessel can be changed over to use combustible crucibles C 14 burns during measuring (page 163, 165)

Ident. No. 8800900 230 V 50/60 Hz 115 V 50/60 Hz 8800901

ldent. No.	
8801400	230 V 50/60 Hz
8801401	115 V 50/60 Hz

C 7000 basic equipment set 1

Consisting of:

C 48 Oxygen station

C 7000 Measurement cell C 7010 Decomposition vessel, standard C 7002 Cooling system C 48 Oxygen station

C 7000 basic equipment set 2

Consisting of: C 7000 Measurement cell C 7012 Decomposition vessel, halogen resistant C 7002 Cooling system

Technical data	
Input power max.	0,1 kW
Power ON-time	continuous operation
Range of measurement	30.000 J
Reproducibility	
based on analysis of 1 g	
benzoic acid NBS 39i NBS 39i	0,2 % RSD
Working modes	patented
	double dry
Measurement time	3 – 7 min
Operating oxygen pressure	30 bar
Cooling medium (C 7002)	tap water
Flow rate (C 7002)	2 – 3 l/h
Temperature	12 – 30 °C
	(cooling water)
Operated at firmly installed water conne	ection
Max. pressure at the tap	9 bar
General Data	
Dimensions (WxDxH)	310 x 490 x 460 mm
Weight	43 kg
Ambient temperature	18 – 30 °C (constant)
Ambient humidity	80 %
Protection class according to DIN EN 6052	29 IP 21

C 5040 CalWin

- of calibrations
- of measurements
- protocols suitable for certification
- Library functions

- CEL (configarable by the user)

C 5041 CalWin plus

different type.

Consisting of: C 5040 CalWin, PCI 8.2 PC Plug-in card (internal), Interface box (8x)

IKA[®] Analytical line

Calorimeters software

CalWin



- Control, monitor and view working procedures - Print and save measurement protocols - Identify and record samples - Administration of sample racks - Flexible administration and evaluation

- Flexible aministration and grouping

- Printing and saving calibration and result

- Data transmission via RS232 interface to Microsoft® EXCEL and Microsoft® Access applications - Preprocessed work sheets for Microsoft® EX-

To control up to 8 calorimeters of the same or







Ident. No. 3166000

IKA® Analytical line

162 Calorimeters accessories

Calorimeters accessories

fas	- ^	2	nn	
101	ີບ	2	JU	

for C 200	Ident. No.
C 5010 Decomposition vessel, standard	7114000
C 5010.4 Attachment for combustible crucible C 14 (for C 5010 / C 5012)	3016900
C 5010.5 Crucible holder, big (for C 5010 / C 5012)	3055900
C 5030 Venting station (for C 5010 / C 5012) with gas wash bottle acc. to DIN 12596 (for gas absorption) $_$	7198000
C 5041.10 Connection cable	3036000
C 21 Pelleting press	1605300
C 29 Pressure gauge, oxygen	0750200
C 248 Oxygen station	3520000
C 200.1 Measuring cup 2.000 ml	3548900

for C 2000	Ident. No.
C 5010 Decomposition vessel, standard	7114000
C 5012 Decomposition vessel, halogen resistant	7215000
C 62 Decomposition vessel, "high pressure"	3265000
C 60 Conversion set for C 62	3187400
C 5010.4 Attachment for combustible crucible C 14 (for C 5010 / C 5012)	3016900
C 5010.5 Crucible holder, big (for C 5010 / C 5012)	3055900
C 5030 Venting station (for C 5010 / C 5012) with gas wash bottle acc. to DIN 12596 (for gas absorption)	7198000
C 5020 Sample rack	7145000
KV 600 Cooling water supply (230 V)	_ 3410500
KV 600 Cooling water supply (115 V)	3410501
C 25 Pressure regulating valve to operate with firmly installed water connection	3197200
C 5041.10 Connection cable	3036000
C 21 Pelleting press	1605300
C 29 Pressure gauge, oxygen	0750200

for	- n	Б	٦N	n
101	U.	J	JU	u

10r C 5000	Ident. No.
C 5010 Decomposition vessel, standard	7114000
C 5012 Decomposition vessel, halogen resistant	7215000
C 5010.4 Attachment for combustible crucible C 14 (for C 5010 / C 5012)	3016900
C 5010.5 Crucible holder, big (for C 5010 / C 5012)	3055900
C 5030 Venting station (for C 5010 / C 5012) with gas wash bottle acc. to DIN 12596 (for gas absorption)	7198000
C 5020 Sample rack	7145000
KV 600 Cooling water supply (230 V)	3410500
KV 600 Cooling water supply (115 V)	3410501
C 5041.10 Connection cable	3036000
C 21 Pelleting press	1605300
C 29 Pressure gauge, oxygen	0750200

for C 7000	Ident. No.
C 7010 Decomposition vessel, standard	3015000
C 7012 Decomposition vessel, halogen resistant	3017000
C 7010.8 Venting handle (for C 7010 / C 7012)	7095000
C 7030 Venting station (for C 7010 / C 7012) with gas wash bottle acc. to DIN 12596 (for gas absorption) $_$	3013300
C 5041.10 Connection cable	3036000
C 7002 Cooling system (230 V)	7011000
C 7002 Cooling system (115 V)	7011001
KV 600 Cooling water supply (230 V)	3410500
KV 600 Cooling water supply (115 V)	3410501

Calorimeters accessories

for C 7000	
C 21 Pelleting press	
C 29 Pressure gauge, oxygen	
C 48 Oxygen station	

Instructions on the IKA® calorimeter systems available on request.

Consumables for all Calorimeters

C 5003.1 Aqua Pro	stabilizing agent (20 ml)
C 710.4 Cotton thre	ead, cut to length (500 pieces)
C 5010.3 Ignition w	/ire, spare (5 pieces)
C 5012.3 Ignition w	vire, platinum (2 pieces)
C 4 Quartz dish	
C 5 Set of VA comb	oustion crucibles (25 pieces)
C 6 Quartz crucible	, big
C 710.2 Set of VA of	combustion crucibles (25 pieces)
C 9 Gelatine capsul	les (100 pieces)
C 10 Acetobutyrate	e capsules (100 pieces)
C 12 Combusion ba	ags 40 x 35 mm (100 pieces)
C 12 A Combusion	bags 70 x 40 mm (100 pieces)
C 14 Combustible of	crucible (100 pieces)
C 15 Paraffin strips	(600 pieces)
C 43 Benzoic acid N	VIST 39i (30 g)
C 723 Benzoic acid	, blister package (50 pieces)
AOD 1.11 Control s	standard for sulfur and chlorine (50 ml)
AOD 1.12 Control s	standard for fluorine and bromine (50 ml)
C 62.1 Set of wear	ing parts (for C 2000 high pressure)
C 62.2 Combustion	crucibles for C 62 (for C 2000 high pressure
	ion wire (for C 2000 bigh pressure) (200 m goil

IKA® Analytical line

Calorimeters accessories / consumables

ldent. No.
1605300
0750200
1560000

Ident. I	٧o.
72077	00
14837	'00
71228	00
29949	00
16955	00
17495	00
03551	00
14835	00
07499	00
07500	00
22014	-00
22015	00
72245	00
31311	00
07506	00
32430	00
30440	00
30802	00
32963	00
32660	00
07496	00

IKA® Analytical line

164 Decomposition system



Protective device AOD 1.3 As per Pressure Vessel Directive 97 / 23 / EC (not included with the delivery), page 165, Ident. No. 3308000

Oxygen filling station C 48 For filling decomposition vessel with oxygen, 30 bar, page 163 Ident. No. 1560000

Venting station C 7030

With gas with DIN 12596 gas wash bottle, for gas absorption (not included with the delivery), page 165 Ident. No. 3013300

Control standard AOD 1.11 (without fig.)

For sulfur and chlorine, page 165 Ident. No. 30440000

Decomposition vessel AOD 1.1

High-alloy, halogen-resistant stainless steel, page 165 Ident. No. 3303000

External ignition unit AOD 1.2

Ignition triggered by pressing the Ignite button Cable length: 5 m, page 165 Ident. No. 3348000

AOD 1 Decomposition system

Consisting of:

- AOD 1.1 Decomposition vessel,
- C 48 Oxygen station,
- AOD 1.2 External ignition unit, AOD 1.11 Control standard (50 ml)
- Oxidative decomposition of solid and liquid organic samples under pressure in a closed system
- Quantitative decomposition of all halogens,
- sulfur, as well as volatile metals, e.g. As and Hg - Absorption of the combustion products in an aqueous medium
- Catalytic support of the oxidation process with

Ident. No. 8801300

auto-regenerating catalytic surfaces - Pressure vessel of high-grade stainless steel

- Decomposition temperature up to 1.200 °C
- Max. operating pressure during decomposition 195 bar
- Decomposition time < 3 min
- The decomposition vessel can be changed over to use combustible crucibles C 14 (page 162 / 163)
- Control standards for Cl, S, F and Br
- Introduction of the combustion gases into the absorption solution via venting station C 7030

Technical data	
Decomposition time	< 3 min
Core temperature	> 1.200 °C
Max. operating temperature	50 °C
Max. operating pressure	195 bar
Volumen of decomposition vessel	210 ml
Oxygen pressure	30 bar

AOD 1.3 Protective device

For use with decomposition vessel AOD 1.1 operated in accordance with Pressure Vessel Directive 97/23/ EC. If the unit is used improperly (e.g. use of unknown explosive substances or high energy overloads) or if the decomposition vessel is worn, bursting can not totally excluded. In this case the protective device protects the user from inquiry.

C 7030 Venting station

vessels AOD 1.1, C 7010 and C 7012.

Decomposition system accessories

AOD 1.1 Control standard
AOD 1.2 External ignition unit
AOD 1.13 Remote ignition head (required where AOD 1.3 is not used)
C 21 Pelleting press
C 29 Pressure gauge, oxygen
C 5010.4 Attachment for combustible crucible. C 14

Decomposition system consumables

IKA[®] Analytical line

Decomposition system accessories and consumables



Ident. No. 3308000

The controls venting of the combustion gases after decomposition. Complete with DIN 12596 gas wash bottle. For use with decomposition



Ident No 3013300

ldent. No.
3303000
3348000
3348100
1605300
0750200
3016900

ldent. No.
1695500
0749900
0750000
2201400
2201500
7224500
3131100
2994900
1483700
3044000
3080200
3243000



IKA[®] General

IKA®-Werke GmbH & Co. KG		Name				Dev
Janke & Kunkel-Str.10		Company				All Ik
D-79219 Staufen		Department				acco to th
		Street				mark appli
Fax: +49 7633 831-98		City / State / Zip				Envir
Plassa condivis a fay or mail in window onvolono						
						Pat
		Phone				Fale
		Fax				Certa right
		E-Mail				withi infor
Type of processing	Mixing	Dissolving	Emulsifying Wet crushing			Gua
	Homogenizing	ouspending	Wet clushing			They
Volume / Quantity	Discontinuous	l/batch				for o
	Continuous	I/h				1 yea
Viscosity		mPas (20 °C)				_
						Сор
Flow behaviour similar to	Water	Motor oil	Honey			Copy refer
Composition	Liquid	%	Material			form
composition	Solid	%	Material			
	Particle size initial	mm	After end of process		µm	
	pH range		Temperature range		°C	
	Vacuum range	mbar	Pressure range		mbar	
Container dimensions	Diameter mm	Total height	mm Filling h	eight	mm	
N. 1. (5						
Voltage / Frequency	V		Hz			
Ex-proof	no	yes, Ex-class				
Bemarks						

Device safety, environment

All IKA® laboratory devices satisfy the international legal regulations according to DIN EN IEC 61010. Any instrument is safety tested according to this norm before it leaves IKA®. Instruments designed for the European market are labeled with the CE mark, to state that they satisfy the applicable EU regulations and norms.

Environmental factors were especially taken into consideration when materials were selected (CFC-free and cadmium-free products).

Patents

Certain products featured in the catalog have been assigned property rights such as patents, trademarks, etc. These property rights only apply within the Federal Republic of Germany. On request, we will gladly provide information with regard to their validity in other countries.

Guarantee, Warranty

The warranty satisfies the relevant legal regulations. The guarantee period for our products is 2 years, for analyzing technology products the period is 1 year.

Copyright

Copying for commercial purposes is expressly permitted. We refer to the copyright with regard to tables, catalog design and formulations. Documentary evidence of used catalog pages is desired.



Illustrations

The glass vessels and containers shown in the photos together with the instruments are generally not included in the product package.

Voltage / Frequency / Plugs

The instruments featured in this catalog require a voltage of 230 V (50/60 Hz), 115 V (50/60 Hz). Please contact us if you have queries concerning different connected loads.

Service

Please contact your specialist dealer or IKA® direct in case of service queries. For spare parts replacement, please indicate the serial number and instrument type.

Certification



DIN EN ISO 9001 Reg. Nr. 4343

AISI steel designation

Refers to the American steel standard.

IKA[®] General 170 Terms and Conditions of Sale

The following terms and conditions shall apply to all sales. unless specifically agreed otherwise:

1. General

All agreements must be made in writing. Any terms and conditions of the buyer in his/its enquiries or orders which deviate from the present Terms and Conditions of Sale shall only apply if the supplier has specifically declared its agreement herewith. Any agreements deviating from the present Terms and Conditions of Sale shall only apply to the business for which they were agreed unless they are specifically prolonged

2. Quotations

The supplier shall be bound to all quoted prices for three months unless otherwise agreed. The right of prior sale shall be reserved. The documents pertaining to the offer, such as illustrations, drawings, weight and dimension details, etc. shall only be approximate unless they are specifically designated as binding. The supplier shall retain the ownership and copyright of cost estimates, drawings and any other documents; they may not be made available to any third parties. Plans received from the buyer and designated as confidential shall only be made available to third parties by the supplier with the consent of the buyer.

3. Conditions of delivery

The written order acknowledgement of the supplier shall be relevant for the scope of delivery. All ancillary agreements and modifications shall require written confirmation by the supplier.

4. Prices and payments

- a) Unless otherwise agreed, prices are ex-works, excluding packaging. INCOTERMS 2000 apply. Unless otherwise agreed, all prices shall apply ex works excluding packing. All prices shall be subject to the statutory rate of valueadded tax. Confirmed prices shall be based on prevailing material prices and wages. The supplier shall reserve the right to charge the material prices and wages prevailing at the time of delivery.
- b) Unless otherwise agreed, all payments shall be made to the cash office of the supplier without deductions or charges, with 2% cash discount for payment within 14 days or net within 30 days. If payments are deferred or not made as agreed, default interest at eight percent above the basic discount rate of the EZB shall be charged. Special payment conditions shall apply to export deliveries.
- c) No withholding of payments, nor any offsetting of counter claims disputed by the supplier, shall be permitted.

5. Deliveries - Delivery period

- a) Unless otherwise agreed in writing, deliveries are ex-works. INCOTERMS 2000 apply.
- b) The delivery period shall commence with the dispatch of the order acknowledgement but not before receipt of the documents, licenses and approvals to be acquired by the
- buyer and not before receipt of the agreed down-payment. c) The delivery period shall be deemed to have been upheld if the object of delivery has left the works of the supplier before the end of the delivery period or if readiness to supply has been notified.
- d) The delivery period shall be reasonable prolonged in the event of labor disputes, particularly strikes or lock-outs, or in the event of unforeseen impediments can be shown to have had a material effect on the production or delivery of the object of supply. This shall also apply if the aforesaid circumstances occur at sub-contractors of the supplier.
- e) If dispatch is delayed at the request of the buyer, the buyer shall be charged with the storage costs incurred commencing one month after the notification of readiness to deliver but not less than 1/2% of the invoice amount for each month if the goods are stored in the works of the supplier.
- f) In case of delayed acceptance by the buyer, and after setting and fruitless course of a reasonable period of time. the supplier has the right of further disposal of the goods.

6. Call-up of goods

Goods ordered on call shall be called up within a reasonable period with special agreement, but no later than 12 months from the date of the order acknowledgement. If ordered goods are not called up on time, the supplier shall be entitled to store the goods which are ready for dispatch, such storage being at the risk of the buyer, and to invoice the goods with all the storage costs incurred as if they had been delivered or to dispatch the goods without having received a dispatch request from the buyer.

7. Transfer of risk and acceptance of goods

- a) Risk shall pass to the buyer no later than the dispatch of goods, also if part-shipments are made or if the supplier has assumed other performances, e.g. dispatch costs or transportation and installation
- b) If specific instructions for the dispatch of goods are not included in the order, goods shall be dispatched at the discretion of the supplier, without any obligation for the cheapest mode of transport
- c) In the interests of the buyer, the supplier shall insure shipments against theft, breakage, transport, fire and water damage and against any other reasonable risks at the cost of the buyer. Only on the specific request of the buyer transport insurance of the aforesaid type shall not be concluded.

Unless otherwise agreed, the supplier shall charge 0.5% of the invoice value for transport insurance and 2% of the invoice value for fragile accessories. Any transport damages shall be notified to the supplier within 8 days, together with the damage report of the transport establishment; such transport damages shall otherwise not be accepted. Any incomplete deliveries shall likewise be notified to the supplier within 8 days; notifications of missing deliveries shall otherwise not be accepted. Shipments destined for export shall only be insured on the specific instructions of the buyer and at the cost of the buyer.

d) If dispatch is delayed for reasons attributable to the buyer, risk shall pass to the buyer on the date of readiness to supply; the supplier shall; however, be obliged to insure the goods at the request of the buyer and at the cost of the buyer. e) Part-shipments shall be admissible

8. Reservation of title

a) The supplier shall reserve title to the goods delivered until all claims of the supplier against the buyer arising from the business relationship have been settled in full, including all future claims arising from simultaneous or subsequent contracts. This shall also apply if individual or all claims of the supplier are placed on a current account and if a balance is drawn and recognized. In the event of any non-contractual conduct by the buyer. in particular payment delay on the part of the buyer, the supplier shall be entitled to demand the return of the reserved goods with prior notification and the buyer shall be obliged to return such goods. The return of goods or the pledging of goods by the supplier shall only constitute withdrawal from the contract if such withdrawal is specifically notified by the supplier in writing unless the German Hire Purchase Law applies. The buyer shall be obliged to notify the supplier immediately in writing if reserved goods are pledged or seized in any other way by a third party.

b) The buyer shall be entitled to sell the delivered goods in the ordinary course of business. The buyer shall, however, hereby assign to the supplier all his/its claims against his/its customers or third parties arising from such resale, irrespective of whether the resreved goods are resold without having been processed or not. The buyer shall also be entitled to collect the aforesaid claims after the aforesaid assignment to the supplier. This shall not prejudice the right of the supplier to collect such claims as long as the buyer discharges his/its payment commitments in an orderly and proper manner. The supplier shall be entitled to demand that the buyer notifies the sassigned claims and the names of the liable parties to the supplier, that all the details required for collection are provided, that the relevant documents are submitted to the supplier and that the liable parties are informed of the assignment. If the reserved goods are sold together with other goods to which the supplier has no title, the claim of the buyer against his/its customer shall be deemed as assigned to the supplier in the amount of the delivery price agreed by the supplier and the buyer.

c) Any processing or transformation of reserved goods by the buyer shall always on behalf of the supplier. If reserved goods are processed with other goods to which the supplier has no title. the supplier shall acquire co-ownership in the new chattel in the ratio of the value of the reserved goods to the value of the new processed chattel at the time of processing. The processed chattel shall also be governed by the provisions relating to the reserved goods. The supplier shall be obliged to release any securities to which he is entitled only if such security exceeds the secured claims by more than 25% provided such claims of the supplier have not already been settled by the buyer. d) The supplier shall, at the cost of the buyer, be entitled to insure

the reserved goods against theft, breakage, fire, water and any other damages unless the buyer is able to prove that he/it has taken out such insurances

e) Any intervention costs incurred by the supplier shall be borne by the buyer.

9. Liability for defects

Notwithstanding Section 11, the supplier shall be liable for defective supplies as follows, to the exclusion of all further claims: a) All those parts which prove unusable or the usability of which is severely impaired within 12 months of putting into service due to circumstances prevailing prior to the transfer of risk shall be rectified or replaced by the supplier without charge and at the reasonable discretion and option of the supplier. The identification of any such defects shall be notified to the supplier in writing immediately. Any replaced parts shall become the property of the supplier. If dispatch, installation or putting into service are delayed for reasons not attributable to the supplier, the aforesaid liability shall lapse no later than 15 months from the transfer of risk.

b) The right of the buyer to enforce claims for defects shall in all cases become statute-barred 6 months from the date of the due complaint by the buyer but no later than the end of the warranty period.

c) No liability shall be assumed for damages arising for the following reasons: improper or incorrect use, defective installation or putting into service by the buyer or third parties, natural wear and tear, incorrect or negligent handling and the use of unsuitable materials, replacement materials, defective construction work, unsuitable foundations, chemical, electrochemical or electrical influences unless they are attributable to negligence or intent on the part of the supplier.

d) The buyer shall, after consultation with the supplier, grant the supplier the necessary time and opportunity to carry out all the rectifications and replacements which the supplier considers necessary at its resonable discretion, otherwise the supplier shall be exempt from its liability for the aforesaid defects. Only in cases of emergency endangering operational safety and to avert disproportionately high damages - were by the supplier is to be informed immediately - or if the supplier is in delay with the rectification of the defect the buyer shall be entitled to rectify the

defect himself/itself, or the have the defect rectified by a third party and to demand reimbursement of the necessary costs from the supplier e) Of the direct costs directly incurred as a result of the

- costs of the replacement parts, including dispatch costs, be borne by the buyer.
- of the supplier
- g) Additional claims of the buyer, particularly compensation goods themselves, shall be excluded if permitted by law.

10. Liability for ancillary obligations

If, for reasons attributable to the supplier, the delivered goods cannot be used by the buyer as specified in the contract due to an omited or defective execution of recommendations and advice given prior to or after the conclusion of the contract - in particular usage a maintenance instructions for the delivered goods - the provisions of Sections 9 and 11 shall apply correspondingly, to the exclusion of any additional claims by the buver.

11. Right of withdraw by the buyer

a) The buyer shall be entitled to withdraw from the contract if the supplier is finally and conclusively unable to perform

prior to the transfer of risk.

b) The buyer shall be entitled to withdraw from the contract if delivery is delayed within the meaning of Section 5 and if the buyer grants the supplier a reasonable period of grace with a specific declaration that he/it will reject acceptance of the goods after such period of grace and if the period of grace is not upheld by the supplier c) If delivery of the goods is not possible during a period of acceptance delay or for reasons attributable to the buyer, the buyer shall be obliged to meet his/its contractual obligations

respond to a period of grace granted for the rectification or replacement of a defect attributable to the supplier within the meaning of the present Terms and Conditions of Sale. Such right of withdrawal by the buyer shall also apply in the event of impossibility to supply or the inability of the supplier to rectify or replace the aforesaid defect

IKA[®] General Terms and Conditions of Sale 171

rectification or replacements - provided the complaints of the buyer prove to be justified - the supplier shall bear the and reasonable dismantling and installation costs and the costs of providing any technicians and auxiliary staff of the buyer if the reimbursement of such costs can be equitably demanded in the specific circumstances. Other costs shall

f) The liability of the supplier shall lapse for the consequences of any improper modification or maintenance work undertaken by the buyer or a third party without the prior consent

claims and claims for damages not sustained by the delivered

d) The buyer shall also have a right of withdrawal from the contract if, through negligence or intent, the supplier fails to e) All other further claims of the buyer shall be excluded, if permitted by law.

12. Rights of withdrawal by the supplier

The contract shall be reasonably modified in case of unforeseen events within the meaning of Section 5 of the present Terms and Conditions of Sale, if such events materially change the financial and substantive implications of the performance of the supplier or if they materially affect the operations of the supplier and if it later transpires that the supplier is unable to perform its contractual obligations. If this is not economically possible, the supplier shall be entitled to withdraw from the whole or part of the contract. Any compensation claims by the buyer due to the exercise of such right of withdraw shall be excluded, if permitted by law. If the supplier makes use of its right to withdraw from the contract, it shall be obliged to notify the buyer immediately after having become aware of the implications of the aforesaid event.

13. Competent court and legal venue

- a))For all disputes arising from the contractual relationship, legal action shall be taken at the competent court for the registered office of the supplier or the branch of the supplier effecting delivery if the buyer is a registered trader, a legal entity under public law or a public-law fund. The supplier shall also be entitled to bring action at the principal place of business of the buyer.
- b) For legal relations in connection with this contract German material law is applicable, whereas the agreement of the United Nations regarding contracts ruling the international purchase of goods (CISG) is excluded.

Issue 04/2007

IKA®-WERKE GmbH & Co. KG D-79219 Staufen



IKA[®] General 172 HANDS for children

HANDS for children

HANDS for children is a nonprofit project of IKA®-Werke in Staufen, Germany with the goal to help and support the needy children of the Third World.

Experienced retirees from the IKA® team volunteer their time to manufacture the laboratory equipment for this program. HANDS for children combines the power of an independent company with the knowledge of experienced retired workers.

The profit gained by these activities is donated, in full, to institutions that help needy children or is used directly to help needy children. The recipients are chosen by the employees of HANDS for children and the donations are closely monitored.



The project »HANDS for Children« is supported by the following products:



EH 4 basic Immersion thermostat

For temperature control of liquids (NFL/I) up to 100 °C in open baths (min. bath depth 160 mm, min. usable depth 75 mm). Page 97



IKA®-PET

High-precision air cushion piston stroke pipette for science, research and routine work in the field of liquid handling. Page 102



VORTEX Genius 3

New vortex shaker suitable for short-time operation (touch function), activated by pres-sing shaker attachment or continuous operation. Page 46

VISION 2020: Christoffel Blind Mission in Peru

According to estimates there are 12.000 blind children under the age of 14 in Peru. Approximately 2.000 live in the capital, Lima, with eight million inhabitants. CBM and its campaign, VISION 2020, made it possible for thousands of Indios to find help for their eye problems. The exemplary work of the ophthalmic clinic for the poor is gaining widespread recognition in the region. Twelve projects for the prevention of blindness have already been implemented by the CBM.

HANDS for children was able to support the VISION 2020 campaign with more than \$ 50.000 in 2001. The money was used for the purchase of needed medication, setup of operating rooms, special microscopes, and anesthetic instruments.

Further information about CBM and its projects can be requested at: Christoffel-Blindenmission e.V., Niebelungenstr. 124, 64625 Bensheim, Germany.



People for People, building schools in Ethiopia project

In his anger over the unjust, inhumane gap between wealth and poverty on our planet the actor, Karlheinz Böhm, founded the association "Stiftung Menschen für Menschen e.V." ("People for People Foundation") in 1981. With this association he has given aid, independent of political, economic, or religious interests, to help Ethiopia. More than 1,7 million people, the majority being children, were helped in the last 20 years through the association.

However, so many still live under inhumane conditions. With profits from the sale of the EH 4 basic 2002 to 2006, HANDS for children will support the project in order to build schools in Babile Woreda - the extremely arid region of eastern Ethiopia.

The major activities of the project will be new buildings, furnishings and equipment as well as the repair of older schools and the building of housing for the teaching staff. At the same time primary education for the rural population will be implemented in order to minimize the illiterate rate (average in Ethiopia: males 60%, females 73%). This should set the stage for better learning and teaching. The long-term goal is to raise the educational level and with that increase the development of Ethiopia.

Further information on "People for People" and about the sponsored project will be available at: www.menschenfuermenschen.de or at: Menschen für Menschen, Brienner Str. 46, 80333 München, Germany.

IKA[®] General HANDS for children 173







Seamless transfer from the laboratory to the production line

Over the last few decades IKA® has built a reputation as a leading innovator for dispersers, kneaders and laboratory systems. Our high-quality stirrers, mixers and kneaders are widely used in the production of pharmaceuticals, chemicals, processed foods, paints, cosmetics and plastics.

Market leaders choose IKA® equipment and systems.







MP 10

1. PF

MP 25

Scale Up

IKA® offers full scaleability from single units to modular series right through to complete systems. IKA® products are renowned worldwide for their high quality. Our products meet the highest specifications for functionality and technology. IKA® products are good to look at, too.

The IKA® Technical Assistance Service works closely with customers to optimise production processes and to develop customized designs.

IKA® is a one-stop supplier offering everything you need from design, prototype development, production through to after-sales service.

For more information, visit us at: www.ikaprocess.com









MP 2000

MP 4000

IKA[®] General

176 Index

Short name	Page	AS 501.3	56	C 5040	161	EUROSTAR power	36	IKA®-PET	103
		AS 501.4	55	C 5041	161	control-visc 6000		Rack 1 und 3	
А		AS 501.5	56	C 5041.10	162			IKA®-PET soft	103
A 11 basic	87	AS 501.6	56	C 6	163	F			
A 11.1	87			C 60	162	FK 1	42	К	
A 11.2	87	В		C 62	162			KMO 2 basic	18
A 11.3	87	BIG SQUID	21, 22	C 62.1	163	Н		KS 130 basic	47
A 11.4	88	BMT-20 G/S	64	C 62.2	163	H 11	29	KS 130 control	47
A 11.5	88			C 7000	160	H 12/135	29	KS 260 basic	48
A 11.6	88	С		C 7000 Set 1	160	H 15	28	KS 260 control	48
A 11.7	88	C 08	163	C 7000 Set 2	160	H 159	98	KS 501 digital	48
AD 25	78	C 10	163, 165	C 7002	162	H 16 V	27	KV 600	162
AK 2.1	148	C 12	163, 165	C 7010	162	H 16.1	27		
AK 2.2	148	C 12 A	163, 165	C 7010.8	162	H 240	98	L	
AK 2.3	148	C 14	163, 165	C 7012	162	H 28	28	LABOR-PILOT 2000	0/4 79
AK 2.4	148	C 15	163, 165	C 7030	162, 165	H 29	28	LABOR-PILOT Con	troller 80
AK 2.5	148	C 200	155	C 710.4	163, 165	H 30	28	lab dancer	45
AK 2.6	148	C 200.1	162	C 723	163, 165	H 38	27	lab disc	18 – 19
AK 2.7	148	C 2000 basic,	156	C 9	163, 165	H 44	27, 116	labworldsoft®	143
AK 2.8	148	Version 1, Version 2		CC3-308B vpc	97	H 52	25	LR 2.1	133
AM 1	120	C 2000 control,	157	COLOR SQUID	20, 22	H 62.51	25	LR 2000.1	133
AOD 1	164	Version 1, Version 2				H 66.51	25	LR 2000.10	132
AOD 1.11	163, 165	C 2000 basic / control	156, 157	D		H 70	25	LR 2000.11	132
AOD 1.12	163, 165	high pressure		DA 2	147	H 99	28	LR 2000.2	133
AOD 1.13	165	C 21	162, 165	DC 2	147	H 100	28	LR 2000.20	132
AOD 1.3	165	C 248	162, 165	DK 25.11	78	H 101	28	LR 2000.21	132
AS 1.10	57	C 25	162	DK 50.11	78	HB 4 basic	95	LR 2000.3	133
AS 1.11	57	C 29	162, 165	DT-20	64	HBR 4 digital	95	LR 2000.30	134
AS 1.12	58	C 4	163, 165	DTM 12	118	HCT basic	96	LR 2000.4	133
AS 1.13	58	C 43	163	DTM 12.10	118, 148	HKD 06.2	59	LR 2000.40	134
AS 1.30	57	C 48	163	DZM control.o	119	HKD-T 06 D	59	LR 2000.52	134
AS 1.31	57	C 5	163	DZM-K	119	HKD 06.10	59	LR 2000.53	133
AS 1.5	57	C 5000 control	158	DZM-M	119	HP 30 digital	96	LR 2000.54	132
AS 1.6	57	Pack. 1/10		DZM-S.o	119	HP 30.1	98	LR 2000.57	132
AS 1.60	57	C 5000 control	158			HS 260 basic	49	LR 2000.60	134
AS 1.61	57	Pack. 1/12		E		HS 260 control	49	LR 2000.70	130
AS 1.7	58	C 5000 control	159	EH 4 basic	97	HS 501 digital	49	LR 2000.75	128
AS 1.8	58	Pack. 2/10		EH 4.1	98			LR 2000.80	132
AS 130.1	54	C 5000 control	159	EH 4.2	98	1		LR 2000.85	132
AS 130.2	54	Pack. 2/12		EH 4.3	98	10 2	148	LR 2000.90	134
AS 130.3	54	C 5000 duocontrol	159	ETS-D 5	24, 117	IO 2.1	148	LR 2000 P	128, 129
AS 130.4	54	Pack. 3/10		ETS-D 6	24, 117	IKAFLON [®] 10	29	LR 2000 V	130, 131
AS 2.1	58	C 5000 duocontrol	159	EUROSTAR digital	33, 34	IKAFLON® 15	29	LR 2000.VK	135
AS 2.2	58	Pack. 3/12		EUROSTAR power	35	IKAFLON [®] 20	29	LR-2.ST	125 – 127
AS 2.3	58	C 5003.1	163	basic		IKAFLON® 25	29	LR-2.SP	135
AS 2.4	58	C 5010	162	EUROSTAR power	30, 31, 35	IKAFLON® 30	29	LT 5.20	99
AS 2.5	58	C 5010.3	163	control-visc		IKAFLON [®] 40	29	LT 5.22	99
AS 260.1	54	C 5010.4	162, 165	EUROSTAR power	36	IKAFLON® 50	29	LT 5.23	99
AS 260.2	55	C 5010.5	162	control-visc P1		IKAFLON® 80	29	LT 5.24	99
AS 260.3	55	C 5012	162	EUROSTAR power	37	IKAFLON® 110	29		
AS 260.5	55	C 5012.3	163, 165	control-visc P4		IKAFLON® 155	29	М	
AS 501.1	55	C 5020	162	EUROSTAR power	37	IKA®-PET Pipetten	102	M 20	89
AS 501.2	56	C 5030	162	control-visc P7		IKA [®] Pipettenspitzen	102	M 20.1	89

M 21	89	R		REO basic	23	S 25 D – 10 G – KS	75	V	
M 22	89	R 1001	42	REO basic C	23	S 25 D – 14 G – KS	75	VC 1.1	121
M 23	89	R 1002	42	RET basic s. c.	13	S 25 KR – 18 G	71	VC 1.3	121
magicLAB®	82 - 85	R 103	114	RET basic C s. c.	13	S 25 KR – 25 F	71	VC 2	120
Maxi MR 1 digital	23	R 104	114	RET control-visc s. c.	14	S 25 KR – 25 G	70	VC 2.4	121
MF 0.25	90	R 1300	40	RET control-visc C s. c.	14	S 25 KV – 18 G	71	VG 3.1	51
MF 0.5	90	R 1302	40	RH 3	116	S 25 KV – 25 F	71	VG 3.2	51
MF 1.0	90	R 1303	40	RH 5	116	S 25 KV – 25 F – IL	71	VG 3.3	51
MF 10 basic	90	R 1311	40	RH basic 2	15	S 25 KV – 25 G	71	VG 3.31	51
MF 10.1	90	R 1312	40	RH basic KT/C s. c.	15	S 25 KV – 25 G – IL	71	VG 3.32	51
MF 10.2	90	R 1313	40	RH digital KT/C s. c.	15	S 25 N – 8 G	71	VG 3.33	51
MF 2.0	90	R 1330	40	RO 5 power	17	S 25 N – 10 G	71	VG 3.34	51
MF 3.0	90	R 1331	40	RO 10 power	17	S 25 N – 18 G	71	VG 3.35	51
Midi MR 1 digital	23	R 1332	40	RO 15 power	17	S 25 N – 25 F	71	VG 3.36	51
Mini MR standard	18	R 1333	40	RS 1	29	S 25 N – 25 G	70	VG 3.37	51
MS 1.21	50	R 1335	40	RSE	29	S 25 N – 10 G – VS	71	VK 60/01	139
MS 1.31	50	R 1342	40	RT 5 power	16	S 25 NK – 19 G	70	VK 600 control	139
MS 1.32	50	R 1345	40	RT 10 power	16	S 50 KB – G 45 F	73	VM 600 basic	135
MS 1.33	50	R 1352	40	BT 15 power	16	S 50 KB – G 45 G	72	VORTEX Genius 3	46
MS 1.34	50	R 1355	40	BV 05.3	109 115	S 50 KB – G 45 M	73	VXR basic	46
MTS 2/4 digital	47	R 1373	40	RV 05 basic 1-B	108	S 50 KV - G 45 G - II	73	VX 1	52
MS 3 basic	45	R 1375	40	RV 05 basic 2-B	108	S 50 N - G 45 F	73	VX 2	52
MS 3 digital	45	R 1376	40	RV 05 basic 2 D	100	S 50 N - G 45 G	72	VX 2 F	52
MS 3 1	-50 50	R 1381	40	RV 06 1	100	S 50 N - G 45 M	72	VX 7	52
MS 3.3	50	R 1382	40	RV 06 11	100	S 50 N - W 65 SK	76	VX 8	52
MS 3 A	50	R 1385	40	RV 06.13	110	S 50 N - W 80 SMK	76	VX 8 1	52
1010 0.4	50	R 1388	40	RV 06.15	110	S 50 KB - W 80 SMK	76	VX 0.1	52
P		R 1389	40	RV 06.2	110	S 65 KG-HH - G 65 F	73	VX 11 1	53
	149	P 1401	40	DV 06.4	110		73	VX 11.1	50
PC 1 2	110 1/10	R 1407	42	RV 06.5	111	S 65 KG-HH - G 65 M	73	VX 11.2	53
PC 1 4	1/10, 140	R 1405	42,70	RV 06.6	111	SI 400	/3	VX 11.3	53
PC 1 5	140	R 192	42,70	RV 06 7	111	SI 472	43	VA 11.4	55
PC 2 1	140	R 1925	11/	RV 06-ML 1-R	107	SI 472	43		
PC 2.2	140	R 1826	114	RV 06-ML 2-B	107	ST-20	43		
PC 2.2	00 1/12	R 1927	114	RW/11 basic	32	SW/ 18	70		
PC 4 1	140	P 2202	40	PW/ 14 basic	22	37710	70		
PC 4.1	147	n 2302	40	RW 14 Dasic	22.24	т.			
PC 0.1	140	N 2303	40	PW/ 20 digital	33, 34	T 10 basia	65		
PCI 0.2 Dragona Dilat 2000/4	70	R 200	40	NVV 20 digital	30	T 10 Dasic	60		
Process-Pilot 2000/4	. 00	R 270	110	DW/ 20 D	20	T 25 digital	66		
Process-Pliot Controller	110	n 270 D 271	110	NW 20 D	39	T 25 digital	60		
PT 100.23	110	n 2/1	110	NVV 47 D	39	T SU DASIC	60		
PT 100.24	110	n 2/22	114	0		T 05 D	69		
PT 100.25	118	R 2/23	115	S	74	1 653	115		
PT 100.27	118	R 301	42	S 10 D - 7 G - KS - 65	74	TRIKA® 25	29		
PT 100.5	99	R 301.1	42	S 10D-7G-KS-110	74	TRIKA® 40	29		
PT 100.50	26	R 302	43	5 10 N - 5 G	/2				
PT 100.51	26	R 350	98	S 10 N - 8 G	72				
PT 100.52	26	R 380	27	S 10 N - 10 G	72	UTL 25 digital	77		
PI 100.7	99	R 4/2	115	S 18 / 25-ET50	75	UTL 50 basic	77		
PT 1000.50	26	K 4/4	115	S 18 D – 10 G – KS	75	ULTRA-TURRAX®	63		
PT 1000.51	26	R 50	76	S 18 D – 14 G – KS	75	Tube Drive			
PT 1000.60	12, 25	RC 1	43	S 18 N – 10 G	70	ULTRA-TURRAX®	63		
		RCT basic s. c.	10 - 12	S 18 N – 19 G	70	Workstation			

IKA[®] General

Index 177

IKA[®] General Notes		

Hotlines

You can reach your contact person directly via IKA® Hotlines.

Germany	Phone:	+49 7633 831-140
Sales	Fax:	+49 7633 831-901
Europe / Middle East / Africa	Phone:	+49 7633 831-250
UK, Scandinavia, Middle East	Fax:	+49 7633 831-98
Belgium, France,	Phone:	+49 7633 831-123
Spain, Africa	Fax:	+49 7633 831-98
Netherlands, Italy, Austria, Switzerland, Eastern Europe	Phone: Phone: Fax:	+49 7633 831-152 +49 7633 831-566 +49 7633 831-98
Russian Federation	Phone: Fax:	+49 7633 831-215 +49 7633 831-98
Service	Dhanai	. 40 7622 021 120

Phone: +49 7633 831-139 Phone: +49 7633 831-202 Fax: +49 7633 500212 e-mail: service@ika.de



Designed to work perfect

IKA®-Werke GmbH & Co. KG Janke & Kunkel-Str. 10 D-79219 Staufen Germany

Phone: +49 7633 831-0 Fax: +49 7633 831-98 e-mail: sales@ika.de www.ika.net

Subject to technical changes. Indications not binding for delivery. EN 04 / 07

