

Measures concentration of  
aqueous solution simply!

**New product!**

Cat.No.3403

# palette Digital Refractometer **PR-32**

- Soluble solid%
- Control of dilution ratio%
- Sugar content (Brix)%



The PR-32 is originally the refractometer developed to measure concentration of various solutions with application of refractive principle.

## Features

### Excellent functions and response

Automatic concentration measurement with the optical sensor, automatic correction of temperature error with the built-in microprocessor, and other many devices realize remarkable operability.

### Highly reliable digital display

Large-sized LCD shows a concentration of measurement result in digits, which are easily and quickly readable.

### Anticorrosive sample stage

The sample stage that is made of stainless steel for dealing with many kinds of solution always keeps cleanliness being free from corrosion.

## Main applications

### In the food and beverage industries

At orchard, fruit shop and beverage factories for measuring sugar content of fruit juice and soft drinks.

### In the industrial fields

Control of dilution ratio(multiples in %)

When water soluble cutting oil, aqueous cleaning agent, surfactant, liquid medicine, liquid chemicals, etc., are watered down for use, their dilution ratios(multiples in %) can be managed by the PR-32 digital refractometer.

## MEASURING METHOD



1. Drip a liquid sample onto the prism surface.



2. Press the START/OFF switch.



3. After the arrow indicator blinks two times, the concentration of the sample is displayed in a Brix value(%).

## Specifications

Cat. No.	3403
Measuring range	Brix 0.0 to 32.0%
Minimum indication	Brix 0.1%
Measuring accuracy	Brix $\pm 0.2\%$
Measuring temperature	5 to 40°C
Measuring time	2 seconds approx.
Sample quantity	0.1ml or more
Zero setting method	By distilled water and ZERO setting switch
Power supply	006P battery(9V)
Dimensions and Weight	17(W) × 9(D) × 4(H)cm, 300g

### In regard of displayed value

The PR-32 digital refractometer measures the refractive index of a liquid sample and displays the measurement result in a Brix value(%) that the refractive index is converted into a concentration(%) of sugar water. Relation between refractive index and concentration(%) depends on the kind of liquid. Therefore, the PR-32 digital refractometer requires to make a conversion table for between Brix(%) and concentration(%) preliminary to measurement.



\*Specifications and design of outer appearance are subject to change without notice.

# ATAGO CO., LTD.

32-10 HONCHO, ITABASHI-KU, TOKYO 173-0001, JAPAN  
 PHONE: 81-3-3964-6156, FAX: 81-3-3964-6137  
<http://www.atago.net/>  
[export@atago.net](mailto:export@atago.net)