

25 years
Tradition :: Quality :: Value



HI 221 & HI 223 Calibration Check™ pH Meters


ISO 9001 CERTIFIED

 **HANNA**®
instruments
hannainst.com

Calibration Check™ pH Meters

A pH electrode, when it is properly manufactured and is kept clean, will retain its physical characteristics for a long time. A common problem associated with pH measurements is the use of a pH electrode that has not been properly cleaned. This is very important because during calibration the instrument assumes that the electrode used is clean, and the standardization curve created during the calibration process will remain the same until the next calibration. pH meters on the market today will allow an offset of approximately ± 60 mV. The deviation from 0 mV is not unusual; in fact, it represents the true characteristics of the pH electrode. The deviation from 0 mV becomes a problem if it is the result of calibrating a dirty electrode. Figure A below shows that the pH measured by a dirty electrode changes over a period of time. This results from the residue on the pH electrode bulb dissolving into the solution and the electrode gradually returning close to its true characteristics. The resulting pH measurements, based upon the calibration of a dirty electrode, will then be incorrect. Hanna's HI 221 and HI 223 compare the characteristics of the pH electrode from one calibration to the next. In cases of large variances in electrode condition that can only be the result of a dirty electrode, these meters alert the user that the electrode needs to be properly cleaned prior to the calibration process.

Figure A



Electrode 1 was properly cleaned prior to calibration.
Electrodes 2 and 3 were not.



GLP

In addition to storing calibration data, the instruments can be programmed to give an alarm when a new calibration is required



Log-On-Demand

Both instruments have a "Log-On-Demand" function to record 100 (HI 221) or 500 (HI 223) readings.



Recall

Readings logged with the HI 221 & HI 223 can be recalled at a later time.

How Calibration Check™ Works

HI 221 and HI 223 are able to detect if the calibration buffer solution is contaminated.



When a pH electrode is dirty, it has a slow response time & unstable reading. HI 221 and HI 223 warn the user in the event a dirty electrode has been used.



Calibration Check™ Features

- :: **Enhanced Calibration Messages**
During calibration, the user is warned if one or more parameters are not suitable to perform an accurate calibration.
- :: **Electrode Condition on LCD Display**
Determined from the electrode offset and slope.
- :: **Electrode Response Time on LCD Display**
Determined from electrode performance during calibration.
- :: **Calibration Alarm Time Out**
Can be programmed from 1 to 7 days or can be disabled.

Other HI 221 & HI 223 features:

- :: Logging up to 500 Samples
- :: Last Calibration Date & Data
- :: Instrument ID Number
- :: Real Time Clock
- :: PC Interface



Specifications		HI 221	HI 223
Range	pH	-2.00 to 16.00	-2.00 to 16.00 or -2.000 to 16.000
	mV	±699.9; ±2000	±999.9; ±2000
	°C	-20.0 to 120.0	-20.0 to 120.0
	°F	-20.0 to 120.0	-20.0 to 120.0
Accuracy	pH	±0.01	±0.01; ±0.002
	mV	±0.2 (±699.9); ±1 (±2000)	±0.2 (±699.9); ±0.5 (±999.9); ±1 (±2000)
	°C	±0.4	±0.4
	°F	±0.4	±0.4
Resolution	pH	0.01	0.01; 0.001
	mV	0.1 (±699.9); 1 (±2000)	0.1 (±999.9 mV); 1 (±2000 mV)
	°C	0.1	0.1
	°F	0.1	0.1
Calibration Check		Yes	Yes
pH Calibration		Automatic 1 or 2 points (pH 1.68, 4.01, 6.86, 7.01, 9.18, 10.01, 12.45)	
Temperature Comp.		Manual (MTC) or Automatic (ATC) temperature compensation from -20.0 to 120.0°C	
pH electrode		HI 1131P glass body, single junction refillable cell, BNC + pin (included)	
Temperature probe		HI 7669/2W stainless steel probe (included)	
PC Interface		Opto-isolated RS232	
Logging		100 points	500 points
Input impedance		10 ¹² ohm	
Power		12VDC adapter (included)	
Environment		0 to 50°C (32 to 122°F) 95% RH	
Dimensions		240 x 182 x 74 mm (9.4 x 7.1 x 2.9")	
Weight		1.1 Kg (2.5 lb.)	
Warranty		5 Years	

Ordering Information

HI 221 and HI 223 are supplied complete with HI 1131P glass-body combination pH electrode with BNC connector + PIN and 1 m cable, HI 7669/2W stainless steel temperature probe, HI 76404 electrode holder, pH 4.01 and pH 7.01 calibration solutions, electrolyte solution, 12V DC adapter and instruction manual.

Authorized Distributor www.clarksonlab.com

Electrodes

All electrodes part numbers ending in P are supplied with a BNC and PIN connector & 1 m (3.3') cable, as shown below:

- HI 1043P** Use: Strong acid/Alkali
Glass-body, double junction, refillable, combination pH electrode
- HI 1053P** Use: Emulsions
Glass-body, triple ceramic, refillable, combination pH electrode
- HI 1083P** Use: Biotechnology
Glass-body, open junction, refillable, combination pH electrode
- HI 1131P** Use: General Purpose
Glass-body, single junction, refillable, combination pH electrode

Solutions

pH Calibration Solutions

- HI 70004P pH 4.01 Buffer Sachets @ 25°C, 20 mL, 25 pcs
- HI 70007P pH 7.01 Buffer Sachets @ 25°C, 20 mL, 25 pcs
- HI 70010P pH 10.01 Buffer Sachets @ 25°C, 20 mL, 25 pcs
- HI 7001L pH 1.68 Buffer Solution @ 25°C, 500 mL
- HI 7004L pH 4.01 Buffer Solution @ 25°C, 500 mL
- HI 7004/1L pH 4.01 Buffer Solution @ 25°C, 1 L
- HI 7004/1G pH 4.01 Buffer Solution @ 25°C, 1 US Gallon
- HI 7006L pH 6.86 Buffer Solution @ 25°C, 500 mL
- HI 7006/1L pH 6.86 Buffer Solution @ 25°C, 1 L
- HI 7006/1G pH 6.86 Buffer Solution @ 25°C, 1 US Gallon
- HI 7007L pH 7.01 Buffer Solution @ 25°C, 500 mL
- HI 7007/1L pH 7.01 Buffer Solution @ 25°C, 1 L
- HI 7007/1G pH 7.01 Buffer Solution @ 25°C, 1 US Gallon
- HI 7009L pH 9.18 Buffer Solution @ 25°C, 500 mL
- HI 7009/1L pH 9.18 Buffer Solution @ 25°C, 1 L
- HI 7009/1G pH 9.18 Buffer Solution @ 25°C, 1 US Gallon
- HI 7010L pH 10.01 Buffer Solution @ 25°C, 500 mL
- HI 7010/1L pH 10.01 Buffer Solution @ 25°C, 1 L
- HI 7010/1G pH 10.01 Buffer Solution @ 25°C, 1 US Gallon

Electrode Storage Solutions

- HI 70300L Storage Solution, 460 mL

Electrode Cleaning Solutions

- HI 70000P Electrode Rinse Sachets, 20 mL, 25 pcs
- HI 7061L General Cleaning Solution, 460 mL
- HI 7073L Protein Cleaning Solution, 460 mL
- HI 7074L Inorganic Cleaning Solution, 460 mL
- HI 7077L Oil & Fat Cleaning Solution, 460 mL

Electrolyte Solution for Electrode Refill

- HI 7071 3.5M KCl + AgCl Electrolyte, 4 x 30 mL, for single junction electrodes
- HI 7082 3.5M KCl Electrolyte, 4 x 30 mL, for double junction electrodes

Connectivity

Software

- HI 92000 Windows® compatible software

Hardware

- HI 920010 9-pin serial cable for PC

HANNA®
instruments
hannainst.com