



Glass Joint REPIPET® Dispensers



Model 3100GR Dispenser

All Glass Dispensers!

- Ground glass joints
- Use with an Erlenmeyer flask on your hot plate
- Completely autoclavable
- ☆ Glassware carries a 50% lifetime guarantee against breakage

Product Description

- A Standard REPIPET dispenser with standard taper, ground glass joints.
- Full range adjustability, 100-Division Scale.
- 1% accuracy and 0.1% reproducibility at full scale guaranteed.
- Includes both magnifying indicator and metal indicator.

ORDERING INFORMATION

Model #	Description	Size (ml)	T Joint Size	US List Price
3005GR	Dispenser with 500 ml Reagent Bottle, glass intake tube	5.0 ml	24/40	\$372.00
3010GR	Dispenser with 1000 ml Reagent Bottle, glass intake tube	10.0 ml	29/42	\$398.00
3020GR	Dispenser with 1000 ml Reagent Bottle, glass intake tube	20.0 ml	29/42	\$450.00
3050GR	Dispenser with 2000 ml Reagent Bottle, glass intake tube	50.0 ml	29/42	\$595.00
3100GR	Dispenser with 2000 ml Reagent Bottle, glass intake tube	100.0 ml	29/42	\$637.00
3001/2GV	Dispenser with 100 ml Conical Vessel, glass intake tube	0.5 ml	24/40	\$405.00
3001GV	Dispenser with 100 ml Conical Vessel, glass intake tube	1.0 ml	24/40	\$405.00
3001/2GE	Dispenser with 125 ml Erlenmeyer Flask, glass intake tube	0.5 ml	24/40	\$352.00
3001GE	Dispenser with 125 ml Erlenmeyer Flask, glass intake tube	1.0 ml	24/40	\$352.00
3001/2G	Dispenser-only, Teflon intake tube, use with catalog #5000125 or 1170L	0.5 ml	24/40	\$316.00
3001G	Dispenser-only, Teflon intake tube, use with catalog #5000125 or 1170L	1.0 ml	24/40	\$316.00
3005G	Dispenser-only, Teflon intake tube, use with catalog #1500500	5.0 ml	24/40	\$342.00
3010G	Dispenser-only, Teflon intake tube, use with catalog #15001000	10.0 ml	29/42	\$358.00
3020G	Dispenser-only, Teflon intake tube, use with catalog #15001000	20.0 ml	29/42	\$409.00
3050G	Dispenser-only, Teflon intake tube, use with catalog #15002000	50.0 ml	29/42	\$501.00
3100G	Dispenser-only, Teflon intake tube, use with catalog #15002000	100.0 ml	29/42	\$546.00

ACCESSORIES

Catalog #	Description	US List Price
1170L	100 ml Conical Vessel-24/40 T joint size	\$82.80
5000125	125 ml Erlenmeyer Flask-24/40 T joint size	\$49.50
1500500	500 ml Corning Bottle-24/40 T joint size	\$40.00
15001000	1000 ml Corning Bottle-29-42 T joint size	\$53.80
15002000	2000 ml Corning Bottle-29/42 T joint size	\$110.60