

Electrothermal Histology



Slide Drying Bench



Higher Temperature Slide-Warmer With Drying Bars!

- Temperature range from slightly above ambient to 100°C
- Accepts up to 50 slides
- Dry slides across the drying bars, angled on the bars or lying flat without bars

CE

Model MH6616

Product Description

- Painted aluminum top reduces glare, offering a high-contrast background for viewing specimen.
- Aluminum housing, black and white exterior is durable and easy to clean.
- Slide capacity: 50 (1" x 3") slides.
- Easily replaceable heating element.

Operation

- Clear light indicates when power is on; amber light indicates when the heater is on.
- Temperature control is provided by a built-in energy regulation with a temperature range to 100°C.
- Reference dial and knob allow the temperature to be closely adjusted and regulated.

ORDERING INFORMATION

Model #	Electrical (50/60 Hertz)		Overall Dimensions Inches (cm)			Shipping Weight Lb. (kg)	US List Price
	Volts	Watts	W	H	D		
MH6616X1'	115	150	15.5 (39)	3.5 (9)	7.12 (18)	5.3 (2.4)	\$829.00
MH6616'	230	150	15.5 (39)	3.5 (9)	7.12 (18)	5.3 (2.4)	\$829.00

'CE