

Thermolyne Histology



RH View Box



Type LL16100

Our Compact Economical Slide-Warmer!

- Soft, fluorescent, glare-free light
- Built-in temperature indicator

Product Description

- Built-in temperature indicator for easy, accurate monitoring of viewing area.
- Viewing area temperature adjusts easily to compensate for ambient temperature changes. (Range: 45°C to 50°C).
- Rh slide reaches temperature of 37° to 39°C in two minutes, eliminating wasted heat up time.
- Soft, fluorescent, glare-free light provides uniform illumination of slide contents.
- View box can be gently rocked using large knob on support cradle.
- Sturdy, easy-to-clean extruded aluminum case is compact; utilizes minimal bench space.
- Meets Rh slide test standards set by the AABB.
- Unit conveniently includes three-wire cord and plug.

Applications

- Blood typing
- Rh determinations
- Warming slides for gram-stain
- Tissue typing
- Reading enzyme antibody screens in microplate

PRODUCT SPECIFICATIONS

Type	Viewing Area Dimensions Inches (cm)		Overall Dimensions Inches (cm)			Shipping Weight Lb. (kg)
	W	D	W	H	D	
LL16100	9.50 (24.1)	3.06 (7.77)	13 (33)	4 (10.1)	4 (10.1)	4.25 (1.92)

ORDERING INFORMATION

Model #	Volts	Electrical (50/60 Hertz) Amps	Watts	Viewing Area Operating Temp Range °C	US List Price
LL16115	120	0.24	30	45-50°C	\$386.00
LL16110-26	220	0.12	25	45-50°C	\$386.00
LL16110	240	0.13	25	45-50°C	\$386.00

-26 Models Supplied with European cord set.