

Thermolyne Ovens



Oven-Incubator/Hot Plate-Oven



Type OV19200

Our Most Compact, Multipurpose Ovens!

- Gravity convection heating
- Excellent temperature uniformity
- Economical

Product Description

- Aluminum chamber provides excellent temperature uniformity.
- Double-pane glass window located in door permits easy viewing of samples.

- Magnetic latch assures positive door closure.
- Supplied with one removable metal shelf for maximum sample capacity, three-wire cord, and plug.

Operation

- Bimetallic thermostat provides time-saving rapid heat-up (56°C in nine minutes).
- Thermometer port at top of oven accommodates any partial immersion thermometer for temperature monitoring.
- Chamber is vented (1" dia. vent) through the top to allow removal of moisture and other exhaust vapors.

Applications

- General lab drying and baking applications
- Culture incubations

PRODUCT SPECIFICATIONS

Type	Chamber Dimensions Inches (cm)			Overall Dimensions Inches (cm)			Chamber Volume cu. ft. (l)	Shipping Weight Lb. (kg)
	W	H	D	W	H	D		
OV19200	12.3 (31.1)	11.7 (29.7)	12.0 (30.5)	13.3 (33.7)	15.5 (39.3)	13.0(33.0)	0.9 (25.5)	30(13.6)

ORDERING INFORMATION

Model #	Electrical (50/60 Hz)		Operating Temp Range °C	Temp Control @ 37°C	Temp Uniformity @ 37°C	US List Price
	Volts	Amps				
OV19225	120	2.9	350	Ambient +5°C-140°C	±1.5°C	\$863.00
OV19220-26	240	1.5	350	Ambient +5°C-140°C	±1.5°C	\$863.00

-26 Models Supplied with European cord set.



Product Description

- Turns your hot plate into a basic gravity convection oven.
- Includes one removable shelf to maximize sample capacity.
- Double-wall construction features baked-on enamel steel exterior and aluminum inner chamber, making the hot plate oven-durable and easy to clean.
- Thermometer port on top of unit accommodates standard laboratory thermometers.

Applications

- Ideal for general drying, baking, preheating and curing applications.

Type OV10600

ORDERING INFORMATION

Model #	Chamber Dimensions Inches (cm)			Overall Dimensions Inches (cm)			Maximum Temperature °C	Chamber Volume cu. ft. (l)	Shipping Weight Lb. (kg)	US List Price
	W	H	D	W	H	D				
OV10600	6.0 (15.2)	6.0 (15.2)	6.0 (15.2)	8.5 (21.5)	7.8 (19.6)	9.3 (23.4)	250°C	0.1 (2.8)	6.8 (3.0)	\$245.00

Fits the following hot plate models: HPA1910M, HPA1915B, HPA1915B-13, HP11515B, RC1910, RC1915, (Sold separately. See Hot plate section of this catalog).