

Lab-Line  
Refrigerators/Freezers



## Refrigerated Incubator



Model 3554-37



Model 3550DT

### All Shapes and Sizes!

- Appropriate for B.O.D. testing
- Available in three chamber sizes and two temperature ranges
- Many options available



### Product Description

- Model 3550 features economical under counter design with high/low safety thermostat back-ups and trim, top-mounted control panel.
- Model 3550DT has easy-access side mounted control panel.
- Models 3550DT and 3550-1DT include a dual set point microprocessor control with a seven day timer that allows two temperatures to be set and switched automatically.
- Models 3554-35 and 3554-37 feature microprocessor-based temperature control with digital LED readout.
- Forced air circulation and quiet-running fan provide temperature control and uniformity.
- Durable enameled steel exterior.
- Many options available—see page 235 for details.
- Door key lock is standard.

### PRODUCT SPECIFICATIONS

Model #	Chamber Volume cu. ft. (ℓ)	Chamber Dimensions Inches (cm)			Overall Dimensions Inches (cm)			Number of Shelves	Shelf Area sq. ft. (sq. m)	B.O.D. Bottle Capacity	Shipping Weight Lb. (kg)
		W	H	D	W	H	D				
3550	5.5 (156)	19.5 (50)	22 (56)	15 (38)	23.5 (60)	34 (86)	23.75 (60)	3 adjustable	6.5(0.6)	78	115 (52)
3550-1	5.5 (156)	19.5 (50)	22 (56)	15 (38)	23.5 (60)	34 (86)	23.75 (60)	3 adjustable	6.5(0.6)	78	115 (52)
3550DT	5.5 (156)	19.5 (50)	22 (56)	15 (38)	27 (69)	34 (86)	23.75 (60)	3 adjustable	6.5(0.6)	78	155 (70)
3550-1DT	5.5 (156)	19.5 (50)	22 (56)	15 (38)	27 (69)	34 (86)	23.75 (60)	3 adjustable	6.5(0.6)	78	155 (70)
3554-35 <sup>1</sup>	19 (538)	27 (69)	56 (142)	18 (46)	32 (81)	73.5 (187)	31 (79)	5 fixed	15.3 (39)	355	420 (191)
3554-36 <sup>1</sup>	19 (538)	27 (69)	56 (142)	18 (46)	32 (81)	73.5 (187)	31 (79)	5 fixed	15.3 (39)	355	420 (191)
3554-37	28 (793)	32 (81)	65.5 (166)	22 (56)	35 (89)	76 (193)	33 (84)	4 adjustable	19.6 (51)	495	570 (259)
3554-38	28 (793)	32 (81)	65.5 (166)	22 (56)	35 (89)	76 (193)	33 (84)	4 adjustable	19.6 (51)	495	570 (259)

<sup>1</sup> ETL listed.

### ORDERING INFORMATION

Model #	Volts	Electrical (50/60 Hz)		Operating Temp Range °C	Temperature Control °C	US List Price	Model #	Volts	Electrical (50/60 Hz)		Operating Temp Range °C	Temperature Control °C	US List Price
		Amps	Watts						Amps	Watts			
3550	120	4.0	350	+5°-50°	±0.25°	\$2,942.00	3554-35	120	7.0	810	-2°-60°	±0.5°	\$4,407.00
3550-1	220	2.2	350	+5°-50°	±0.25°	\$2,942.00	3554-36	220	3.8	810	-2°-60°	±0.5°	\$4,722.00
3550DT	120	4.0	350	+5°-50°	±0.25°	\$3,887.00	3554-37	120	11.0	1320	-2°-60°	±0.5°	\$7,532.00
3550-1DT	220	2.2	350	+5°-50°	±0.25°	\$3,887.00	3554-38	220	6.0	1320	-2°-60°	±0.5°	\$7,242.00