



## Variable Speed Platform Angle Rockers



Model 4638

### Rocking Motion Provides Gentle to Vigorous Mixing!

- Variable speed control of 0-100 cycles per minute
- Adjustable rocking angle from 10-15 degrees from horizontal
- ☆ 24 Month Warranty



#### Product Description

- Provides gentle, adjustable-angle rocking.
- Ideal for applications such as hybridization blotting or staining and destaining gels.
- 15.4" x 15.4" (39 x 39 cm) non-skid platform accommodates many types of glassware.
- White platform mat allows easy viewing of glassware contents.
- Durable steel body encloses motor to prevent hazardous spills onto drive mechanism.

#### Operation

- Electronic speed control from 0 to 100 cycles per minute.
- Adjustable rocking angle from 10 to 15° from horizontal.
- Operates in temperatures from 4° C to 40°C and humidity conditions from 20 to 80% non-condensing.

#### Applications

- Hybridizations
- Blotting
- Staining and destaining gels

#### PRODUCT SPECIFICATIONS

Type	Maximum Load	Platform Dimensions W x D Inches (cm)	Overall Dimensions Inches (cm)			Shipping Weight
			W	H	D	
4638	5 Lbs (2.3 kg)	15.4(39) x 15.4(39)	15.4(39)	7(15)	15.4(39)	19 Lbs (9 kg)
4638-1CE	5 Lbs (2.3 kg)	15.4(39) x 15.4(39)	15.4(39)	7(15)	15.4(39)	19 Lbs (9 kg)

#### ORDERING INFORMATION

Model #	Volts	Amps	Watts	Hz	Speed Cycles per minute	Rocking Angle Degrees	Ambient Conditions	US List Price
4638 <sup>1</sup>	120	0.4	50	50/60	0-100	10-15 from horizontal	4 to 40°C, 20 to 80% RH non-condensing	\$977.00
4638-1CE <sup>1,2</sup>	240	0.2	50	50/60	0-100	10-15 from horizontal	4 to 40°C, 20 to 80% RH non-condensing	\$977.00

<sup>1</sup> ETL, <sup>2</sup> CE