



## Fixed Speed Platform Angle Rockers



Model 4635

### When Your Application Requires Safe, Gentle Rocking!

- Rocking speed of 10 cycles/min
- Platform moves through a 10° angle from horizontal
- Mechanical timer for time dependent studies
- ☆ 24 Month Warranty

#### Product Description

- Ideal for applications requiring smooth, gentle rocking motion; gel staining, mixing vacutainer tubes or aeration of cultures in petri dishes, culture bottles or flasks.
- Rocking motion is a constant 10 oscillations per minute.
- Versatile—Can be used in incubators and cold rooms between 4°C and 40°C.

- Durable steel body totally encloses motor to prevent hazardous spills onto drive mechanism.
- 11.5" x 13.5" (29 x 34 cm) platform holds samples securely on a removable grooved vinyl pad.
- Platform moves through a 10° angle from horizontal.
- Mechanical timer can be set for 1 to 60 minutes or for continuous operation.

#### PRODUCT SPECIFICATIONS

Type	Maximum Load	Platform Dimensions Inches (cm)	Overall Dimensions Inches (cm)			Shipping Weight
			W	H	D	
4635	5 Lbs (2.3 kg)	14.5 (37) x 11.5 (29)	14.5 (37)	6 (15)	11.5 (29)	21 Lbs (10 kg)
4635-1	5 Lbs (2.3 kg)	14.5 (37) x 11.5 (29)	14.5 (37)	6 (15)	11.5 (29)	21 Lbs (10 kg)

#### ORDERING INFORMATION

Model #	Volts	Amps	Watts	Hz	Speed (cycles per minute)	Rocking Angle Degrees	Ambient Conditions	Timer (min)	US List Price
4635	120	1.0	12	60	10	10 from horizontal	4° to 40°C, 20 to 80% RH non-condensing	1-60	\$1,115.00
4635-1	240	1.0	24	50	10	10 from horizontal	4° to 40°C, 20 to 80% RH non-condensing	1-60	\$1,115.00