




Harvey Sterilizers



Temperature Indicators



Chemitest™



Chemitest Pouches

Simple Daily Monitor for Your Sterilizer!

- Verify your sterilizer's temperature with every cycle

Product Description

Chemitest See-Thru Pouches

- Available in six sizes to meet all needs of the customer.
- Self sealing edges.
- Perforation on adhesive strip ensures a stronger seal.
- Solid rails on sides for strong seal—easy to open.
- Used for steam, ethylene oxide, or chemical vapor sterilization.
- Medical grade paper.

Daily Temperature Verification Accessories

- Available in strips or tape.
- Color change indicates that your sterilizer has reached the proper temperature.
- Indicators assure that items being processed have been run through a cycle, preventing the accidental use of items which have not been processed.
- Alternate the use of Harvey Temperature Indicators with the Harvey sterilization process validation procedure for guaranteed assurance of your sterilizer's performance.

ORDERING INFORMATION

Instrument Packaging for Clean Process and Temperature Verification

Model #	Description	US List Price
230300001	Self-Sealing Chemitest Indicator Paper Bags, Tray Size Box of 150 - 5" W (12.7 cm) x 2" H (5.1 cm) x 15" L (38.1 cm)	\$44.00
230301001	Self-Sealing Chemitest Indicator Paper Bags, Instrument Size Box of 250 - 2.5" W (6.4 cm) x 1.5" H (3.8 cm) x 10.5" L (26.7 cm)	\$40.00
ZSX152	Self-Sealing Chemitest See-Thru Indicator Pouches Box of 200 - 2.25" x 4" (5.7 cm x 10.2 cm)	\$14.10
ZSX153	Self-Sealing Chemitest See-Thru Indicator Pouches Box of 200 - 2.75" x 9" (7 cm x 22.9 cm)	\$15.20
ZSX154	Self-Sealing Chemitest See-Thru Indicator Pouches Box of 200 - 3.5" x 5.25" (8.9 cm x 13.3 cm)	\$14.60
ZSX155	Self-Sealing Chemitest See-Thru Indicator Pouches Box of 200 - 3.5" x 9" (8.9 cm x 22.9 cm)	\$16.00
ZSX156	Self-Sealing Chemitest See-Thru Indicator Pouches Box of 200 - 5.25" x 10" (13.3 cm x 25.4 cm)	\$25.20
ZSX157	Self-Sealing Chemitest See-Thru Indicator Pouches Box of 200 - 7.5" x 13" (19.1 cm x 33 cm)	\$45.30

* See-thru pouches are a combination of mylar and paper. May be sealed with Chemitest indicator tape.

ACCESSORIES

Temperature Verification (Include in every sterilizer load for daily confidence.)

Model #	Description	US List Price
AY759X2	Chemical Indicator Strips - Box of 250, Steam Sterilizers only	\$27.00
260577	Chemitest Indicator Tape-Closure system for packs, bags or other containers, 1.25" (3.2 cm) x 360" L (914.4 mm) Roll	\$29.00
230305001	Chemitest Indicator Strips - Box of 250, Chemiclave Sterilizers only	\$25.00
230306001	Chemitest Indicator Tabs, Roll of 250 self-adhering tabs 1.25" (3.2 cm) sq.	\$25.00