


Single Scale Graduated Cylinders

Roll-Proof and Easy to Read

These graduated cylinders have a single scale with raised numerals for easy reading. Flared top rims and large spouts facilitate pouring, and hexagon-shaped bases provide stability. The "non-wetting" polypropylene surface eliminates concave meniscus and makes reading easier and more accurate. Chemically resistant and autoclavable at 121°C (250°F). **Per each, cases listed below.**



CATALOG NO.	CAPACITY	SUBDIVISION	CASE	PRICE/EACH
F28452-0000	10ml	0.2ml	24/each	\$2.25
F28453-0000	25ml	0.5ml	24/each	3.00
F28454-0000	50ml	1.0ml	12/each	4.00
F28455-0000	100ml	1.0ml	12/each	5.00
F28456-0000	250ml	2.0ml	12/each	5.75
F28457-0000	500ml	5.0ml	6/each	8.00
F28458-0000	1000ml	10.0ml	6/each	11.00
F28459-0000	2000ml	20.0ml	4/each	20.00



Clear Graduated Cylinders

Shatter-Proof and Easy to Read

Clear TPX® (polymethylpentene) cylinders create no concave meniscus, making the raised numeral, single scale graduations easy to read. They are shatter-proof with a wide octagonal base to prevent tipping and rolling. Superior chemical resistance and heat resistant to 177°C (350°F) make them a great value. Steam autoclavable at 121°C (250°F). **Per each, cases listed below.**



CATALOG NO.	CAPACITY	SUBDIVISION	CASE	PRICE/EACH
F28690-0000	10ml	0.2ml	24/each	\$3.90
F28691-0000	25ml	0.5ml	24/each	4.40
F28692-0000	50ml	1.0ml	12/each	5.40
F28693-0000	100ml	1.0ml	12/each	6.45
F28694-0000	250ml	2.0ml	12/each	7.50
F28695-0000	500ml	5.0ml	6/each	9.80
F28696-0000	1000ml	10.0ml	6/each	13.00
F28697-0000	2000ml	20.0ml	4/each	30.50

Made of TPX® Brand PMP for Easy Viewing



Holdfast™ Graduated Cylinders

Grip Handle for Steady Pouring

The handle on a Holdfast™ Cylinder provides a secure, solid grip for superior control. The result is safer, more convenient pouring. The one-piece polypropylene cylinder will not create a meniscus. The "Total Contained, Total Delivered" double scale is accurate and easy to read, and the hexagonal base prevents tipping and rolling. The hand opening is 89 x 29 mm (3 x 1 1/8"). Steam autoclavable at 121°C (250°F).



CATALOG NO.	CAPACITY	PRICE/EACH
F28461-1000	1000ml	\$16.50
F28461-2000	2000ml	32.50

See Beaker Section for Additional Containers, Page 6

For Graduated Pitchers, See Pages 101

