



### Poxygrid® Wire Test Tube Racks

Choice of Four Colors and Multiple Sizes

Chemical and corrosion-resistant steel wire racks are coated with heat cured epoxy. Epoxy coatings are unaffected by organic solvents, salts, weak acids and weak alkalis, will withstand a temperature range of -223°C (-369°F) to +121°C (250°F) and can be autoclaved.

All racks have square holes, three decks, a slotted grid bottom and reinforced edges. Orange racks can be used for biohazard warning, yellow for radioactive warning. **Per each, 24 per case.**

**To Order:**

1. Select a tube size.
2. Determine number of places required.
3. Find corresponding catalog number.
4. Add catalog number suffix for desired color.

**For Example: F18757-0003 is an orange, 72 place rack for 15-16mm tubes.**

Custom Racks Made to Order Call 1 800-4BEL-ART.



CATALOG NO	GREEN	BLUE	YELLOW	ORANGE	HOLDS TUBES	PLACES	ROWS	LENGTH X WIDTH X HEIGHT	PRICE/EACH
F18750-	0000	0001	0002	0003	10-13mm	72	6 x 12	205 x 113 x 64mm (8 <sup>1</sup> / <sub>16</sub> x 4 <sup>3</sup> / <sub>16</sub> x 2 <sup>1</sup> / <sub>2</sub> "")	\$22.70
F18753-	0000	0001	not available	0003	10-13mm	108	9 x 12	205 x 164 x 64mm (8 <sup>1</sup> / <sub>16</sub> x 6 <sup>1</sup> / <sub>16</sub> x 2 <sup>1</sup> / <sub>2</sub> "")	24.75
F18754-	0000	0001	not available	not available	15-16mm	40	4 x 10	203 x 92 x 64mm (8 x 3 <sup>3</sup> / <sub>16</sub> x 2 <sup>1</sup> / <sub>2</sub> "")	21.90
F18755-	0000	0001	0002	0003	15-16mm	48	4 x 12	242 x 92 x 64mm (9 <sup>1</sup> / <sub>2</sub> x 3 <sup>3</sup> / <sub>16</sub> x 2 <sup>1</sup> / <sub>2</sub> "")	21.90
F18757-	0000	0001	0002	0003	15-16mm	72	6 x 12	242 x 130 x 64mm (9 <sup>1</sup> / <sub>2</sub> x 5 <sup>1</sup> / <sub>16</sub> x 2 <sup>1</sup> / <sub>2</sub> "")	21.90
F18759-	0000	0001	0002	0003	15-16mm	108	9 x 12	242 x 189 x 64mm (9 <sup>1</sup> / <sub>2</sub> x 7 <sup>1</sup> / <sub>16</sub> x 2 <sup>1</sup> / <sub>2</sub> "")	24.75
F18762-	0000	0001	0002	0003	18-20mm	40	4 x 10	243 x 110 x 83mm (9 <sup>5</sup> / <sub>16</sub> x 4 <sup>1</sup> / <sub>4</sub> x 3 <sup>1</sup> / <sub>4</sub> "")	21.90
F18764-	0000	0001	not available	not available	18-20mm	80	8 x 10	243 x 203 x 83mm (9 <sup>5</sup> / <sub>16</sub> x 8 x 3 <sup>1</sup> / <sub>4</sub> "")	24.75
F18767-	0000	0001	not available	not available	20-24mm	40	4 x 10	263 x 116 x 83mm (10 <sup>1</sup> / <sub>2</sub> x 4 <sup>3</sup> / <sub>16</sub> x 3 <sup>1</sup> / <sub>4</sub> "")	24.75
F18772-	0000	0001	0002	0003	22-25mm	40	4 x 10	293 x 128 x 83mm (11 <sup>1</sup> / <sub>16</sub> x 5 <sup>1</sup> / <sub>16</sub> x 3 <sup>1</sup> / <sub>4</sub> "")	22.70
F18773-	0000	0001	0002	0003	22-25mm	80	8 x 10	293 x 242 x 83mm (11 <sup>1</sup> / <sub>16</sub> x 9 <sup>5</sup> / <sub>16</sub> x 3 <sup>1</sup> / <sub>4</sub> "")	24.75



### Poxygrid® "Half - Size" Test Tube Racks

When Space is Limited... Use Half-Rack

Excellent space-savers for use in water baths, under fume hoods, and any place where space is limited. Epoxy coated and available in three colors: orange for biohazard warning, green and blue for general lab use.

Per each, 24 per case.



CATALOG NO.	COLOR	LENGTH X WIDTH X HEIGHT	PLACES	ROWS	HOLDS TUBES	PRICE EACH
F18788-1300	Green	113 x 107 x 64mm (4 <sup>1</sup> / <sub>16</sub> x 4 <sup>3</sup> / <sub>16</sub> x 2 <sup>1</sup> / <sub>2</sub> "")	36	6 x 6	10-13mm	\$16.50
F18788-1301	Blue	113 x 107 x 64mm (4 <sup>1</sup> / <sub>16</sub> x 4 <sup>3</sup> / <sub>16</sub> x 2 <sup>1</sup> / <sub>2</sub> "")	36	6 x 6	10-13mm	16.50
F18788-1600	Green	125 x 99 x 64mm (4 <sup>1</sup> / <sub>16</sub> x 3 <sup>3</sup> / <sub>16</sub> x 2 <sup>1</sup> / <sub>2</sub> "")	24	6 x 4	15-16mm	16.50
F18788-1601	Blue	125 x 99 x 64mm (4 <sup>1</sup> / <sub>16</sub> x 3 <sup>3</sup> / <sub>16</sub> x 2 <sup>1</sup> / <sub>2</sub> "")	24	6 x 4	15-16mm	16.50
F18788-1603	Orange	125 x 99 x 64mm (4 <sup>1</sup> / <sub>16</sub> x 3 <sup>3</sup> / <sub>16</sub> x 2 <sup>1</sup> / <sub>2</sub> "")	24	6 x 4	15-16mm	16.50
F18788-2000	Green	125 x 108 x 64mm (4 <sup>1</sup> / <sub>16</sub> x 4 <sup>1</sup> / <sub>4</sub> x 2 <sup>1</sup> / <sub>2</sub> "")	20	5 x 4	18-20mm	16.50
F18788-2001	Blue	125 x 108 x 64mm (4 <sup>1</sup> / <sub>16</sub> x 4 <sup>1</sup> / <sub>4</sub> x 2 <sup>1</sup> / <sub>2</sub> "")	20	5 x 4	18-20mm	16.50
F18788-2003	Orange	125 x 108 x 64mm (4 <sup>1</sup> / <sub>16</sub> x 4 <sup>1</sup> / <sub>4</sub> x 2 <sup>1</sup> / <sub>2</sub> "")	20	5 x 4	18-20mm	16.50

### Poxygrid® "Rack And A Half"

For Volume Work

The same quality rack as the Scienceware® F18750 through F18788 series in a convenient extra-long size for volume test tube work. **Per each, 12 per case.**



CATALOG NO.	LENGTH X WIDTH X HEIGHT	PLACES	ROWS	HOLDS TUBES	PRICE/EACH
F18786-0750	331 x 110 x 60mm (13 <sup>1</sup> / <sub>2</sub> x 4 <sup>3</sup> / <sub>16</sub> x 2 <sup>1</sup> / <sub>2</sub> "")	120	6x20	10-13mm	\$24.75
F18786-0770	386 x 183 x 64mm (15 <sup>1</sup> / <sub>16</sub> x 7 <sup>1</sup> / <sub>16</sub> x 2 <sup>1</sup> / <sub>2</sub> "")	180	9x20	15-16mm	24.75
F18786-0780	386 x 108 x 64mm (16 <sup>1</sup> / <sub>16</sub> x 4 <sup>1</sup> / <sub>4</sub> x 2 <sup>1</sup> / <sub>2</sub> "")	100	5x20	15-16mm	24.75

### Poxygrid® 96-Place Test Tube Rack



Now you can maintain the 96-place format of high throughput applications with large, 15ml test tube capacity.

The Poxygrid® 96-Place Test Tube Rack is ideal for mimicking the 96-place format used in microplates, deep-well plates and other tools configured for drug discovery and robotic

applications. Holds 96 tubes, 15ml in capacity, in an 8 x 12 configuration. Flow-through design allows rapid temperature changes of samples by varying the water bath temperature. Rigid wire construction with blue epoxy finish. Steam autoclavable at 121°C (250°F).



CATALOG NO.	LENGTH X WIDTH X HEIGHT	PRICE/EACH
F18765-0001	242 x 170 x 64mm (9 <sup>1</sup> / <sub>2</sub> x 6 <sup>1</sup> / <sub>16</sub> x 2 <sup>1</sup> / <sub>2</sub> "")	\$28.35

Refer to Tube Sizing Chart on Page 173 to Select the Ideal Scienceware® Rack for Your Application