

Polygon Spinbar® Magnetic Stirring Bars




Teflon® PTFE

Per each, 50 per case.

CATALOG NO.	SIZE	PRICE/EACH
F37120-0012	12.7 x 8mm (½ x ⅝")	\$3.50
F37120-0058	15.9 x 8mm (⅝ x ⅝")	3.25
F37120-0078	19.1 x 9.5mm (¾ x ⅜")	4.10
F37120-7836	22.2 x 4.7mm (¾ x ⅜")	3.55
F37120-0038	25.4 x 9.5mm (1 x ⅜")	4.30
F37120-0112	38.1 x 9.5mm (1½ x ⅜")	5.10
F37120-0002	50.8 x 9.5mm (2 x ⅜")	6.85
F37120-0212	63.5 x 12.7mm (2½ x ½")	14.00
F37120-0003	76.2 x 12.7mm (3 x ½")	15.35

Polygon Spinbar® Magnetic Stirring Bars

F37120-0010	10 x 6mm (⅜ x ¼")	\$3.80
F37120-0125	12 x 5mm (½ x ⅜")	3.50
F37120-0015	15 x 6mm (⅝ x ¼")	4.00
F37120-0020	20 x 8mm (¾ x ⅝")	4.25
F37120-0025	25 x 8mm (1 x ⅝")	4.35
F37120-0030	30 x 8mm (1⅞ x ⅝")	4.50
F37120-0040	40 x 8mm (1⅞ x ⅝")	5.10
F37120-0050	50 x 8mm (2 x ⅝")	5.60
F37120-0060	60 x 7mm (2⅞ x ¼")	5.70
F37120-0070	70 x 10mm (2⅞ x ⅜")	9.55
F37120-0080	80 x 10mm (3⅞ x ⅜")	9.90




Giant Polygon Spinbar® Magnetic Stirring Bars

Teflon® PTFE Bars for Large Volume Stirring

Heavy duty Spinbar® Magnetic Stirring Bars for use in stirring large volumes of fluid in vessels on large magnetic stirring machines such as the Scienceware® H37028 or other large machines. Will fit into minimum vessel neck diameter of 32 mm. **Per each, 10 per case.**

CATALOG NO.	SIZE	PRICE/EACH
F37118-0002	57 x 27mm (2¼ x 1⅞")	\$80.00
F37118-0004	108 x 27mm (4¼ x 1⅞")	88.00
F37118-0006	159 x 27mm (6¼ x 1⅞")	106.00



Circulus™ Magnetic Stirring Bars

Teflon® PTFE

The Circulus™ Magnetic Stirring Bar offers fast mixing action, high turbulence, greater stability, and less drifting at higher speeds. Ideal for all "problem" liquids, highly viscous mixtures, metal particles in fluids, and containers with round or convex bottoms. All sizes have solid Teflon® PTFE end disks 20mm (.79") in diameter and a bar diameter of 8mm (.315). **Per each, 50 per case.**



CATALOG NO.	COLOR	SIZE	PRICE/EACH
F37170-0000	white	32mm long (1¼")	\$9.30
F37170-0001	red	32mm long (1¼")	9.30
F37170-0002	yellow	32mm long (1¼")	9.30
F37170-0003	blue	32mm long (1¼")	9.30
F37172-0000	white	52mm long (2")	10.25
F37172-0001	red	52mm long (2")	10.25
F37172-0002	yellow	52mm long (2")	10.25
F37172-0003	blue	52mm long (2")	10.25

Round Spinbar® Magnetic Stirring Bars

Pyrex® Glass

Glass stirring bars are utilized when high temperature ranges, up to 274°C (525°F) are necessary and when "zero absorption" of the stirred solution is preferred. The glass encapsulation of the magnet is completely smooth, making the bars easy to clean off. **Per each, 50 per case.**



CATALOG NO.	SIZE	PRICE/EACH
F37101-0012	12.7 x 9.5mm (½ x ⅜")	\$11.00
F37101-7814	22 x 6.4mm (⅝ x ¼")	10.30
F37101-0001	25.4 x 9.5mm (1 x ⅜")	12.65
F37101-0112	38 x 9.5mm (1½ x ⅜")	14.20

Spinning® Magnetic Stirring Bars

Teflon® PTFE

Octagon Spinbar® Magnetic Stirring Bar with a removable stabilizing ring to keep it centered in the vessel. Creates a deep vortex without splashing or chattering and resists spinout. **Per each, 50 per case.**



CATALOG NO.	BAR SIZE	RING I.D.	PRICE/EACH
F37140-0005	12.7 x 8mm (½ x ⅝")	19mm (¾")	\$22.00
F37140-0010	25.4 x 8mm (1 x ⅝")	31.8mm (1¼")	25.00
F37140-0015	38.1 x 8mm (1½ x ⅝")	44.5mm (1⅞")	25.75
F37140-0020	50.8 x 8mm (2 x ⅝")	57.2mm (2¼")	28.00

Disposable Spinbar® Magnetic Stirring Bars

Economically Avoid Cross Contamination

Disposable Spinbar® Magnetic Stirring Bars are ideal for single use applications or when cross contamination might be a concern. All sizes are offered in economical bags of 100 bars to ensure a large supply for most users. The multifaceted stirring bar surface adds turbulence relative to similar size smooth cylindrical stirring bars. **100 per bag.**

CATALOG NO.	SIZE	PRICE/BAG
F37112-0012	12.7 x 3.2 mm (½ x ⅛")	\$77.00
F37112-0001	25.4 x 8 mm (1 x ⅝")	108.00
F37112-0016	41.3 x 8 mm (1⅞ x ⅝")	130.00
F37112-0002	50.8 x 8 mm (2 x ⅝")	180.00