

Safe-Lab® Stirrer Bearings

Patent Number D294, 969

Shown with Stopper-Loc® Nut

Thread size for use with Stopper-Loc® Nut. Allows stirrer bearing to be locked to outer ground joint when stirring viscous materials for long periods

of time. Sold complete with Stopper-Loc® Nut and two Teflon® PTFE inserts. Not autoclavable.



CATALOG NO.	JOINT SIZE	PRICE/EACH
F20320-1661	24/40	\$77.50
F20320-1662	29/42	89.00



Stopper-Loc® Nut

A dual-action device that locks and unlocks threaded products to and from outer ground joints. Autoclavable polypropylene nut fits Safe-Lab® inner ground joints, glass adapters, condensers, funnels, stoppers and any other products used without sleeves. **6 per pack.**



CATALOG NO.	LENGTH	STL THREAD SIZE	USE WITH	PRICE/PACK
F20202-9205	short	19	19/22, 19/38	\$13.25



Safe-Lab® Joint Vacuum Adapters

Teflon® PTFE body with Viton® O-ring. Teflon® hose connection for vacuum work or introduction and removal of gases. Rotate threaded polypropylene extracting nut to remove adapter from flask. Not autoclavable.



CATALOG NO.	JOINT SIZE	PRICE/EACH
F20314-1690	14/20	\$53.75
F20314-1720	19/22	56.50
F20314-1729	24/40	63.50

Safe-Lab® Joint Stirrer Bearing

Teflon® PTFE body with Viton® O-ring for use with 10 mm stirring rods. Teflon® PTFE insert has inside O-ring for fine adjustment of stirring rod tension. Threaded polypropylene locking compression cap on top. Extracting nut above ground joint allows easy and quick removal from flask. Comes complete with two Teflon® PTFE inserts. Not autoclavable.



CATALOG NO.	JOINT SIZE	PRICE/EACH
F20320-1660	24/40	\$65.00

Safe-Lab® Teflon® PTFE Bearing Insert

Replacement Teflon® PTFE insert with Viton® O-ring for 10 mm Teflon® PTFE stirrer bearings. Not autoclavable. **2 per pack.**



CATALOG NO.	PRICE/PACK
F20321-1665	\$22.95

Cup-Type Stoppers

Double as Containers for Micro Work


These polyethylene cup-type stoppers fit tightly, yet remove easily. All are 25mm high. **24 per bag, 6 bags per case.**



CATALOG NO.	SIZE NO.	TOP	BOTTOM	PRICE/BAG
H20193-0000	00	13mm	7mm	\$21.00
H20193-0001	1	17mm	11mm	21.00
H20193-0002	2	20mm	14mm	21.00
H20193-0003	3	23mm	17mm	21.00
H20193-0004	4	26mm	20mm	21.00
H20193-0006	6	32mm	26mm	21.00
H20193-0008	8	38mm	32mm	21.00

Flaskup™ Flask Holders

Patent Number D307,869



Holds Round Bottom Flasks

Unique flask holders are molded from glass-filled polypropylene, and are durable as well as easy to clean. The inside wall has a slight taper to the open bottom, allowing a variety of flask shapes to fit into one holder. The wide base keeps flasks stable.



Side slots can be used as sight windows. Autoclavable at 121°C (250°F). **3 per package, 8 packages per case, except assortment.**



CATALOG NO.	FITS FLASK SIZE	PRICE/PACKAGE
F38951-2000	10ml	\$ 11.85/3
F38951-2002	25ml	11.85/3
F38951-2004	50ml	12.00/3
F38951-2006	100ml	13.95/3
F38951-2010	250ml	15.20/3
F38951-2012	500ml	17.00/3
F38951-2008	Assortment Pack; 3 each 10, 25, 50, and 100ml sizes	48.95/12