


Precision Products for Making the Best Beer and Wine



HANNA
instruments

Manufacturers Since 1978

pH Meters



pHep[®] 4 & pHep[®] 5

Revolutionary NEW high accuracy pH testers with replaceable electrodes and temperature in °C and °F.

- High Resolution, Large LCD
- Waterproof Casing
- Percent Battery Power Illumination
- Replaceable pH Cartridge
- Temperature in °C and °F
- Tri-level LCD
- Automatic Calibration

HI 98127..... w/0.1 pH Resolution

HI 98128..... w/0.01 pH Resolution



Checker[®]

The Checker family of pH meters offers a reliable measurement at an affordable price.

- pH Range 0.00 to 14.00
- Replaceable pH Electrode

HI 98103Checker[®] 1

HI 98104Checker[®] 2

HI 98105Checker[®] 3

Temperature

Checktemp®

Professional high accuracy pocket thermometers with exclusive calibration check device.

- Celsius (-50.0 to 150.0°) or Fahrenheit (-58.0 to 302.0°)
- Stainless Steel Probe

HI 98501 °C

HI 98502 °F

Checktemp®-Dip

- Celsius (-20 to 100°) or Fahrenheit (-4 to 212°)
- 3000 Hours of Battery Life from One AAA Battery
- Stainless Steel Probe

HI 98509-01 °C w/10 ft. Cable

HI 98510-01 °F w/10 ft. Cable



HI 190M

Compact magnetic mini-stirrers with electronic controls that allow the user to regulate speed with greater precision.

- 1 Liter (1/4 Gallon) Capacity
- Minimum 100 rpm/Maximum 800 to 1000 rpm
- Spill Resistant ABS Plastic Cover

HI 190M



Mixing

pH Meter Accessories

HIGH ACCURACY CHEMICAL TEST KITS



PROFESSIONALLY MADE

RANGE

HI 3820 Acidity Test Kit	0-500 ppm
HI 3841 Hardness Test Kit	40-500 ppm
HI 3835 Salinity Test Kit	0-40 g/Kg
HI 3874 Nitrate Test Kit	0-50 ppm
HI 38001 Sulfate Test Kit	100-1000ppm
HI 3834 Iron Test Kit	0-5 ppm

CALIBRATION, STORAGE, AND CLEANING



CONVENIENT SINGLE USE 20 ML SACHETS

HI 70004P	Pack of 25 Sachets, pH 4.01
HI 70007P	Pack of 25 Sachets, pH 7.01
HI 70010P	Pack of 25 Sachets, pH 10.01
HI 7061L	Cleaning Solution for pH Electrodes
HI 70300L ..	Storage Solution for pH Electrodes

Authorized Distributor www.clarksonlab.com

HANNA
instruments
Manufacturers Since 1978

Certified
ISO 9001