



Product	Grade Certifications	(CH <sub>3</sub> ) <sub>2</sub> CO	CAS 67-64-1	FW 58.08
---------	----------------------	------------------------------------	-------------	----------

### Acetone - Lab Grade - Hospital Clinical Grade

Catalog Number	Package Sizes	Product Specifications	Limits
<b>329000LAB</b>	4x1 gallon poly	Assay (GC), min	95.5%
	5 gallon Cube w/ sp.	Water, max	0.5%
	5 gallon Shelf w/sp.	Appearance	Clear
	5 gallon Metal	Color (APHA), max	10
	5 gallon Poly	Residue After Evaporation, max	10 ppm
		Aldehydes & Alcohols	To Pass Test
	55 gallon New Lined	Density (g/ml) @ 25C, max	0.7857
	55 gal Poly	Micron Filtered	
		6.63 lbs/gal	360 lbs. net
		<b>* Note: Not intended for use as a medical device -</b>	

**When this product is ordered the ACS grade and label is substituted!**

***Approved by Paul DiMarco, Vice President Quality Management***

58 Vale Road, Brookfield, CT 06804, USA  
 Telephone: (203) 740-3471 ext. 134, 800-243-5360 ext. 134, Fax (203) 740-3481  
 E-Mail: [paul@pharmco-prod.com](mailto:paul@pharmco-prod.com)  
 Website: [www.pharmco-prod.com](http://www.pharmco-prod.com)