



PRODUCT SPECIFICATION SHEET

Sec-Butyl Alcohol

Reagent Grade

Catalog No: 305000REA

Available in the following sizes:

4x4 liter Amber Glass, 5 gallon Metal w/Rieke Spout, 55 gallon New Metal Drum

PRODUCT SPECIFICATIONS	LIMITS	Typical Lot Analysis
Assay	NLT 99.0, wt %	99.98 wt %
Color, Pt-Co	NMT 10	<5
Water Content	0.20% max	0.015%
Specific Gravity @ 20°C	0.807-0.809	0.808
Odor	Characteristic	Pass
Appearance	Clear	Pass

Form Sec-Butyl Alcohol, Reagent, #101, Rev. 2.0, 10/01

Approved by: Paul DiMarco, Vice President Quality Management

FOR INDUSTRIAL OR LAB USE ONLY. Not intended as an active substance in Food or Drug. Not to be considered a Medical Device. Not intended for use as a Disinfectant as defined by the EPA. The expiration date of this product is three years from the date of manufacture. (Rev. # Disclaimer only: Rev 2.0 11/19/01 PJM)

This document has been electronically generated and is therefore not signed. A signed copy can be obtained from Pharmco's Quality Control office.

58 Vale Road, Brookfield, CT 06804, USA
Telephone: (203) 740-3471 Fax (203) 740-3481
E-Mail: paul@pharmco-prod.com
Web Site: www.pharmco-prod.com

DOT INFORMATION
 Butanols
 Hazard Class: 3
 UN 1120
 PG: II
 Butyl Alcohol (sec)
 CAS 78-92-2

LOT NO:
 BATCH NO:

 **pharmco**
 PRODUCTS INC
 58 Vale Road, Brookfield, Connecticut 06804 USA
 Phone (203) 740-3471 FAX (203) 740-3481 DSP-CT-18

**BUTYL ALCOHOL (Sec)
 REAGENT GRADE**

CATALOG NUMBER: 305000REA 5 Gallons

FLAMMABLE LIQUID
 CHEMTREC NO. 1-800-424-9300 3

Product Specifications	Limits	Typical Lot Analysis
Assay	NLT 99.9, wt%	99.99 wt %
Color, Pt-Co	NMT 10	5
Water Content	0.20% max.	0.000%
Specific Gravity @ 20C	0.807-0.809	0.8077
Odor	Characteristic	Pass
Appearance	Clear	Pass

  

DANGER: FLAMMABLE LIQUID & VAPOR. Keep away from all sources of ignition. Improper Exposure Can Result in Injury or Death! Wear all recommended protective equipment. Refer to MSDS for additional information. *

STORAGE AND HANDLING: Keep this seal, dry place. Do not expose to direct sunlight. Use only with proper personal protective equipment and ventilation. Refer to MSDS for additional information. *

FIRST AID: SEEK IMMEDIATE MEDICAL ATTENTION. REFER TO MSDS FOR APPROPRIATE FIRST AID FOR THIS CHEMICAL.

REFUSE TO USE PRIOR TO USE. Use Only Under Appropriate Conditions. Use Only Under Trained Supervision. *

FOR INDUSTRIAL OR LAB USE ONLY. Not intended as an active substance in Food or Drug. Not to be considered a Medical Device. Not intended for use as a Disinfectant as defined by the EPA. *

800-424-9300

Approved by: Paul DiMarco, Vice President Quality Management

FOR INDUSTRIAL OR LAB USE ONLY. Not intended as an active substance in Food or Drug. Not to be considered a Medical Device. Not intended for use as a Disinfectant as defined by the EPA. The expiration date of this product is three years from the date of manufacture. (Rev. # Disclaimer only: Rev 2.0 11/19/01 PJM)

This document has been electronically generated and is therefore not signed. A signed copy can be obtained from Pharmco's Quality Control office.

58 Vale Road, Brookfield, CT 06804, USA
 Telephone: (203) 740-3471 Fax (203) 740-3481
 E-Mail: paul@pharmco-prod.com
 Web Site: www.pharmco-prod.com