



810012



810013R



810018



810022



810029AR



810029R



810033




810035



810036



810037R



810050

Timers SPER SCIENTIFIC

Sper Scientific professional laboratory timers use only highly accurate quartz crystal movements guaranteed to be within ± 3 seconds per day. Certificates of traceability to NIST are available for all models. To order certified add "C" to the end of the part number. At Sper Scientific even our least expensive timers undergo extensive quality control and come with a 5 year warranty.

Catalog Number	Retail	Description	Count Up Chnl.	Count Down Chnl.	Prev. Setting Memory	Clock	Alarm Clock	Range	Resolution
810044	\$9.00	Digital Hourglass Timer	0	1	Yes	Yes	No	100 Mn	1 Sec
810017	\$11.00	Project Timer	1	2	No	Yes	No	24Hrs/Clock Up to 8 years	1/100 Sec
810040	\$11.00	Visual Alarm Timer	0	1	Yes	No	No	99 Mn	1 Sec
810041	\$11.00	Visual Alarm Timer	0	1	Yes	No	No	99 Hrs	1 Sec
810042	\$11.00	Visual Alarm Timer	1	1	Yes	No	No	9999 Mn	1 Mn/1 Sec
810043	\$11.00	Visual Alarm Timer	1	1	Yes	No	No	9999 Hrs	1 Hr/1 Sec
810021	\$12.00	Mini Timer	0	1	Yes	No	No	16 Hrs	1 Min
810019	\$14.00	Mini Timer	0	1	Yes	No	No	99 Mn	1 Sec
810026R	\$16.00	Count Up/Count Down	1	1	Yes	Yes	Yes	100 Hrs	1 Sec
810024R	\$18.00	Count Up/Count Down	1	1	Yes	No	No	99 Mn	1 Sec
810050	\$19.00	TIME-IT Software	8	8	Yes	Yes	Yes	99 Hrs	1/1000 Sec
810015	\$22.00	5 Channel Timer	4	4	Yes	Yes	No	100 Hrs	1 Sec
810027	\$22.00	Simultaneous Start Timer	4 1	4 0	Yes	Yes	Yes	100 Hrs 10 Hrs	1 Sec 1/100 Sec
810025R	\$24.00	Stopwatch Timer	1	1	Yes	No	No	100 Hrs	1/100 Sec
810005	\$25.00	White Large Display Timer	1	2	Yes	Yes	No	24 Hrs	1 Sec
810006	\$25.00	Black Large Display Timer	1	2	Yes	Yes	No	24 Hrs	1 Sec
810038	\$39.00	Talking Timer	1	1	Yes	Yes	No	24 Hrs	1 Sec



810005

810006

810015

810017

810019/21



810024/25/26R

810027

810038

810040/41/42/43

810044

810050

Stopwatches

SPER SCIENTIFIC

Sper Scientific professional laboratory stopwatches use only highly accurate quartz crystal movements guaranteed to be within ± 3 seconds per day. Certificates of traceability to NIST are available for all models. To order certified add C to the end of the part number. At Sper Scientific even our least expensive stopwatches undergo extensive quality control and come with a 5 year warranty.

Catalog Number	Retail	Description	Count Up Chnl.	Count Down Chnl.	Clock	Alarm Clock	Range	Resolution
810012	\$8.00	Water Resistant Stopwatch	1	0	Yes	Yes	24 Hrs	1/100 Sec
810018	\$17.00	Stopwatch Recorder	1	0	Yes	Yes	24 Hrs	1/100 Sec
810050	\$19.00	TIME-IT Software	8	8	Yes	Yes	99 Hrs	1/100 Sec
810013R	\$19.95	Lightweight Stopwatch	1	0	Yes	Yes	24 Hrs	1/100 Sec
810033	\$24.00	60 Memory Stopwatch	2	1	Yes	Yes	10 Hrs	1/100 Sec
810022	\$25.00	Large Display Stopwatch	1	0	Yes	No	24 Hrs	1/100 Sec
810029R	\$27.00	8 Memory Stopwatch	2	1	No	No	10 Hrs	1/100 Sec
810035	\$37.00	Large Display Water Resistant Stopwatch	1	0	Yes	Yes	24 Hrs	1/100 Sec
810036	\$38.00	100 Memory Water Resistant Stopwatch	1	0	Yes	No	10 Hrs	1/100 Sec
810037	\$43.00	Electro-Luminescent Stopwatch	1	0	Yes	Yes	24 Hrs	1/100 Sec
810029AR	\$65.00	Observational Research Stopwatch	2	1	No	No	10 Hrs	1/100 Sec