



STOPCLOCK 810032

Ruggedly designed for prolonged use as both a bench-top count-up timer and a stopclock. Non-skid feet keep the unit firmly in place during usage, yet it is small and light enough to be easily moved around the lab as needed. A large LCD and buttons make it easy to operate. Counts-up to 10 hours with 1/100th second resolution in selectable cumulative or interval split modes.

Comes complete with one AA battery and two leads for connection to an external trigger such as a photo gate. Dim: 4½" x 4" x 2" (120 x 100 x 50 mm). Weight: 7 oz (210 g).



With NIST traceable certificate of compliance 810032C