

FREEDOM™

Freedom™ Microscope Features

- Free from cords
- Free from burned out bulbs
- Free from the risk of A/C current
- Free from charging downtime
- Free to use anywhere, totally portable
- perfect for field use
- Free from heat

Its affordability frees your science budget.

Specifications

- Life-lite illumination. Lamp life over 100,000 hours (Brilliant white light LED for years of trouble free illumination).
- Battery powered, just three simple AA batteries for one full school year of power (over 240 hours with a single battery change). A/C adapter is available as an option.
- Safer!! No need for risky power cords and A/C current. All student-proof design.
- Portable: Use at any desk, workstation or in the field.
- Locked battery compartment so that little fingers are locked out.
- The highest quality optics and construction for which UNICO has become known for over a decade.
- Triple nosepiece with click stop.
- Widefield 10X eyepiece, rotatable with safety lock and built-in pointer.
- 4X-10X-40XR DIN objectives.
- Monocular head inclined at 45 degrees, rotates 360 degrees for ease of use.
- Large stage, acid and reagent resistant with lock-on stage clips and safety stop.
- Stage drilled to accept optional attachable mechanical stage.
- Separate coarse and fine adjustments.
- Built-in N.A. 0.65 condenser with 5 - position disc diaphragm.
- Optional A/C Adapter

New Generation

- now with new revolutionary illumination



Actual Size: 8 1/2" x 5" x 11 1/2"
Shipping Dimensions : 12" x 8" x 16"
Shipping Weight: 9 lbs.



United Products & Instruments, Inc.

182 Ridge Road, Suite E, Dayton, NJ 08810 USA Tel: 800-588-9776 / 732-274-1155 Fax: 732-274-1151
E-MAIL: sales@unico1.com Visit us at www.unico1.com