



L-TTR200 Test Tube Rocker



The UNICO® model TTR200 Test Tube Rocker is an almost silent heavy-duty unit designed for years of everyday use.

Features

- Capacity of 15 tubes or mix of tubes and vials
- 25 rocking cycles per minute
- Reversible and auto-clavable tube pad

Size	13.7 in x 4.7 in x 5.0 in (347 mm x 119 mm x 126 mm)
Shipping Size	15.7 in x 6.7 in x 6.9 in (400 mm x 170 mm x 175 mm)
Net Weight	4.8 lbs (2.2 kg)
Shipping Weight	5.8 lbs (2.6 kg)
Limited Warranty	One year

United Products & Instruments, Inc. 182 Ridge Road, Suite E, Dayton, New Jersey, 08810, USA
Tel: 1-800-588-9776 Fax: 732-274-1151 E-mail: sales@unicol.com Website: <http://www.unicol.com>