

SELECTION GUIDE

Refrigerators/Freezers



**HAZARDOUS ENVIRONMENT?
SAFETY FIRST!**



SAFE COLD STORAGE!



SIZES FOR YOUR EVERY NEED!



**OUR LARGEST COLD STORAGE
UNITS!**

REFRIGERATORS	Product	Total Volume cu. ft.	Cabinet Style	Page
HAZARDOUS ENVIRONMENT? SAFETY FIRST!	Explosion Proof	5.6 20	Undercounter Upright	266
SAFE COLD STORAGE!	Flammable Material Storage	5.6 20	Undercounter Upright	268
SIZES FOR YOUR EVERY NEED!	General Purpose	5.6 20	Undercounter Upright	270
READY TO USE!	Chromatography	23, 47, 69	Upright	273
OUR LARGEST COLD STORAGE UNITS!	Large Capacity	23, 49, 72	Upright	275
FREEZERS				
HAZARDOUS ENVIRONMENT? SAFETY FIRST!	Explosion Proof	5.6 6, 18 20	Undercounter Chest Upright	266
SAFE COLD STORAGE!	Flammable Material Storage	5.6 20	Undercounter Upright	268
SIZES FOR YOUR EVERY NEED!	General Purpose	1.9 5.6 6, 18 20	Undercounter Chest Upright Benchtop	270
REACH IN FOR EASY ACCESS!	Large Capacity	23, 49, 72	Upright	274
REFRIGERATOR/FREEZER COMBINATIONS				
HAZARDOUS ENVIRONMENT? SAFETY FIRST!	Explosion Proof	5.6 10, 11	Upright Undercounter	266
SAFE COLD STORAGE!	Flammable Material Storage	5.6 10, 11	Upright Undercounter	268
SIZES FOR YOUR EVERY NEED!	General Purpose	5.6 10, 14.4	Undercounter Upright	270

Refrigerators / Freezers

EXPLOSION - PROOF

Lab-Line

Model 3560



HAZARDOUS ENVIRONMENT? SAFETY FIRST!

- — Engineered to prevent triggering an explosion inside or outside the unit
- — Suitable for use in Class I, Group C and D hazardous environments
- — Meets OSHA and National Fire Protection Association safety requirements



PRODUCT DESCRIPTION

- Available in a wide range of sizes and styles, including uprights, combination refrigerator/freezers, undercounter and horizontal chest units.
- Thermostat and compressor components are isolated within a vapor-proof enclosure.
- Electrical connections are hardwired to prevent sparking and accidental disconnection.
- Magnetic vinyl gaskets for positive door seal.
- CFC-free insulation and refrigerant.

Model 3560

- 20 cubic ft. corrosion-resistant model designed to store caustic materials and acids.
- Five adjustable stainless steel shelves.
- Enameled steel interior and exterior with stainless steel door liner for greater corrosion resistance.
- Built-in key lock.
- UL listed, CSA certified.

Model 3551

- 11 cubic ft. combination refrigerator/freezer with separate doors.
- Three adjustable shelves plus four door shelves for maximum storage space.
- Hydraulic thermostat; adjustable natural air flow vent.
- Manual defrost.
- UL listed, CSA certified.

Model 3552

- 20 cubic ft. upright freezer.
- Built-in evaporators in shelves for even cooling.
- Four permanently positioned wire shelves.
- Door equipped with key lock.
- Manual defrost.
- UL listed, CSA certified.

Model 3557

- 5.6 cubic ft. undercounter refrigerator.
- Hydraulic thermostat.
- Three adjustable wire shelves plus three door shelves.
- Off cycle self defrost.
- UL listed.

Model 3557-4

- 5.6 cubic ft. undercounter refrigerator/freezer.
- Hydraulic thermostat.
- Two adjustable wire shelves, three door shelves.
- ABS plastic interiors with enameled steel exterior.
- Manual defrost.

Model 3557-2

- 5.6 cubic ft. undercounter freezer.
- Hydraulic thermostat.
- Three fixed shelves, three door shelves.
- Aluminum evaporator coil.
- Off cycle defrost.

Model 3555

- 18 cubic ft. horizontal chest freezer.
- Storage basket.
- Built-in key lock in lid.
- Manual defrost.

Model 3555-6

- 7 cubic ft. horizontal chest freezer.
- Storage basket.
- Cover with key lock.
- UL listed, CSA certified.

Model 3566

- 20 cubic ft. one-door refrigerator.
- Copper evaporator coil for efficient thermal transfer.
- Five adjustable shelves.
- Door with built-in key lock.
- UL listed, CSA certified.

Model 3559

- 10 cubic ft. single door refrigerator/freezer with interior freezer.
- Three adjustable shelves plus four door shelves.
- Hydraulic thermostat.
- Manual defrost.
- UL listed, CSA certified.

Model 3551



Refrigerators / Freezers

EXPLOSION-PROOF

Lab-Line

PRODUCT SPECIFICATIONS

Model #	Chamber Dimensions Inches (cm) Refrigerator Section			Chamber Dimensions Inches (cm) Freezer Section			Total Volume Cubic ft.	Exterior Dimensions Inches (cm)			Shipping Weight Lb. (kg)
	W	H	D	W	H	D		W	H	D	
3551, 3551-1	20.75 (53)	41.5 (105)	14 (36)	19 (48)	12.75 (32)	15 (38)	11	23.25 (59)	63 (160)	24.5 (62)	150(68)
3552, 3552-1		N/A		27 (69)	57 (145)	18 (46)	20	32 (81)	71.5 (182)	33 (84)	240 (109)
3555, 3555-1		N/A		50.25 (128)	28 (71)	22.25 (57)	18	60.25 (153)	35 (89)	29.5 (75)	200 (91)
3555-6, 3555-7		N/A		26.25 (67)	28 (71)	16.5 (42)	6	31.5 (80)	35 (89)	26.25 (67)	125 (57)
3557	19 (48)	27.5 (70)	13 (33)		N/A		5.6	23.5 (59)	33.5 (85)	23.75 (60)	120 (55)
3557-2		N/A		18 (46)	27.25 (70)	15 (38)	5.6	21.25 (54)	33.5 (85)	25.25 (64)	120 (55)
3557-4	18 (46)	21 (53)	15.25 (39)	15.25 (39)	5.25 (13)	11 (28)	5.6	21.25 (54)	33.5 (85)	25.25 (64)	120 (55)
3559, 3559-1	19.75 (50)	32 (81)	17.5 (44)	17.25 (44)	8 (20)	12 (30)	10	23 (58)	58 (147)	23.5 (60)	150 (68)
3560, 3560-1	27 (69)	55 (140)	20 (51)		N/A		20	32 (81)	71.5 (182)	28.5 (72)	280 (127)
3566, 3566-1	27 (69)	55 (140)	20 (51)		N/A		20	32 (81)	71.5 (182)	28.5 (72)	280 (127)

ORDERING INFORMATION

Model #	Electrical			Temperature Range °C (°F)		Freezer Section	US List Price
	Volts	Hertz	Watts	Refrigerator Section	Refrigerator Section		
3551	120	60	360	0° to 10° (32° to 50°)	-20° to 0° (-4° to 32°)		\$3274.00
3551-1	240	50	360	0° to 10° (32° to 50°)	-20° to 0° (-4° to 32°)		\$3932.00
3552	120	60	516	N/A	-23° to -5° (-9° to 23°)		\$3769.00
3552-1	240	50	516	N/A	-23° to -5° (-9° to 23°)		\$4447.00
3555	120	60	516	N/A	-28° to -18° (-18° to 0°)		\$4162.00
3555-1	240	50	516	N/A	-28° to -18° (-18° to 0°)		\$4753.00
3555-6	120	60	276	N/A	-21° to -4° (-5° to 24°)		\$3233.00
3555-7	240	50	276	N/A	-21° to -4° (-5° to 24°)		\$3907.00
3557	120	60	160	2° to 7° (36° to 45°)	N/A		\$2295.00
3557-2	120	60	200	N/A	-26° to -15° (-15° to 5°)		\$2295.00
3557-4	120	60	160	2° to 7° (36° to 45°)	-26° to -6° (-15° to 21°)		\$2295.00
3559	120	60	360	1° to 10° (34° to 50°)	-17° to -6° (-1° to 21°)		\$3101.00
3559-1	240	50	360	1° to 10° (34° to 50°)	-17° to -6° (-1° to 21°)		\$3749.00
3560	120	60	480	-5° to 15° (23° to 59°)	N/A		\$5212.00
3560-1	240	50	480	-5° to 15° (23° to 59°)	N/A		\$5885.00
3566	120	60	480	-5° to 15° (23° to 59°)	N/A		\$2963.00
3566-1	240	50	480	-5° to 15° (23° to 59°)	N/A		\$3198.00

Model 3555



Refrigerators / Freezers

FLAMMABLE MATERIAL STORAGE

Lab-Line

Model 3552-10



SAFE COLD STORAGE!

- — Designed for safe cold storage of volatile materials
- — Engineered with no internal electrical components which could trigger an explosion
- — Meets OSHA and National Fire Protection Association safety requirements



PRODUCT DESCRIPTION

- Provide a wide variety of units to keep your lab safe from potentially destructive explosions.
- Available in many sizes and styles, including uprights, combination refrigerator/freezers, and undercounter models.
- Thermostat located on exterior of cabinet to eliminate interaction with stored materials.
- Magnetic vinyl gaskets for positive door seal.
- CFC-free insulation and coolant.

Model 3551-10

- 11 cubic ft. refrigerator/freezer with separate doors.
- Three adjustable shelves and four door shelves.
- Hydraulic thermostat; adjustable natural air flow vent.
- ABS plastic interior with enameled steel exterior.
- Manual defrost.
- UL listed.

Model 3552-10

- 20 cubic ft. upright freezer.
- Built-in evaporators in shelves for even cooling.
- Four permanently positioned wire shelves.
- Enameled steel exterior with ABS plastic interior.
- Door equipped with key lock.
- Manual defrost.
- UL listed.

Model 3556-4

- 5.6 cubic ft. undercounter freezer.
- Hydraulic thermostat.
- Three fixed wire shelves, three door shelves.
- Aluminum evaporator coil.
- Off cycle defrost.

Model 3556-6

- 5.6 cubic ft. undercounter refrigerator/freezer.
- Hydraulic thermostat.
- Two adjustable wire shelves, three door shelves.
- Manual defrost.

Model 3556

- 5.6 cubic ft. undercounter refrigerator.
- Hydraulic thermostat.
- Three adjustable wire shelves plus three door shelves.
- Aluminum evaporator coil.
- Off cycle self defrost.
- Key locking door.
- UL listed.

Model 3559-10

- 10 cubic ft. single door refrigerator with interior freezer.
- Three adjustable shelves plus four door shelves.
- Hydraulic thermostat.
- ABS plastic interior with enameled steel exterior.
- Manual defrost.
- UL listed.

Model 3556-4



Model 3566-10

- 20 cubic ft. one-door refrigerator.
- Five adjustable shelves.
- Enameled steel exterior.
- Manual defrost.
- UL listed.

Model 3560-10

- 20 cubic ft. corrosion-resistant model designed to store caustic materials and strong acids.
- Stainless steel evaporators.
- Five adjustable stainless steel shelves.
- Enameled steel interior and exterior with stainless steel door liner for greater corrosion resistance.
- Door comes with built-in key lock.
- Manual defrost.
- UL listed.

Refrigerators/Freezers

FLAMMABLE MATERIAL STORAGE

Lab-Line

PRODUCT SPECIFICATIONS

Model #	Chamber Dimensions Inches (cm) Refrigerator Section			Chamber Dimensions Inches (cm) Freezer Section			Total Volume Cubic ft.	Exterior Dimensions Inches (cm)			Shipping Weight Lb. (kg)
	W	H	D	W	H	D		W	H	D	
3551-10, 3551-11	20.75 (53)	41.5 (105)	14 (36)	19 (48)	12.75 (32)	15 (38)	11	23.25 (59)	63 (160)	24.5 (62)	150 (68)
3552-10, 3552-11		N/A		27 (69)	57 (145)	18 (46)	20	32 (81)	71.5 (182)	33 (84)	240 (109)
3556	19 (48)	27.5 (70)	13 (33)		N/A		5.6	23.5 (59)	33.5 (85)	23.75 (60)	120 (55)
3556-4		N/A		18 (46)	27.5 (70)	15 (38)	5.6	21.25 (54)	33.5 (85)	25.25 (64)	120 (55)
3556-6	18 (46)	21 (53)	15.25 (39)	15.25 (39)	5.25 (13)	11 (28)	5.6	21.25 (54)	33.5 (85)	25.25 (64)	120 (55)
3559-10, 3559-11	19.75 (50)	32 (81)	17.5 (45)	17.25 (44)	8 (20)	12 (30)	10	23 (58)	58 (147)	23.5 (60)	150 (68)
3560-10, 3560-11	27 (69)	55 (140)	20 (51)		N/A		20	32 (81)	71.5 (182)	28.5 (72)	280 (127)
3566-10, 3566-11	27 (69)	55 (140)	20 (51)		N/A		20	32 (81)	71.5 (182)	28.5 (72)	280 (127)

ORDERING INFORMATION

Model #	Electrical			Temperature Range °C (°F)		Freezer Section	US List Price
	Volts	Hertz	Watts	Refrigerator Section	Refrigerator Section		
3551-10	120	60	360	0° to 10° (32° to 50°)	-20° to 0° (-4° to 32°)		\$2815.00
3551-11	240	50	360	0° to 10° (32° to 50°)	-20° to 0° (-4° to 32°)		\$3463.00
3552-10	120	60	516	N/A	-23° to -5° (-9° to 23°)		\$3203.00
3552-11	240	50	516	N/A	-23° to -5° (-9° to 23°)		\$3896.00
3556	120	60	160	2° to 7° (36° to 45°)		N/A	\$2397.00
3556-4	120	60	200	N/A	-26° to -15° (-15° to 5°)		\$2397.00
3556-6	120	60	160	2° to 7° (36° to 45°)	-26° to -6° (-15° to 21°)		\$2397.00
3559-10	120	60	360	1° to 10° (34° to 50°)	-17° to -6° (-1° to 21°)		\$2621.00
3559-11	240	50	360	1° to 10° (34° to 50°)	-17° to -6° (-1° to 21°)		\$3269.00
3560-10	120	60	480	-5° to 15° (23° to 59°)		N/A	\$5059.00
3560-11	240	50	480	-5° to 15° (23° to 59°)		N/A	\$5738.00
3566-10	120	60	480	-5° to 15° (23° to 59°)		N/A	\$2866.00
3566-11	240	50	480	-5° to 15° (23° to 59°)		N/A	\$3101.00

Model 3559-10



Refrigerators/Freezers

GENERAL PURPOSE

Lab-Line

Model 3767



SIZES FOR YOUR EVERY NEED!

- — Available as all-refrigerator, combo refrigerator/freezer, upright freezer, and chest freezer models in many popular sizes
- — Adjustable shelves or baskets included
- — CFC-free insulation and coolant



PRODUCT DESCRIPTION

Model 3750

- 5.6 cubic ft. undercounter refrigerator/freezer.
- Adjustable thermostatic control.
- Two adjustable wire shelves, three door shelves.
- ABS plastic interior with enameled steel exterior.
- Manual defrost.
- UL listed.

Model 3751

- 5.6 cubic ft. undercounter refrigerator.
- Off-cycle defrost.
- Adjustable hydraulic thermostat control.
- Three adjustable wire shelves, three door shelves.
- Key lock.
- UL listed.

Model 3752

- 5.6 cubic ft. undercounter freezer.
- Hydraulic thermostat.
- Three fixed shelves.
- Off cycle defrost.

Model 3761

- 7 cubic ft. chest-style freezer.
- Includes storage basket.
- Textured enameled steel exterior.
- Key lock in lid.
- Manual defrost.

Model 3762

- 18 cubic ft. chest-style freezer.
- Largest chest-style model.
- Storage basket included.
- Textured enameled steel exterior.
- Key lock in lid.

Model 3763

- 14 cubic ft. upright two-door combination refrigerator/freezer.
- Hydraulic thermostat control with separate "cold air flow" control for freezer.
- Three adjustable shelves in refrigerator.
- Magnetic vinyl gaskets for positive door seal.
- Automatic defrost.
- UL listed.

Model 3764

- 10 cubic ft. upright combination refrigerator/freezer.
- Single-door model with internal freezer compartment.
- Three adjustable wire shelves.
- Manual defrost.

Model 3766

- 20 cubic ft. frost-free refrigerator.
- Five permanently mounted shelves plus storage basket.
- Hydraulic thermostat.
- Automatic defrost.
- Locking door comes standard.


Model 3767

- 20 cubic ft. upright freezer.
- Four fixed wire shelves and six door shelves.
- Door comes with built-in key lock.
- Manual defrost.
- Enameled steel interior and exterior.

Model 3770

- 1.9 cubic ft. bench top combination refrigerator/freezer.
- Single door with interior freezer compartment.
- One removable shelf included.
- Manual defrost.
- UL listed.

Model 3751



48 Hour or Sooner
Express Shipping
Guaranteed!*



*see inside front cover for details

Refrigerators/Freezers

GENERAL PURPOSE

Lab-Line

PRODUCT SPECIFICATIONS

Model #	Chamber Dimensions Inches (cm) Refrigerator Section			Chamber Dimensions Inches (cm) Freezer Section			Total Volume Cubic ft.	Exterior Dimensions* Inches (cm)			Shipping Weight Lb. (kg)
	W	H	D	W	H	D		W	H	D	
3750, 3750-1	18 (46)	21 (53)	15.25 (39)	15.25 (39)	5.25 (13)	11 (28)	5.6	21.25 (54)	33.5 (85)	25.25 (64)	120 (55)
3751, 3751-1	19 (48)	27.5 (70)	13 (33)		N/A		5.6	23 1/2(59)	33.3 (85)	23.75 (60)	120 (55)
3752, 3752-1		N/A		18 (46)	27.5 (70)	15 (38)	5.60	21.25 (54)	33.5 (85)	25.25 (64)	120 (55)
3761, 3761-1		N/A		26.25 (67)	28 (71)	16.5 (42)	7	31 1/2(80)	35 (89)	26.25 (67)	125 (57)
3762, 3762-1		N/A		50.25 (128)	28 (71)	22.25 (57)	18	60.25 (153)	35 (89)	29.5 (75)	200 (91)
3763, 3763-1	25 (64)	36.5 (93)	15.5 (39)	24 (61)	15.25 (39)	16 (41)	14.4	28 (71)	61.5 (156)	28 (71)	190 (86)
3764, 3764-1	19.75 (50)	32 (81)	17.5 (45)	17.25 (44)	8 (20)	12 (30)	10	23 (58)	58 (147)	23.5 (60)	150 (68)
3766, 3766-1	27 (69)	57 (145)	17.75 (45)		N/A		20	32 (81)	71 (180)	28.5 (72)	280 (129)
3767, 3767-1		N/A		27 (69)	57 (145)	18 (46)	20	32 (81)	71.5 (182)	33 (84)	240 (109)
3770	16.25 (41)	15.25 (39)	11 (28)	6.25 (16)	6 (15)	8.25 (21)	1.9	18 1/2(47)	19 (48)	19 (48)	52 (24)

* Upon installation, we recommend clearance of 3" in rear, 2" on each side and 12" above unit.

Model 3766



Refrigerators / Freezers

GENERAL PURPOSE

Lab-Line

Model 3750



ORDERING INFORMATION

Model #	Electrical			Temperature Range °C(°F)		US List Price
	Volts	Hertz	Watts	Refrigerator Section	Freezer Section	
3750	120	60	160	2° to 7° (36° to 45°)	-26° to -6° (-15° to 21°)	\$1209.00
3750-1	240	50	160	2° to 7° (36° to 45°)	-26° to -6° (-15° to 21°)	\$1255.00
3751	120	60	160	2° to 7° (36° to 45°)	N/A	\$1209.00
3751-1	240	50	160	2° to 7° (36° to 45°)	N/A	\$1255.00
3752	120	60	200	N/A	-26° to -15° (-15° to 5°)	\$1209.00
3752-1	240	50	200	N/A	-26° to -15° (-15° to 5°)	\$1255.00
3761	120	60	276	N/A	-21° to -4° (-5° to 24°)	\$1148.00
3761-1	240	50	276	N/A	-21° to -4° (-5° to 24°)	\$1499.00
3762	120	60	516	N/A	-28° to -18° (-18° to 0°)	\$1969.00
3762-1	240	50	516	N/A	-28° to -18° (-18° to 0°)	\$2198.00
3763	120	60	420	0° to 10° (32° to 50°)	-20° to 0° (-4° to 32°)	\$1464.00
3763-1	240	50	420	0° to 10° (32° to 50°)	-20° to 0° (-4° to 32°)	\$1800.00
3764	120	60	360	1° to 10° (34° to 50°)	-17° to -6° (-1° to 21°)	\$1321.00
3764-1	240	50	360	1° to 10° (34° to 50°)	-17° to -6° (-1° to 21°)	\$1658.00
3766	120	60	480	-5° to 15° (23° to 59°)	N/A	\$2387.00
3766-1	240	50	480	-5° to 15° (23° to 59°)	N/A	\$2632.00
3767	120	60	516	N/A	-23° to -5° (-9° to 23°)	\$1504.00
3767-1	240	50	516	N/A	-23° to -5° (-9° to 23°)	\$1805.00
3770	120	60	132	0° to 7° (32° to 45°)	-5° to 0° (23° to 32°)	\$393.00

Model 3763



Refrigerators/Freezers

CHROMATOGRAPHY REFRIGERATORS

Lab-Line

Model 3796



READY TO USE!

- — Three models available—23, 47, and 69 cu. ft.
- — CFC-free insulation and coolant
- ★ — 60 Month Warranty on Compressor



PRODUCT SPECIFICATIONS

Model #	Capacity Cu. Ft.	Access Ports 2" Diameter	Interior Dimensions Inches (cm)			Exterior Dimensions* Inches (cm)			Shipping Weight Lb. (kg)
			W	H	D	W	H	D	
3790, 3790-1	23	1	24 (61)	61.5 (157)	27 (69)	27.25 (70)	78.25 (199)	29.5 (75)	311 (141)
3793, 3793-1	47	2	50.75 (129)	61.25 (156)	24.5 (63)	54 (138)	78.25 (199)	29.5 (75)	495 (225)
3796, 3796-1	69	2	75 (191)	61.5 (157)	27 (69)	78 (199)	79.25 (202)	29.5 (75)	678 (308)

*Exterior height does not include leveling feet. Add 6" for feet.

PRODUCT DESCRIPTION

- Temperature uniformity of $\pm 0.5^\circ$ at 4°C .
- Adjustable thermostat control from 1° to 10°C .
- Triple tempered glass doors.
- 2" foamed-in-place polyurethane insulation.
- Two 2" diameter PVC access ports for gassing applications.
- White vinyl exterior coating is durable and chemically resistant.
- Interior construction is .040 aluminum with baked-on enamel finish.
- Vapor-proof outlets.
- Off-cycle defrost.
- Optional chromatography support rods and casters available.

ORDERING INFORMATION

Model #	Electrical		Duplex Outlets	Chromatography			US List Compressor Price
	Volts	Hertz		Rods	Shelves	Doors	
3790	120	60	1	N/A	3	1 (swinging)	1/3 hp \$3583.00
3790-1	240	50	1	N/A	3	1 (swinging)	1/3 hp \$3843.00
3793	120	60	2	1	6	2 (sliding)	1/2 hp \$5585.00
3793-1	240	50	2	1	6	2 (sliding)	1/2 hp \$5760.00
3796	120	60	2	2	9	3 (sliding)	1/2 hp \$7020.00
3796-1	240	50	2	2	9	3 (sliding)	1/2 hp \$7342.00

ACCESSORIES

Model #	Description
C	Casters
CR	Chromatography Rods
DT	Digital Temperature Readout
P	Access Ports
V	Vapor-proof Duplex Outlet

When ordering, add accessory letter to end of model number, i.e.: 3790C

Refrigerators / Freezers

LARGE CAPACITY FREEZERS

Lab-Line

Model 3785



REACH IN FOR EASY ACCESS!

- — Storage of biological materials and other nonflammable substances
- — Three models available—from 23, 49 and 72 cu. ft.
- — CFC-free insulation and coolant
- ★ — 60 Month Warranty on Compressor



PRODUCT DESCRIPTION

- Automatic defrost.
- Key locks on all doors.
- Adjustable thermostat control from -17°C to -23°C.
- Temperature uniformity of $\pm 0.5^{\circ}\text{C}$ at -23°C.
- .040 polished aluminum exterior.
- 22 gauge stainless steel floor.
- Covered interior corners for ease of cleaning.
- Foamed-in-place insulation.
- PVC coated adjustable shelves.
- Compressor bottom-mounted for easy access.

Model 3788



PRODUCT SPECIFICATIONS

Model #	W	Interior Dimensions Inches (cm)			Exterior Dimensions Inches (cm)			Shipping Weight Lb. (kg)
		H	D	W	H	D		
3782 3782-1	23.75 (60)	61.5 (156)	26 (66)	27 (69)	78.25 (199)	29.5 (75)	286 (130)	
3785 3785-1	50.5 (128)	61.5 (156)	26 (66)	54 (137)	78.25 (199)	29.5 (75)	425 (193)	
3788 3788-1	74.75 (190)	61.5 (156)	26 (66)	78 (198)	79.25 (201)	29.5 (75)	598 (271)	

ORDERING INFORMATION

Model #	Electrical		Capacity (cu. ft.)	Number of Shelves (sq. ft.)	Solid Doors	Compressor	US List Price
	Volts	Hertz					
3782	120	60	23	3 (13)	1	1/2 hp	\$4621.00
3782-1	240	50	23	3 (13)	1	1/2 hp	\$4621.00
3785	120	60	49	6 (29)	2	3/4 hp	\$6217.00
3785-1	240	50	49	6 (29)	2	3/4 hp	\$6217.00
3788	240	50	72	9 (43)	3	3/4 hp	\$8619.00
3788-1	240	50	72	9 (43)	3	3/4 hp	\$8619.00

Refrigerators / Freezers

LARGE CAPACITY REFRIGERATORS

Lab-Line

MODEL 3784



OUR LARGEST COLD STORAGE UNITS!

- — Storage of biological materials and other nonflammable substances
- — Six models available—from 23, 49, and 72 cu. ft.
- — CFC-free insulation and coolant
- ★ — 60 Month Warranty on Compressor

PRODUCT DESCRIPTION

- Available in glass or solid door styles.
- Key locks on all doors.
- Adjustable thermostatic control from 1° to 10°C.
- Temperature uniformity of ±0.5°C at 4°C.
- .040 polished aluminum exterior.
- 22 gauge stainless steel floor.
- Foamed-in-place insulation.
- PVC coated adjustable shelves.
- Compressor bottom-mounted for easy access.
- UL listed.

ACCESSORIES

Model #	Description
B	Digital Alarm with Remote Contacts
L	6" Adjustable Stainless Steel Legs
DT	Digital Temperature Readout
S	Stainless Steel Interior and Exterior

When ordering, add accessory letter to end of model number, i.e.:3782DT

PRODUCT SPECIFICATIONS



Model #	Interior Dimensions Inches (cm)			Exterior Dimensions Inches (cm)			Shipping Weight Lb. (kg)
	W	H	D	W	H	D	
3780 3780-1	23.75 (60)	61.5 (157)	26 (66)	27 (69)	78.25 (198)	31.75 (81)	245 (111)
3781 3781-1	23.75 (60)	61.5 (157)	26 (66)	27 (69)	78.25 (198)	31.75 (81)	277 (126)
3783 3783-1	50.5 (129)	61.5 (157)	26 (66)	54 (138)	78.25 (198)	31.75 (81)	415 (188)
3784 3784-1	50.5 (129)	61.5 (157)	26 (66)	54 (138)	78.25 (198)	31.75 (81)	480 (218)
3786 3785-1	74.75 (190)	61.5 (157)	26 (66)	78 (198)	78.25 (198)	31.75 (81)	585 (265)
3787 3787-1	74.75 (190)	61.5 (157)	26 (66)	78 (198)	78.25 (198)	31.75 (81)	655 (297)

*Exterior height does not include castors. Add 6" for castor height.

ORDERING INFORMATION

Model #	Electrical		Capacity (cu.ft.)	Number of Shelves (sq.ft.)	Solid Doors	Compressor	US List Price
	Volts	Hertz					
3780	120	60	23	3 (13)	Solid	1/3 hp	\$3667.00
3780-1	240	50	23	3 (13)	Solid	1/3 hp	\$3667.00
3781	120	60	23	3 (13)	Glass	1/3 hp	\$4167.00
3781-1	240	50	23	3 (13)	Glass	1/3 hp	\$4167.00
3783	120	60	49	6 (29)	Solid	1/2 hp	\$4891.00
3783-1	240	50	49	6 (29)	Solid	1/2 hp	\$4891.00
3784	120	60	49	6 (29)	Glass	1/2 hp	\$5748.00
3784-1	240	50	49	6 (29)	Glass	1/2 hp	\$5748.00
3786	120	60	72	9 (43)	Solid	1/2 hp	\$6579.00
3786-1	240	50	72	9 (43)	Solid	1/2 hp	\$6579.00
3787	120	60	72	9 (43)	Glass	1/2 hp	\$7864.00
3787-1	240	50	72	9 (43)	Glass	1/2 hp	\$7864.00