

ADVANTEC[®] Method & Application Quick Reference Guide

STANDARD METHODS

Number	Application	ADVANTEC Products	
9222	Fecal Coliform	A045H047A(W)	Membrane Filter
	Total Coliform	501030(KP47U)	Filter Holder
2540D	TSS	GA-55, GC-50	Glass Fiber Filter
5520	Oil and Grease	No.5B	Quantitative Filter
9213	Pseudomonas sp.	A045H047A(W)	Membrane Filter
9215	Total Bacteria	A045H047A(W)	Membrane Filter
9230	Fecal Streptococcus	A045H047A(W)	Membrane Filter
		A080R047A(W)	Membrane Filter
9610	Yeast & Mold	A080H047A(W)	Membrane Filter
		A080R047A(W)	Membrane Filter

AIR SAMPLING

Method Number	Application	ADVANTEC Products	
EPA Method 5	Air Sampling & Analysis	No.88R 25 x 90mm No.88R 19 x 90mm	Quartz Thimble Filter
EPA Method 17	Stack Gas Emissions Analysis	PD-2	Petri Dishes
EPA Method 201a, b	PM10	QR-100	Quartz Fiber Filter
NIOSH 7400	Asbestos Sampling(PCM)	A045A025A A080A025A, A080X025A	Membrane Filter
NIOSH 7402	Asbestos Sampling(TEM)	M-085	Absorbent Pad

OIL & SDI TESTING

Method Number	Application	ADVANTEC Products	
SAE ARP 958	Aviation Fuel Testing	A080A047A	Membrane Filter
ASTM D2276	Aviation Fuel Testing	(A080H047A)	
ASTM D4189-07	SDI Test	A045A047A	Membrane Filter