

WESCO

SERIES

enjoy the VU

Bio VU™ 2000 Series

Bio VU™ 2000 series biological microscopes incorporate economical and durable design with quality optics and mechanics to meet the standards of medical students as well as those features needed by clinical technicians.

There are currently nine models in the Bio VU™ 2000 Series with many additional accessories available to further customize your viewing needs.



Bio VU 2000 Series

SPECIFICATIONS	MODELS								
	2200	2300	2320	2330	2335	2500	2520	2530	2535
Monocular observation body	•								
Binocular observation body		•	•	•	•				
Trinocular observation body						•	•	•	•
10X widefield eyepieces 18mm field of view	•	•	•	•	•	•	•	•	•
Quintuple revolving nosepiece	•	•	•	•	•	•	•	•	•
DIN 4X, 10X, 40XR & 100XR oil semi-plan achromatic objective	•	•				•			
DIN 4X, 10X, 40XR & 10XR oil plan achromatic objective			•				•		
DIN plan phase contrast PH10X, PH20X, PH40XR & PH100XR oil objectives with turret phase condenser, telescope eyepiece and G.I. filter				•				•	
DIN 4X, PH10X, PH40XR & 100XR oil objectives with slider phase condenser, telescope eyepiece and G.I. filter					•				•
120mm X 128mm graduated mechanical stage with low position 76mm X 50mm travel X-Y controls	•	•	•	•	•	•	•	•	•
Coaxial coarse and fine focusing with tension control	•	•	•	•	•	•	•	•	•
Focusable N.A. 1.25 Abbe condenser with filter holder	•	•	•	•	•	•	•	•	•
6 volt 20 watt halogen illumination with variable light control	•	•	•	•	•	•	•	•	•
Options									
Dual viewing attachment	•	•	•	•	•	•	•	•	•
Attachable mirror	•	•	•	•	•	•	•	•	•

WESCO



Distributed By: