


Abbe Refractometer
ATAGO'S STANDARD MODEL
Cat.No.1210

NAR-1T

The 1T is an Abbe refractometer of standard type, measures refractive indices of liquids and solids quickly and accurately.

It is designed to give refractive indices of D line without using a sodium lamp.



Specifications

Measurement range

:Refractive index(nD) 1.3000 to 1.7000, Brix(%) 0.0 to 95.0%

Minimum scale

:Refractive index(nD)0.001, Brix(%)0.5%

Measurement accuracy

:Refractive index(nD) ± 0.0002 , Brix(%) $\pm 0.1\%$

Average dispersion value

:nF-nC(to be calculated according to conversion table)

Digital thermometer

:Indication range 0.0 to 50.0°C

(Accuracy $\pm 0.2^\circ\text{C}$, Minimum indication 0.1°C)

Measurement temperature range :0 to 50 °C

Light source :8V, 0.15A tungsten lamp

Power supply :AC 100 to 240 V, 50/60 Hz

Power consumption : 6 VA

Dimensions and weight : Main unit 13(W) × 28(D) × 23(H) cm, 3.4kg

Accessories

| | | |
|--|-------|-------|
| Test piece (Refractive index standard) | | 1pc |
| Monobromonaphthalene (4 ml) | | 1pc |
| Thermometer (Digital Type) | | 1pc |
| Screwdriver for scale adjustment | | 1pc |
| Lamp cord | | 1pc |
| Illuminator lamp (8V, 0.15A) | | 5pcs |
| Tube band | | 10pcs |
| Instruction manual | | 1pc |

*Specifications and appearance are subject to change without notice.

ATAGO CO., LTD.

Main Office:32-10, Honcho, Itabashi-ku, Tokyo 173-0001, Japan

Phone:81-3-3964-6156, Fax:81-3-3964-6137

<http://www.atago.net/> export@atago.net

