

Barnstead Deionization



Cartridge Kits



Cartridge Kits for NANOpure® Diamond™, NANOpure Infinity™, NANOpure, NANOpure II, E-pure™ Ultrapure Water Systems

Barnstead ultrapure water systems require expendable cartridges and filters not supplied with the unit. Cartridge kits include everything needed to allow the Barnstead ultrapure water systems to supply 18.2 megohm-cm resistivity Type I reagent grade water.

ORDERING INFORMATION

| Model # | Description | Content | Price |
|--|---|--|----------|
| Cartridge Kits for QuadStyle NANOpure Systems (Final Filters sold separately) - Includes canister o-rings | | | |
| D4801 | Type I Cartridge Kit; Tap Feed | 1-D0835 Pretreatment Cartridge, 1-D0803 High Capacity Cartridge, 2-D5027 Ultrapure Cartridges | \$347.00 |
| D4802 | Type I ORGANICfree Cartridge Kit; Tap Feed | 1-D0836 MACROpure Cartridge, 1-D0803 High Capacity Cartridge, 1-D5027 Ultrapure Cartridge, 1-D5021 ORGANICfree Cartridge | \$415.00 |
| D5025 | Type I Cartridge Kit; Pretreated Feed | 1-D0835 Pretreatment Cartridge, 2-D5027 Ultrapure Cartridges, 1-D0809 Ultrapure Cartridge | \$363.00 |
| D5026 | Type I ORGANICfree Cartridge Kit; Pretreated Feed | 1-D0836 MACROpure Cartridge, 1-D0809 Ultrapure Cartridge, 1-D5027 Ultrapure Cartridges, 1-D5021 ORGANICfree Cartridge, | \$426.00 |
| Cartridge Kits for E-pure Systems (Final Filters sold separately) - Includes canister o-rings | | | |
| D5029 | 3 Module Type I Cartridge Kit | 1-D0835 Pretreatment Cartridge, 2-D5027 Ultrapure Cartridges | \$273.00 |
| D5028 | 4 Module Type I Cartridge Kit | 1-D0835 Pretreatment Cartridge, 1-D0803 High Capacity Cartridge, 2-D5027 Ultrapure Cartridges | \$354.00 |
| D50227 | 4 Module Type I Cartridge Kit; Pretreated Feed | 1-D0835 Pretreatment Cartridge, 1-D0809 Ultrapure Cartridge, 2-D5027 Ultrapure Cartridges | \$360.00 |
| D5022 | 3 Module ORGANICfree Type I Cartridge Kit | 1-D0836 MACROpure Cartridge, 1-D5027 Ultrapure Cartridge, 1-D5021 ORGANICfree Cartridge, | \$415.00 |
| D5023 | 4 Module ORGANICfree Cartridge Kit | 1-D0836 MACROpure Cartridge, 1-D0803 High Capacity Cartridge, 1-D5027 Ultrapure Cartridge, 1-D5021 ORGANICfree Cartridge | \$418.00 |
| D50228 | 4 Module ORGANICfree Cartridge Kit; Pretreated Feed | 1-D0836 MACROpure, 1-D5027 Ultrapure Cartridges, 1-D5021 ORGANICfree Cartridge 1-D0809 Ultrapure Cartridge | \$422.00 |
| Cartridge Kits for NANOpure and NANOpure II Systems (Final Filters sold separately) | | | |
| D3805 | 3 Module Type I Cartridge Kit | 1-D0835 Pretreatment Cartridge, 2-D0809 Ultrapure Cartridges | \$273.00 |
| D3804 | 4 Module Type I Cartridge Kit | 1-D0835 Pretreatment Cartridge, 1-D0803 High Capacity Cartridge, 2-D0809 Ultrapure Cartridges | \$354.00 |
| D3806 | 3 Module ORGANICfree Type I Cartridge Kit | 1-D0836 MACROpure Cartridge, 1-D0809 Ultrapure Cartridge, 1-D0820 ORGANICfree | \$327.00 |
| D3807 | 4 Module ORGANICfree Type I Cartridge Kit | 1-D0836 MACROpure Cartridge, 1-D0803 High Capacity Cartridge, 1-D0809 Ultrapure Cartridge, 1-D0820 ORGANICfree | \$407.00 |
| Cartridge Kits for NANOpure Infinity | | | |
| D50256 | Type I Cartridge Kit TAP or DI Feed | 1 -D50251 Pretreatment Cartridge, 3 -D50253 Ultrapure Cartridges, 2 -FL703X2 Final Filters | \$516.00 |
| D50257 | Type I Cartridge Kit RO or Dist. Feed | 1 -D50250 Pretreatment Cartridge, 3 -D50253 Ultrapure Cartridges, 2 -FL703X2 Final Filters | \$479.00 |
| D50254 | ORGANICfree Cartridge Kit TAP or DI Feed | 1 -D50251 Pretreatment Cartridge, 2 -D50253 Ultrapure Cartridges, 1 -D50252 ORGANICfree Cartridge, 2-FL703X2 Final Filters | \$501.00 |
| D50255 | ORGANICfree Cartridge Kit RO or Dist. Feed | 1-D50250 Pretreatment Cartridge, 2-D50253 Ultrapure Cartridges, 1-D50252 ORGANICfree Cartridge, 2-FL703X2 Final Filters | \$492.00 |
| Cartridge Kits for NANOpure Diamond | | | |
| D50280 | Diamond Kit Organic Free RO and Dist. Feed | Ultra-Low Organics Type I Water Reverse Osmosis or Distilled Water Feed | \$473.00 |
| D50281 | Diamond Kit Organic Free Deionized Feed | Ultra-Low Organics Type I Water Deionized Water Feed | \$493.00 |
| D50282 | Diamond Kit Type I Reverse Osmosis or Dist. Feed | Low Organics Type I Water Reverse Osmosis or Distilled Water Feed | \$464.00 |
| D50283 | Diamond Kit Type I Deionized Feed | Low Organics Type I Water Deionized Water Feed | \$480.00 |