



PRODUCTS



APPLICATIONS

Marine Biology

Substance Abuse Testing

Water Testing

Dri-Baths™

SELECTION GUIDE

MODULAR BLOCK

	Type	Temperature Range	Capacity	Removable Blocks	Temperature Control	Page
<i>High Temperature Models!</i>	DB104200	AMB+5 to 160°C	One Block	Yes	±0.1°C	29
	DB104300	AMB+5 to 160°C	Two Block	Yes	±0.1°C	29
	DB104400	AMB+5 to 160°C	Three Block	Yes	±0.1°C	29
<i>Analog Standard Temperature!</i>	DB17600	AMB+5 to 115°C	One Block	Yes	±0.5°C	30
	DB28100	AMB+5 to 140°C	Three Block	Yes	±0.5°C	30
	DB16500	2°C-110°C	Three Block	Yes	±0.1°C	30
<i>Standard and High Temperature Models!</i>	5401	AMB+5 to 200°C	One Block	Yes	±0.1°C	32
	5111	AMB+5 to 110°C	Four Block	Yes	±0.1°C	32
<i>Analog and Digital Standard Temperature!</i>	2000	AMB+5 to 130°C	One to Six Block	Yes	±0.15°C	33
	2050	AMB+5 to 130°C	One to Six Block	Yes	±0.1°C	33
<i>Digital High Temperature!</i>	2005	AMB+5 to 200°C	One to Six Block	Yes	±0.2°C	34
	2006					34
	2007					34

FIXED BLOCK

<i>Simple Operation!</i>	DB5900	AMB+5 to 60°C	10 or 20 Well Block	No	±0.5°C	36
	DB12200	AMB+5 to 60°C	40 Well Block	No	±0.5°C	36

HEATING AND COOLING

<i>Heating and Cooling in One Unit!</i>	DB82200	-10°C to 70°C	Two Block	Yes	±0.2°C	37
---	---------	---------------	-----------	-----	--------	----