

Thermolyne
Dri-Baths™



Modular Block Accessories



A Block for Any Application!

- Versatile
- Anodized coating resists the harshest environments



Modular Blocks

*Optional
Thermometer:
Optional 6" (15.2 cm)
long thermometer for
exact temperature
measurement of alu-
minum block. 30-
100°C range in 1°
divisions. Accuracy of
±1.0°C; ±0.5°C at
50°C. Part # ME384X1.*

Product Description

- Solid anodized aluminum blocks resist the harshest environments and provide maximum heat retention.
- All blocks include thermometer well for easily measuring block temperature.
- Fit all Thermolyne brand modular block dri-baths.

ORDERING INFORMATION

Block Dimensions: Inches (cm)—3 (7.6)W x 2 (5)H x 4 (10.1)D
Block Well Depth: 48 mm

Model #	Tube Size (mm)	# of Wells	Max Tube Size (mm)	Well Volume (ml)	US List Price
BK165X3A	6	30	6.5	.62	\$89.00
BK165X4A	10	20	10.75	2.4	\$89.00
BK165X12A	12	20	12.75	4.1	\$89.00
BK165X5A	13	20	13.75	4.5	\$89.00
BK165X6A	15	12	15.75	6.3	\$89.00
BK165X7A	16	12	16.75	7.4	\$89.00
BK165X13A	18	12	18.75	9.8	\$89.00
BK165X8A	20	12	20.75	12.2	\$89.00
BK165X9A	25	6	25.75	19.9	\$89.00
BK165X17A	1.5 ml micro tube block, 20 tube capacity				\$93.00
BK165X20A	0.5 ml micro tube block, 30 tube capacity				\$124.00
BK165X2A	Solid aluminum block to be machined by user				\$89.00
BK165X21A	1 Cuvette block, 12 cuvette capacity				\$101.00
ME384X1	Optional 6" (15.2 cm) Thermometer				\$91.00