

## Lab-Line Incubators



# Air-Jacketed Automatic CO<sub>2</sub>



Model 495

## Exclusive Decontam Cycle!

- Choose from IR or TC sensors
- Exclusive SMARTGAS® system for accurate CO<sub>2</sub> control
- Exclusive DECONTAM cycle automatically decontaminates - TC sensor models only
- Security password protection of set points for unintended changes

☆ 24 Month Warranty

- Fast 7 minute CO<sub>2</sub> recovery to 5% after 30 second door opening protects cultures.
- Dual chamber models come with separate, independent control systems.
- Corrosion-resistant polished stainless steel chamber features rounded corners for easy cleaning.
- Temperature adjustable heated door keeps inner glass door free of condensation for improved viewing of cultures.
- Keypad programming for fast, easy adjustment of parameters.
- Five shelves and supports can be removed and reinstalled without tools for complete decontamination and easier cleaning. Twelve shelf positions per chamber.
- Accessory HEPA filter (99.7% efficient at 0.5 microns) can be easily installed to ensure aseptic cell culture.



## Product Description

- Built-in RS232 port for transferring temperature and CO<sub>2</sub> data to a printer.
- Easy-to-read LED display scans between actual temperature and % CO<sub>2</sub>.

## Operation

- CO<sub>2</sub> concentration from 0-20% with ±0.1% control.
- Audible and visual alarms indicate temperature deviations of ±1°C or more, CO<sub>2</sub> concentration deviations of 1% or greater, or when incubator door has been left open for more than five minutes.
- Independent hydraulic high limit thermostat protects cultures in case of microprocessor system failure.

## PRODUCT SPECIFICATIONS

Model #	Sensor Type	Operating Temp Range °C	Temperature Control °C	Temperature Uniformity °C	Chamber Volume cu. ft. (l)	Chamber Dimensions Inches (cm)			Overall Dimensions Inches (cm)			Shipping Weight Lb. (kg)
						W	H	D	W	H	D	
490, 490-1CE 495, 495-1	TC IR	Ambient +5° to 60°C Ambient +5° to 55°C	±0.1°C	±0.25°C	single chamber 6.2 (176)	22 (56)	22 (56)	22 (56)	27.3 (69)	34 (86)	26.8 (68)	245 (111)
492, 492-1CE 497, 497-1	TC IR	Ambient +5° to 60°C Ambient +5° to 55°C	±0.1°C	±0.25°C	dual chamber (2) 12.4 (352)	44 (112)	side-by-side 44 (112)	22 (56)	54.6 (139)	34 (86)	26.8 (68)	490 (222)
494, 494-1CE 499, 499-1	TC IR	Ambient +5° to 60°C Ambient +5° to 55°C	±0.1°C	±0.25°C	dual chamber (2) 12.4 (352)	22 (56)	stacked 44 (112)	22 (56)	27.3 (69)	68 (173)	26.8 (68)	490 (222)

## ORDERING INFORMATION

Model #	Volts	Electrical (50/60 Hz)		US List Price	Model #	Volts	Electrical (50/60 Hz)		US List Price
		Amps	Watts				Amps	Watts	
<b>TC Sensor</b>					<b>IR Sensor</b>				
490	120	7.1	850	\$5,028.00	495 <sup>5</sup>	120	7.1	850	\$6,222.00
490-1CE	230/240	3.5	850	\$5,028.00	495-1 <sup>5</sup>	230/240	3.5	850	\$8,110.00
492	120	14.2	1700	\$9,975.00	497 <sup>5</sup>	120	14.2	1700	\$12,362.00
492-1CE	230/240	7.1	1700	\$9,975.00	497-1 <sup>5</sup>	230/240	7.1	1700	\$12,362.00
494 <sup>5</sup>	120	14.2	1700	\$9,975.00	499 <sup>5</sup>	120	14.2	1700	\$12,362.00
494-1CE	230/240	7.1	1700	\$9,975.00	499-1 <sup>5</sup>	230/240	7.1	1700	\$12,362.00

## ACCESSORIES

Model #	Description	Dimensions, Inches (cm)			US List Price	Model #	Description	Dimensions, Inches (cm)			US List Price
		W	H	D				W	H	D	
<b>Factory Installed Options for 495, 497, 499 Series (IR models only)</b>					<b>Factory installed Options for 495, 497, 499 Series (IR models only)</b>						
370-24	Printer				\$711.00	490-4	Stacking Stand for 460/490	34.2 (86.8)	41.7 (105.9)	26.9 (68.3)	\$572.00
370-25	Cable for printer				\$46.00	490-9	Hepa Filter				\$185.00
490-5	Mobile cart for 460/490	27 (68.6)	3.25 (8.3)	23.9 (60.7)	\$450.00	RH	RH Display				\$594.83

If a 495 with RH display is required, order 495RH