

## Lab-Line Incubators



# Thermal Rocker® Incubator



## Ideal for Hybridization of Transfer Membranes!

- Accommodates different size containers and heat sealed plastic bags
- ☆ 24 Month Warranty on Components



Model 4637

### Product Description

- Electronic proportional temperature controller maintains temperatures from ambient  $\pm 5^\circ$  to  $70^\circ\text{C}$  with  $\pm 0.5^\circ$  accuracy.
- Automatic resetting thermal cutoff shuts down power to heaters if temperature set point is exceeded.
- Easy-to-read LED temperature display.
- Rocking speed is adjustable from 0 to 100 cycles per minute.
- 14.19" x 14.19" (36 cm x 36 cm) non-skid rubber mesh pad holds utility trays in place.
- Convenient front-mounted switch for operation with or without heat.
- Angle of rotation adjusts from  $10^\circ$  to  $15^\circ$ .
- Load capacity: 10 pounds.

### Applications

- Blood procedures
- Cell cultures
- Hybridizations
- Southern blotting techniques
- Northern blotting techniques

### PRODUCT SPECIFICATIONS

Model #	Platform Dimensions Inches (cm)		Overall Dimensions* Inches (cm)			Shipping Weight Lb. (kg)
	W	D	W	H	D	
4637, 4637-1CE	14.19 (36)	14.19 (36)	16.5 (42)	8.5 (22)	16.0 (41) (without cover)	25 (11)

\* Clear acrylic cover adds 5" (12.7 cm) to height.

### ORDERING INFORMATION

Model #	Volts	Electrical (50/60 Hz)		Watts	US List Price
		Amps			
4637 <sup>1</sup>	120	3.3		400	\$1,899.00
4637-1CE <sup>2</sup>	240	1.7		400	\$1,899.00

<sup>1</sup> ETL, <sup>2</sup> CE