

Electrothermal Pyrometers



Digital



Model 1006

Quicker and Easier Than a Thermometer!

- Rugged, compact and easy to use
- Accompanied by a range of temperature probes
- Use with Mel-Temp for quick melting point determinations

Product Description

- Provides 0.1°C resolution on an easy to read digital display.
- HOLD key allows uninterrupted observation of the sample rather than dividing attention between sample and thermometer.
- Choose from different accessory probes.
- Includes both J and K thermocouple inputs.

Applications

- Product temperature measurement testing.
- Temperature measurement for research & development.
- Frozen product quality control.
- Temperature measurement in production processes.
- Melting point determination with the Mel-Temp.

ORDERING INFORMATION

Model #	Read-out Scale	Range °C	US List Price
1006	°C	- 200 to 1372°C	\$345.00