


# Harvey Sterilizers



## MC8/10 Steam



MC10



*Non-recirculating option keeps internal components and filters free of debris and contamination*

### Easy-To-Program Tabletop Sterilizers!

- Self cleaning
- Compact design
- Choose from two chamber sizes
- One button, preset programs cycle the unit from start to finish
- ☆ 24 Month Warranty
- MC10 chamber holds three 1 liter flasks.
- 4 preset cycles; wrapped, unwrapped, packs, and liquids and an automatic drying cycle.
- Closed-door drying protects against airborne contaminants.
- Special venting for liquids to reduce boil over; allows cycle to run under 45 minutes for EPA Standards.
- ASME coded chamber.
- Has received 510 (K) clearance.
- Optional printer for hard copy records.

### Product Description

- Easy to set up; simple to operate one-touch controls.
- Easily accessible reservoir and removable racks and trays make cleaning easy.
- Built and tested to rigorous industry standards.

### PRODUCT SPECIFICATIONS

	MC8	MC10		MC8	MC10
Overall Dimensions (W x H x D)	13.5" x 16" x 22" (34.3 cm x 40.6 cm x 55.9 cm)	16" x 18" x 23" (40.6 cm x 45.7 cm x 58.4 cm)	Fixed Cycle Parameters	Liquids: 30 min. @ 121°C Unwrapped: 3 min. @ 135°C Wrapped: 5 min. @ 135°C Packs: 30 min. @ 121°C	
Chamber Dimensions	8" diameter x 14.5" usable depth (20.3 cm diameter x 36.8 cm usable depth)	10" diameter x 15.5" usable depth (25.4 cm diameter x 39.3 cm usable depth)	Chamber Capacity	.42 cubic ft., 11.9 l	.7 cubic ft., 19.8 l
Tray Size	Large  Small	15.25" x 9" x 1.25" (38.7 cm x 22.8 cm x 3.2 cm 14" x 6.25" x 0.88" (35.6 cm x 15.8 cm x 2.2 cm)	Operational Weight	55 lb (25 kg)	70.5 lb (32 kg)
			Shipping Weight	73 lb (33 kg)	92 lb (42 kg)
			Heater Wattage	1350 Watts	1425 Watts
			Reservoir Capacity	3.0 l (approximately 3 qts.)	4.0 l (approximately 4 qts.)
			Accessories Included:	Rack, 1 large tray, 1 small tray	

### ORDERING INFORMATION

Model #	Description	Electrical (50/60 Hz)		US List Price	Model #	Description	Electrical (50/60 Hz)		US List Price
		Volts	Amps				Volts	Amps	
265560	Harvey MC10	120	12	\$5,616.00	265595	Harvey MC8	120	12	\$4,779.00
265565*	Harvey MC10	230	6	\$5,616.00	265600*	Harvey MC8	230	6	\$4,779.00
265588	MC10 with non-recirculating feature <sup>2</sup>	120	12	\$5,849.00	265610	MC8 with non-recirculating feature	120	12	\$5,018.00
265589*	MC10 with non-recirculating feature <sup>2</sup>	230	6	\$5,849.00	265615*	MC8 with non-recirculating feature	230	6	\$5,018.00
265591	MC10D with built-in data logger	120	12	\$6,334.00	265585	MC8 with built-in data logger	120	12	\$5,496.00
265840*	MC10D with built-in data logger	230	6	\$6,334.00	265617*	MC8 with built-in data logger	230	6	\$5,496.00
265575	MC10D with non-recirculating feature and built-in data logger	120	12	\$6,567.00	265608	MC8 with non-recirculating feature and built-in data logger	120	12	\$5,736.00
265580*	MC10D with non-recirculating feature and built-in data logger	230	6	\$6,567.00	265612*	MC8 with non-recirculating feature and built-in data logger	230	6	\$5,736.00

<sup>1</sup>Non-recirculating feature patent # 5,480,623.  
<sup>2</sup>European Cord Set. Not intended for use in US.

### ACCESSORIES

Model #	Description	US List Price	Model #	Description	US List Price
265656	Non-Recirculating Feature for MC8 or MC10, retrofittable and easily installed.	\$289.00	265001	Tray, small 6.25" x 14" (15.9 x 35.6), for MC8	\$45.30
260284	Datalogger Paper for MC units (Box of 6 rolls)	\$13.00	232010202	Metal Polish, pint bottle with spray top (Case of 6)	\$96.10
265593	Plug-in Data Logger for MC8	\$795.20	232010102	Chamber Cleaner, pint bottle with dispenser top (Case of 6)	\$77.70
2-513025	Standard Tray Rack for MC10	\$120.60	2-69904	Door Gasket for Harvey MC8 (Replace once per year to assure best seal)	\$36.00
265002	Standard Tray Rack for MC8	\$74.30	2-513005	Door Gasket for Harvey MC10 (Replace once per year to assure best seal)	\$57.90
2-67811	Tray, large 9" x 15.25" (22.9 x 38.7), for MC10	\$101.90	HRX5	Pouch Rack	\$81.10