

## Thermolyne Stirrers



# Aluminum Top and Stir Light



Models S7200 and SL7200

**Includes FREE Stirbar**

## Illuminated Top For Observing Titration Endpoints!

- Stir Light is ideal for spotting titration endpoints
- Aluminum top stirrer accommodates 4 l vessels



### Product Description

#### Stir Light (SL7200)

- Acrylic top plate is illuminated with cool fluorescent bulbs.
- Light and stirring functions can be used independently.
- Rugged cast aluminum body includes an integral ring stand holder.
- Perfect for titration.

#### Stirrer (S7200)

- Aluminum top plate remains cool while stirring.
- Rugged cast aluminum body includes an integral ring stand holder.
- 7" x 7" top plate accommodates 4 l flasks.
- Top side drip edge channels spilled liquid away from components.

### Applications

- Standard reagent mixing
- Illuminate hard-to-see endpoints in titrations
- Stir bacteriological cultures
- Stir up to four liters of aqueous solutions

### PRODUCT SPECIFICATIONS

Type	Top Plate Dimensions Inches (cm)		Overall Dimensions Inches (cm)		
	W	D	W	H	D
S7200	7 (17.7)	7 (17.7)	7.68 (19.5)	4.5 (11.4)	8.5 (21.5)
SL7200	7 (17.7)	7 (17.7)	7.68 (19.5)	4.5 (11.4)	8.5 (21.5)

### ORDERING INFORMATION

Model #	Volts	Electrical (50/60 Hz)		Stirring Speed Range (rpm)	Shipping Weight Lb. (kg)	US List Price
		Amps	Watts			
S7225 <sup>1</sup>	120	0.3	18	100-1000	6 (2.8)	\$305.00
S7220 <sup>1</sup>	230	0.1	11	100-1000	6 (2.8)	\$305.00
SL7225 <sup>1</sup>	120	0.6	34	100-1000	9 (4.1)	\$486.00
SL7220-33 <sup>1,2</sup>	230	0.3	34	100-1000	9 (4.1)	\$486.00

<sup>1</sup> CSA, <sup>2</sup> CE, Non-CE -26 Units with European cord set are available.