


Thermolyne Stirrers



Maxi Stirrer



Type S25500

Powerful Large Volume Stirring!

- 1/8 hp motor provides maximum torque
- Powerful drive magnet accommodates 6" stir bars

Product Description

- 196 square inches (1260 sq. cm) of stainless steel stirring surface accommodates up to 300 lb. (136 kg) loads.
- Powerful 1/8 hp permanent split capacitor motor provides maximum torque, even when mixing viscous solutions (> 3,000 centipoise).
- Powerful 7.5" (19.1 cm) drive magnet accommodates 3" to 6" (7.6 to 15.2 cm) stir bars for stirring up to 35 gallons (132 liters) of aqueous solution.

Applications

- Stirring viscous gels
- Stirring latex mixtures
- Mixing polymers and starches

PRODUCT SPECIFICATIONS

Type	Top Plate Dimensions Inches (cm)		W	Overall Dimensions Inches (cm)	
	W	D		H	D
S25500	14 (35.5)	14 (35.5)	14.75 (37.4)	8 (20.3)	14.75 (37.4)

ORDERING INFORMATION

Model #	Volts	Electrical		Stirring Speed Range (rpm)	Shipping Weight Lb. (kg)	US List Price
		Amps	Watts	Hz		
S25535	120	1.4	160	60	48 (21.8)	\$1,709.00
S25530	240	0.7	160	50/60	48 (21.8)	\$1,709.00
S25530-26	230	0.7	160	50/60	48 (21.8)	\$1,709.00

ACCESSORIES

Model #	Description	Length Inches (cm)	Diameter Inches (cm)	US List Price
MGX7	Teflon-coated stir magnet	4.25 (10.8)	1.06 (2.7)	\$221.20