



Table-Top Buchner Funnel

With Fritware® Porous Filter or HDPE Perforated Plate

Polyethylene funnels are one piece design with a FIXED or REMOVABLE porous filter plate

Custom sizes available upon request.

supported by a multiple ring grid. Filter plates are made of high density polyethylene 6.4mm (1/4") thick with a non-porous ring at the periphery of the plate which seals filter paper when used. Vacuum connector accepts 12.7mm (1/2") I.D. tube. Available in two porosities (medium and coarse). Coarse porosity plate is 90-130 microns, medium porosity plate is 45-90 microns. Use below 52°C (125°F). **Special sizes and features made to order.**

With Perforated Filter Plate Section

Similarly, these buchner funnels are also available with a FIXED or REMOVABLE perforated filter plate. The plate is made of 4.8mm (3/16") thick high density polyethylene and has 4.8mm (3/16") perforations on 11mm (7/16") centers. The FIXED version is for coarse filtration or use with a cloth or paper filter. A REMOVABLE filter plate is available for use in precious metals recovery. Use below 52°C (125°F). 24" and 36" sizes are crated for shipment.



Filter Paper Discs

Sized to Fit the Table-Top Buchner Funnel

Filter Paper Discs to fit table top Buchner funnels are E-D 615 type for fast filtering. They have a rough crepe surface with a thickness of 0.0105mm, basic weight 20 lb. G/Sq. meter, flow rate 235 ml/min., retention of 25 microns. **100 discs per box, 6 boxes per case.**

CATALOG NO.	DESCRIPTION	PRICE/BOX
H14632-0010	For 260mm (10.25") Inside Diameter Funnel	\$28.80
H14632-0018	For 457mm (18") Inside Diameter Funnel	68.25
H14632-0024	For 610mm (24") Inside Diameter Funnel	83.70
H14632-0036	For 914mm (36") Inside Diameter Funnel	204.00



Buchner Funnel

Economical, Lightweight and Durable

Easy to clean two-piece, polypropylene Buchner funnels with perforated HDPE filter plate. Lightweight, but not top-heavy. Steam autoclavable at 121°C (250°F). **Per each, except F14600-0043 which is 6 per bag.**



CATALOG NO.	TAKES FILTER PAPER DIA.	RIM TO PLATE	CAPACITY	CASE	PRICE/EACH
F14600-0043	42.5mm	24mm	40ml	2 Bags	\$31.95/6
H14601-0000	55mm	30mm	75ml	12 Each	7.95
H14602-0000	70mm	40mm	165ml	6 Each	9.75
H14603-0000	90mm	47mm	315ml	6 Each	12.99
H14604-0000	110mm	67mm	875ml	6 Each	38.25



Buchner Funnel

260MM (10.25") INSIDE DIAMETER, 203MM (8") OVERALL HEIGHT, 127MM (5") RIM TO PLATE

CATALOG NO.	DESCRIPTION	PRICE/EACH
H14620-0000	Funnel with Coarse Porosity Fixed Plate	\$216.00
H14620-1260	Funnel with Coarse Porosity Removable Plate	216.00
H14625-3510	Funnel with Medium Porosity Fixed Plate	216.00
H14626-3510	Funnel with Medium Porosity Removable Plate	216.00
H14627-0000	Funnel with Perforated Fixed Plate	216.00
H14627-1260	Funnel with Perforated Removable Plate	216.00

457MM (18") INSIDE DIAMETER, 292MM (11.5") OVERALL HEIGHT, 203MM (8") RIM TO PLATE

CATALOG NO.	DESCRIPTION	PRICE/EACH
H14621-0000	Funnel with Coarse Porosity Fixed Plate	\$540.00
H14621-1457	Funnel with Coarse Porosity Removable Plate	540.00
H14625-3518	Funnel with Medium Porosity Fixed Plate	540.00
H14626-3518	Funnel with Medium Porosity Removable Plate	540.00
H14628-0000	Funnel with Perforated Fixed Plate	465.00
H14627-1457	Funnel with Perforated Removable Plate	465.00

610MM (24") INSIDE DIAMETER, 330MM (13") OVERALL HEIGHT, 267MM (10.25") RIM TO PLATE

CATALOG NO.	DESCRIPTION	PRICE/EACH
H14622-0000	Funnel with Coarse Porosity Fixed Plate	\$1125.00
H14622-1610	Funnel with Coarse Porosity Removable Plate	1125.00
H14625-3524	Funnel with Medium Porosity Fixed Plate	1125.00
H14626-3524	Funnel with Medium Porosity Removable Plate	1125.00
H14629-0000	Funnel with Perforated Fixed Plate	741.60
H14627-1610	Funnel with Perforated Removable Plate	741.60

914MM (36") INSIDE DIAMETER, 375MM (14.75") OVERALL HEIGHT, 305MM (12") RIM TO PLATE

H14623-0000	Funnel with Coarse Porosity Fixed Plate	\$1770.00
H14623-1914	Funnel with Coarse Porosity Removable Plate	1770.00
H14625-3536	Funnel with Medium Porosity Fixed Plate	1770.00
H14626-3536	Funnel with Medium Porosity Removable Plate	1770.00
H14630-0000	Funnel with Perforated Fixed Plate	1335.00
H14627-1914	Funnel with Perforated Removable Plate	1335.00

Replacement Discs

FOR 260MM (10.25") INSIDE DIAMETER FUNNEL

H91462-4026	Replacement Coarse Porous Plate	\$159.00
H91462-4326	Replacement Medium Porous Plate	156.00
H91463-1260	Replacement Perforated Plate	144.00

FOR 457MM (18") INSIDE DIAMETER FUNNEL

H91462-4457	Replacement Coarse Porous Plate	306.00
H91462-4345	Replacement Medium Porous Plate	270.00
H91463-1457	Replacement Perforated Plate	147.00

FOR 610MM (24") INSIDE DIAMETER FUNNEL

H91462-4610	Replacement Coarse Porous Plate	603.00
H91462-4361	Replacement Medium Porous Plate	600.00
H91463-1610	Replacement Perforated Plate	261.00

FOR 914MM (36") INSIDE DIAMETER FUNNEL

H91462-4091	Replacement Coarse Porous Plate	1089.00
H91462-4914	Replacement Medium Porous Plate	861.00
H91463-1914	Replacement Perforated Plate	357.00

Filter Paper Retaining Rings Are Available, Call for Current Pricing

Buchner Funnel Single Piece

Fast Flowing with Fritware® Porous Plates

One piece HDPE funnels with sealed-in, polyethylene filter plate of coarse porosity (90-130 microns) offer good retention and fast flow under vacuum. Useable with or without filter paper. Use below 52°C (125°F).



CATALOG NO.	CAPACITY	PLATE DIAMETER	RIM TO PLATE	CASE	PRICE/EACH
F14608-0000	50ml	36mm	35mm	12 Each	\$8.80
F14609-0000	150ml	70mm	51mm	6 Each	10.30
F14610-0000	250ml	104mm	50mm	6 Each	14.70
F14611-0000	1000ml	122mm	104mm	4 Each	20.60