



SOLAR THERMOMETER 800018 — NEW

Powered forever by its own solar panel this unit perpetually monitors ambient temperature without batteries. Large $\frac{5}{8}$ " digital display easily readable from a distance. Switches from $-5/50^{\circ}\text{C}$ to $23/122^{\circ}\text{F}$. Stands alone or wall mounts. Dim: $2\frac{7}{8}$ " x $2\frac{3}{4}$ " x $\frac{3}{4}$ ". Weight: 1.7 oz.



**800019 AMBIENT TEMPERATURE THERMOMETER
(PREVIOUSLY 840019)**

The large LCD continuously displays ambient temperature in laboratories, inside incubators, or anywhere. Compact, measuring only $3\frac{3}{8}$ x $1\frac{3}{4}$ x $\frac{1}{2}$ ". Switches from $^{\circ}\text{C}$, ($-50/70^{\circ}\text{C}$), to $^{\circ}\text{F}$, ($-58/160^{\circ}\text{F}$), with a resolution of 0.1° . Can also display the time. Designed to stand alone or mount on the wall. Comes with a 1.5 Volt button cell battery. Weight: 1.7 oz.



**800020 INDOOR/OUTDOOR THERMOMETER
(PREVIOUSLY 840020)**

Same as 840019 but with the addition of an external sensor with a 9 ft. cord. A switch changes the readout from indoor to outdoor temperature. Ideal for use as an indoor/outdoor thermometer, in incubators, or any place where simultaneous internal and external readings are required.