



## RED LED CLOCK 810031

This bench-top alarm clock features a bright red LED (3" x 1") easily seen from anywhere in the lab. Especially useful for work in dark rooms, HPLC, TLC, and electrophoresis. In the event of a power failure a 9V battery back up (not included) maintains the correct time internally for up to 12 hours while the LED remains off to preserve battery life. UL listed. Works on AC120V and comes with a 60" electric power cord. Weight: 9 oz (260 g).

Dim: 4½" x 4" x 2" (101 x 50 x 99 mm).