




WATER RESISTANT STOPWATCH 810012 — NEW

Counts up to 24 hours in 1/100th second resolution for the first 30 minutes with a 1 second resolution thereafter. Split interval function. Also has a 12/14 hour alarm clock, calendar, and an hourly chime option. Dim: 2 1/2" x 2" x 1/2". Weight: 1.2 oz.



With N.I.S.T. traceable certificate of compliance: 810012C



LARGE DISPLAY STOPWATCH 810022

Click stop controls give the definitive feel of a mechanical stopwatch but with hundreds of times the accuracy. Large LCD for good visibility. Rugged ABS plastic case. Counts to 24 hours in 1/100th second resolution for the first 30 minutes with 1 second resolution thereafter. Split interval, clock, and calendar functions. Comes with leatherette case. Dim: 3 3/4" x 2 1/4" x 5/8". Weight: 2.6 oz.



With N.I.S.T. traceable certificate of compliance: 810022C