



## NEW - LARGE DISPLAY TIMER 810005

Features the largest display available on a timer; 3" X 2 1/2". Can be easily read from across a room. Simultaneously displays either one count up channel and a clock, or two count down channels. All channels have a range of 24 hours with 1 second resolution. Also indicates how much time has elapsed after a count down alarm sounds and conveniently recalls the last timer setting for reuse. White case has extra large start/stop button and a loud alarm. Wall mount or fold for bench top use. Runs for approximately 2 years on one AAA battery, (included). Dim: 5" X 3 3/4" X 1 1/2". Weight: 6.7 oz.



With N.I.S.T. traceable certificate of compliance: 810005C

## NEW - LARGE DISPLAY TIMER 810006

Same as 810005 in black.



With N.I.S.T. traceable certificate of compliance: 810006C

## WATER RESISTANT WALL CLOCK 810034

A quartz crystal movement in a sealed water and dust resistant case ensures accuracy and reliability even in conditions of 100% humidity. Use in areas of high dust and humidity such as in factories, greenhouses, outdoors, or in the lab. The attractive contemporary case features a green face and bright red second hand and is 7" in diameter. Weight: 3.8 oz. Powered by one AA battery which is included.



With N.I.S.T. traceable certificate of compliance: 810034C