



5 CHANNEL TIMER 810015

4 independent channels each of which counts up or down for any interval from 1 second to 100 hours. Extra large display indicates channels in use. Each one has a distinctive 1 minute alarm. Automatically begins to count up to indicate time elapsed after count down reaches zero. Memory recalls last timer setting for reuse. Also has a clock mode. Has a stand, pocket clip, and a magnet. Dim: 2⁷/₈" x 2¹/₂" x ⁹/₁₆". Weight: 2 oz.



With N.I.S.T. traceable certificate of compliance: 810015C



MINI TIMER 99 MIN. 810019

Compact and easy to hold. Large, easy-to-read LCD and controls. Counts down from 99 minutes in intervals of 1 second. 30 second alarm when time is up. Memory recalls last timer setting for reuse. Touch Tone buttons. Magnetic back. Dim: 3¹/₂" x 1⁷/₈" x ¹/₂". Weight: 1 oz.



With N.I.S.T. traceable certificate of compliance: 810019C



MINI TIMER 16 HOUR 810021

Same ergonomic design as 810019. Counts up to 16 hours with resolution in minutes. Warning tones at 10 min. and 5 min. 20 second alarm when time is up.



With N.I.S.T. traceable certificate of compliance: 810021C