

KILOTECH



February 28, 2010 to September 1, 2010

PORTION CONTROL



Item #	Model	Capacity & Graduation	Platter	PACK	LB	List (e)
851165	KPC 160	5kg x 1g / 180oz. x 0.1oz. +/-3d	Pastic bowl 7.25" diam.	12	2.52	\$60.00
851167	KPC2000-03 (kg, lb, g, oz, lb-oz)	3kg x 1g / 6lb. x 0.002lb.	9.5 x 8.5	6	8.1	\$135.00
851168	KPC2000-06 (kg, lb, g, oz, lb-oz)	6kg x 2g / 12lb. x 0.005lb.				
851286	KPC 1000SS-0.5	500g x 0.1g / 1 lb x 0.0002 lb	5 7/8 x 4 3/4"	4	3.75	\$99.00
851287	KPC 1000SS-2.5	2500g x 0.5g / 5.5 lb x 0.001 lb				
851288	KPC 1000SS-2	2000g x 1g / 4.4 lb x 0.002 lb				
851289	KPC 1000SS-5	5000g x 1g / 11 lb x 0.002 lb				
851146	KPC 5000-2	2kg x 0.5g / 4 lb 6 oz x 1/8 oz	6 3/4" x 7	4	3	\$99.00
851147	KPC 5000-5	5kg x 1g / 11 lb. x 1/4 oz.				
851154	KD-200-110 (g, oz.)	1kg x 1g / 35 oz. x 0.05 oz.	7.2 x 7.2	6	4.3	\$168.00
851155	KD-200-210 (g, oz.)	2kg x 2g / 70 oz. x 0.1 oz.				
851156	KD-200-510 (g, oz.)	5kg x 5g / 11 lb. x 0.2 oz.				
851157	KPC 1144	6kg x 1g / 13 lb x 0.002 lb. +/- 3d	6.75 x 6.75"	4	2.9	\$65.00

Accessories

Item #	Description	List (e)
861355	KPC 2000/KWS 152 clear dust and water cover	\$4.00
861523	KPC 5000 clear dust and water cover	\$4.00

Mix and Match any scale on this page to receive Quantity Discount Pricing. Discount only applies to single orders that are shipped to one location.

ELECTRONIC BENCH SCALES



Item #	Model	Capacity & Graduation	Platter	PACK	LB	List	Inspection
851169	KWS 152 (kg, lb, g, oz, lb-oz)	15kg x 2g / 30lb. x 0.005lb.	9.5 x 8.5	6	8.1	\$135.00	
851280	KPS-68MS Bench/Parcel scale	68kg x 100g / 150 lb. x 0.2 lb.	12 1/2 x 11 1/4	4	11.46	\$150.00	
851281	KPS-180MS Bench/Parcel scale	180kg x 200g / 400lb. x 0.5lb.					
851282	KPS 25-2 Bench/Parcel scale w/RS232	25kg x 10g / 55 lb. x 0.02 lb. (DNR)	12 1/2 x 11 1/4	4	11.46	\$180.00	
851283	KPS 50-2 Bench/Parcel scale w/RS232	50kg x 20g / 110 lb. x 0.05 lb.					
851284	KPS 100-2 Bench/Parcel scale w/RS232	100kg x 50g / 220 lb. x 0.1 lb.					
851285	KPS 150-2 Bench/Parcel scale w/RS232	150kg x 100g / 330lb. x 0.2 lb.					
851290	KPS-060SS Bench/Parcel scale	60kg x 20g / 130 lb. x 0.05 lb. / 2116 oz x 1 oz.	14" x 16" S/S	4	14.4	\$230.00	
851291	KPS-150SS Bench/Parcel scale	150kg x 50g / 330 lb. x 0.1 lb. / 5290 oz x 2 oz.					
853182	KWD500-05 (front & rear display)	2.5kg x 0.001kg / 5lb. x 0.002lb. 2500g x 1g / 80oz. x 0.05oz.	8 1/2 x 7	4	9.1	\$220.00	NTEP
853183	KWD500-10 (front & rear display)	5kg x 0.002kg / 10 lb. x 0.005 lb. 5000g x 2g / 160oz. x 0.1oz.					
853184	KWD500-20 (front & rear display)	10kg x 0.005kg / 20 lb. x 0.01 lb. 10000g x 5g / 400oz. x 0.2oz.					
853185	KWD500-50 (front & rear display)	20kg x 0.01kg / 50 lb. x 0.02 lb. 20000g x 10g / 800oz. X 0.5oz.					
851836	KSB-4KR weighing scale x 1gram	4kg x 1g / 8 lb x 0.002 lb.	11 1/4 x 9 1/4	6	12.75	\$200.00	
851837	KSB-6KR weighing scale x 1gram	6kg x 1g / 13 lb x 0.002 lb.				\$225.00	
851838	KSB-10KR weighing scale x 1gram	10kg x 1g / 22 lb x 0.002 lb.				\$250.00	
851253	KWS301 with Back light	15kg x 0.5g / 30 lb. x 0.005 lb.	9 1/2 x 11" s/s	2	11.75	\$340.00	
851253	KWS301 with Back light	30kg x 1g / 60 lb. x 0.002 lb.	9 1/2 x 11" s/s	2	11.75	\$340.00	

ACCESSORIES

Item #	Description	List (e)
861355	KPC 2000/KWS 152 clear dust and water cover	\$4.00
861253	KCS/KWS 301 clear dust and water cover	\$4.00

Prices subject to change Les prix sont sujet à changements

KILOTECH



February 28, 2010 to September 1, 2010

PRICE COMPUTING SCALE

COUNTING SCALES

 KRS 1000 Price computing 	 Printer for KRS 1000 & KCS 1000 	 Labels for KRS 1000 2 9/32" x 1 5/8" 58 x 40mm Food grade glue Pre-printed	 Labels for KCS 1000 2 1/6" x 1 5/8" 56 x 40mm Food grade glue Blank	 KCS 1000 Counting 250 item memory 	 KCS 301 Counting
--	--	---	--	--	------------------------------------

Item #	Model	Capacity & Graduation	Platter	PACK	LB	List		
851302	KRS 1000 Scale with 261 PLU and RS232 to connect to Item #853500 label printer	0-15 lb x 0.005 lb / 15-30 lb x 0.01 lb 0-6kg x 2g / 6-15kg x 5g	11 1/4" x 8 1/4"	1	9	\$400.00		
853500	Printer for KRS 1000 and KCS 1000 with RTC card and label separator kit	Uses Kilotech Item #862049 58 x 40mm labels	5 x 9 x 6.5"	1	3	\$490.00		
862049	KRS 1000 labels lb. / kg	58 x 40 mm (700 / roll ~ 30 rolls / case)		1 cs	22	\$185.00		
862051	KCS 1000 labels lb. / kg	56 x 40 mm (600 / roll ~ 30 rolls / case)		1 cs	22	\$185.00		
851310	KCS 1000-6 256 item memory	6kg x 1g / 12 lb. x 0.001 lb.	11 1/4" x 8 1/4"	4	9	\$300.00		
851311	KCS 1000-15 256 item memory	15kg x 1g / 30 lb. x 0.002 lb.						
851312	KCS 1000-30 256 item memory	30kg x 2g / 60 lb. x 0.005 lb.						
851201	KCS301-3 with Back light	3kg x 0.1g / 6.6 lb. x 0.0002 lb.	11 " x 9 1/4" s/s	2	11.8	\$380.00		
851202	KCS301-6 with Back light	6kg x 0.2g / 13.2 lb. x 0.0004 lb.						
851203	KCS301-15 with Back light	15kg x 0.5g / 33 lb. x 0.001 lb.						
851204	KCS301-30 with Back light	30kg x 1g / 66 lb. x 0.002 lb.						

Mix and Match any scale or printer on this page to receive Quantity Discount Pricing. Discount only applies to single orders that are shipped to one location.


PLATFORMS, INDICATORS, CRANE SCALES

 KBS 	 KWS 1000 	 KIN 	 VHS 300
--------------------	-------------------------	--------------------	------------------------

Item #	Model	Capacity & Graduation	Platform	PACK	LB	List		Inspection
851270	KBS 30	30kg x 10g / 60 lb. x 0.02 lb.	12 x 16" s/s	2	30.2	\$690.00		
851271	KBS 60	60kg x 20g / 150 lb. x 0.05 lb.						
851272	KBS 150	150kg x 50g / 300 lb. x 0.1 lb.						
851684	KWD 1000-60 IP64 base	60kg x 20g / 120 lb x 0.05lb	13 1/2 x 17" s/s	1	37.4	\$750.00		NTEP
851685	KWD 1000-150 IP64 base	150kg x 50g / 300lb x 0.1lb	17 1/4 x 22 1/4 s/s		62.7	\$900.00		NTEP
851688	KWD 1000-300 IP64 base	300kg x 100g / 600lb x 0.2lb						NTEP
						List (e)		
851704	KIN 201 6 button programmable	ABS housing (rechargeable battery in option)		8	3.1	\$450.00		NTEP
851705	KIN 202 18 button programmable				3.1	\$510.00		
851706	KIN 203 18 button programmable	Stainless steel housing rechargeable battery in option			5	\$700.00		
						List (e)		
851711	VHS 300-1000 kg/lb., IP54	1000kg x 0.2kg / 2204lb. x 0.5lb.	N/A	1	30.8	\$2,600.00		
851714	VHS 300-2000 kg/lb., IP54	2000kg x 0.5kg / 4000lb. x 1lb.	N/A	1	30.8			
851713	VHS 300-5000 kg/lb., IP54	5000kg x 1kg / 11023lb. x 2lb.	N/A	1	55.0	\$3,000.00		
851715	VHS 300-10 000 kg/lb., IP54	10 000kg x 2kg / 22 046lb. x 5lb.	N/A	1	79.2	\$3,600.00		
851716	VHS 300-20 000 kg/lb., IP54	20 000kg x 5kg / 44 092lb. x 10lb.	N/A	1	165.0	\$4,100.00		

Prices subject to change Les prix sont sujet à changements

KILOTECH



February 28, 2010 to September 1, 2010

KILOTECH TOPLOADER / BENCH ENTRY LEVEL LAB & EDUCATIONAL BALANCES



Item #	Model	Capacity & Graduation	Platter / Plateau	PACK	LB	List (e)		
852925	KLE-200 Econo	200g x 0.1g (g, oz, lb., ozt, dwt, cnt)	6½ x 5" s/s	12	3.1	\$99.00		
852926	KLE-2000 Econo	2000g x 1g (g, oz, lb., ozt, dwt, cnt)						
851940	KHA200R Draft Shield	200g x 0.01g	120mm round		3.93	\$240.00		
851941	KHA600R Draft Shield	600g x 0.01g			4.18	\$320.00		
851942	KHA1000R Draft Shield	1000g x 0.01g			3.68	\$380.00		
851943	KHA201R Draft Shield	200g x 0.1g			3.87	\$180.00		
851944	KHA601R Draft Shield	600g x 0.1g			4.1	\$210.00		
851945	KHA1201R	1200g x 0.1g			6¾ x 6¾" s/s	4.63	\$240.00	
851946	KHA2401R	2400g x 0.1g	4.86			\$270.00		
851947	KHA5001R	5000g x 0.1g	4.81			\$300.00		

Note: All KHA have CSA/UL adapter and g, oz, lb., ozt, dwt and ct weighing modes and CNT feature. The KHA has a mono-directional RS232.
 ozt = Troy Onces dwt = Pennyweight ct = Carats cnt = Counting

Mix and Match any scale on this page to receive Quantity Discount Pricing. Discount only applies to single orders that are shipped to one location.


INDIVIDUAL TEST WEIGHTS							
Item #	Model	Capacity		Tolerance	PACK	LB.	List
853555	Class M1 (OIML R111) with statement of accuracy; non-traceable	100g	Cylindrical, stainless steel, packed in plastic screw cap container with cushioned interior	+/- 5mg	1	0.22	\$20.00
853554		200g		+/- 10mg		0.44	\$25.00
853553		500g		+/- 25mg		1.1	\$38.00
853552		1kg		+/- 50mg		2.2	\$60.00
853551		2kg		+/- 100mg		4.4	\$88.00
853556		20kg		+/- 1g		44	\$100.00
853557		50 lb.	+/- 1.25g	50		\$110.00	
853429	No class and no statement of accuracy	20g	Polished alloy	N/A		0.08	\$3.00
853431		100g		N/A		0.44	\$6.00

KILOTECH



February 28, 2010 to September 1, 2010

DIAL SCALES



See below or web site www.kilotech.com for more details. Some models with "RAPID STOP" (KDP) and/or N.S.F. approved (NS).

Item #	Model	Capacity & Graduation	Platter	PACK	LB	List (e)				
852408	KDP 6-1000SPLRD ROTATING	1kg x 5g / 2 lb. x 1/4 oz.	8" x 8"	8	4.4	\$60.00				
852409	KDP 6-2000SPLRD ROTATING	2kg x 10g / 5 lb. x 1/2 oz.	8" x 8"	8	4.4	\$60.00				
852404	KAM 6-2000SPLRD ROTATING	2kg x 10g / 5 lb. x 1/2 oz.	8" x 8"	8	4.4	\$60.00				
852405	KDP 6-1000SPL RAPID STOP	1kg x 5g / 35 oz. x 1/4 oz.	8" x 8"	8	5	\$60.00				
852406	KDP 6-2000SPL RAPID STOP	2kg x 10g / 4 lb. 6 oz. x 1/2 oz.	8" x 8"	8	5	\$60.00				
852407	KDP 6-5000SPL RAPID STOP	5kg x 25g / 11 lb. x 1 oz.	8" x 8"	8	5	\$60.00				
852281	KAM 3218PL	1kg x 5g / 34 oz. x 1/8 oz.	7½ x 7½	4	6.78	\$90.00				
852285	KAM 512PL	2kg x 5g / 5 lb. x 1/2 oz.	7½ x 7½	4	6.78	\$90.00				
852291	KAM 1101PL	5kg x 25g / 11 lb. x 1 oz.	7½ x 7½	4	6.78	\$90.00				
852295	KAM 2501PL	12kg x 25g / 26 lb. x 1 oz.	7½ x 7½	4	6.78	\$90.00				
852299	KAM 5002PL	22kg x 50g / 50 lb. x 2 oz.	7½ x 7½	4	6.78	\$90.00				
852300	KAM 1004PL	45kg x 100g / 100 lb. x 4 oz.	13 x 13"	2	24.0	\$200.00				
852350	KAM 2008PL	90kg x 200g / 200 lb. x 8 oz.	15½ x 15½"	2	27.6	\$220.00				
852351	KBI 64 FL. OZ. & OZ.	64 fl oz. x 1 fl oz. / 160 oz. X 1/2 oz.	7½ x 7½	4	7.75	\$180.00				


Mix and Match any scale on this page to receive Quantity Discount Pricing. Discount only applies to single orders that are shipped to one location.

KILOTECH



February 28, 2010 to September 1, 2010

HANGING DIAL, TRIPPLE BEAM, BAKER'S, BEAM



Item #	Model	Capacity & Graduation	Platter Plateau	PACK	LB	List (e)			
852410	KHS 1022	10kg x 25g / 22 lb. x 1 oz. 7" dial	12 1/4" x 2"	4	7.2	\$129.00			
851681	KHS-C360 with large "J" hook	60kg x 200g / 125lbs. x 0.5lbs. 8" dial	N/A	8	7.44	\$140.00			
851682	KHS-C3120 with large "J" hook	120kg x 500g / 250lbs. x 1lb. 8" dial	N/A	8	7.87	\$140.00			
851683	KHS-C3250 with large "J" hook	250kg x 1kg / 540lbs. x 2 lb. 8" dial	N/A	8	7.95	\$140.00			
853180	KW 2610 Includes weight set	2610g x 0.1g	6" diameter	4	6.4	\$90.00			
852820	BK-8	8kg x 10g	9" chrome	2	30.64	\$260.00			
852825	BK-16	16lb. x 1/4oz.	9" chrome	2	30.64	\$260.00			

Accessories

852850	BK stainless steel scoop		12 x 18½ x 5"	1	3	\$75.00			
852855	BK counterweight for BK s/s scoop		N/A	1	2	\$25.00			
862851	BK plastic scoop and counter weight		N/A	1	2	\$65.00			
862852	BK plastic scoop only		11 x 18¾ x 4¾"	1	2	\$45.00			

Item #	Model	Capacity & Graduation	Platter Plateau	PACK	LB	List (e)			Inspection
850100	KT-SP50KG - Bench (kg in front)	50kg x 20g / 100 lb. x 1 oz.	12 x 16"	1	57.0	\$280.00			NTEP
850101	KT-SP100KG - Bench (kg in front)	100kg x 50g / 200 lb. x 1 oz.	12 x 16"	1	59.0	\$320.00			NTEP
850102	KT-SP300KG - Platform (kg in front)	300kg x 100g / 600 lb. x 4 oz.	16 x 24"	1	178.0	\$540.00			NTEP
850103	KT-SP500KG - Platform (kg in front)	500kg x 100g / 1200 lb. x 4 oz.	24 x 25"	1	231.5	\$600.00			NTEP
850104	KT-SP1000KG - Platform (kg in front)	1000kg x 500g / 2300 lb. x 8 oz.	29 x 29"	1	540.0	\$1,000.00			NTEP
850110	KT-SP100LB - Bench (lb in front)	100 lb. x 1 oz. / 50kg x 20g	12 x 16"	1	57.0	\$280.00			NTEP
850111	KT-SP200LB - Bench (lb in front)	200 lb. x 1 oz. / 100kg x 50g	12 x 16"	1	59.0	\$320.00			NTEP
850112	KT-SP600LB - Platform (lb in front)	600 lb. x 4 oz. / 300kg x 100g	16 x 24"	1	178.0	\$540.00			NTEP
850113	KT-SP1200LB - Platform (lb in front)	1200 lb. x 4 oz. / 500kg x 100g	24 x 25"	1	231.5	\$600.00			NTEP
850114	KT-SP2300LB - Platform (lb in front)	2300 lb. x 8 oz. / 1000kg x 500g	29 x 29"	1	540.0	\$1,000.00			NTEP