



# STOCK AND STANDARD PRODUCT CATALOG

Comar, Inc.  
One Comar Place  
Buena, NJ 08310

Toll Free: 800-962-6627  
Direct: 856-692-6100  
Fax: 856-692-9251

[www.comar.com](http://www.comar.com)  
[information@comar.com](mailto:information@comar.com)



# Comar Stock and Standard Products Catalog

Thank you for using the Comar Stock and Standard Products Catalog!

This catalog represents Comar's Stock and Standard Product Offerings found in Comar's new Online Catalog. To view an interactive version of this catalog, with additional information, detailed product pictures, product drawings and the ability to request samples or quotes, please visit, [www.comar.com/ecatalog](http://www.comar.com/ecatalog).

## Product Categories

Products are separated into the following categories and may be represented in this printed version:

### Plastics Products

Oral Dispensers  
Plastic Dropper Assemblies  
3cc, 6cc, 10cc & 15cc Bottle Packages  
Plastic Bottles  
Closures  
Injection Molded Products

### Glass Products

Glass Dropper Assemblies  
Media Tubes  
Screw Thread Vials  
Serum Vials  
Shell Vials  
Sure Thread Vials  
EPA/Water Vials

## Product Information

Information on each product reflects a combination of individual attribute information, including: Size, Length, Diameter, Fill Volumes, Capacity, Neck Finish, Style, Color, Plastic Material, and Glass Type.

Each product is sorted in its category by Size, Capacity, Neck Finish or Fill Volumes in Ascending Order.

Products represented in this catalog are Comar's Stock and Standard Product Offerings.

## Stock Products

Comar's "Stock Program" is designed to maintain constant availability of selected products to smooth out delivery patterns caused by volume fluctuations and its immediate customer requirement for product. Comar will continue to maintain a safety stock for selected "stock items" based on typical order patterns and quantities. For your convenience, Comar will replenish the stock as necessary. Stock Items may be added or removed from our list at any time.

## Standard Products

Comar's Standard Products are manufactured with production molds owned by Comar. Standard Products may or may not be available in specific quantities at the time of request. In this case, the requested product will be scheduled for production and manufactured for your order. Comar does keep Standard Products in inventory for sampling purposes.

## Custom Searches with the Comar Online Catalog

To create custom searches by product attribute, please visit the Comar Online Catalog at [www.comar.com/ecatalog](http://www.comar.com/ecatalog), log-in, click the Custom Search link on the Left Navigation Bar and begin your search.

## Custom Products

If you don't see a product here that matches your requirements, Comar specializes in developing custom packaging for its customers. Please visit [www.comar.com](http://www.comar.com), and click the Comar [ASR] "A Solutions Resource" logo on the Home Page. Here, we outline our flexible custom packaging process and walk you through each step of the custom packaging process with a product case study. You may also request a no-charge consultation with Comar's Packaging Development Team.

## Thank You

The Comar Team would like to thank you for using the Comar Stock and Standard Products Catalog. We hope that you will visit and use the Comar Online Catalog as a resource for all of your packaging needs. [www.comar.com/ecatalog](http://www.comar.com/ecatalog).

Please feel free to provide us with any feedback that you feel might help Comar better serve your needs.

# About Comar

---



## **Mission**

We will embrace innovation and customer responsiveness as the principal tools for delivering long-term value to customers and Comar. We will accomplish this by combining product value with affordable quality and reliability to minimize the total cost of ownership for customers who value what we sell.

## **Vision**

We will combine innovation, leading-edge technology, uncompromising quality and exemplary customer responsiveness to create and deliver packaging solutions that delight our customers throughout every aspect of our relationship.

For over 55 years, Comar, Inc. has offered the Pharmaceutical, Diagnostic, Health and Personal Care Markets with Innovative Packaging Solutions with an emphasis on Quality and Service.

“Quality on Time” is a way of life at Comar, and the way things are done. The many years of Innovative Design Experience continue to make Comar the leader in customized glass and plastic packaging.

Comar’s Stock and Standard Products Catalog provides you with a detailed view of what Comar offers to the Market.

Contact Customer Service or your Sales Representative for additional information or to place an order.

Thank you for using the Comar Stock and Standard Products Catalog.

# Oral Dispensers




Part Number	Description	Size	Tip Style	Capacity	Color	Type
14-0209-001	Cover Cap				Natural	STOCK
14-0209-004	Cover Cap				Natural	STOCK
14-0223-001	Cover Cap Stand				White	STOCK
41-0001-002	Oral Dispenser		Short Tip	3mL		STOCK
41-0002-002	Oral Dispenser		Long Tip	3mL		STOCK
41-0003-002	Oral Dispenser		Short Tip	5mL		STOCK
41-0004-002	Oral Dispenser		Long Tip	5mL		STOCK
41-0005-002	Oral Dispenser		Short Tip	10mL		STOCK
41-0006-002	Oral Dispenser		Long Tip	10mL		STOCK
41-0008-001	Oral Dispenser		Short Tip	1mL		STOCK
22-0218-001	Press-In Bottle Adapter	24mm			Natural	STOCK
22-0198-001	Press-In Bottle Adapter	28mm			Natural	STOCK
22-0231-002	Press-In Bottle Adapter	33mm			Natural	STOCK

Comar, Inc. One Comar Place, Buena, NJ 08310

Toll Free: 800-962-6627 • Fax: 856-692-9251 • Distributor Hotline: 866-807-9070

www.comar.com • information@comar.com

# Plastic Dropper Assemblies



Part Number	Closure Color	Bulb Color	Pipette Style	Bulb Material	Closure Size	Pipette Diameter	Pipette Length	Type
20-0627-001	Black	Black	Stub	Rubber	13-425	4.5mm	45mm	STANDARD
20-0678-001	White	White	Stub	Rubber	13-425	4.5mm	45mm	STANDARD
20-0821-001	Black	Black	Blunt	Rubber	15-415	5mm	66mm	STANDARD
20-0822-001	White	White	Blunt	Rubber	15-415	5mm	66mm	STANDARD
20-0029-002	Black	Black	Straight	Rubber	18-400	7mm	62mm	STOCK
20-0527-002	Black	Black	Straight	Monprene	18-400	7mm	62mm	STOCK
20-0152-003	Black	Black	Straight	Rubber	20-400	7mm	89mm	STOCK
20-0425-003	Black	Black	Straight	Rubber	20-400	7mm	72mm	STOCK
20-0483-003	Black	Black	Straight	Monprene	20-400	7mm	72mm	STOCK
20-0626-003	Black	Black	Straight	Monprene	20-400	7mm	89mm	STOCK
20-0825-001	Black	Black	Blunt	Rubber	22-400	7mm	59mm	STANDARD
20-0826-001	White	White	Blunt	Rubber	22-400	7mm	59mm	STANDARD

Comar, Inc. One Comar Place, Buena, NJ 08310

Toll Free: 800-962-6627 • Fax: 856-692-9251 • Distributor Hotline: 866-807-9070

www.comar.com • information@comar.com

# 3CC Bottle Package



Part Number	Description	Shape	Neck Finish	Capacity	Color	Finish Size	Style	Drop Size	Type
12-0004-001	3cc Bottle	Cylinder	8-425	3cc	Natural				STOCK
12-0004-002	3cc Bottle	Cylinder	8-425	3cc	White				STOCK
12-0004-003	3cc Bottle	Cylinder	8-425	3cc	Black				STOCK
14-0165-013	Dome Closure				White	8-425	Dome		STOCK
14-0165-014	Dome Closure				White	8-425	Dome		STOCK
14-0165-017	Dome Closure				Natural	8-425	Dome		STOCK
14-0165-021	Dome Closure				Black	8-425	Dome		STOCK
22-0057-001	Controlled Drop Tip		8mm		Natural			42uL	STOCK
22-0057-007	Controlled Drop Tip		8mm		White			42uL	STOCK
22-0057-008	Controlled Drop Tip		8mm		Black			42uL	STOCK

Comar, Inc. One Comar Place, Buena, NJ 08310

Toll Free: 800-962-6627 • Fax: 856-692-9251 • Distributor Hotline: 866-807-9070

www.comar.com • information@comar.com

# 6CC Bottle Package



Part Number	Description	Shape	Neck Finish	Capacity	Color	Finish Size	Style	Drop Size	Type
12-0093-001	6cc Bottle	Cylinder	13-425	6cc	Natural				STOCK
12-0093-002	6cc Bottle	Cylinder	13-425	6cc	White				STOCK
12-0093-003	6cc Bottle	Cylinder	13-425	6cc	Black				STOCK
14-0166-012	Dome Closure				Black	13-425	Dome		STOCK
14-0166-015	Dome Closure				White	13-425	Dome		STOCK
14-0166-018	Dome Closure				Natural	13-425	Dome		STOCK
22-0200-001	Controlled Drop Tip		13mm		Natural			40uL	STOCK
22-0200-002	Controlled Drop Tip		13mm		Natural			40uL	STOCK
22-0200-006	Controlled Drop Tip		13mm		White			40uL	STOCK
22-0200-007	Controlled Drop Tip		13mm		White			40uL	STOCK
22-0200-010	Controlled Drop Tip		13mm		Black			40uL	STOCK
22-0200-011	Controlled Drop Tip		13mm		Black			40uL	STOCK

Comar, Inc. One Comar Place, Buena, NJ 08310

Toll Free: 800-962-6627 • Fax: 856-692-9251 • Distributor Hotline: 866-807-9070

www.comar.com • information@comar.com

# 10CC Bottle Package



Part Number	Description	Shape	Neck Finish	Capacity	Color	Finish Size	Style	Drop Size	Type
12-0016-001	10cc Bottle	Boston Round	13-425	10cc	Natural				STOCK
12-0016-003	10cc Bottle	Boston Round	13-425	10cc	White				STOCK
12-0016-008	10cc Bottle	Boston Round	13-425	10cc	Black				STOCK
14-0166-012	Dome Closure				Black	13-425	Dome		STOCK
14-0166-015	Dome Closure				White	13-425	Dome		STOCK
14-0166-018	Dome Closure				Natural	13-425	Dome		STOCK
22-0200-001	Controlled Drop Tip		13mm		Natural			40uL	STOCK
22-0200-002	Controlled Drop Tip		13mm		Natural			40uL	STOCK
22-0200-006	Controlled Drop Tip		13mm		White			40uL	STOCK
22-0200-007	Controlled Drop Tip		13mm		White			40uL	STOCK
22-0200-010	Controlled Drop Tip		13mm		Black			40uL	STOCK

Comar, Inc. One Comar Place, Buena, NJ 08310

Toll Free: 800-962-6627 • Fax: 856-692-9251 • Distributor Hotline: 866-807-9070

www.comar.com • information@comar.com



# 15CC Bottle Package (13mm Tip)



Part Number	Description	Shape	Neck Finish	Capacity	Color	Finish Size	Style	Drop Size	Type
12-0094-001	15cc Bottle	Boston Round	15-415	15cc	Natural				STOCK
12-0094-002	15cc Bottle	Boston Round	15-415	15cc	White				STOCK
12-0094-003	15cc Bottle	Boston Round	15-415	15cc	Black				STOCK
14-0175-020	Dome Closure				Black	15-415	Dome		STOCK
14-0175-021	Dome Closure				Natural	15-415	Dome		STOCK
14-0175-022	Dome Closure				White	15-415	Dome		STOCK
22-0200-001	Controlled Drop Tip		13mm		Natural			40uL	STOCK
22-0200-002	Controlled Drop Tip		13mm		Natural			40uL	STOCK
22-0200-006	Controlled Drop Tip		13mm		White			40uL	STOCK
22-0200-007	Controlled Drop Tip		13mm		White			40uL	STOCK
22-0200-010	Controlled Drop Tip		13mm		Black			40uL	STOCK

Comar, Inc. One Comar Place, Buena, NJ 08310

Toll Free: 800-962-6627 • Fax: 856-692-9251 • Distributor Hotline: 866-807-9070

www.comar.com • information@comar.com

# 15CC Bottle Package (15mm Tip)



Part Number	Description	Shape	Neck Finish	Capacity	Color	Finish Size	Style	Drop Size	Type
12-0008-001	15cc Bottle	Boston Round	15-415	15cc	Natural				STOCK
12-0008-002	15cc Bottle	Boston Round	15-415	15cc	White				STOCK
12-0008-008	15cc Bottle	Boston Round	15-415	15cc	Black				STOCK
14-0175-020	Dome Closure				Black	15-415	Dome		STOCK
14-0175-021	Dome Closure				Natural	15-415	Dome		STOCK
14-0175-022	Dome Closure				White	15-415	Dome		STOCK
22-0215-001	Controlled Drop Tip		15mm		Natural			40uL	STOCK
22-0215-002	Controlled Drop Tip		15mm		Natural			40uL	STOCK
22-0215-004	Controlled Drop Tip		15mm		White			40uL	STOCK
22-0215-005	Controlled Drop Tip		15mm		White			40uL	STOCK
22-0215-006	Controlled Drop Tip		15mm		Black			40uL	STOCK
22-0215-007	Controlled Drop Tip		15mm		Black			40uL	STOCK

Comar, Inc. One Comar Place, Buena, NJ 08310

Toll Free: 800-962-6627 • Fax: 856-692-9251 • Distributor Hotline: 866-807-9070

www.comar.com • information@comar.com

# Plastic Bottles



Part Number	Shape	Neck Finish	Capacity	Material	Color	Type
12-0042-004	Cylinder	15-415	1/4oz	High Density Polyethylene	Natural	STANDARD
12-0042-005	Cylinder	15-415	1/4oz	High Density Polyethylene	White	STANDARD
12-0050-001	Cylinder	15-415	1/4oz	Low Density Polyethylene	Natural	STOCK
12-0104-001	Cylinder	13-425	6cc	High Density Polyethylene	Natural	STANDARD
12-0104-002	Cylinder	13-425	6cc	High Density Polyethylene	White	STANDARD
12-0104-003	Cylinder	13-425	6cc	High Density Polyethylene	Black	STANDARD
12-0014-002	Modified Cylinder	20-400	20mL	High Density Polyethylene	Natural	STANDARD
12-0052-002	Cylinder	22-400	20mL	High Density Polyethylene	Natural	STANDARD
12-0012-002	Boston Round	28-400	25cc	High Density Polyethylene	Natural	STOCK
12-0012-003	Boston Round	28-400	25cc	High Density Polyethylene	White	STOCK
12-0048-001	Boston Round	28-480	100cc	High Density Polyethylene	Natural	STANDARD
12-0084-003	Blake	33-400	75mL	High Density Polyethylene	Natural	STANDARD
12-0028-002	Square	33-400	95mL	High Density Polyethylene	White	STANDARD
12-0049-001	Boston Round	28-480	150cc	High Density Polyethylene	Natural	STANDARD
12-0026-001	Square	38-400	165mL	High Density Polyethylene	White	STANDARD

Comar, Inc. One Comar Place, Buena, NJ 08310

Toll Free: 800-962-6627 • Fax: 856-692-9251 • Distributor Hotline: 866-807-9070

www.comar.com • information@comar.com

# Closures




Part Number	Description	Color	Finish Size	Style	Type
14-0168-004	Closure	Natural	13-415	Linerless	STANDARD
14-0168-006	Closure	Black	13-415	Linerless	STANDARD
14-0168-007	Closure	White	13-415	Linerless	STANDARD
14-0269-025	Closure	White	13-415	Lined	STOCK
14-0269-026	Closure	Black	13-415	Lined	STOCK
14-0438-102	Closure	Black	13-425	Lined	STOCK
14-0438-119	Closure	White	13-425	Lined	STOCK
14-0170-020	Closure	Natural	15-400	None	STANDARD
14-0252-005	Closure	Black	15-415	Lined	STOCK
14-0252-006	Closure	White	15-415	Lined	STOCK
14-0252-007	Closure	Black	15-415	Lined	STOCK
14-0176-072	Closure	White	15-425	Lined	STOCK
14-0572-001	Combo Closure	White	15-425	Combo/Grey Butyl Stopper	STANDARD
14-0572-002	Combo Closure	Black	15-425	Combo/Grey Butyl Stopper	STANDARD
14-0248-005	Closure	Natural	15mm	Plug	STANDARD
14-0174-039	Closure	Black	18-400	Lined	STOCK
14-0177-003	Closure	Black	18-415	Linerless	STANDARD
14-0177-004	Closure	Natural	18-415	Linerless	STANDARD
14-0177-005	Closure	White	18-415	Linerless	STANDARD
14-0251-018	Closure	Black	18-415	Lined	STOCK
14-0251-019	Closure	White	18-415	Lined	STOCK
14-0441-004	Closure	Natural	18-415	Linerless	STANDARD
14-0441-005	Closure	White	18-415	Linerless	STANDARD
14-0441-006	Closure	Black	18-415	Linerless	STANDARD
14-0135-001	Combo Closure	Black	20-400	Combo/Grey Butyl Stopper	STANDARD
14-0135-002	Combo Closure	White	20-400	Combo/Grey Butyl Stopper	STANDARD
14-0214-003	Closure	Black	20-400	Lined	STOCK

Comar, Inc. One Comar Place, Buena, NJ 08310

Toll Free: 800-962-6627 • Fax: 856-692-9251 • Distributor Hotline: 866-807-9070

www.comar.com • information@comar.com

# Injection Molded Products



Part Number	Description	Capacity	Color	Style	Drop Size	Type
14-0142-010	Dispette Cover		White			STANDARD
14-0574-006	Dispette Cover		White			STANDARD
22-0020-002	Dosage Cup		Natural	With Ribs		STANDARD
22-0021-003	Dosage Cup		Natural	Without Ribs		STANDARD
22-0026-002	Dispette	5cc	Natural			STANDARD
22-0034-008	Dispette	1cc	Natural			STANDARD
22-0196-001	Controlled Drop Tip	13mm	Natural		60uL	STANDARD
22-0211-001	Controlled Drop Tip	13mm	Natural		50uL	STANDARD
22-0330-006	Dispette	2cc	Natural			STANDARD

Comar, Inc. One Comar Place, Buena, NJ 08310

Toll Free: 800-962-6627 • Fax: 856-692-9251 • Distributor Hotline: 866-807-9070

www.comar.com • information@comar.com

# Glass Dropper Assemblies



Part Number	Closure Color	Bulb Color	Pipette Style	Bulb Material	Closure Size	Pipette Diameter	Pipette Length	Type
19-0005-002	Black	Black	Straight	Rubber	13-425	5mm	46mm	STANDARD
19-0631-001	White	White	Straight	Rubber	13-425	5mm	46mm	STANDARD
19-0022-002	Black	Black	Straight	Rubber	15-425	5mm	66mm	STANDARD
19-0629-001	Black	Black	Straight	Rubber	15-425	5mm	61mm	STANDARD
19-0630-001	White	White	Straight	Rubber	15-425	5mm	66mm	STANDARD
19-0636-001	Black	Black	Straight	Rubber	15-415	5mm	38mm	STANDARD
19-0637-001	White	White	Straight	Rubber	15-415	5mm	38mm	STANDARD
19-0403-002	Black	Black	Straight	Rubber	18-400	7mm	66mm	STOCK
19-0542-001	White	White	Straight	Monprene	18-400	7mm	66mm	STOCK
19-0557-002	Black	Black	Straight	Monprene	18-400	7mm	66mm	STOCK
19-0632-001	Black	Black	Straight	Rubber	18-400	7mm	70mm	STANDARD
19-0633-001	White	White	Straight	Rubber	18-400	7mm	70mm	STANDARD
19-0062-003	Black	Black	Straight	Rubber	20-400	7mm	75mm	STOCK
19-0067-003	Black	Black	Straight	Rubber	20-400	7mm	89mm	STOCK
19-0469-002	White	White	Straight	Rubber	20-400	7mm	89mm	STANDARD
19-0523-002	White	White	Straight	Monprene	20-400	7mm	75mm	STOCK
19-0529-002	White	White	Straight	Rubber	20-400	7mm	75mm	STOCK
19-0763-001	Black	Black	Straight	Monprene	20-400	7mm	75mm	STOCK
19-0764-001	Black	Black	Straight	Monprene	20-400	7mm	89mm	STOCK
19-0560-001	Black	Black	Straight	Monprene	22-400	7mm	108mm	STOCK
19-0799-001	Black	Black	Straight	Rubber	22-400	7mm	108mm	STOCK
19-0634-001	Black	Black	Straight	Rubber	24-400	7mm	58mm	STANDARD
19-0635-001	Black	Black	Straight	Rubber	24-400	7mm	95mm	STANDARD
19-0020-002	White	White	Straight	Rubber	115-415	5mm	61mm	STANDARD

Comar, Inc. One Comar Place, Buena, NJ 08310

Toll Free: 800-962-6627 • Fax: 856-692-9251 • Distributor Hotline: 866-807-9070

www.comar.com • information@comar.com

# Media Tubes



Part Number	Neck Finish	Base of Neck Fill	Overflow	Vial Diameter	Vial Length	Glass Type	Style	Type
24-0024-003	13-415	7mL	8.59mL	13mm	100mm	Type I Clear	Round Bottom	STANDARD
24-0034-001	13-415	8mL	8.90mL	13mm	100mm	Type I Clear	Flat Bottom	STOCK
24-0004-003	15-415	12mL	13.36mL	16mm	100mm	Type I Clear	Round Bottom	STOCK
24-0025-001	15-415	12mL	13.71mL	16mm	100mm	Type I Clear	Flat Bottom	STOCK
24-0023-001	15-415	15mL	17.10mL	16mm	125mm	Type I Clear	Flat Bottom	STOCK
24-0026-001	15-415	15mL	16.81mL	16mm	125mm	Type I Clear	Round Bottom	STOCK
24-0018-003	15-415	18mL	20.14mL	16mm	150mm	Type I Clear	Round Bottom	STANDARD
24-0035-001	15-415	18mL	20.64mL	16mm	150mm	Type I Clear	Flat Bottom	STANDARD
24-0027-003	18-415	23mL	26mL	20mm	125mm	Type I Clear	Round Bottom	STANDARD
24-0036-001	18-415	24mL	27.11mL	20mm	125mm	Type I Clear	Flat Bottom	STOCK
24-0028-001	18-415	29mL	32.88mL	20mm	150mm	Type I Clear	Flat Bottom	STANDARD
24-0031-003	18-415	29mL	32.66mL	20mm	150mm	Type I Clear	Round Bottom	STANDARD
24-0070-001	18-415	30mL	33.38mL	20mm	150mm	Type I Clear	Flat Bottom	STANDARD

Comar, Inc. One Comar Place, Buena, NJ 08310

Toll Free: 800-962-6627 • Fax: 856-692-9251 • Distributor Hotline: 866-807-9070

www.comar.com • information@comar.com

# Screw Thread Vials



Part Number	Neck Finish	Base of Neck Fill	Overflow	Vial Diameter	Vial Length	Glass Type	Type
24-0120-001	8-425	1.8mL	2.07mL	12mm	32mm	Type I Clear	STANDARD
24-0124-008	8-425	1.8mL	2mL	12mm	32mm	Amber	STANDARD
24-0249-001	13-425	1.8mL	2.11mL	16mm	50mm	Type I Clear	STANDARD
24-0046-001	13-425	2.3mL	2.56mL	15mm	26.5mm	Amber	STANDARD
24-0050-001	13-425	2.3mL	2.61mL	15mm	26.5mm	Type I Clear	STANDARD
24-0016-001	13-425	4.4mL	4.97mL	15mm	45mm	Type I Clear	STOCK
24-0055-001	13-425	4.3mL	4.85mL	15mm	45mm	Amber	STOCK
24-0066-001	13-425	4.1mL	4.65mL	17mm	35mm	Type I Clear	STANDARD
24-0012-001	15-415	8mL	9.73mL	16mm	75mm	Type I Clear	STANDARD
24-0048-002	15-425	8mL	9.01mL	17mm	60mm	Type I Clear	STOCK
24-0080-001	15-425	8.mL	8.96mL	17mm	60mm	Amber	STANDARD
24-0231-001	15-425	10mL	12.10mL	19mm	65mm	Amber	STANDARD
24-0054-004	15-425	11mL	12.25mL	19mm	65mm	Type I Clear	STANDARD
24-0032-001	18-400	14mL	16.63mL	21mm	70mm	Type I Clear	STOCK
24-0194-001	18-400	14mL	16.43mL	21mm	70mm	Amber	STANDARD
24-0006-002	24-400	22mL	25.37mL	27mm	57.5mm	Type I Clear	STOCK

Comar, Inc. One Comar Place, Buena, NJ 08310

Toll Free: 800-962-6627 • Fax: 856-692-9251 • Distributor Hotline: 866-807-9070

www.comar.com • information@comar.com



# Serum Vials



Part Number	Neck Finish	Base of Neck Fill	Overflow	Vial Diameter	Vial Length	Glass Type	Type
24-0123-001	11mm	1.7mL	1.91mL	12mm	32mm	Amber	STANDARD
24-0125-001	11mm	1.8mL	2.03mL	12mm	32mm	33 Exp. Clear	STANDARD
24-0085-002	13mm	2.8mL	3.16mL	15mm	32mm	Type I Clear	STOCK
24-0658-001	13mm	2.8mL	3.18mL	15mm	32mm	Amber	STANDARD
24-0063-001	13mm	4.1mL	4.65mL	17mm	37mm	Amber	STOCK
24-0365-002	13mm	4.4mL	4.98mL	17mm	37.7mm	Type I Clear	STOCK
24-0038-001	20mm	8mL	9.74mL	22mm	40mm	Type I Clear	STANDARD
24-0002-001	20mm	9mL	9.97mL	22mm	40mm	Amber	STOCK
24-0003-001	20mm	13mL	15.30mL	24mm	50mm	Type I Clear	STOCK
24-0007-001	20mm	13mL	15.07mL	24mm	50mm	Amber	STANDARD
24-0103-002	20mm	25mL	28.23mL	29mm	62mm	Amber	STANDARD
24-0289-001	20mm	25mL	27.79mL	29mm	62mm	Type I Clear	STANDARD

Comar, Inc. One Comar Place, Buena, NJ 08310

Toll Free: 800-962-6627 • Fax: 856-692-9251 • Distributor Hotline: 866-807-9070

www.comar.com • information@comar.com

# Shell Vials



Part Number	Neck Finish	Base of Neck Fill	Overflow	Vial Diameter	Vial Length	Glass Type	Type
24-0466-006	None	5mL	6.04mL	15mm	45mm	Type I Clear	STANDARD
24-0558-001	None	9mL	10.75mL	17mm	60mm	Type I Clear	STANDARD
24-0803-001	None	48mL	53.82mL	30mm	93.1mm	Type I Clear	STANDARD

Comar, Inc. One Comar Place, Buena, NJ 08310

Toll Free: 800-962-6627 • Fax: 856-692-9251 • Distributor Hotline: 866-807-9070

www.comar.com • information@comar.com

# Sure Thread Vials



Part Number	Neck Finish	Base of Neck Fill	Overflow	Vial Diameter	Vial Length	Glass Type	Type
24-0259-001	18-400	9mL	10.86mL	24mm	40mm	33 Exp. Clear	STANDARD
24-0092-001	18-400	13mL	15.30mL	22mm	60mm	Type I Clear	STANDARD
24-0258-001	18-400	14mL	15.77mL	24mm	54mm	33 Exp. Clear	STANDARD

Comar, Inc. One Comar Place, Buena, NJ 08310

Toll Free: 800-962-6627 • Fax: 856-692-9251 • Distributor Hotline: 866-807-9070

www.comar.com • information@comar.com

# EPA/ Water Vials



Part Number	Neck Finish	Base of Neck Fill	Overflow	Vial Diameter	Vial Length	Glass Type	Type
24-0163-004	24-400	39mL	43.75mL	27mm	95mm	Type I Clear	STOCK
24-0165-001	24-400	60mL	66.93mL	27mm	140mm	33 Exp. Clear	STOCK