

Qorpak®



value
service
integrity

Laboratory Container Solutions



Qorpak® Delivers The Containers You Want – When You Want Them – How You Need Them



Qorpak Cartons Carry Convenience to Your Door.

- Full height partitions and top pads minimize breakage.
- Pressure Sensitive Labels are included for fast and easy product identification.
- Cartons are bar-coded with a catalog number and product description to streamline re-ordering.
- Bottles are packed in small quantities... the way you actually use them.
- Packaging is lightweight, easy to handle and store.

Vials	3
Glass	4-7
Safety Coated Glass	8-9
Plastic	10-14
Caps & Liners	15-17
Cleaning Services	18-21
Cans, Paper, Weigh Dishes & Bags	22
UN Compliant Shippers	23



We Are Committed to Helping You.

Our friendly Customer Service Representatives are available to help you with all of your container needs. They can also assist you with product information, and help you locate non-stock containers.



In Stock.

Most Qorpak products are readily available. We can drop-ship directly to you.



Free Samples: Try It Before You Buy.

Don't try to guess if a container will work. Free samples available on most products.

Qorpak is an ISO Compliant Company.

1-800-922-7558

Storage & Sampling Vials



CLOSURES SOLD SEPARATELY.
R.A.M.™ VIAL IS A TRADEMARK OF
J. G. FINNERAN ASSOCIATES.

LARGE OPENING R.A.M.™ VIALS

TYPE I BOROSILICATE GLASS, 2.0 ML VOLUME,
9 MM NECK FINISH, 12 X 32 MM

	CLEAR	AMBER
Plain	GLD-09527	GLD-09528
White Marking Spot	GLD-09527WM	GLD-09528WM
White Graduation Spot	GLD-09527WG	GLD-09528WG
Blue Graduation Spot	GLD-09527BG	GLD-09528BG
Green Graduation Spot	GLD-09527GG	GLD-09528GG
Rust Graduation Spot	GLD-09527RG	GLD-09528RG
Yellow Graduation Spot	GLD-09527YG	GLD-09528YG
Wt. Case (lbs.)	8	8
No. Per Case	1000	1000



The R.A.M.™ vial features a 40% larger opening and is designed to work in robotic arm samplers. Incorporates the unique StepVial design that precisely centers a limited volume insert in the vial neck. Available with graduated marking spots which can be written on with a # 2 pencil.

SCREW THREAD CLOSURES/ SEPTAS FOR R.A.M.™ VIALS

POLYPROPYLENE HOLE CAPS (9MM)	
Black	CPD-09529
Blue	CPD-09530
Green	CPD-09531
Natural	CPD-09532
Red	CPD-09533
Yellow	CPD-09534
Wt. Case (lbs.)	2
No. Per Case	1000

POLYPROPYLENE CAPS (9MM)	
Black, unlined	CPD-09535
Black w/F217 liner	CPD-06930
Natural w/molded 10mil Septa	CPD-09536
Wt. Case (lbs.)	2
No. Per Case	1000

SEPTAS (9MM)	
White PTFE Disk	0.010" CPD-09537
PTFE/Red Rubber Septa	0.040" CPD-09538
PTFE/Silicone Septa	0.040" CPD-09539
Red PTFE/Silicone	0.040" CPD-09540
Wt. Case (lbs.)	2
No. Per Case	1000



AMBER BOROSILICATE GLASS THREADED VIALS

Capacity (ml)	1/2dr (1.875)	1dr (3.75)	2dr (7.5)	4dr (15)
Neck Size (MM)	8-425	13-425	15-425	18-400

CAPS ATTACHED				
SOLID GREEN TEFLON® caps	GLC-06786	GLC-06787	GLC-06788	GLC-06789
Wt. Case (lbs.)	2	2	3	4
No. Per Case	144	144	144	144

GLASS ONLY				
Cat. No.	2524	2525	2522	2526
Wt. Case (lbs.)	2	2	3	4
No. Per Case	144	144	144	144



BOROSILICATE COMPOUND VIALS

Capacity	1/4 dr	1/3 dr	1/2 dr	1/4 dr	1/3 dr	1/2 dr
Neck Size (MM)	13-425	13-425	13-425	13-425	13-425	13-425
Dimensions	14.65 x 19	14.65 x 22	14.65 x 26	14.75 x 19	14.75 x 22	14.75 x 26

	CLEAR			AMBER		
Cat. No.	GLA-07794	GLA-07795	GLA-07796	GLA-07797	GLA-07798	GLA-07799
Wt. Case (lbs.)	2	2	2	2	2	2
No. Per Case	144	144	144	144	144	144



CLEAR BOROSILICATE GLASS THREADED VIALS

Capacity (ml)	1/2dr (1.875)	1dr (3.75)	2dr (7.5)	3dr (11.25)	4dr (15)	6dr (22.5)	8dr (30)	5dr (20)	6dr (22.5)
Neck Size (MM)	8-425	13-425	15-425	15-425	18-400	20-400	22-400	24-400	24-400

CAPS ATTACHED - CONVENIENCE PACK									
TEFLON® PTFE LINED caps	2501T	2502T	2503T	2504T	2505T	2506T	2507T	2508T	2509T
PULP/VINYL LINED caps	2501Q	2502Q	2503Q	2504Q	2505Q	2506Q	2507Q	2508Q	2509Q
POLYSEAL caps	—	2502B	2503B	2504B	2505B	2506B	2507B	2508B	2509B
14B WHITE RUBBER caps	—	2502D	2503D	2504D	2505D	2506D	2507D	2508D	2509D
Wt. Case (lbs.)	2	3	4	4	6	6	9	3	4
No. Per Case	144	144	144	144	144	144	144	72	72

GLASS ONLY									
Cat. No.	2501	2502	2503	2504	2505	2506	2507	2508	2509
Wt. Case (lbs.)	2	3	4	4	6	6	9	3	4
No. Per Case	144	144	144	144	144	144	144	72	72



20, 40 & 60 ML BOROSILICATE GLASS THREADED VIALS

	CLEAR			NEW!		AMBER		NEW!
Capacity (ml)	2/3 oz (20)	1 1/3oz (40)	1 1/3oz (40)	2 oz (60)	2/3 oz (20)	1 1/3oz (40)	2 oz (60)	
Neck Size (MM)	24-400	24-400	24-400	24-410	24-400	24-400	24-410	

	TALL		SQUAT		TALL		TALL	
TEFLON® PTFE SEPTA caps	2V20MHSC	2517T/144	2515T/144	GLC-07863	2V22MHSC	2518T/144	GLC-07867	
SOLID GREEN TEFLON® caps	GLC-06064	7588T/144	7587T/144	GLC-07876	2V22MGTC	7589T/144	GLC-07881	
PULP/VINYL LINED caps	—	2517Q	2515Q	GLC-07877	2516Q	2518Q	GLC-07882	
POLYSEAL caps	—	2517B	2515B	GLC-07878	2516B	2518B	GLC-07883	
14B WHITE RUBBER caps	—	2517D	2515D	GLC-07879	2516D	2518D	GLC-07884	
BLACK P/P - .015 TEFLON disc	2V20QEDC	2V40QEDC	2V41QEDC	GLC-07880	2V22QEDC	2V42QEDC	GLC-07885	
Wt. Case (lbs.)	10	11	11	8	10	11	8	
No. Per Case	144	144	144	72	144	144	72	

CAPS ATTACHED - BULK PACK								
TEFLON® PTFE SEPTA caps	—	2517T	2515T	—	—	2518T	—	—
Wt. Case (lbs.)	—	40	40	—	—	40	—	—
No. Per Case	—	576	576	—	—	576	—	—

GREEN P/P CLOSURE / .125 TEFLON®
SILICONE SEPTUM COMMONLY USED FOR
VOA ANALYSIS. ALSO AVAILABLE
PRE-CLEANED, SEE PAGE 18.



CLEAR STRAIGHT SIDED ROUNDS

Capacity (ml)	1oz (30)	2oz (60)	4oz (120)	6oz (180)	9oz (270)	16oz (480)	32oz (960)
Neck Size (MM)	43-400	53-400	58-400	63-400	70-400	89-400	89-400

CAPS ATTACHED - CONVENIENCE PACK

TEFLON® PTFE LINED caps ●	7939	7940	7943	7945	7946	7947	7950
PULP/VINYL LINED caps	7739	7740	7743	7745	7746	7747	7750
14B WHITE RUBBER caps	—	7740D	7743D	—	7746D	—	—
Wt. Case (lbs.)	7	6	8	10	13	10	15
No. Per Case	48	24	24	24	24	12	12

CAPS ATTACHED - BULK CASE

PULP/VINYL LINED caps	7529	7530	7531	7553	7533	7534	7535
Wt. Case (lbs.)	46	29	8	10	13	10	15
No. Per Case	384	144	24	24	24	12	12

GLASS ONLY - BULK CASE

Cat. No.	2039	2040	2043	2045	2046	2047	2050
Wt. Case (lbs.)	45	28	7	9	12	9	14
No. Per Case	384	144	24	24	24	12	12



AMBER STRAIGHT SIDED ROUNDS

Capacity (ml)	4oz (120)	8oz (240)	32oz (960)
Neck Size (MM)	58-400	70-400	89-400

CAPS ATTACHED

SOLID GREEN TEFLON® caps ●	GLC-06777	GLC-06778	GLC-06779
-----------------------------------	------------------	------------------	------------------

GLASS ONLY

No Closure	GLA-06759	GLA-06760	GLA-06761
Wt. Case (lbs.)	7	6	14
No. Per Case	24	12	12



AMBER BOSTON ROUNDS

Capacity (ml)	1/2oz (15)	1oz (30)	2oz (60)	4oz (120)	8oz (240)	16oz (480)	32oz (960)	33oz (1000)
Neck Size (MM)	18-400	20-400	20-400	22-400	24-400	28-400	33-400	33-430

CAPS ATTACHED - CONVENIENCE PACK

TEFLON® PTFE LINED caps ●	—	—	7918	7919	7920	7922	7923	7724T*
PULP/VINYL LINED caps	—	—	7718	7719	7720	7722	7723	—
POLYSEAL caps	—	—	7718B	7719B	7720B	7722B	7723B	—
14B WHITE RUBBER caps	—	—	7718D	7719D	7720D	7722D	7723D	—
Wt. Case (lbs.)	—	—	4	6	10	8	14	15
No. Per Case	—	—	24	24	24	12	12	12

CAPS ATTACHED - BULK CASE

PULP/VINYL LINED caps	7558	7559	7515	7516	7517	7518	7519	—
POLYSEAL caps	7558B	7559B	7515B	7516B	7517B	7518B	7519B	—
Wt. Case (lbs.)	46	38	39	38	40	39	33	—
No. Per Case	720	432	288	160	96	60	30	—

GLASS ONLY - BULK CASE

Cat. No.	2025	2021	2018	2019	2020	2022	2023	2024
Wt. Case (lbs.)	45	37	38	37	39	38	32	15
No. Per Case	720	432	288	160	96	60	30	12

*BLACK PHENOLIC TEFLON LINED.



AMBER WIDE MOUTH PACKERS

Capacity (ml)	1/5 oz (6)	1/4 oz (8.5)	1/2 oz (15)	1oz (30)	2oz (60)	4oz (120)	8oz (250)	10oz (300)
Neck Size (MM)	20-400	20-400	28-400	28-400	33-400	38-400	45-400	53-400

CAPS ATTACHED - CONVENIENCE PACK

TEFLON® PTFE LINED caps ●	—	—	—	7930	7931	7933	7934	7932
PULP/VINYL LINED caps	—	—	—	7730	7731	7733	7734	7732
POLYSEAL caps	—	—	—	7730B	7731B	7733B	—	—
14B WHITE RUBBER caps	—	—	—	7730D	7731D	7733D	—	7732D
Wt. Case (lbs.)	—	—	—	4	5	7	10	54
No. Per Case	—	—	—	24	24	24	24	96

CAPS ATTACHED - BULK CASE

PULP/VINYL LINED caps	7561	7563	7562	7522	7523	7525	7526	7532
Wt. Case (lbs.)	38	52	47	55	38	47	36	54
No. Per Case	720	912	624	432	216	180	84	96

GLASS ONLY - BULK CASE

Cat. No.	2091	2098	2092	2030	2031	2033	2034	2095
Wt. Case (lbs.)	31	46	41	53	36	46	35	53
No. Per Case	720	912	624	432	216	180	84	96

Ultimate Long-Term Storage Container



CLEAR BOSTON ROUNDS

Capacity (ml)	1oz (30)	2oz (60)	4oz (120)	8oz (240)	16oz (480)	32oz (960)
Neck Size (MM)	20-400	20-400	22-400	24-400	28-400	33-400

CAPS ATTACHED - CONVENIENCE PACK

TEFLON® PTFE LINED caps ●	7909	7910	7912	7913	7914	7917
PULP/VINYL LINED caps	7709	7710	7712	7713	7714	7717
POLYSEAL caps	7709B	7710B	7712B	7713B	7714B	7717B
14B WHITE RUBBER caps	7709D	7710D	7712D	7713D	7714D	7717D
Wt. Case (lbs.)	5	4	6	10	8	14
No. Per Case	48	24	24	24	12	12

CAPS ATTACHED - BULK CASE

PULP/VINYL LINED caps	7508	7509	7511	7512	7513	7514
POLYSEAL caps	7508B	7509B	7511B	7512B	7513B	7514B
Wt. Case (lbs.)	38	39	38	40	39	15
No. Per Case	432	288	108	96	60	12

GLASS ONLY - BULK CASE

Cat. No.	2009	2010	2012	2013	2015	2017
Wt. Case (lbs.)	37	38	37	39	38	14
No. Per Case	432	288	108	96	60	12



CLEAR STANDARD WIDE MOUTHS

Capacity (ml)	4oz (120)	8oz (240)	16oz (480)	32oz (960)	64oz (1920)	130oz (3810)
Neck Size (MM)	48-400	58-400	63-400	70-400	83-400	89-400

CAPS ATTACHED - CONVENIENCE PACK

TEFLON® PTFE LINED caps ●	7967	7968	7970	7972	7973	7974
PULP/VINYL LINED caps	7767	7768	7770	7772	7773	7774
Wt. Case (lbs.)	6	10	14	13	10	11
No. Per Case	24	24	24	12	6	4

GLASS ONLY - BULK CASE

Cat. No.	2067	2068	2070	2072	2073	2074
Wt. Case (lbs.)	6	9	13	12	10	11
No. Per Case	24	24	24	12	6	4



CLEAR TALL STRAIGHT SIDED JARS

Capacity (ml)	8oz (240)	16oz (480)
Neck Size (MM)	58-400	63-400

CAPS ATTACHED - CONVENIENCE PACK

TEFLON® PTFE LINED caps ●	7578T	GLC-05728
PULP/VINYL LINED caps	7578	7580
14B WHITE RUBBER caps	GLC-05721	—
Wt. Case (lbs.)	14	10
No. Per Case	24	12

GLASS ONLY - CONVENIENCE PACK

Cat. No.	2078	2080
Wt. Case (lbs.)	13	9
No. Per Case	24	12

16oz (480) 53-400	32oz (950) 53-400	42oz (1250) 70-400	84oz (2500) 70-400
----------------------	----------------------	-----------------------	-----------------------

7935	7938	7937	7936
7735	7738	7737	7736
—	—	—	—
7735D	7738D	7737D	7736D
11	12	8	28
12	12	6	12

7527	7528	7560	7536
38	36	8	28
60	36	6	12

2035	2038	2029	2036
37	35	6	27
60	36	6	12

Amber Glass Provides
Protection from Harmful
Light and UV Rays.

TEFLON® IS A REGISTERED TRADEMARK OF
E.I. DUPONT DE NEMOURS AND COMPANY.
● THE COLOR GREEN FOR CAPS IS A TRADEMARK OF QORPAK,
REGISTERED IN THE U.S. PATENT AND TRADEMARK OFFICE.

1-800-922-7558



GRADUATED

QORPAK BOTTLE BEAKERS® - CLEAR MEDIUM ROUNDS

Capacity (ml)	1oz (30)	2oz (60)	4oz (120)	8oz (240)	16oz (480)
Neck Size (MM)	33-400	38-400	48-400	58-400	70-400

CAPS ATTACHED - CONVENIENCE PACK

TEFLON® PTFE LINED caps ●	7981	7982	7983	7984	7985
PULP/VINYL LINED caps	7781	7782	7783	7784	7785
TINFOIL LINED caps	7781A	7782A	7783A	7784A	7785A
POLYSEAL caps	7781B	7782B	—	—	—
14B WHITE RUBBER caps	7781D	7782D	7783D	7784D	7785D
Wt. Case (lbs.)	7	10	8	12	18
No. Per Case	48	48	24	24	24

CAPS ATTACHED - BULK CASE

PULP/VINYL LINED caps	7547	7548	7549	7550	7551
POLYSEAL caps	7547B	7548B	—	—	—
Wt. Case (lbs.)	55	54	44	43	33
No. Per Case	432	288	144	96	48

GLASS ONLY - BULK CASE

Cat. No.	2081	2082	2083	2084	2085
Wt. Case (lbs.)	54	53	43	42	33
No. Per Case	432	288	144	96	48



SPACE SAVING DESIGN
32 OZ. GRADUATED IN MILLILITERS AND OUNCES.

WIDE MOUTH FRENCH SQUARES

Capacity (ml)	1/2oz (15)	1oz (30)	2oz (60)	4oz (120)	8oz (240)	16oz (480)	32oz (960)
Neck Size (MM)	20-400	24-400	28-400	33-400	43-400	48-400	58-400

CAPS ATTACHED - CONVENIENCE PACK

TEFLON® PTFE LINED caps ●	7901	7902	7903	7904	7905	7906	7907
PULP/VINYL LINED caps	7701	7702	7703	7704	7705	7706	7707
POLYSEAL caps	7701B	7702B	7703B	7704B	—	—	—
14B WHITE RUBBER caps	7701D	7702D	7703D	7704D	7705D	7706D	7707D
Wt. Case (lbs.)	4	6	10	8	13	18	14
No. Per Case	48	48	48	24	24	24	12

CAPS ATTACHED - BULK CASE

PULP/VINYL LINED caps	7500	7501	7502	7504	7505	7506	—
Wt. Case (lbs.)	40	31	46	32	45	29	—
No. Per Case	576	280	240	96	84	40	—

GLASS ONLY - BULK CASE

Cat. No.	2001	2002	2003	2004	2005	2006	2007
Wt. Case (lbs.)	39	30	45	31	44	28	26
No. Per Case	576	280	240	96	84	40	24



CLEAR COMPOSITE TEST JARS

Capacity (ml)	2oz (60)
Neck Size (MM)	43-400

CAPS ATTACHED - CONVENIENCE PACK

TEFLON® PTFE LINED caps ●	7537T
PULP/VINYL LINED caps	7537
Wt. Case (lbs.)	39
No. Per Case	180

GLASS ONLY - CONVENIENCE PACK

Cat. No.	2057
Wt. Case (lbs.)	35
No. Per Case	180



SPECIMEN MAILERS

COMPLETE WITH FIBER MAILER AND METAL SCREW THREAD LID

Capacity (ml)	1oz (30)	2oz (60)	4oz (120)	8oz (240)
----------------------	----------	----------	-----------	-----------

CAPS ATTACHED

GLASS NM BOSTON ROUNDS	2300	2301	2302	2303
GLASS WM ROUNDS	2304	2305	2306	2307
HDPE NM CYLINDERS	—	3300	3301	3302
HDPE WM ROUNDS	3303	3304	3305	—
Specimen Mailer Only	2910	2911	2912	2913
No. Per Case	24	24	24	24



SQUARE TABLET BOTTLES

Capacity (ml)	1oz (30)	2oz (60)
Neck Size (MM)	33-400	33-400

CAPS ATTACHED - CONVENIENCE PACK		
TEFLON® PTFE LINED caps ●	7520T	7521T
PULP/VINYL LINED caps	7520	7521
POLYSEAL caps	7520B	GLC-04869
Wt. Case (lbs.)	40	29
No. Per Case	288	144

GLASS ONLY - CONVENIENCE PACK		
Cat. No.	2026	2027
Wt. Case (lbs.)	36	26
No. Per Case	288	144



CLEAR OIL SAMPLE BOTTLE

Capacity (ml)	4oz (120)
Cork Size	# 6 Cork

CAPS ATTACHED	
Cat. No.	2090
Wt. Case (lbs.)	35
No. Per Case	144



AMBER OVAL DROPPER BOTTLES

Capacity (ml)	1oz (30)	2oz (60)
Neck Size (MM)	20-400	20-400

CAPS ATTACHED - CONVENIENCE PACK		
PLASTIC DROPPER ASSEMBLY	2088	2089
Wt. Case (lbs.)	10	7
No. Per Case	72	36



AMBER SQUARE WIDE MOUTHS

Capacity (ml)	8oz (240)
Neck Size (MM)	45-400

WHITE PP CAPS ATTACHED	
F217 P/E FOAM caps	2538
Wt. Case (lbs.)	23
No. Per Case	48

SPACE SAVING DESIGN



ROUND DROPPER BOTTLES

Capacity (ml)	1oz (30)	2oz (60)	4oz (120)	2oz (60)	4oz (120)
Neck Size (MM)	20-400	20-400	22-400	20-400	22-400

CAPS ATTACHED - CONV. PACK	CLEAR			AMBER	
	GLASS DROPPER ASSEMBLY	GLC-05723	GLC-05724	GLC-05726	GLC-05725
Wt. Case (lbs.)	6	5	7	5	7
No. Per Case	48	24	24	24	24



CLEAR JUGS

Capacity (ml)	64oz (1920)	128oz (3840)
Neck Size (MM)	38-400	38-400

CAPS ATTACHED		
TEFLON® PTFE LINED caps ●	7963	7964
PULP/VINYL LINED caps	7763	7764
POLYSEAL caps	7763B	7764B
14B WHITE RUBBER caps	7763D	7764D
Wt. Case (lbs.)	14	13
No. Per Case	6	4

GLASS ONLY		
Cat. No.	2063	2064
Wt. Case (lbs.)	14	13
No. Per Case	6	4



AMBER JUGS

Capacity (ml)	80oz (2400)	128oz (3785)
Neck Size (MM)	38-430	38-400

CAPS ATTACHED		
TEFLON® PTFE LINED caps ●	7562TA*	7965
PULP/VINYL LINED caps	—	7765
POLYSEAL caps	—	7765B
14B WHITE RUBBER caps	—	7765D
Wt. Case (lbs.)	16	13
No. Per Case	6	4

GLASS ONLY		
Cat. No.	2062	2065
Wt. Case (lbs.)	16	13
No. Per Case	6	4

* BLACK PHENOLIC TEFLON LINED CLOSURE

The Glass Barrier.

Chemicals, oils and gases cannot penetrate glass. Which means your samples and stored materials are safe from permeation.



The Coating that Prevents Accidents and Provides Peace of Mind.

- Prevent lab mishaps. Bottles are slip resistant — which gives technicians improved handling stability.
- Reduce breakage. Impact-resistant coating means fewer broken bottles in the lab, or in transit.

- Quick clean-up. If a bottle should break, the clear plastisol coating contains the acids, solvents, alcohols, surfactants, esters, ethers, and glass particles long enough to allow proper disposal.
- Helps ensure compliance with OSHA, D.O.T., and UPS requirements.
- Protects employees. Safety glass removes some of the risk in dealing with certain chemicals and compounds.
- Protects against pressure build up. Overfilling or atmo-

- spheric changes can cause pressure within a container. The safety coating protects the container from being affected by pressure build up.
- Eliminates permeation. Glass is inert and provides complete protection against gas and product permeation.
- Reduces product liability.
- Saves you time. Safety glass containers are supplied with lined Thermoset caps already attached — which makes the container ready to use right out of the box.

SAFETY COATED



CLEAR BOSTON ROUNDS

Capacity (ml)	8oz (240)	16oz (480)	32oz (960)
Neck Size (MM)	24-400	28-400	33-400

CAPS ATTACHED - CONVENIENCE PACK

TEFLON® PTFE LINED caps ●	7713TW	7714TW	7717TW
PULP/VINYL LINED caps	7713W	7714W	7717W
POLYSEAL caps	7713BW	7714BW	7717BW
14B WHITE RUBBER caps	7713DW	7714DW	7717DW
Wt. Case (lbs.)	14	9	17
No. Per Case	24	12	12

GLASS ONLY - BULK CASE

Cat. No.	2013W	2015W	2017W
Wt. Case (lbs.)	43	48	40
No. Per Case	96	60	30

We Offer Additional Shapes & Sizes With 1,000 Piece Minimum Order. Call Us for Details.

SAFETY COATED



* BLACK PHENOLIC TEFLON LINED CLOSURE

AMBER BOSTON ROUNDS

Capacity (ml)	4oz (120)	8oz (240)	16oz (480)	32oz (960)	33oz (1000)
Neck Size (MM)	22-400	24-400	28-400	33-400	33-430*

CAPS ATTACHED - CONVENIENCE PACK

TEFLON® PTFE LINED caps ●	7719TW	7720TW	7722TW	7723TW	7724TW
PULP/VINYL LINED caps	7719W	7720W	7722W	7723W	—
POLYSEAL caps	7719BW	7720BW	7722BW	7723BW	—
14B WHITE RUBBER caps	7719DW	7720DW	7722DW	7723DW	—
Wt. Case (lbs.)	8	13	10	17	18
No. Per Case	24	24	12	12	12

GLASS ONLY - BULK CASE

Cat. No.	2019W	2020W	2022W	2023W	2024W
Wt. Case (lbs.)	35	50	50	40	43
No. Per Case	126	108	60	30	30

SAFETY COATED



AMBER WIDE MOUTH PACKERS

Capacity (ml)	8oz (250)	16oz (500)	32oz (950)	83oz (2500)
Neck Size (MM)	45-400	53-400	53-400	70-400

CAPS ATTACHED - CONVENIENCE PACK

TEFLON® PTFE LINED caps ●	7526TW	7527TW	7528TW	4P83QGTC
PULP/VINYL LINED caps	7526W	7527W	7528W	4P83QPVC
Wt. Case (lbs.)	14	12	15	10
No. Per Case	24	12	12	4

GLASS ONLY - BULK CASE

Cat. No.	2034W	2035W	2038W	2036W
Wt. Case (lbs.)	13	12	14	30
No. Per Case	24	12	12	12

TEFLON® IS A REGISTERED TRADEMARK OF E.I. DUPONT DE NEMOURS AND COMPANY. ● THE COLOR GREEN FOR CAPS IS A TRADEMARK OF QORPAK, REGISTERED IN THE U.S. PATENT AND TRADEMARK OFFICE.

Helps to Prevent Accidents

SAFETY COATED



CLEAR WIDE MOUTHS

Capacity (ml)	8oz (270)	16oz (480)	32oz (960)
Neck Size (MM)	70-400	70-400	70-400

CAPS ATTACHED - CONVENIENCE PACK

TEFLON® PTFE LINED caps ●	7533TW	7785TW	7772TW
PULP/VINYL LINED caps	7533W	7785W	7772W
Wt. Case (lbs.)	15	21	14
No. Per Case	24	24	12

GLASS ONLY - BULK CASE

Cat. No.	2046W	2085W	2072W
Wt. Case (lbs.)	13	31	33
No. Per Case	24	40	30

SAFETY COATED



CLEAR WIDE MOUTHS- LARGER CAPACITY

Capacity (ml)	64oz (1920)	128oz (3840)
Neck Size (MM)	83-400	89-400

CAPS ATTACHED - CONVENIENCE PACK

TEFLON® PTFE LINED caps ●	GLC-05729	GLC-05731
PULP/VINYL LINED caps	GLC-05730	GLC-05732
Wt. Case (lbs.)	12	13
No. Per Case	6	4

GLASS ONLY - BULK CASE

Cat. No.	GLA-00972	GLA-05580
Wt. Case (lbs.)	11	12
No. Per Case	6	4

SAFETY COATED



CLEAR ACID JUGS

Capacity (ml)	80oz (2400)
Neck Size (MM)	38-439

CAPS ATTACHED - CONVENIENCE PACK

STURDEE SEAL® F-422 LINED caps	7761GW
Wt. Case (lbs.)	18
No. Per Case	6

GLASS ONLY - BULK CASE

Cat. No.	2061W
Wt. Case (lbs.)	34
No. Per Case	12

SAFETY COATED



AMBER JUGS

Capacity (ml)	133oz (3990)
Neck Size (MM)	38-430

CAPS ATTACHED - CONVENIENCE PACK

TEFLON® PTFE LINED caps*	7765TW
Wt. Case (lbs.)	16
No. Per Case	4

GLASS ONLY - BULK CASE

Cat. No.	2065W
Wt. Case (lbs.)	23
No. Per Case	6

* BLACK PHENOLIC TEFLON LINED CLOSURE



Get Your
FREE Qorpak
Measurer.

**HANDY
MEASURER
MAKES
CAP
ORDERING
A SNAP.**

This two-sided card helps you measure cap sizes and bottle neck finishes with confidence. It's also loaded with lots of useful information like a conversion chart for commonly used measurements, a benefits list for glass and plastic to help you to decide which container to order, a six-inch ruler, and answers to a number of commonly asked ordering questions.

Call customer service at **1-800-922-7558** to request your measurer today!

1-800-922-7558



NARROW MOUTH COMMERCIAL CYLINDERS

WITH WHITE PP STURDEESEAL P/E FOAM LINED CAPS ATTACHED

Capacity (ml)	2oz (60)	4oz (120)	8oz (250)	12oz (360)	16oz (500)	32oz (960)
Neck Size (MM)	24-410	24-410	24-410	38-400	24-410	28-400

CAPS ATTACHED - CONVENIENCE PACK

Cat. No.	7336G	7337G	7338G	7330G	7339G	7333G
Wt. Case (lbs.)	2	3	5	3	4	3
No. Per Pack	12	12	12	—	6	—
No. Per Case	48	48	48	24	24	12

CAPS ATTACHED - BULK CASE

Cat. No.	7136G	7137G	7138G	7130G	7139G	7133G
-----------------	--------------	--------------	--------------	--------------	--------------	--------------

BOTTLES ONLY - BULK CASE

Cat. No.	3036	3037	3038	3030	3039	3033
Wt. Case (lbs.)	20	17	14	13	15	12
No. Per Case	800	425	230	216	180	77



NARROW MOUTH CLEAR CYLINDER ROUNDS

WITH WHITE PP STURDEESEAL P/E FOAM LINED CAPS ATTACHED

Capacity (ml)	4oz (120)	8oz (240)	16oz (500)	32oz (960)
Neck Size (MM)	20-410	24-410	28-410	28-410

CAPS ATTACHED - BULK CASE

Cat. No.	PLC-05771	PLC-05772	PLC-05773	PLC-05774
-----------------	------------------	------------------	------------------	------------------

BOTTLES ONLY - BULK CASE

Cat. No.	PLA-05756	PLA-05757	PLA-05758	PLA-05759
Wt. Case (lbs.)	26	31	18	11
No. Per Case	559	489	243	128



NARROW MOUTH WHITE CYLINDER ROUNDS

WHITE PP STURDEESEAL P/E FOAM LINED CAPS

Capacity (ml)	2oz (60)	6oz (180)	8oz (240)	12oz (360)	16oz (500)	32oz (960)
Neck Size (MM)	20-410	24-410	24-410	24-210	28-400	28-410

CAPS ATTACHED - BULK CASE

Cat. No.	PLC-05764	PLC-05766	PLC-05767	PLC-05768	PLC-05769	PLC-05770
-----------------	------------------	------------------	------------------	------------------	------------------	------------------

BOTTLES ONLY - BULK CASE

Cat. No.	PLA-05749	PLA-05751	PLA-05752	PLA-05753	PLA-05754	PLA-05755
Wt. Case (lbs.)	22	18	14	19	14	11
No. Per Case	1000	300	230	220	162	77



NARROW MOUTH BOSTON ROUNDS

WITH BLACK PHENOLIC POLYSEAL LINED CAPS ATTACHED

Capacity (ml)	1/2oz (15)	1oz (30)	2oz (60)	4oz (120)	8oz (250)	16oz (500)	32oz (960)
Neck Size (MM)	15-415	18-410	18-410	24-410	24-410	28-400	38mm Buttress*

CAPS ATTACHED - CONVENIENCE PACK

Cat. No.	7367B	7311B	7312B	7314B	7315B	7316B	7117G
Wt. Case (lbs.)	2	2	2	4	8	5	18
No. Per Pack	12	12	12	12	12	6	—
No. Per Case	48	48	48	48	48	24	65

BOTTLES ONLY - BULK CASE

Cat. No.	3067	3011	3012	3014	3015	3016	3017
Wt. Case (lbs.)	32	29	21	18	19	19	15
No. Per Case	3100	1550	850	430	225	140	65

*POLYETHYLENE LINERLESS BUTTRESS CAP ATTACHED



WIDE MOUTH ROUNDS

WHITE PP STURDEESEAL P/E FOAM LINED CAPS

Capacity (ml)	1oz (30)	2oz (60)	4oz (120)	8oz (250)	16oz (500)	32oz (950)	66oz(2000)	100oz(3000)	128oz(3840)
Neck Size (MM)	28-400	33-400	38-400	53-400	53-400	53-400	89-400	89-400	89-400

CAPS ATTACHED - CONVENIENCE PACK

Cat. No.	7352G	7307G	7354G	7355G	7356G	7357G	7358G	7359G	7320G
Wt. Case (lbs.)	2	2	3	3	2	3	3	2	3
No. Per Pack	12	12	12	12	6	—	—	—	—
No. Per Case	48	48	48	48	24	12	6	4	4

CAPS ATTACHED - BULK CASE

Cat. No.	7152G	7107G	7154G	7155G	7156G	7157G	7158G	7159G	—
-----------------	--------------	--------------	--------------	--------------	--------------	--------------	--------------	--------------	---

PLASTIC BOTTLE ONLY - BULK CASE

Cat. No.	3052	3007	3054	3055	3056	3057	3058	3059	3020
Wt. Case (lbs.)	25	21	18	15	14	13	19	11	3
No. Per Case	1650	950	500	260	137	72	75	33	4

Lightweight, Shatterproof Alternative

Select the Plastic Resin That's Right for You.

Plastic containers are unbreakable, lightweight, functional and compact. If your product can be packed in plastic, you need to select the plastic resin that suits your application.



POLYETHYLENE TEREPHTHALATE

- Clear and transparent
- Good alcohol and solvent barrier
- Good stress and impact resistance



HIGH DENSITY POLYETHYLENE

- Excellent moisture barrier
- Impact resistant
- Translucent
- Ideal for high-speed filling, demanding pharmaceutical specifications, and the rough handling of interstate shipping.



VINYL/POLYVINYL CHLORIDE

- Transparent, rigid plastic.
- Excellent resistance to oxygen permeation, oils and good water barrier characteristics.



LOW DENSITY POLYETHYLENE

- Used for applications requiring squeezability and resilience.



POLYPROPYLENE

- Greater tolerance to higher temperatures than HDPE and better contact clarity.



POLYSTYRENE

- Offers great rigidity and glass-like clarity.



WIDE MOUTH OBLONG BOTTLES

WHITE PP STURDEESEAL P/E FOAM LINED CAPS

Capacity (ml)	2oz (60)	4oz (125)	5oz (160)	8oz (250)	16oz (500)	33oz(1000)	40oz (1280)
Neck Size (MM)	33-400	38-400	38-400	43-400	43-400	53-400	53-400

CAPS ATTACHED - CONVENIENCE PACK

Cat. No.	3M02QWFC	3M04QWFC	3M05QWFC	3M08QWFC	3M17QWFC	3M33QWFC	3M41QWFC
Wt. Case (lbs.)	3	4	4	5	5	4	3
No. Per Pack	12	12	12	12	6	—	—
No. Per Case	48	48	48	48	24	12	6

CAPS ATTACHED - BULK CASE

Cat. No.	3M02BWFC	3M04BWFC	3M05BWFC	3M08BWFC	3M17BWFC	3M33BWFC	3M41BWFC
----------	----------	----------	----------	----------	----------	----------	----------

BOTTLES ONLY - BULK CASE

Cat. No.	3M02B	3054Z	3010	3M08B	3066	3M33B	3M41B
Wt. Case (lbs.)	23	20	17	15	15	21	16
No. Per Case	950	500	400	250	160	116	84



CLEAR PET BOTTLES

Capacity (ml)	60cc	120cc	200cc	250cc	500cc
Neck Size (MM)	33-400	38-400	38-400	45-400	45-400

CONVENIENCE PACK

WHITE PP/F217 caps	PLC-06589	PLC-06591	—	PLC-06594	PLC-06596
WHITE PP/F217 caps - ULTRA CLEAN	—	PLC-06597	—	—	PLC-06806
GREEN TEFLON PTFE caps ●	—	PLC-06794	—	—	—
No. Per Case	48	48	—	48	24

BULK CASE

WHITE PP/F217 caps	PLC-06764	PLC-06762	PLC-06884	—	PLC-06877
GREEN TEFLON PTFE caps ●	—	PLC-06795	—	—	—
No. Per Case	1150	576	302	270	140

BOTTLES ONLY - BULK CASE

Cat. No.	PLA-06565	PLA-06572	PLA-06883	PLA-06568	PLA-06570
No. Per Case	1150	576	302	270	140

HEAVY WEIGHT PET BOTTLES

24 GRAM THICKNESS

Capacity (ml)	120cc
Neck Size (MM)	38-400

CONVENIENCE PACK

WHITE PP/F217 caps - ULTRA CLEAN	PLC-06796
No. Per Case	48

BULK CASE

WHITE PP/F217 caps	PLC-06793
WHITE PP/F217 caps - SUPER CLEAN	PLC-06996
No. Per Case	450

BOTTLES ONLY - BULK CASE

Cat. No.	PLA-06792
No. Per Case	450

AMBER PET BOTTLES

Capacity (ml)	60cc	120cc	250cc	500cc
Neck Size (MM)	33-400	38-400	45-400	45-400

CONVENIENCE PACK

WHITE PP/F217 caps	PLC-06590	PLC-06592	PLC-06593	PLC-06595
No. Per Case	48	48	48	24

BOTTLES ONLY - BULK PACK

Cat. No.	PLA-06566	PLA-06567	PLA-06569	PLA-06571
No. Per Case	1150	576	270	140

PET BOTTLES ARE IDEAL FOR OIL ANALYSIS. AVAILABLE IN YOUR CHOICE OF CASE PACKS. BOTTLES CAN BE SUPER CLEANED TO ISO SPECIFICATIONS. PLEASE CALL FOR FURTHER INFORMATION.

TEFLON® IS A REGISTERED TRADEMARK OF E.I. DUPONT DE NEMOURS AND COMPANY. ● THE COLOR GREEN FOR CAPS IS A TRADEMARK OF QORPAK, REGISTERED IN THE U.S. PATENT AND TRADEMARK OFFICE.




QORPAK BOTTLES RETAIN THEIR CONTAINER SHAPE.

The Plastic that Performs Like Glass.

Until now, you had to choose between a container that was lightweight and impact resistant or one that could stop chemical permeation. Qorpak Plastic Barrier Bottles give you the benefits of plastic and the containment qualities necessary for handling most hydrocarbon-based chemicals.

Qorpak Barrier Plastic containers provide:

- A guard against product weight loss and permeation.
- Improved product odor control.
- A guard against discoloration of the container.
- Container shape retention.



Qorpak Barrier Plastic provides a barrier to permeation of hydro-carbon solvents 100 times better than high density polyethylene!



HDPE FLUORINATED BARRIER PLASTIC BOTTLES

WITH TEFLON PTFE LINED CAPS

Capacity (ml)	2oz (60)	4oz (120)	8oz (250)	16oz (500)	32oz (1000)
Neck Size (MM)	24-410	24-410	24-410	28-400	28-400

CAPS ATTACHED - CONVENIENCE PACK

Cat. No.	7890T	7891T	7892T	7893T	7894T
Wt. Case (lbs.)	2	3	4	4	3
No. Per Case	48	48	48	24	12



NARROW MOUTH HDPE LABORATORY STYLE BOTTLES

WITH SPECIAL FINISH PP LINERLESS CAPS

Capacity (ml)	1oz (30)	2oz (60)	4oz (120)	8oz (250)	16oz (500)	32oz (1000)
Neck Size (MM)	20	20	24	24	28	38

CAPS ATTACHED - CONVENIENCE PACK

Cat. No.	3512	3513	3514	3515	3516	3517
Wt. Case (lbs.)	3	3	3	6	7	7
No. Per Case	72	72	72	72	48	24



WIDE MOUTH HDPE LABORATORY STYLE BOTTLES

WITH SPECIAL FINISH PP LINERLESS CAPS

Capacity (ml)	1oz (30)	2oz (60)	4oz (120)	8oz (250)	16oz (500)	32oz (1000)
Neck Size (MM)	28	28	38	45	53	63

CAPS ATTACHED - CONVENIENCE PACK

Cat. No.	3522	3523	3524	3525	3526	3527
Wt. Case (lbs.)	3	3	4	7	8	7
No. Per Case	72	72	72	72	48	24



WIDE MOUTH POLYPROPYLENE JARS

AUTOCLAVABLE WITH WHITE PP LINERLESS CAPS ATTACHED

Capacity (ml)	1oz (30)	2oz (60)	4oz (120)	8oz (250)	16oz (500)	32oz (960)
Neck Size (MM)	43-400	53-400	58-400	89-400	89-400	120-400

CAPS ATTACHED - CONVENIENCE PACK

Cat. No.	3820	3821	3822	3823	3824	3825
Wt. Case (lbs.)	3	3	3	5	5	8
No. Per Case	72	48	36	36	24	24

JAR ONLY - BULK CASE

Cat. No.	3220	3221	3222	3223	3224	3225
Wt. Case (lbs.)	22	15	18	13	14	12
No. Per Case	1188	546	432	168	112	60



WIDE MOUTH WHITE POLYPROPYLENE JARS

AUTOCLAVABLE WITH WHITE PP LINERLESS CAPS ATTACHED

Capacity (ml)	1oz (30)	2oz (60)	4oz (120)	8oz (240)	16oz (480)	32oz (960)
Neck Size (MM)	43-400	53-400	58-400	70-400	89-400	120-400

CAPS ATTACHED - CONVENIENCE PACK

Cat. No.	PLC-07184	PLC-07185	PLC-07186	PLC-07187	PLC-07188	PLC-07189
Wt. Case (lbs.)	20	3	4	6	6	10
No. Per Case	72	48	36	36	24	24

JAR ONLY - BULK CASE

Cat. No.	3Z01B	3229	PLA-06007	3223Z1	3S17B	3Z33B
Wt. Case (lbs.)	20	18	18	27	14	12
No. Per Case	980	588	432	336	112	60



WIDE MOUTH POLYSTYRENE JARS

WITH WHITE PP STURDEESEAL P/E FOAM LINED CAPS ATTACHED

Capacity (ml)	1oz (30)	2oz (60)	4oz (120)	8oz (240)	16oz (480)	32oz (960)
Neck Size (MM)	43-400	53-400	58-400	89-400	89-400	120-400

CAPS ATTACHED - CONVENIENCE PACK

Cat. No.	PLC-07572	3811	3891	3812	3892	3813
Wt. Case (lbs.)	4	3	4	6	6	10
No. Per Case	72	48	36	36	24	24

JAR ONLY - BULK CASE

Cat. No.	PLA-07565	3211	3091	3212	3092	3213
Wt. Case (lbs.)	26	20	19	17	15	16
No. Per Case	1188	546	324	168	112	60

With or Without Attached Caps



OTHER STYLES AVAILABLE. PLEASE INQUIRE.

LABORATORY GRADE POLYPROPYLENE HINGED VIALS

	CAPACITY	NO./PACK	CAT. NO.
Shorty Cornell Vial	11ml	500	PLD-08849
Tall Graduated ML Vet Vial - specify plain or lock tab	13ml	500	PLD-08850
Hi-Pro Vial	2oz (60)	600	PLD-05721
Hi-Pro Custody w/Sodium Thiosulfate w/string & lock tab	4oz (120)	200	PLD-04452C
Hi-Pro Vial - marked in oz and every 10ml	4oz (120)	250	PLD-06446
Hi-Pro Custody - w/string and lock tab	4oz (120)	200	PLD-06446C
Cup Graduated ml & oz	8oz (240)	150	PLD-05722
Short Vial w/Locking String Graduated ml & oz	10oz (300)	90	PLD-71028
Hi-Pro - specify w/custody string or lock tab	10oz (300)	100	PLD-71028C



POLYSTYRENE SNAP-CAP VIALS & CAPS

Capacity (ml)	5/8oz (18)	7/8oz (26)	1 5/8oz (48)	1 7/8oz (55)
Neck Size (MM)	7/8 inch	1 inch	1 5/16 inch	1 1/4 inch

CAPS ATTACHED - CONVENIENCE PACK				
Cat. No.	3075	3076	3078	3079
Wt. Case (lbs.)	3	3	6	5
No. Per Case	144	144	144	144



POLYETHYLENE JUGS

Capacity (ml)	64oz (1920)	128oz (3840)	128oz (3840)
Neck Size (MM)	38-400	38-400	38-400

CAPS ATTACHED - CONVENIENCE PACK			
Cat. No.	7119G	7121G	7322G
Wt. Case (lbs.)	16	23	3
No. Per Case	40	60	4



SQUARE HANDLED POLYETHYLENE JUGS

Capacity (ml)	128oz (3840)
Neck Size (MM)	38-400

CAPS ATTACHED - CONV. PACK	
Cat. No.	7161G
Wt. Case (lbs.)	3
No. Per Case	4



NATURAL HDPE LIGHTWEIGHT DAIRY JUGS

Capacity (ml)	8oz (240)	16oz (480)	32oz (960)	128oz (3785)
Neck Size (MM)	38-400	38-400	38-400	38-400

BOTTLES ONLY - BULK CASE				
Cat. No.	PLA-06816	PLA-06817	PLA-06818	PLA-05565
Wt. Case (lbs.)	35	24	22	10
No. Per Case	500	400	224	36



NATURAL HDPE HEAVYWEIGHT JUGS

Capacity (ml)	4 liter
Neck Size (MM)	38-439

BOTTLES ONLY - BULK CASE	
Cat. No.	3042
Wt. Case (lbs.)	32
No. Per Case	55



WHITE HDPE F-STYLE JUGS

Capacity (ml)	16oz (480)	32oz (960)	64oz (1920)	128oz (3785)
Neck Size (MM)	33-400	33-400	38-400	38-400

BOTTLES ONLY - BULK CASE				
Cat. No.	3070	3071	3072	3073
Wt. Case (lbs.)	26	17	16	4
No. Per Case	240	120	60	6



TALL WIDE MOUTH NATURAL HDPE ROUND

Capacity (ml)	32 (960)
Neck Size (MM)	70-400

BOTTLES ONLY - BULK CASE	
Cat. No.	3003
Wt. Case (lbs.)	16
No. Per Case	106



FLUORINATED NATURAL HDPE JUG WITH GREEN TEFLON CLOSURE

Capacity (ml)	128 (3840)
Neck Size (MM)	38-400

CAPS ATTACHED - CONVENIENCE PACK	
Cat. No.	PLC-07567
Wt. Case (lbs.)	3
No. Per Case	4



NATURAL HDPE WIDE MOUTH HANDLE ROUND

Capacity (ml)	128 (3840)
Neck Size (MM)	89-400

BOTTLES ONLY - BULK CASE	
Cat. No.	3023
Wt. Case (lbs.)	24
No. Per Case	60



UN RATINGS ON PAIL AND LID COMBINATIONS ARE AVAILABLE. VARIOUS COLOR PAILS AND LIDS AVAILABLE IN PALLET QUANTITIES. *ALSO AVAILABLE IN BULK. PLEASE CALL FOR MORE INFORMATION.

WHITE HIGH-DENSITY POLYETHYLENE PAILS

Capacity 128 oz 2 Gallon 3.5 Gallon 5 Gallon 6 Gallon

PAILS					
Cat. No.	7011	7001	3J03B	7044/4*	PLA-07351
Wt. Case (lbs.)	58	56	64	11	40
No. Per Case	100	60	36	4	18



LIDS

FOR WHITE HIGH-DENSITY POLYETHYLENE PAILS

Pail Capacity 128 oz 2 gallon 3.5, 5, 6 gal 3.5, 5, 6 gal 3.5, 5 gal 3.5, 5, 6 gal

LID STYLE:	PLAIN	PLAIN	PLAIN	EZ TEAR	JR REIKE SPOUT	SR REIKE SPOUT
Cat. No.	7012	7010	7046/4*	PLA-07597	7045/4*	7047
Wt. Case (lbs.)	50	62	4	37	4	48
No. Per Case	200	200	4	50	4	60



CUBITAINER IN ASSEMBLED CARTON

Capacity (ml) 32oz (960)* 128oz (3840) 5 Gallon (20,000)
Neck Size (MM) 38-400 38-400 38-400

WITH SCREW CAP			
Cat. No.	7201	3H00Q	7204
Wt. Case (lbs.)	16	5	3
No. Per Case	144	8	1

*NOT ASSEMBLED. CARTON NOT INCLUDED.

QUICK SERVE TAP

FOR CUBITAINER
Neck Size (MM) 38-400

Cat. No.	7206
No. Per Case	1



CARBOY

Capacity 5 gallon
Neck Size (MM) 53-400

Natural HDPE	3043
UN No. (Ⓢ)	3H1/Y1.8/250
No. Per Case	8

SPIGOTS

FOR CARBOYS
Neck Size (MM) 3/4"

Cat. No.	1Y03BY1
No. Per Case	25

CAPS

FOR CARBOYS
Neck Size (MM) 53-400 w/3/4" insert

Cat. No.	3051
No. Per Case	50



WINPAK

Capacity 5 gallon
Neck Size (MM) 70-400

WITH SCREW CAP	
Natural HDPE	7225
UN No. (Ⓢ)	1H1/Y1.8/100
No. Per Case	1



HEDPAK

Capacity 5 gallon
Neck Size (MM) 70-400

WITH SCREW CAP	
Natural HDPE	7221
No. Per Case	1

UN PACKING GROUP 2 1Y03BY1 SPIGOT CAN BE USED WITH THIS CAP.



INDUSTRIAL SPRAYER*

Capacity (ml) 16oz (480)
Neck Size (MM) 28-400

Cat. No.	PLC-06934
Wt. Case (lbs.)	2
No. Per Case	6

*NYLON / POLYETHYLENE BOTTLE HELPS REDUCE RISK OF INCOMPATIBLE PRODUCT. POLYESTER / POLYPROPYLENE SPRAYER.



TRIGGER SPRAYER BOTTLES

Capacity (ml) 16oz (480) 22oz (660) 32oz (960)
Neck Size (MM) 28-400 28-400 28-400

Cat. No.	7331X	7340X	7333X
Wt. Case (lbs.)	2	2	4
No. Per Case	6	4	12



YORKER TOP DISPENSING BOTTLES

Capacity (ml) 2oz (60) 4oz (120) 8oz (250)
Neck Size (MM) 24-410 24-410 24-410

CAPS ATTACHED - CONVENIENCE PACK			
Cat. No.	7336H	7337H	7338H
Wt. Case (lbs.)	3	3	4
No. Per Case	48	48	48



FLIP TOP DISPENSING BOTTLES

Capacity (ml) 2oz (60) 4oz (120) 8oz (240) 16oz (500)
Neck Size (MM) 24-410 24-410 24-410 24-410

CAPS ATTACHED - CONVENIENCE PACK				
Cat. No.	7336J	7337J	7338J	7339J
Wt. Case (lbs.)	3	3	4	4
No. Per Case	48	48	48	24

1-800-922-7558

Specialty Caps & Liner Systems



WHITE POLYPROPYLENE FLIP-TOP CAPS WITH 3.0MM ORIFICE

Neck Size (MM)	24-410	28-400
Cat. No.	CAP-06047	CAP-06048
No. Per Case	144	144



BLACK PHENOLIC DROPPER ASSEMBLIES*

Neck Size (MM)	20-400	20-400	22-400
Dropper Length (MM)	74-76	86-88	100-102
Cat. No.	CAP-05740	CAP-05742	CAP-05744
No. Per Case	12	12	12

* FIT BOTTLES WITH SAME NECK SIZE AND HEIGHT



NATURAL LDPE YORKER SPOUTS WITH RED CAP AND .030 ORIFICE

Neck Size (MM)	18-400	24-400	28-400	38-400
Cat. No.	CAP-06046	CAP-06053	CAP-06049	CAP-06050
No. Per Case	144	144	144	144



WHITE POLYPROPYLENE HOLE CLOSURES / SONICALLY WELDED SEPTA

Neck Size (MM)	28-400	33-400	38-400	48-400	53-400	58-400
AUTOCLAVABLE TEFLON® PTFE SILICONE						
Convenience (120/Case)	CAP-06680	CAP-06681	CAP-06682	CAP-06683	CAP-06684	CAP-06685
Bulk (1000/Case)	5Y28BHM	5Y33BHM	5Y38BHM	5Y48BHM	5Y53BHM	5Y58BHM



TEFLON PTFE DISCS

Neck Size (MM)	22-400	24-400	28-400	33-400	33-430	38-400	38-430	43-400
Cat. No.	5Y22QYD	5Y24QYD	5Y28QYD	5Y33QYD	5Y34QYD	5Y38QYD	5Y30QYD	5Y43QYD
AUTOCLAVABLE								
Neck Size (MM)	45-400	48-400	53-400	58-400	63-400	70-400	83-400	89-400
Cat. No.	5Y45QYD	5Y48QYD	5Y53QYD	5Y58QYD	5Y63QYD	5Y70QYD	5Y83QYD	5Y89QYD



HOLE CAPS AND SEPTAS

Neck Size (MM)	20-400	22-400	24-410
OPEN TOP HOLE CAPS			
Convenience Pack	CAP-07004	CAP-07850	5Y24QH Y
No. Per Case	500	200	100
Bulk Pack	CAP-00330	—	5240
No. Per Case	5300	—	3500
GREEN PP			
Neck Size (MM)	20MM	22MM	24MM
.125 TEFLON SILICONE SEPTA			
Convenience Pack	—	CAP-05188	5Y24QYS
No. Per Case	—	100	100
Bulk Pack	CAP-06864	—	5240A
No. Per Case	1000	—	4000

.125 TEFLON/SILICONE SEPTA FITS 24MM HOLE CAP. IDEAL FOR VOA ANALYSIS.



Why You Should Order Caps Attached to Bottles.

- Saves you time. No more searching for the proper sized cap.
- Prevents mistakes. The right liner with the right cap, every time.
- Eliminates washing prior to use. Carton dust can't get inside the bottle or behind the liner.
- Simplifies ordering. Consolidates two inventory items into one.



SPECIALTY CAPS

Specialty caps include spray and pump dispensers, and dispensing applicator caps. Specialty caps allow for proper application of a product, dispensing measured amounts and one-handed application. Please call for catalog numbers and specific information.

CAPS & LINERS



GREEN TEFLON® PTFE LINED THERMOJET CAPS

No need to worry about permeation of the product or an attack of the liner system by the product. These caps provide a wide range of chemical compatibility and temperature tolerance. The foam backed Teflon PTFE resists attack from virtually all chemicals at room temperature and provides resilience for a tight seal. Not autoclavable.

BLACK PHENOLIC THERMOJET CAPS

These caps are known for providing the widest range of chemical compatibility and for being the most temperature tolerant of plastic closures. They also provide the most consistent adherence to close dimensional tolerances. Seven different liner systems are available from stock. (See chart.)

* POLYPRO SHELL
 ** AUTOCLAVABLE — ACCEPTABLE FOR STEAM HEATING OR STERILIZING WITH HEAT.



White polypropylene caps are available with a SturdeeSeal® Polyethylene foam liner or linerless. The SturdeeSeal® liner has excellent sealing characteristics and good chemical resistance. Linerless caps are suitable for many dry products.



For use in the widest range of temperatures, metal caps are also very resistant to fracture from impact. Available in a full range of sizes with several liner systems including Plastisol.

GREEN TEFLON® PTFE LINED THERMOJET CAPS

Neck Size (MM)	8-425	13-425	15-425	18-400	18-410	20-400	22-400	24-400	24-410	24-414
Convenience Pack	5218/100	5200/100	5201/100	5202/100	—	5203/100	5204/100	5205/100	—	—
Bulk Case	5218	5200	5201	5202	—	5203	5204	5205	—	—
Bulk Case Quantity	24000	13500	11000	6500	—	4800	4600	4000	—	—

BLACK PHENOLIC THERMOJET CAPS

Neck Size (MM)	8-425	13-425	15-425	18-400	18-410	20-400	22-400	24-400	24-410	24-414
PTFE LINED										
Convenience Pack	—	—	—	—	—	CAP-06689	CAP-06690	CAP-06691	—	—
Bulk Case	—	—	—	—	—	5203A	5204A	5205A	—	—
Bulk Case Quantity	—	—	—	—	—	5500	4400	3600	—	—

PULP/VINYL LINER

Convenience Pack	5001/12	5002/12	5003/12	5004/12	5005/12	5007/12	5009/12	5011/12	5012/12	—
Bulk Case	5001	5002	5003	5004	5005	5007	5009	5011	5012	—
Bulk Case Quantity	24000	13500	11000	6500	4900	4800	4600	4000	2200	—

POLYSEAL CONE LINER

Convenience Pack	—	5045/12	5083/12	5069/12	—	5084/12	5070/12	5026/12	5027/12	—
Bulk Case	—	5045	5083	5069	—	5084	5070	5026	5027	—
Bulk Case Quantity	—	13000	10000	6500	—	5500	4700	4200	2900	—

.030 SOLID PE LINER

Convenience Pack	—	—	—	—	—	—	5071/12	5072/12	—	—
Bulk Case	—	—	—	—	—	—	5071	5072	—	—
Bulk Case Quantity	—	—	—	—	—	—	4600	4000	—	—

PULP/TINFOIL LINER

Convenience Pack	5116/12	5117/12	5118/12	5119/12	—	5120/12	5121/12	5122/12	—	—
Bulk Case	5116	5117	5118	5119	—	5120	5121	5122	—	—
Bulk Case Quantity	24000	13500	11000	6500	—	4800	4600	4000	—	—

14B WHITE RUBBER LINER**

Convenience Pack	—	5086/12	5087/12	5088/12	—	5089/12	5090/12	5091/12	—	5097/12
Bulk Case	—	5086	5087	5088	—	5089	5090	5091	—	5097
Bulk Case Quantity	—	13500	11000	6500	—	4800	4600	4000	—	3000

PULP/ALUMINUM FOIL LINER

Convenience Pack	—	5150/12	5151/12	—	—	—	5153/12	5154/12	—	—
Bulk Case	—	5150	5151	—	—	—	5153	5154	—	—
Bulk Case Quantity	—	13500	11000	—	—	—	4600	4000	—	—

WHITE POLYPROPYLENE CAPS

Neck Size (MM)	24-410	28-400	33-400	38-400	43-400	45-400	48-400	53-400	58-400	63-400
FOAMED POLYETHYLENE STURDEESEAL® LINER										
Convenience Pack	CAP-06679	CAP-06624	CAP-06623	CAP-06625	CAP-06626	CAP-06627	CAP-06628	CAP-06629	CAP-06630	CAP-06631
Bulk Case	5401	5402	5403	5404	5412	CAP-05019	5416	5405	5413	5414
Bulk Case Quantity	3500	4300	3000	2250	1800	2000	1700	1300	1000	900

LINERLESS AUTOCLAVABLE

Convenience Pack	—	—	CAP-06633	CAP-06634	CAP-06635	—	CAP-06636	CAP-06637	CAP-06638	CAP-06639
Bulk Case	5Y25BWU	—	5507	5508	5509	—	5510	5511	5512	5515
Bulk Case Quantity	4000	—	4000	2900	2200	—	1700	1300	1100	900

BLACK POLYPROPYLENE CAPS

Neck Size (MM)	22-400	24-400	28-400	33-400	33-430	38-400	38-430	43-400	45-400	48-400
LINERLESS AUTOCLAVABLE										
Convenience Pack	5Y22QEU	5Y24QEU	5Y28QEU	5Y33QEU	5Y34QEU	5Y38QEU	5Y30QEU	5Y43QEU	5Y45QEU	5Y48QEU
Bulk Case	5Y22BEU	5Y24BEU	5Y28BEU	5Y33BEU	5Y34BEU	5Y38BEU	5Y30BEU	5Y43BEU	5Y45BEU	5Y48BEU
Bulk Case Quantity	7000	6500	5000	4000	2200	2900	1000	2200	2000	1700

WHITE METAL CAPS

Neck Size (MM)	24-400	28-400	33-400	38-400	43-400	48-400	53-400	58-400	63-400	70-400
PLASTISOL LINER										
Convenience Pack	—	—	—	CAP-06642	—	CAP-06643	CAP-06644	CAP-06645	4050/12	4057/12
Bulk Case	—	—	—	4072	—	4048	4058	4049*	4050	4057

PULP/ALUMINUM FOIL LINER

Convenience Pack	CAP-06646	CAP-06647	CAP-06648	CAP-06649	CAP-06650	CAP-06651	CAP-06652	CAP-06653	CAP-06654	CAP-06655
Bulk Case	5Y24BMA	5Y28BMA	5Y33BMA	5Y38BMA	5Y43BMA	5Y48BMA	5Y53BMA	5Y58BMA	5Y63BMA	5Y70BMA

PULP/POLYETHYLENE LINER

Convenience Pack	4033/12	4034/12	4035/12	4036/12	CAP-06658	CAP-06659	4039/12	CAP-06660	CAP-06661	CAP-06662
Bulk Case	4033	4034	4035	4036	4037	4038	4039	4040	4041	4042
Bulk Case Quantity	6500	4600	3400	3800	3100	2300	1900	2000	1700	1400

*BULK CASE 1,900 / CASE.

tion of In-Stock Closure Sizes & Liner Systems

CAPS & LINERS

28-400 33-400 33-430 38-400 38-430 43-400 45-400 48-400 53-400 58-400 63-400 70-400 83-400 89-400 100-400

5206/100	5207/100	—	5208/100	—	5209/100	5210/100	5211/100	5212/100	5213/100	5214/12	5215/12	5216/12	5217/12	—
5206	5207	—	5208	—	5209	5210	5211	5212	5213	5214	5215	5216	5217	—
3000	2100	—	2160	—	1760	1520	1400	1080	980	740	640	340	304	—

28-400 33-400 33-430 38-400 38-430 43-400 45-400 48-400 53-400 58-400 63-400 70-400 83-400 89-400 100-400

CAP-06692	CAP-06693	CAP-06700	CAP-06694	CAP-06701	—	CAP-06695	CAP-06696	CAP-06697	CAP-06698	—	CAP-06699	—	—	—
5206A	5207A	5220A	5208A	5221A	—	5210A	5211A	5212A	5213A	—	5215A	—	—	—
2800	2150	1100	2160	900	—	1520	1400	1080	980	—	640	—	—	—

5013/12	5014/12	—	5015/12	—	5016/12	5017/12	5018/12	5019/12	5020/12	5021/12	5022/12	5023/12	5024/12	5025/12*
5013	5014	—	5015	—	5016	5017	5018	5019	5020	5021	5022	5023	5024	5025
3000	2100	—	2160	—	1760	1520	1400	1080	980	740	640	340	304	256

5028/12	5029/12	—	5030/12	—	—	—	—	—	—	—	—	—	—	—
5028	5029	—	5030	—	—	—	—	—	—	—	—	—	—	—
3100	2300	—	1600	—	—	—	—	—	—	—	—	—	—	—

5073/12	5074/12	—	5075/12	—	5076/12	—	5077/12	5078/12	5079/12	5080/12	5081/12	5067/12	5082/12	—
5073	5074	—	5075	—	5076	—	5077	5078	5079	5080	5081	5067	5082	—
3000	2100	—	2160	—	1720	—	1400	1080	980	740	640	340	304	—

5123/12	5124/12	—	5125/12	—	5127/12	5128/12	5129/12	5130/12	5131/12	5132/12	5133/12	—	5134/12	—
5123	5124	—	5125	—	5127	5128	5129	5130	5131	5132	5133	—	5134	—
3000	2100	—	2160	—	1760	1520	1400	1080	980	740	640	—	304	—

5092/12	5093/12	—	5094/12	—	5095/12	—	5098/12	5099/12	5100/12	—	5101/12	—	—	—
5092	5093	—	5094	—	5095	—	CAP-05781	CAP-05782	CAP-05783	—	CAP-05784	—	—	—
3000	2800	—	2160	—	1720	—	1400	1080	980	—	640	—	—	—

5152/12	5156/12	—	5157/12	—	5158/12	—	5159/12	5160/12	5161/12	5162/12	5163/12	5165/12	5164/12	—
5152	5156	—	5157	—	5158	—	5159	5160	5161	5162	5163	5165	5164	—
3000	2100	—	2160	—	1760	—	1400	1080	980	740	640	340	304	—

70-400 83-400 89-400 120-400

5417/12	CAP-06632	5406/12	5415/6
5417	5407	5406	5415
900	500	370	210

CAP-05088	—	CAP-06640	CAP-06641
—	—	5513	5514
760	—	456	210

53-400 58-400 63-400 70-400 83-400 89-400

5Y53QEU	5Y58QEU	5Y63QEU	5Y70QEU	5Y83QEU	5Y89QEU
5Y53BEU	5Y58BEU	5Y63BEU	5Y70BEU	5Y83BEU	5Y89BEU
1300	1100	900	760	500	456

83-400 89-400 100-400

4067/12	4069/12	—
4067	4069	—


CAP-06656	CAP-06657	—
5Y83BMA	5Y89BMA	—

CAP-06663	CAP-06664	CAP-06665
4043	4044	4045
875	700	500

TEFLON® IS A REGISTERED TRADEMARK OF E.I. DUPONT DE NEMOURS AND COMPANY.
 ● THE COLOR GREEN FOR CAPS IS A TRADEMARK OF QORPAK, REGISTERED IN THE U.S. PATENT AND TRADEMARK OFFICE.

Caps & Neck Finishes

Screw thread cap sizes are expressed in two numbers, which refer to the cap diameter (in millimeters) and the GPI Finish (the top of the container). For example, a size 28-430 cap has a diameter of 28 mm, and fits a container with a 430 GPI thread finish.



Kapclean

Environmental Containers – Available Options:

KAPT Catalog Numbers:

Bottle is fully assembled with polypropylene closure and .015” Teflon PTFE disc.

CLEAN Catalog Numbers :

Bottles and closure systems have been pre-cleaned and certified to meet designated standards. Bottle is fully assembled with polypropylene closure and .015” Teflon PTFE disc. Certificate of analysis is included in each case.

Save Time.

No need to wash bottles before using.

Provide A Clean Final Appearance.

There won't be dust particles floating in your samples.

Sealed To Stay Clean.

Since each container is individually sealed, it won't collect dust or other debris in storage.

Extremely Cost Effective.

You get outstanding cleanliness at minimal cost.

Pre-cleaned for Volatile Organic Analysis.

VIALS

AVAILABLE IN CLEAR OR AMBER BOROSILICATE GLASS WITH GREEN OPEN TOP CAP AND TEFLON® PTFE/SILICONE SEPTA

Capacity	Neck Size (MM)	No./Case	Lbs./Case	Kapt Cat.#	Clean Cat.#
CLEAR (PHOTO PAGE 3)					
2/3 oz (20)	24-400	144	10	0200-0005	0200-0017
1 1/3oz (40) Squat	24-400	144	11	0200-0007	0200-0019
1 1/3oz (40) Tall	24-400	144	11	0200-0008	0200-0020
NEW! 2 oz (60) Tall	24-410	72	8	GLC-07863	GLC-07865
AMBER (PHOTO PAGE 3)					
2/3 oz (20)	24-400	144	10	0200-0006	0200-0018
1 1/3oz (40)	24-400	144	11	0200-0009	0200-0021
NEW! 2 oz (60) Tall	24-410	72	8	GLC-07867	GLC-07869

PRE-PRESERVED VIALS – .5ML 1:1 HCL

Capacity	Neck Size (MM)	No./Case	Lbs./Case	Kapt Cat.#	Clean Cat.#
TALL CLEAR (PHOTO PAGE 3)					
1 1/3oz (40)	24-400	144*	11	—	GLC-05606
TALL AMBER (PHOTO PAGE 3)					
1 1/3oz (40)	24-400	144*	11	—	GLC-05543

* VIALS PACK 2 TRAYS OF 72

NARROW MOUTH GLASS BOTTLES

WITH GREEN OPEN TOP CAP WITH TEFLON® PTFE/SILICONE SEPTA

Capacity	Neck Size (MM)	No./Case	Lbs./Case	Kapt Cat.#	Clean Cat.#
CLEAR BOSTON ROUND (PHOTO PAGE 5)					
8oz (240)	24-400	12	6	0200-0010	0200-0022
AMBER BOSTON ROUND (PHOTO PAGE 4)					
8oz (240)	24-400	12	6	0200-0011	0200-0023

WIDE MOUTH GLASS BOTTLES

WITH ULTRASONICALLY WELDED SEPTA CAP

Capacity	Neck Size (MM)	No./Case	Lbs./Case	Kapt Cat.#	Clean Cat.#
CLEAR MEDIUM ROUND (PHOTO PAGE 6)					
4oz (120)	48-400	24	8	0200-0012	0200-0024
CLEAR STRAIGHT SIDED ROUNDS (PHOTO PAGE 4)					
2oz (60)	53-400	24	6	0200-0013	0200-0025
4oz (120)	58-400	24	8	0200-0014	0200-0026
AMBER WIDE MOUTH PACKERS (PHOTO PAGE 4)					
16oz (480)	53-400	12	11	0200-0015	0200-0027
32oz (950)	58-400	12	12	0200-0016	0200-0028

COLIFORM TESTING - STERILE VIAL

PRE-PRESERVED WITH SODIUM THIOSULFATE

Capacity	Neck Size (MM)	No./Case	Lbs./Case	Kapt Cat.#	Clean Cat.#
(PHOTO PAGE 13)					
120ml*	Snap	200	15	—	PLD-04452C

*GRADUATED

TEFLON® IS A REGISTERED TRADEMARK OF E.I. DUPONT DE NEMOURS AND COMPANY. THE COLOR GREEN FOR CAPS IS A TRADEMARK OF QORPAK, REGISTERED IN THE U.S. PATENT AND TRADEMARK OFFICE.

Pre-cleaned for Semivolatiles, Pesticides, PCBs and Metals.

Pre-cleaned for Metals Analysis.

High Density Polyethylene bottles are available in cylinder round style, wide-mouth packer style or straight sided cylindrical style with white polypropylene closure with foamed polyethylene liners.

NARROW MOUTH GLASS BOTTLES

Capacity	Neck Size (MM)	No./Case	Lbs./Case	Kapt Cat.#	Clean Cat.#
----------	----------------	----------	-----------	------------	-------------

CLEAR BOSTON ROUNDS (PHOTO PAGE 5)

8oz (240)	24-400	12	5	0100-0001	0100-0025
16oz (480)	28-400	12	8	0100-0002	0100-0026
32oz (960)	33-400	12	14	0100-0003	0100-0027

CLEAR JUG (PHOTO PAGE 7)

128oz (3785)	38-400	4	13	0100-0004	0100-0028
--------------	--------	---	----	-----------	-----------

AMBER BOSTON ROUNDS (PHOTO PAGE 4)

4oz (120)	22-400	24	6	0100-0005	0100-0029
8oz (240)	24-400	12	6	0100-0006	0100-0030
16oz (480)	28-400	12	8	0100-0007	0100-0031
1 Liter	33-430	12	14	0100-0008	0100-0032
32oz	33-400	30	12	0100-0068	0100-0073

AMBER JUGS (PHOTO PAGE 7)

80oz (2400)	38-430	6	15	0100-0009	0100-0033
4 Liters	38-430	4	13	0100-0010	0100-0034

WIDE MOUTH GLASS BOTTLES

Capacity	Neck Size (MM)	No./Case	Lbs./Case	Kapt Cat.#	Clean Cat.#
----------	----------------	----------	-----------	------------	-------------

CLEAR STRAIGHT SIDED ROUNDS (PHOTO PAGE 5)

2oz (60)	53-400	24	5	0100-0011	0100-0035
4oz (120)	58-400	24	7	0100-0012	0100-0036
9oz (270)	70-400	24	13	2S080EDC	0100-0037A
16oz (480)	89-400	12	10	0100-0014	0100-0038
32oz (950)	89-400	12	15	0100-0015	0100-0039

AMBER STRAIGHT SIDED ROUNDS (PHOTO PAGE 4)

4oz (120)	58-400	24	8	GLC-07259	GLC-06774
8oz (240)	70-400	12	7	GLC-07260	GLC-06775
32oz (960)	89-400	12	15	GLC-07261	GLC-06776

AMBER WIDE MOUTH PACKERS (PHOTO PAGE 4)

2oz (60)	33-400	24	5	0100-0016	0100-0040
4oz (120)	38-400	12	4	0100-0017	0100-0041
8oz (250)	45-400	12	10	0100-0018	0100-0042
16oz (480)	53-400	12	8	0100-0019	0100-0043
32oz (950)	53-400	12	12	0100-0020	0100-0044

CLEAR TALL MEDIUM ROUND (PHOTO PAGE 6)

4oz (120)	48-400	24	8	0100-0021	0100-0045
-----------	--------	----	---	-----------	-----------

CLEAR TALL STRAIGHT SIDED JARS (PHOTO PAGE 5)

8oz (240)	58-400	24	14	0100-0022	0100-0046
16oz (480)	63-400	12	10	0100-0023	0100-0047

CLEAR WIDE MOUTH JARS (PHOTO PAGE 5)

64oz (1920)	83-400	6	10	0100-0024	0100-0048
128oz (3810)	89-400	4	11	0100-0087	0100-0063

CLEAR STANDARD WIDE MOUTH JARS (PHOTO PAGE 5)

4oz (120)	48-400	24	6	2W04QEDC	0100-0084
8oz (250)	58-400	24	10	0100-0077	0100-0083
16oz (480)	63-400	24	14	0100-0078	0100-0082
32oz (950)	70-400	12	13	0100-0079	0100-0081

HDPE NARROW MOUTH PLASTIC BOTTLES

Capacity	Neck Size (MM)	No./Case	Lbs./Case	Kapt Cat.#	Clean Cat.#
----------	----------------	----------	-----------	------------	-------------

HDPE CYLINDER ROUNDS (PHOTO PAGE 10)

120ml	24-410	48	3	PLC-05733	0300-0250
240ml	24-410	24	5	PLC-05734	0300-0021
500ml	28-410	24	3	0300-0001	0300-0006

HDPE WIDE MOUTH PLASTIC BOTTLES

Capacity	Neck Size (MM)	No./Case	Lbs./Case	Kapt Cat.#	Clean Cat.#
----------	----------------	----------	-----------	------------	-------------

HDPE PLASTIC STRAIGHT SIDED JARS

500ml	89-400	24	4	0300-0016	0300-0017
1000ml	89-400	12	3	0300-0019	0300-0020

HDPE PLASTIC WIDE MOUTH ROUNDS (PHOTO PAGE 10)

30cc	28-400	48	2	PLC-05735	0300-0012
60cc	33-400	48	2	PLC-05736	PLC-05058
120cc	38-400	48	3	PLC-05737	PLC-03645
250ml	53-400	24	3	0300-0003	0300-0008
500ml	53-400	12	2	0300-0004	0300-0009
950ml	53-400	12	3	0300-0005	0300-0010

Other Bottle Styles and Cleaning Procedures are Available.

Call for Details.
1-800-922-7558

1-800-922-7558

ULTRA CLEAN



CLEAR JUGS

Capacity (ml)	64oz (1920)	128oz (3785)
Neck Size (MM)	38-400	28-400

CAPS ATTACHED - CONVENIENCE PACK

TEFLON® PTFE LINED caps ●	8663T	8664T
PULP/VINYL LINED caps	8663	8664
Wt. Case (lbs.)	14	13
No. Per Case	6	4
Minimum Case Quantity	8	7

ULTRA CLEAN



AMBER WIDE MOUTH PACKERS

Capacity (ml)	32oz (960)
Neck Size (MM)	53-400

CAPS ATTACHED - CONVENIENCE PACK

TEFLON® PTFE LINED caps ●	8638T
Wt. Case (lbs.)	12
No. Per Case	12
Minimum Case Quantity	6

ULTRA CLEAN



AMBER BOSTON ROUNDS

Capacity (ml)	8oz (240)	16oz (480)	32oz (960)
Neck Size (MM)	24-400	28-400	33-400

CAPS ATTACHED - CONVENIENCE PACK

TEFLON® PTFE LINED caps ●	8620T	—	—
PULP/VINYL LINED caps	—	8622	8623
POLYSEAL caps	—	8622B	—
Wt. Case (lbs.)	10	8	14
No. Per Case	24	12	12
Minimum Case Quantity	6	9	8

ULTRA CLEAN



CLEAR BOSTON ROUNDS

Capacity (ml)	1oz (30)	2oz (60)	4oz (120)	8oz (240)	16oz (480)	32oz (960)
Neck Size (MM)	20-400	20-400	22-400	24-400	28-400	33-400

CAPS ATTACHED - CONVENIENCE PACK

TEFLON® PTFE LINED caps ●	8609T	8610T	8612T	8613T	8614T	8617T
PULP/VINYL LINED caps	8609	8610	8612	8613	8614	8617
POLYSEAL caps	—	—	—	—	8614B	—
Wt. Case (lbs.)	5	4	6	10	8	14
No. Per Case	48	24	24	24	12	12
Minimum Case Quantity	4	6	6	6	9	8

ULTRA CLEAN



CLEAR MEDIUM ROUNDS

Capacity (ml)	1oz (30)	2oz (60)	4oz (120)	8oz (240)	16oz (480)
Neck Size (MM)	33-400	38-400	48-400	58-400	70-400

CAPS ATTACHED - CONVENIENCE PACK

TEFLON® PTFE LINED caps ●	8681T	8682T	8683T	8684T	8685T
PULP/VINYL LINED caps	8681	8682	8683	8684	8685
Wt. Case (lbs.)	6	10	8	11	17
No. Per Case	48	48	24	24	24
Minimum Case Quantity	4	4	6	6	5

ULTRA CLEAN



CLEAR STANDARD WIDE MOUTHS

Capacity (ml)	4oz (120)	8oz (240)	16oz (480)	32oz (960)	64oz (1920)	128oz (3810)
Neck Size (MM)	48-400	58-400	63-400	70-400	83-400	89-400

CAPS ATTACHED - CONVENIENCE PACK

TEFLON® PTFE LINED caps ●	8667T	8668T	8670T	8672T	8673T	8674T
Wt. Case (lbs.)	6	10	14	13	10	11
No. Per Case	24	24	24	12	6	4
Minimum Case Quantity	6	6	5	6	8	7

ULTRA CLEAN



WIDE MOUTH FRENCH SQUARES

Capacity (ml)	1/2oz (15)	1oz (30)	2oz (60)	4oz (120)	8oz (240)	16oz (480)	32oz (960)
Neck Size (MM)	20-400	24-400	28-400	33-400	43-400	48-400	58-400

CAPS ATTACHED - CONVENIENCE PACK

TEFLON® PTFE LINED caps ●	8601T	8602T	8603T	8604T	8605T	8606T	8607T*
PULP/VINYL LINED caps	8601	8602	8603	8604	8605	8606	8607*
Wt. Case (lbs.)	4	6	9	9	12	17	25
No. Per Case	48	48	48	24	24	24	12
Minimum Case Quantity	4	4	4	6	6	5	5

Ultra Clean

This cleaning service is designed for those needing containers free from particulate contaminants.

Ideal for chemical and resin sampling applications. Also well suited to small filling line operations and general lab use.

By using a three-step process, all particles 5 microns or larger are eliminated and traces of bottle release agents and other contaminants are removed.

● Our cleaning procedure includes ultra-sonic cleaning, chemistry, a de-ionized water rinse and double sealing.

● After the bottles pass our rigid quality control tests, they are certified.

● All processes take place under Class 100 Environmental Conditions per Federal Standard 209 on Clean Rooms.

Vacuum & Ionization

– Eliminate Carbon Dust and Visible Dirt Particles

Vacuum & Ionization

When you select vacuum and ionization-treated containers, you don't have to be concerned with dust, dirt, carbon lint or fine glass particles. Contaminants are removed with air that's been compressed in a special area and then filtered.

As the air moves through a condensate separator and coalescing filter, it removes aerosols and fine dust particles. (100% of 3 microns). The air is then run through an absorbent

filter to remove gasses.

A vacuum draws the particulate-laden atmosphere out of the container, while closures are individually cleaned and placed immediately on the treated containers.

Odor free and tasteless.

The air inside the bottle will not interfere with your product in any way.

Photo: Actual accumulation of foreign material removed from containers over an 80 hour period.

Additional Sizes in
Either Glass or
Plastic are Available
for all Cleaning
Processes.

Please Inquire.
1-800-922-7558



AMBER WIDE MOUTH PACKERS

Capacity (ml)	1oz (30)	2oz (60)	4oz (120)	8oz (240)	16oz (500)	32oz (950)
Neck Size (MM)	28-400	33-400	38-400	45-400	53-400	53-400

CAPS ATTACHED - CONVENIENCE PACK

TEFLON® PTFE LINED caps ●	2P01QGTV	2P02QGTV	2P04QGTV	2P08QGTV	2P16QGTV	2P32QGTV
PULP/VINYL LINED caps	2P01QPVV	2P02QPVV	2P04QPVV	2P08QPVV	2P16QPVV	2P32QPVV
Wt. Case (lbs.)	4	5	7	10	11	12
No. Per Case	24	24	24	24	12	12



AMBER BOSTON ROUNDS

Capacity (ml)	2oz (60)	4oz (120)	8oz (240)	16oz (480)	32oz (950)
Neck Size (MM)	20-400	22-400	24-400	28-400	33-400

CAPS ATTACHED - CONVENIENCE PACK

TEFLON® PTFE LINED caps ●	2A02QGTV	2A04QGTV	2A08QGTV	2A16QGTV	2A32QGTV
PULP/VINYL LINED caps	2A02QPVV	2A04QPVV	2A08QPVV	2A16QPVV	2A32QPVV
Wt. Case (lbs.)	4	6	10	8	14
No. Per Case	24	24	24	12	12



CLEAR BOSTON ROUNDS

Capacity (ml)	1oz (30)	2oz (60)	4oz (120)	8oz (240)	16oz (480)	32oz (960)
Neck Size (MM)	20-400	20-400	22-400	24-400	28-400	33-400

CAPS ATTACHED - CONVENIENCE PACK

TEFLON® PTFE LINED caps ●	2B01QGTV	2B02QGTV	2B04QGTV	2B08QGTV	2B16QGTV	2B32QGTV
PULP/VINYL LINED caps	2B01QPVV	2B02QPVV	2B04QPVV	2B08QPVV	2B16QPVV	2B32QPVV
Wt. Case (lbs.)	5	4	6	10	8	14
No. Per Case	48	24	24	24	12	12



CLEAR STRAIGHT SIDED ROUNDS

Capacity (ml)	1.5 oz (45)	2oz (60)	4oz (120)	9oz (240)	16oz (480)	32oz (950)
Neck Size (MM)	43-400	53-400	58-400	70-400	89-400	89-400

CAPS ATTACHED - CONVENIENCE PACK

TEFLON® PTFE LINED caps ●	2S01QGTV	2S02QGTV	2S04QGTV	2S08QGTV	2S16QGTV	2S32QGTV
PULP/VINYL LINED caps	2S01QPVV	2S02QPVV	2S04QPVV	2S08QPVV	2S16QPVV	2S32QPVV
Wt. Case (lbs.)	7	6	8	13	10	15
No. Per Case	48	24	24	24	12	12



CLEAR WIDE MOUTH FRENCH SQUARES

Capacity (ml)	1/2oz (15)	1oz (30)	2oz (60)	4oz (120)	8oz (240)	16oz (500)	32oz (950)
Neck Size (MM)	20-400	24-400	28-400	33-400	43-400	48-400	58-400

CAPS ATTACHED - CONVENIENCE PACK

TEFLON® PTFE LINED caps ●	2F00QGTV	2F01QGTV	2F02QGTV	2F04QGTV	2F08QGTV	2F16QGTV	2F32QGTV
PULP/VINYL LINED caps	2F00QPVV	2F01QPVV	2F02QPVV	2F04QPVV	2F08QPVV	2F16QPVV	2F32QPVV
Wt. Case (lbs.)	4	6	10	8	13	18	14
No. Per Case	48	48	48	24	24	24	12



**Tin Cans...
When Only
an Airtight
Container
Will Do.**



*128 OZ. COMES COMPLETE WITH EARS AND BAIL HANDLES

Tin cans provide a number of excellent performance characteristics, including:

- The ability to withstand physical shock. No need to worry about breakage during shipping and storage since cans are rugged and durable.

- Protection from permeation. Air, oxygen and other gases can't get in or out of the container.

- Shapes and sizes to fit your application. Choose from various sizes of plain tin paint cans and screw cap oblong cans.

.25 PLAIN TIN PAINT CANS

COMPLETE WITH LIDS AND SECURING CLIPS

Capacity (ml) 4oz (120) 8oz (250) 16oz (500) 32oz (960) 64oz (1920) 128oz (3785)

Cat. No.	1M00Q	1M03Q	1M06Q	1M07Q	1M08Q	1M10Q*
Wt. Case (lbs.)	6	10	14	10	9	6
No. Per Case	48	48	24	24	12	6



.25 PLAIN TIN SCREW CAP OBLONG CANS

COMPLETE WITH LINED CAPS AND INNER SEALS

Capacity (ml) 16oz (500) 32oz (960) 128oz (3785)
Cap Size 1" 1" 1-3/4"

Cat. No.	1M11Q	1M14Q	1M16Q
Wt. Case (lbs.)	7	10	6
No. Per Case	24	24	6



PAPER CYLINDERS

Capacity (ml) 8oz (240) 16oz (500) 32oz (960)

Cat. No.	2901	2902	2903
Wt. Case (lbs.)	14	19	28
No. Per Case	250	250	250



DISPOSABLE ALUMINUM WEIGH DISHES

Size (MM) 43 57 70

Cat. No.	4101	4100	4102
No. Per Case	1000	1000	1000



2 MIL POLYETHYLENE CLEAR OPEN END BAGS

Size (Inches) 14 x 15 15 x 20 15 x 24

Cat. No.	BAG-00039	BAG-00040	BAG-00041
Wt. Case (lbs.)	27	21	27
No. Per Case	1000	500	500

WARNING

ALL PACKAGES ARE SUPPLIED COMPLETE AND THE SHIPPER MUST NOT MAKE ANY SUBSTITUTIONS. FAILURE TO HEED THIS WARNING WOULD INVALIDATE PACKAGE COMPLIANCE. THE SELLER MAKES NO CLAIM AS TO THE SUITABILITY OF THESE PACKAGES FOR THE CUSTOMER'S PRODUCT. DETERMINATION OF THE SUITABILITY IS THE RESPONSIBILITY OF THE CUSTOMER.

More and more labs must comply with stringent regulations governing the shipment of hazardous materials. The benefits for choosing Qorpak to provide this type of packaging include:

- United Nations compliant packaging that is acceptable for shipping anywhere in the world.
- Packaging that is certified. No need for you to pay for costly certification since Qorpak's packaging is in compliance with domestic and international regulations.
- Packaging available for liquids or solids.
- Easy to use packaging. Already assembled, the container, closure, corrugated outer packaging, tapes, cushioning material and absorbent are all included.
- All containers have passed pressure testing.
- No additional freight costs. These versatile packages can be shipped via air freight companies like Federal Express.
- Convenient quantities. Case lot quantities are available from your distributor.

Shipping Info Key:

- LIQUIDS
- DRY SUBSTANCES
- AIR SHIPPING
- GROUND SHIPPING

1-800-922-7558

CANS, PAPER, WEIGH DISHES & BAGS

Transport Infectious Substances & Hazardous Materials with Durable Shippers



EXEMPTION PACKAGING DOT-E 9168

(FOR CLASS B POISONS - NO HAZARD LABEL REQUIRED FOR DOMESTIC SHIPMENTS)
INCLUDES SAFETY COATED GLASS CONTAINER AND THERMOSET CAP WITH FOAMED BACK TEFLON® PTFE LINER

Capacity (ml)	4oz (120)	4oz (120)	8oz (240)	8oz (250)	16oz (480)	16oz (500)	32oz (960)	32oz (950)
Neck Size (MM)	22	38	24	45	28	53	33	53

CONVENIENCE PACK

Cat. No.	6654M	6653M	6656M	6651M	6645M	6644M	6642M	6643M
Wt. Case (lbs.)	7	7	7	8	22	23	24	23
No. Per Case	4	4	4	4	4	4	4	4

NOTE: ACCEPTABLE FOR GREAT DANGER HAZARDOUS MATERIALS WHEN GROSS MASS OF 4KG. IS NOT EXCEEDED. PACKING GROUP I.



4-PAK GLASS MULTIPLE PACKAGE

INCLUDES SAFETY COATED GLASS CONTAINER AND THERMOSET CAP WITH FOAM BACKED TEFLON® PTFE LINER

Capacity (ml)	4oz (120)	4oz (120)	16oz (480)	16oz (500)	32oz (960)	32oz (950)
Neck Size (MM)	22	38	28	53	33	53

CONVENIENCE PACK

Cat. No.	6454M	6453M	6445M	6444M	6442M	6443M
Wt. Case (lbs.)	7	8	21	21	21	21
No. Per Case	4	4	4	4	4	4
Max Gross Mass (kg)	1.4	1.6	4.9	5.2	8.9	8.6

NOTE: ACCEPTABLE FOR MEDIUM DANGER HAZARDOUS MATERIALS WHEN GROSS MASS IN THIS CHART IS NOT EXCEEDED. PACKING GROUP II.



4 LITER AMBER JUG

CONVENIENCE PACK

UN MARKING ®	4G/Y 6.6/S
K/A CAPABILITY	250
Cat. No.	6551M
Wt. Case (lbs.)	23
No. Per Case	4

NOTE: ACCEPTABLE FOR MEDIUM DANGER HAZARDOUS MATERIALS WHEN GROSS MASS 6.6 KG. IS NOT EXCEEDED. PACKING GROUP II.



4-PAK PLASTIC MULTIPLE PACKAGE

INCLUDES HIGH DENSITY POLYETHYLENE LABORATORY STYLE BOTTLE AND POLYPROPYLENE CAP

Capacity (ml)	8oz (250)	16oz (480)	16oz (500)	32oz (960)	32oz (1000)
Neck Size (MM)	43	28	53	38	63

CONVENIENCE PACK

Cat. No.	6763M	6764M	6765M	6766M	6767M
Wt. Case (lbs.)	5	7	7	11	12
No. Per Case	4	4	4	4	4
Max Gross Mass (kg)	2.5	4	4.5	8.6	8.6

NOTE: ACCEPTABLE FOR MEDIUM DANGER HAZARDOUS MATERIALS WHEN GROSS MASS IN THIS CHART IS NOT EXCEEDED. PACKING GROUP II.



6.2 DRY ICE INFECTIOUS SUBSTANCE SHIPPERS

CONVENIENCE PACK

UN MARKING ®	4GU 6.2
Cat. No.	HMS-69110
Wt. Case (lbs.)	10
No. Per Case	4

ADD YOUR OWN PRIMARY CONTAINER.



INFECTIOUS SUBSTANCE SHIPPERS

IDEAL FOR LABS AND HOSPITALS TO TRANSPORT BLOOD COLLECTION TUBES, PATHOLOGY SAMPLES AND BIOPSIES.

CONVENIENCE PACK

UN MARKING ®	4GU 6.2
Cat. No.	HMS-69070
Wt. Case (lbs.)	17
No. Per Case	12

ADD YOUR OWN PRIMARY CONTAINER.

NOTE: ACCEPTABLE FOR CLASS 6.2 INFECTIOUS SUBSTANCE.



DRY ICE INFECTIOUS SUBSTANCE SHIPPERS

CONVENIENCE PACK

UN MARKING ®	4G Class 6.2
	4G/X14.5/S
Cat. No.	6905M
Wt. Case (lbs.)	16
No. Per Case	2

VIALS NOT INCLUDED.
NOTE: ACCEPTABLE FOR CLASS 6.2 INFECTIOUS SUBSTANCE.




INFECTIOUS SUBSTANCE SHIPPERS

CONVENIENCE PACK

UN MARKING ®	4G Class 6.2
Cat. No.	6901
Wt. Case (lbs.)	36
No. Per Case	12

VIALS NOT INCLUDED.
NOTE: ACCEPTABLE FOR CLASS 6.2 INFECTIOUS SUBSTANCE.



ONE GALLON STEEL CLOSED HEAD CONTAINER

CONVENIENCE PACK

UN MARKING ®	4G/X 9.5/S
K/A CAPABILITY	250
Cat. No.	6740M*
Cat. No.	6741M**
Wt. Case (lbs.)	21
No. Per Case	4

* WITH EPOXY-PHENOLIC LINED INTERIOR.

** WITH UNLINED INTERIOR.

NOTE: ACCEPTABLE FOR GREAT DANGER HAZARDOUS MATERIALS WHEN GROSS MASS 9.5 KG. IS NOT EXCEEDED. PACKING GROUP I.



ONE GALLON PAINT CAN 4G COMBINATION PACKAGE

CONV. PACK

	1 CAN/PK	2 CANS/PK	1 TALL PAINT
UN MARKING ®	4G/Y 7.7/S	4G/Y 16.6/S	4G/Y 4/S
K/A CAPABILITY	95	95	95
Cat. No.	6768M	6770M	6750M*
Wt. Case (lbs.)	11	17	10
No. Per Case	4	4	4
Max Gross Mass (kg)	7.7	16.6	4

NOTE: ACCEPTABLE FOR MEDIUM DANGER HAZARDOUS MATERIALS WHEN GROSS MASS IN THIS CHART IS NOT EXCEEDED. PACKING GROUP II.

*THIS SPECIALLY DESIGNED PACKAGE IS A COMBINATION PACKAGE FOR PAINT SHIPMENTS AS WELL AS AN OVERPACK FOR CONTAINERS THAT HAVE UNKNOWN INTERNAL PRESSURE TESTS.

new!

- Wide GL45 neck opening is designed for easy filling and pouring.
- Bottles, caps and pouring rings are autoclavable up to 140°C.
- Ideal for storage as well as mixing and sampling.
- Permanent white enamel graduations in ML increments and marking spot.
- Cap and pouring rings available separately. Please call for more information.



LAB STYLE SODA-LIME GLASS MEDIA BOTTLES

Capacity	100 ml	250 ml	500 ml	1000 ml
MEDIA BOTTLE WITH STANDARD POURING RING AND BLUE CAP				
Cat. No.	GLC-08561	GLC-08562	GLC-08563	GLC-08564
Wt. Case (lbs.)	5	7	9	15
No. Per Case	10	10	10	10
MEDIA BOTTLE ONLY				
Cat. No.	GLA-08565	GLA-08566	GLA-08567	GLA-08568
Wt. Case (lbs.)	4	6	8	14
No. Per Case	10	10	10	10

Visit Our
Web Site at:
www.qorpak.com

Or Call Today!
1-800-922-7558



CORPORATE ONE WEST
1195 WASHINGTON PIKE
BRIDGEVILLE, PA 15017-2808

PRESORTED
STANDARD
US POSTAGE
PAID
PITTSBURGH, PA
PERMIT NO. 1715