

# WESCO

## SERIES

*enjoy the VU*

9000

### Uni VU™ 9000 Series

*Uni VU™ 9000 series universal microscope features versatility for reflected light and transmitted light microscopy with quality optics and mechanics for a variety of applications.*

*There are currently six models in the Uni VU™ 9000 Series with many additional accessories available to further customize your viewing needs.*



# Uni VU 9000 Series

SPECIFICATIONS	MODELS					
	9200T	9200R	9200TR	9300T	9300R	9300TR
Transmitted light stand	•			•		
Reflected brightfield, darkfield and polarizing light stand		•			•	
Universal reflected and transmitted light stand			•			•
30 degree binocular observation body	•	•	•			
30 degree trinocular observation body				•	•	•
10X widefield eyepieces 20mm field of view	•	•	•	•	•	•
Reverse angle revolving quintuple nosepiece	•		•	•		•
Front angle revolving quintuple nosepiece		•	•		•	•
DIN 4X, 10X, 25X, 40XR & 100XR oil plan achromatic objectives	•		•	•		•
DIN BF4X, DF10X, DR25X & DR40X & BF 63X plan achromatic objectives		•	•		•	•
Coaxial coarse and fine focus with tension control	•	•	•	•	•	•
150mm X 150mm flat top graduated mechanical stage with slide holder and low position 76mm X 50mm travel X-Y controls	•		•	•		•
145mm X 185mm metallurgical stage with low position 45mm X 55mm travel X-Y controls		•			•	
N.A. 1.25 Abbe condenser	•		•	•		•
12 volt 50 watt halogen illumination with variable light control	•	•	•	•	•	•
Field iris diaphragm	•		•	•		•

# WESCO



Distributed By: